

510(k) Summary

K112173

Submission Date: 29 September 2011
OCT - 5 2011

Submitter: Spacelabs Healthcare
5150 220th Avenue SE
Issaquah, WA 98029

Submitter Contact: Spacelabs Healthcare
5150 220th Avenue SE
Issaquah, WA 98029
Mr. David J. Geraghty
Spacelabs Healthcare
Phone: +1 (425) 657-7200, ext 5889
Fax: +1 (425) 657-7210
Email: david.geraghty@spacelabs.com

Official Contact: Thomas Kroenke
Principal Consultant
Speed To Market, Inc.
PO Box 3018
Nederland, CO 80466 USA
tkroenke@speedtomarket.net
303 956 4232

Manufacturing Site: Spacelabs Healthcare
5150 220th Avenue SE
Issaquah, WA 98029

Trade Name: Spacelabs Multigas Module, Model 92518

Common Name: Multigas analyzer

Classification Name: Carbon dioxide gas analyzer, Enflurane gas analyzer, Desflurane gas analyzer, Sevoflurane gas analyzer, Isoflurane gas analyzer, Halothane gas analyzer, Nitrous oxide gas analyzer, Oxygen gas analyzer

Classification Regulation: 21 CFR §868.1700

Primary Product Code: CBR

Secondary Product Codes: CCK, CBQ, NHO, NHP, NHQ, CBS, CCL

510(k) Summary

**Substantially
Equivalent Devices:**

<i>New Spacelabs Model</i>	<i>Predicate 510(k) Number</i>	<i>Predicate Manufacturer / Model</i>
Spacelabs Multigas Module, Model 92518	K053599	Spacelabs Medical, Inc. / Multigas Analyzer Module 91518 and Accessories

Device Description:

The Spacelabs Multigas Module, Model 92518 (Model 92518) is an easy-to-use modular unit in the Spacelabs Healthcare Ultraview-family of monitors (Ultraview). The Module 92518 is inserted into a bay within an Ultraview monitor. The Ultraview monitor is then used to control the Model 92518, and provide the user interface for the Model 92518.

The Module 92518 is a sidestream analyzer intended to provide a measurement of the following parameters: Carbon Dioxide (CO₂); Oxygen (O₂) and Nitrous Oxide (N₂O); Anesthetic agents including Desflurane; Enflurane; Halothane; Isoflurane; and Sevoflurane; Respiratory rate; and Calculated MAC and age-dependent MAC values.

The Spacelabs Ultraview monitors display information from the Module 92518 and are the user interface for the Model 92518.

Ultraview monitors provide a number display for the anesthetic agent concentrations and respiratory rate, and a capnograph waveform display for O₂ and CO₂. The Module 92518 automatically identifies which anesthetic agent or mixture of anesthetic agents is present, and measures the concentration of the identified agent(s). An alarm is issued if a mixture of more than two anesthetic agents is detected.

The Module 92518 is intended to be used primarily in the operating room environment

510(k) Summary

Intended Use:

The Spacelabs Multigas Module, Model 92518 (Model 92518) is intended to provide a means of monitoring a variety of gas concentrations and alert clinical personnel when the concentration of anesthetic agent(s), oxygen, carbon dioxide or nitrous oxide moves outside of user-defined limits. The multigas module is capable of automatically identifying which anesthetic agent(s) is being administered.

The 92518 Multigas Module is intended to be used with and controlled by a Spacelabs Healthcare Ultraview-family monitor. The 92518 Multigas Module is intended to be used for monitoring all hospitalized patients, under the direction of qualified medical personnel.

Although the 92518 Multigas Module alarms when the duration between breaths exceeds user defined limits, it is not intended to be a primary diagnostic apnea monitor and/or recording device.

Technology Comparison:

The Model 92518 employs the same technological characteristics as the predicate device.

<i>Characteristic</i>	<i>Predicate Device</i>	<i>Proposed Device</i>
<i>Sampling Technique</i>	Sidestream	Same
<i>Carbon-Dioxide (CO₂) Measurement</i>	Infrared	Same
<i>Oxygen (O₂) Measurement</i>	Paramagnetic	Same
<i>Anesthetic Agent Measurement</i>	Infrared	Same

Summary of Performance Testing:

Electrical Safety

The Model 92518 was tested for performance in accordance with the following Standards:

- *IEC 60601-1: 2005, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance; and*
- *UL 60601-1: 2003, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.*

Test results indicated that the Model 92518 complies with the Standards.

510(k) Summary

Electromagnetic Compatibility (EMC) Testing

The Model 92518 was tested for performance in accordance with the following Standard:

- *IEC 60601-1-2: 2007, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests.*

Test results indicated that the Model 92518 complies with the Standards.

Software Testing

Software device modifications made to the Model 92518 were designed and developed according to a robust software development process, and were rigorously verified and validated.

Software information is provided in accordance with internal documentation and the following Standards and guidance documents:

- *FDA guidance: The content of premarket submissions for software contained in medical devices, 11 May 05;*
- *FDA guidance: Off-the-shelf software use in medical devices, 09 Sep 99;*
- *FDA guidance: General principles of software validation; Final guidance for industry and FDA staff, 11 Jan 02;*
- *IEC 60601-1-4: 2000, Medical electrical equipment Medical electrical equipment – Part 1-4: General requirements for safety – Collateral Standard: Programmable electrical medical systems.*

Test results indicate that the Model 92518 complies with its predetermined specification and the Standards and guidance documents.

Performance Testing

The Model 92518 was tested for performance in accordance with internal documentation and the following Standards:

- *IEC 60601-1-8: 2006, Medical electrical equipment – Part 1-8: General requirements for basic safety and essential performance – Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems; and*
- *ISO 21647: 2004, Medical electrical equipment – Particular requirements for the basic safety and essential performance of respiratory gas monitors.*

Test results indicated that the Model 92518 complies with its predetermined specification and with the applicable Standards.

510(k) Summary

Conclusion

Verification and validation activities were conducted to establish the performance and safety characteristics of the software device modifications made to the Model 92518. The results of these activities demonstrate that the Model 92518 is safe and effective when used in accordance with its intended use and labeling.

Therefore, the Model 92518 is considered substantially equivalent to the predicate device.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Spacelabs Healthcare
C/O Mr. Thomas Kroenke
Principal Consultant
Speed To Market, Incorporated
P.O. Box 3018
Nederland, Colorado 80466

OCT - 5 2011

Re: K112173
Trade/Device Name: Spacelabs Multigas Module, Model 92518
Regulation Number: 21 CFR 868.1700
Regulation Name: Nitrous Oxide Gas Analyzer
Regulatory Class: II
Product Code: CBR, CCK, CBQ, NHO, NHP, NHQ, CBS, CCL
Dated: September 6, 2011
Received: September 8, 2011

Dear Mr. Kroenke:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K112173

Device Name: Spacelabs Multigas Module, Model 92518

Indications for Use: The Spacelabs Multigas Module, Model 92518 (Model 92518) is intended to provide a means of monitoring a variety of gas concentrations and alert clinical personnel when the concentration of anesthetic agent(s), oxygen, carbon dioxide or nitrous oxide moves outside of user-defined limits. The multigas module is capable of automatically identifying which anesthetic agent(s) is being administered.

The 92518 Multigas Module is intended to be used with and controlled by a Spacelabs Healthcare Ultraview-family monitor. The 92518 Multigas Module is intended to be used for monitoring all hospitalized patients, under the direction of qualified medical personnel.


Although the 92518 Multigas Module alarms when the duration between breaths exceeds user defined limits, it is not intended to be a primary diagnostic apnea monitor and/or recording device.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K112173