



**510(k) Notification**  
**Audit® MicroLQ™ Spinal Fluid**

JUN 15 2012

K112742

**510(k) Summary**

**A. Submitter**

Aalto Scientific, Ltd.  
1959 Kellogg Ave.  
Carlsbad, CA 92008  
Telephone: (760) 431-7922  
Fax: (760) 431-6824

**B. Contact Person**

Dessi Lyakov  
Regulatory Affairs Manager  
Telephone: (760) 431-7922 Ext. 118  
E-mail: dlyakov@aaltoscientific.com

**C. Date of Summary Preparation**

June 13, 2012

**D. Device Identification**

Product Trade Name: Audit® MicroLQ™ Spinal Fluid  
Common Name: Spinal Fluid  
Classification Name: Assay QC Material  
Device Classification: Class I  
Regulation Number: 21 CFR 862.1660  
Panel: 75  
Product Code: JJY

*Device to Which Substantial Equivalence is Claimed:*

Product Trade Name: Audit® MicroCV™ Protein Linearity Set  
Aalto Scientific, Ltd., Carlsbad, CA  
K101216

**E. Description of the Device**

The Audit® MicroLQ™ Spinal Fluid is a human based, liquid set of QC material. Each level of the set contains Chloride, Glucose, Immunoglobulin A (IgA), Immunoglobulin G (IgG), Immunoglobulin M (IgM), Lactate, Lactate Dehydrogenase (LD), Microalbumin, Microprotein, and Sodium analytes. It is used to confirm the proper calibration of Chloride, Glucose, Immunoglobulin A (IgA), Immunoglobulin G (IgG), Immunoglobulin M (IgM), Lactate, Lactate Dehydrogenase (LD), Microalbumin, Microprotein, and Sodium

**F. Statement of Intended Use**

Audit® MicroLQ™ Spinal Fluid Control is a quality control material intended for monitoring the



## 510(k) Notification Audit® MicroLQ™ Spinal Fluid

precision of laboratory testing procedures.

When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

The Audit® MicroLQ™ Spinal Fluid is for In Vitro Diagnostic use only.

### **G. Summary of Performance Data**

Stability studies have been performed to determine the shelf life for the Audit® MicroLQ™ Spinal Fluid Set. All supporting data is retained on file at Aalto Scientific, Ltd. Product claims are as follows:

*Shelf Life:* Two years, when stored unopened at 2 - 8° C.

*Open Vial:* 30 days, when stored at 2 - 8° C.

*Note:* Real time studies are ongoing to support the shelf life of this product.



## 510(k) Notification Audit<sup>®</sup> MicroLQ<sup>™</sup> Spinal Fluid

### H. Technical Characteristics Compared to Predicate Device

Characteristics	Audit <sup>®</sup> MicroLQ <sup>™</sup> Spinal Fluid Set (New Device)	Audit <sup>™</sup> MicroCV <sup>™</sup> Protein Linearity Set (K101216)
Intended Use	<p>Audit<sup>®</sup> MicroLQ<sup>™</sup> Spinal Fluid Control is a quality control material intended for monitoring the precision of laboratory testing procedures.</p> <p>When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.</p> <p>The Audit<sup>®</sup> MicroLQ<sup>™</sup> Spinal Fluid is for In Vitro Diagnostic use only.</p>	Same
Number of Analytes per vial	10	7
Contents	6 x 3 mL	5 x 5 mL
Matrix	Human Based Serum	Same
Type of Analytes	Chloride, Glucose, Immunoglobulin A (IgA), Immunoglobulin G (IgG), Immunoglobulin M (IgM), Lactate, Lactate Dehydrogenase (LD), Microalbumin, Microprotein, and Sodium	Alpha-1-Antitrypsin, Complement C3, Complement C4, Immunoglobulin G, Immunoglobulin A, Immunoglobulin M and Transferrin
Form	Liquid	Same
Storage	2 to 8° C Until expiration date	Same
Open Bottle Stability	30 days at 2 to 8° C	24 hours at 2 to 8°

### I. Conclusions

Based upon the purpose of the device and the descriptions of the predicate device, the safety and efficacy, and the stability data generated, the product is substantially equivalent to the predicate device.



10903 New Hampshire Avenue  
Silver Spring, MD 20993

Aalto Scientific, Ltd.  
c/o Dessi Lyakov  
1959 Kellogg Avenue  
Carlsbad, CA 92008

JUN 15 2012

Re: k112742  
Trade Name: Audit® MicroLQ™ Spinal Fluid Control  
Regulation Number: 21 CFR §862.1660  
Regulation Name: Quality Control material (assayed and unassayed)  
Regulatory Class: Class I, reserved  
Product Codes: JJY  
Dated: April 23, 2012  
Received: April 25, 2012

Dear Dessi Lyakov:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (301) 796-5760. For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance...

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-5680 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>

Sincerely yours,



Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

**510(k) Number:** K112742

**Device Name:** AUDIT® MicroLQ™ Spinal Fluid Control Set

**Indications For Use:**

Audit® MicroLQ™ Spinal Fluid Control is a quality control material intended for monitoring the precision of laboratory testing procedures.

When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

The Audit® MicroLQ™ Spinal Fluid is for In Vitro Diagnostic use only.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



Division Sign-Off Office of In Vitro  
Diagnostic Device Evaluation and Safety

510(k)           K112742