

NOV - 1 2011

Section 5 - 510(k) Summary

Date of Summary Preparation: 08/25/2011

1. Submitter's Identifications

Submitter's Name: ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD.

Address: Jin'an Road, Minzhong, Zhongshan 528441, Guangdong, China

Contact Person: Lisa Li

Contact Email Address: lishal@transtek.cn

Telephone: +86(760)88282982 ext. 876

Fax: +86(760)85339231

2. Correspondent's Identifications

Correspondent's Name: A03 Lab of BTS

Address: No.1 Fanghua Street, Hi-tech Zone, Chengdu 610041, Sichuan, China

Contact Person: Leo Wang

Contact Email Address: leo.w@hibts.com

Telephone: 086-28-86083300

Fax: 086-4006070050 ext. 8

3. Name of the Device

Device Classification Name: Analyzer, Body Composition (Impedance Plethysmograph)

Product Name: TRANSTEK Glass Body Fat Analyzer

Trade Name: TRANSTEK

Models: GBF-733-W1, GBF-1012, BF-1039, BF-1041-A, GBF-950-D

Classification Panel: Cardiovascular

Common/Usual Name: Body Composition Analyzer/Scales

Product Code: MNW

Device Classification: Class II

Contraindications: Do not use the Analyzer if you have a pacemaker or other internal medical device.

4. The Predicate Devices

TRANSTEK, Glass Body Analyzer, Model GBF-950, K102191

5. Device Description

The TRANSTEK Glass Body Fat Analyzer uses BIA (Bio Impedance Analysis) technology which passes an electrical current through the body to estimate body fat mass, lean mass, total body water

and bone mass. The electrical current is small and may not be felt. Contact with the body is made via glass and stainless steel pads on the platform of the analyzer.

This method simultaneously calculates your personal weight, body fat, total body water, bone mass and muscle mass, giving you a more accurate reading of your overall health and fitness.

This scale stores the personal data of up to 4 or 8 users. As well as being an analyzer, this device can be used as a conventional scale.

6. Intended Use of Device

The TRANSTEK Glass Body Fat Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally healthy adults 18 years of age or older. It is intended for use in the home/domestic setting only.

7. Summary of Substantial Equivalence

Table 1: The difference between TRANSTEK Glass Body Fat Analyzer and the predicate device, Transtek Glass Body Analyzer (Model GBF-950)

Feature	Proposed Device: TRANSTEK Glass Body Fat Analyzer Models: BF-1039, BF-1041-A GBF-1012, GBF-733-W1, GBF-950-D	Predicate Device: Transtek Glass Body Analyzer Model: GBF-950
510(k) Number	—	K102191
Manufacturer	ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD	ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD
Classification	21 CFR 870.2770	21 CFR 870.2770
Product Code	MNW	MNW
Indication for use	The TRANSTEK Glass Body Fat Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally healthy adults 18 years of age or older. It is intended for use in the home/domestic setting only.	The TRANSTEK Glass Body Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally healthy adults 18 years of age or older. It is intended for use in the home/domestic setting only.

Feature	Proposed Device: TRANSTEK Glass Body Fat Analyzer Models: BF-1039, BF-1041-A GBF-1012, GBF-733-W1, GBF-950-D	Predicate Device: Transtek Glass Body Analyzer Model: GBF-950
Device description	TRANSTEK Glass Body Fat Analyzer utilizes a "foot-to-foot" bioelectrical impedance analysis (BIA) technology to determine internal body composition.	TRANSTEK Glass Body Analyzer utilizes a "foot-to-foot" bioelectrical impedance analysis (BIA) technology to determine internal body composition.
Analysis method	BIA (Bioelectrical Impedance Analysis)	BIA (Bioelectrical Impedance Analysis)
Operating parameters	50 KHz	50 KHz
Power source	BF-1039, 2*CR2032 lithium batteries BF-1041-A, 2*CR2032 lithium batteries GBF-1012, 4*AAA batteries GBF-733-W1, 2*CR2032 lithium batteries GBF-950-D, 4*AAA batteries	6V, 4*AAA batteries
Operating keys	4	4
Number of electrodes	4	4

8. Conclusions

The subject devices have all features of the predicate device, Transtek Glass Body Analyzer (Model GBF-950) except the new features and the power source voltage of the device. These differences do not affect the safety and effectiveness of the subject devices. BIA (Bioelectrical Impedance Analysis) technology is same as what is used in Glass Body Analyzer (Model GBF-950). Thus, the subject devices are substantially equivalent to the predicate devices.

--- End of this section ---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
10903 New Hampshire Avenue
Document Mail Center - WO66-G609
Silver Spring, MD 20993-0002

ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD.
% Mr. Leo Wang, Consulting Manager
A03 Lab of BTS
No. 1 Fanghua Street, Hi-tech District
CHENGDU SICHUAN 610041
CHINA

NOV - 1 2011

Re: K112932

Trade/Device Name: Transtek Glass Body Fat Analyzer
Models: GBF-733-W1, GBF-1012, BF-1039, BF-1041-A, GBF-950-D
Regulation Number: 21 CFR§ 870.2770
Regulation Name: Impedance plethysmograph
Regulatory Class: II
Product Code: MNW
Dated: October 3, 2011
Received: October 3, 2011

Dear Mr. Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related

adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Herbert P. Lerner, M.D., Director (Acting)
Division of Reproductive, Gastro-Renal
and Urological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Section 4 - Indications for Use

510(k) Number (if known): K112932

Device Name:

Transtek Glass Body Fat Analyzer

Models: GBF-733-W1, GBF-1012, BF-1039, BF-1041-A, GBF-950-D

Indications for Use:

The Transtek Glass Body Fat Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally healthy adults 18 years of age or older. It is intended for use in the home/domestic setting only.

Prescription Use _____

AND/OR


Over-The-Counter Use X

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Reproductive, Gastro-Renal, and
Urological Devices
510(k) Number K112932