

K 113083

510(k) Summary

JAN 19 2012

Date of Summary Preparation: August.16.2011

1. Submitter's Identifications

Submitter's Name: Truly Instrument Limited
Address: Site 2, Truly Industrial Area, Shanwei City,
Guangdong Province, China
Contact Person: Manager Yang Jian-Hao
Telephone: 86-0660-3380070
Fax: 86-0660-3380377

2. Name of the Device

Device Classification Name: System, Measurement, Blood-Pressure, Non-invasive
Trade Name: Truly Automatic Arm Blood Pressure Monitor
Models: DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93
Classification Panel: cardio-vascular
Common/Usual Name: Automatic Arm Blood Pressure Monitor
Product Code: DXN
Device Classification: Class II
Contraindications : N/A

3. The Predicate Devices

Truly Automatic Arm Blood Pressure Monitor, Model DB61M, K091434

4. Device Description

Truly Automatic Arm Blood Pressure Monitor DB series, Models DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93 are designed to measure the systolic and diastolic blood pressure and pulse rate of an individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm. Our method to define systolic and diastolic pressure is similar to the auscultatory method but uses an electronic pressure sensor rather than a stethoscope and mercury manometer. The sensor converts tiny alterations in cuff pressure to electrical signals, by analyzing those signals to define the systolic and diastolic blood pressure and calculating pulse rate, which is a well-known technique in the market called the "oscillometric method".

The main components of the Truly Automatic Arm Blood Pressure Monitor DB series are the main unit and cuff unit. ABS is used to outer housing of the main unit. The preformed cuff unit, which is applicable to arm circumference approximately between 220 and 340 mm, includes the inflatable bladder and nylon shell. All models of the arm blood pressure monitor use a single size of cuff. The device consists of the microprocessor, the pressure sensor, the operation keys, the pump, the electromagnetic deflation control valve and the LCD. The subject devices are powered by four AA alkaline batteries.

The device also compares the longest and the shortest time intervals of detected pulse waves to mean time interval and displays a warning signal with the reading to indicate the detection

K 113083

TRULY® Truly Instrument Limited

of irregular pulse rhythm when the difference of the time intervals is over 25%.

5. Intended use of device

Truly Automatic Arm Blood Pressure Monitor DB series, Models DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93 are a series devices intended to measure the systolic and diastolic blood pressure and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm.

The devices features include irregular pulse rhythm detection during measurement, and display a warning signal with the reading once the irregular heartbeat is detected.

6. Summary of Substantial Equivalence

Table-1: The comparison table

Parameter	Predicate Devices DB61M	DB 33	DB 81	DB 82	DB 83	DB 85	DB 91	DB 92	DB 93	Result
Intended use	Measuring systolic and diastolic blood pressure and pulse rate of adult individual	Measuring systolic and diastolic blood pressure and pulse rate of adult individual								Same
Indications for use	Measuring systolic and diastolic blood pressure and pulse rate of adult individual, Including irregular pulse rhythm detection . Over-The-Counter Use	Measuring systolic and diastolic blood pressure and pulse rate of adult individual, Including irregular pulse rhythm detection . Over-The-Counter Use								Same
Target Population	Adult	Adult								Same
Anatomical sites	Upper Arm	Upper Arm								Same
Where used (hospital, home, ambulance, etc)	Home	Home								Same
Energy used and / or delivered	4x 1.5V AA Battery	4x 1.5V AA Battery								Same
Human factors	Blood pressure	Blood pressure								Same
design	Refer to Table-2									Same
performance	Measuring systolic and diastolic blood pressure and pulse rate of adult individual, Including irregular pulse rhythm detection	Measuring systolic and diastolic blood pressure and pulse rate of adult individual, Including irregular pulse rhythm detection								Same
materials	Refer to Table-2									Same
biocompatibi	Cuff	Cuff								Same

K 113083

TRULY® Truly Instrument Limited

Parameter	Predicate Devices DB61M	DB 33	DB 81	DB 82	DB 83	DB 85	DB 91	DB 92	DB 93	Result
lity	According to ISO-10993	According to ISO-10993								
Compatibility with the environment and other devices	Operation Environment: 10°C~ 40°C, 15%~90%RH Storage Environment: -20°C~ 60°C, 10%~95%RH	Operation Environment: 10°C~ 40°C, 15%~90%RH Storage Environment: -20°C~ 60°C, 10%~95%RH								Same
sterility	N/A	N/A								Same
Electrical safety	According to IEC60601-1-2 According to IEC60601-1	According to IEC60601-1-2 According to IEC60601-1								Same
Mechanical safety	Same	Same								Same
Chemical safety	N/A	N/A								Same
Thermal safety	N/A	N/A								Same

Table-2: The comparison table in Design and Materials

Parameter	Predicate Devices DB61M	DB 33	DB 81	DB 82	DB 83	DB 85	DB 91	DB 92	DB 93	
Measurement algorithm Method	Oscillometric method	No change ,all same								
Measurement site of body	Arm	No change ,all same								
Pressure Sensor	MSP-2107	No change ,all same								
Cuff		No change ,all same								
Software		No change ,all same								
Irregular heartbeat detection		More than ±25% to the mean interval of pulse intervals. About the more detailed description of the IH detection algorithm, please refer to "Software validation report I-5. Algorithm description 4. Determination method of irregular heartbeat"..								
Memory Size	4 x 99	2X60	1X40	1X40	4X99	4X99	4X99	4X99	4X99	
Measurement Pressure Range	20 ~ 280 mmHg	No change ,all same								
Measurement Pulse Range	40 ~ 195 beats/min	No change ,all same								
Mesasuring resolution	1 mmHg	No change ,all same								
Accuracy	±3mmHg	No change ,all same								

TRULY®

Truly Instrument Limited

Parameter	Predicate Devices DB61M	DB 33	DB 81	DB 82	DB 83	DB 85	DB 91	DB 92	DB 93
Pressure									
Accuracy Pulse	±5%	No change ,all same							
Pressurization Source	Automatic internal pump	No change ,all same							
Cuff Deflation	Automatic deflation	No change ,all same							
Operating Environment	10~40℃ 15~90%RH	No change ,all same							
Power Voltage	4X 1.5V	No change ,all same							
Hardware circuit		No change ,all same							
Electronic element		No change ,all same							
PCB		No change ,all same							
Display Type	Liquid crystal display	Liquid crystal display ,Only difference size							
Cover		Difference							

7. Summary of Clinical study

1). Subjects:

Ninety subjects including 15 hypertensive patients in the hospital were participated in clinical study.

2). Method:

A standard mercury sphygmomanometer was used as a reference standard. Simultaneous and blinded blood pressure determinations were performed by two doctors.

3). Criteria:

The ANSI/AAMI SP10 Standard recommended :

- A. a mean difference of ± 5 mmHg, with standard deviation of differences of ± 8 mmHg between test device and reference method.
- B. Between-observer agreement should be:95% or more of readings made simultaneously by observers agree to within ± 10 mmHg and 85% or more agree to within ± 5 mmHg.

4). Result

Through clinical research, we can convinced that the clinical device is safe and effective. The results of the clinical data refer the follow two tables.

Table-1 Between test device and reference method

Criteria	Test Result	
	Systolic Pressure	Diastolic Pressure

TRULY[®] Truly Instrument Limited

Mean differences	x 5 mmHg	3.3 mmHg	4.5 mmHg
SD differences	x 8 mmHg	5.2 mmHg	4.6 mmHg

Table-2 Between-observer agreement

Criteria	Test Result	
	Systolic Pressure	Diastolic Pressure
At least 95% of readings agree to within x 10mmHg	100%	100%
At least 85% of readings agree to within x 5mmHg	95%	99%

8. Conclusions

The new subject series devices of Truly Automatic Arm Blood Pressure Monitor continue to follow principles of hardware and software design of the predicate device DB61M(K091434), and the feature, safety, effectiveness are also as same as DB61M., just only in case and LCD's size changes. Thus, the subject devices are substantially equivalent to the predicate devices.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Truly Instrument Co., Ltd.
c/o: Mr. Jack Lee
Truly (USA) INC
2620 Concord Avenue, Suite 106
Alhambra, CA 91803

JAN 19 2012

Re: K113083

Truly Automatic Arm Blood Pressure Monitor (DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93) - 8 models

Regulation Number: 21 CFR 870.1130

Regulation Name: Noninvasive blood pressure measurement system

Regulatory Class: II

Product Code: DXN

Dated: January 13, 2012

Received: January 13, 2012

Dear Mr. Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K 113083

Indications for Use

510(K) Number (if Known)

Device Name: Truly Automatic Arm Blood Pressure Monitor DB Series:

Models DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93

Indication For Use:

Truly Automatic Arm Blood Pressure Monitor, Models DB33, DB81, DB82, DB83, DB85, DB91, DB92, DB93 are a series device intended to measure the systolic and diastolic blood pressure and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm.

The devices features include irregular pulse rhythm detection during measurement, and display a warning signal with the reading once the irregular heartbeat is detected.

Prescription Use _____ AND/OR Over-The-Counter Use X

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



Page 1 of 1

(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K 113083