

JUN 12 2012

## 510(k) Summary

**Date Summary Prepared:** June 12, 2012

**510(k) Owner:** Ventus Medical, Inc.  
101 Metro Drive, Suite 620  
San Jose, CA 95110  
Phone: 408-200-5299  
Fax: 408-200-5298  
[info@ventusmedical.com](mailto:info@ventusmedical.com)

**Contact Person:** Cindy Domecus, R.A.C. (US & EU)  
Principal, Domecus Consulting Services LLC  
Consultant to Ventus Medical, Inc.  
(650) 343-4813

**Trade Name:** InVent Snoring Device

**Common Name:** Nasal Dilator

**Classification Name:** 21 CFR 874.3900 Nasal Dilator; Product Code LWF

**Predicate Device:** BR2 Nasal Dilator

### Device Description:

As with the predicate device, the subject device is a nasal dilator intended to reduce or eliminate snoring. While the predicate device is an internal nasal dilator, the subject device is an external nasal dilator in which the device is placed over the patient's nostrils. As with the predicate device, the subject device creates dilating pressure or expiratory positive airway pressure (EPAP) to dilate and open the nasal airway. The device directs expiratory flow through selected pathways, which maintain nasal airway pressure and help maintain nasal passage dilation.

**Intended Use:** The InVent Snoring Device is intended to reduce or eliminate snoring.

### Technological Characteristics Comparison:

The predicate device and subject device have the same mechanism of action – expiratory resistance. Minor changes have been made to the dimensions and materials.

**Performance Data:**

Bench testing of the subject device confirms that it meets pre-defined specifications and is substantially equivalent to the predicate device. Clinical testing also supports that the device is safe and effective for its intended use.

**Conclusion:**

The InVent Snoring Device is substantially equivalent to the BR2 predicate device.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Ventus Medical, Incorporated  
% Ms. Cindy Domecus, R.A.C. (US & EU)  
Principal  
Domecus Consulting Services LLC  
1171 Barroilhet Drive  
Hillsborough, CA 94010

JUN 12 2012

Re: K120665  
Trade/Device Name: InVent Snoring Device  
Regulation Number: 21 CFR 874.3900  
Regulation Name: Nasal dilator  
Regulatory Class: Class I  
Product Code: LWF  
Dated: May 2, 2012  
Received: May 3, 2012

Dear Ms. Domecus:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

We note that your device exceeded the Limitations of exemptions from section 510(k) of the Federal Food, Drug, and Cosmetic Act (21 CFR Part 874.9), and therefore required the submission and clearance of a premarket notification prior to commercial distribution in the United States. Future devices of this same type that meet the exemption criteria and do not exceed the limitations of exemptions found in 21 CFR Part 874.9 will be exempt from the premarket notification requirements of the Act.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Malvina B. Eydelman, M.D.  
Director  
Division of Ophthalmic, Neurological,  
and Ear, Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

# 1. Indications for Use Statement

## Indications for Use Form

### Indications for Use

510(k) Number (if known): \_\_\_\_\_

Device Name: InVent Snoring Device

Indications for Use:

The InVent Snoring Device is intended to reduce or eliminate snoring.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR


Over-The-Counter Use X  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

-----

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page \_\_\_ of \_\_\_

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Ophthalmic, Neurological and Ear,  
Nose and Throat Devices

510(k) Number K120665