

K120708

510(k) Summary
(as required by section 807.92(c).

FEB 7 2013

Sponsor	LifeScan Europe, a Division of Cilag GmbH International Landis and Gyr Strasse 1 Zug, Switzerland 6300
Correspondent	Nadine Nasr, Regulatory Project Manager LifeScan, Inc. 1000 Gibraltar Drive Milpitas, California 95035 Phone: 408-956-4079, 408-942-5906 Email: nnasr@its.jnj.com Alternate 510(k) Contact: Yvonne Middlefell, Director Regulatory Affairs LifeScan Scotland Ltd Beechwood Business Park North Inverness, Scotland IV2 3ED United Kingdom Phone: 44 01463 721250 Mobile: 44 (0) 7900 138 650 Fax: 44 01463 722000 Email: ymiddlef@its.jnj.com
Date Prepared	07th Feb 2013
Device Trade Name	OneTouch Verio Sync Blood Glucose Monitoring System
Common Name	Glucose Test System

Classification	OneTouch Verio Sync Blood Glucose Meters and OneTouch Verio Test Strips are Class II devices (21 CFR § 862.1345), Product Code NBW, LFR
System Description	The OneTouch® Verio™ Sync Blood Glucose Monitoring System consists of the OneTouch® Verio™ Sync Blood Glucose Meter, OneTouch® Verio™ Test Strips, OneTouch® Verio™ Level 3 and Level 4 Control Solutions, Lancing Device and Sterile Lancets. The OneTouch® Verio™ Sync Blood Glucose Monitoring System measures the glucose content of a blood sample by means of an electrical current produced in the test strip and sent to the meter for measurement.
Predicate Device	OneTouch Verio Blood Glucose Monitoring System (K093745, cleared on February 11, 2011)
Intended Use/Indications for Use	<p>The OneTouch Verio Sync Blood Glucose Monitoring System (BGMS) is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertip. The system is intended to be used by a single person and should not be shared.</p> <p>The OneTouch Verio Sync BGMS is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control. The OneTouch Verio Sync BGMS should not be used for the diagnosis of or screening of diabetes or for neonatal use.</p> <p>The OneTouch Verio Test Strips are for use with the OneTouch Verio Sync Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh</p>

	capillary whole blood samples drawn from the fingertips.
Comparison to Predicate Device	<p>The Subject device is different from the predicate device for the following aspects:</p> <ul style="list-style-type: none"> • Meter: shape/size, ergonomic design, simplified user interface, wireless communication with the OneTouch Reveal Diabetes Management Application and a rechargeable battery. <p>There have been no changes to the OneTouch® Verio™ Test Strip principle of operation, materials of construction, design of strip or electrodes, and they are identical to the test strips in the OneTouch® Verio™ Blood Glucose Monitoring System, K093745. The OneTouch® Verio™ Control Solutions composition and functional performance are also the same as those cleared in the OneTouch® Verio™ Blood Glucose Monitoring System, K093745. However, changes were made to the labeling provided with those components of the system at the request of the Food and Drug Administration during review of the OneTouch Verio IQ 510(k) Submission, K110637. There are some differences in the intended use, specifically around alternate site testing and the removal of the indication for forearm or palm testing. The OneTouch Verio Sync Blood Glucose Monitoring System is not validated for alternate site testing. Only</p>

	fingertip testing is indicated for this system.
Technological Characteristics	There has been no change to the fundamental scientific technology, which is amperometric detection. The operating principle remains electrochemical reaction.
Summary of Performance Characteristics	The OneTouch [®] Verio [™] Sync Blood Glucose Monitoring System (meter, strips, and control solutions) was tested in accordance with ISO 15197:2003(E). Analytical performance testing included system accuracy, repeatability, intermediate precision and linearity testing. A user performance evaluation assessed accuracy of results and usability of the device in the hands of intended users. The OneTouch [®] Verio [™] Sync Blood Glucose Monitoring System performed similarly to both the predicate device as well as to a laboratory reference method, the Yellow Springs Instrument (YSI).

System Accuracy

A comparison of system accuracy performance demonstrated that the OneTouch® Verio™ Sync Blood Glucose Monitoring System and the OneTouch® Verio™ Blood Glucose Monitoring System are substantially equivalent.

System Accuracy Results for Glucose Concentrations <75 mg/dL

Number (and percent) of meter results that match the laboratory test

Within ±5 mg/dL	Within ±10 mg/dL	Within ±15 mg/dL
40/57 (70%)	55/57 (96%)	57/57 (100%)

System Accuracy Results for Glucose Concentrations ≥75 mg/dL

Number (and percent) of meter results that match the laboratory test

Within ±5%	Within ±10%	Within ±15%	Within ±20%
174/243 (72%)	235/243 (97%)	242/243 (99.59%)	243/243 (100%)

Total Precision

(600 Control Solution Tests)

Glucose Level Ranges (mg/dL)	Mean Glucose (mg/dL)	Standard Deviation (mg/dL)	Coefficient of Variation (%)
Low (38-62)	38.63	0.88	2.27
Mid (102-138)	117.39	2.00	1.70
High (298-403)	335.91	6.25	1.86

User Performance Evaluation

Subject Fingertip Results for Glucose Concentrations <75 mg/dL

Tester	Within ±5 mg/dL	Within ±10 mg/dL	Within ±15 mg/dL
Subject	5/23 (22%)	15/23 (65%)	22/23 (96%)

Subject Fingertip Results for Glucose Concentrations ≥75 mg/dL

Tester	Within ±5%	Within ±10%	Within ±15%	Within ±20%
Subject	141/250 (56%)	219/250 (88%)	242/250 (97%)	248/250 (99%)

Design verification and validation testing confirmed that the performance, safety, and effectiveness of the OneTouch® Verio™ Sync Blood Glucose Monitoring System was equivalent to that of the predicate device. The OneTouch® Verio™ Sync Meter also met recognized electrical and safety standards.

Conclusions

The OneTouch Verio Sync Blood Glucose Monitoring System is substantially equivalent in its intended use, performance, safety, effectiveness and underlying scientific and operating principles used to the predicate, the OneTouch Verio Blood Glucose Monitoring System (K093745, 2/11/11).



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

February 7, 2013

Cilag GMBH International
c/o Nadine Nasr, Regulatory Project Manager
Lifescan, Inc.
1000 Gibraltar Drive
Milpitas, CA 95035

Re: k120708

Trade/Device Name: OneTouch Verio Sync Blood Glucose Monitoring System
Regulation Number: 21 CFR 862.1345
Regulation Name: Glucose test system
Regulatory Class: II
Product Code: NBW, LFR
Dated: January 09, 2013
Received: January 30, 2013

Dear Ms. Nasr:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Carol C. Benson for

Courtney H. Lias, Ph.D.
Director
Division of Chemistry and Toxicology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use Form

510(k) Number (if known): k120708

Device Name: OneTouch Verio Sync Blood Glucose Monitoring System

Indications for Use:

The OneTouch Verio Sync Blood Glucose Monitoring System (BGMS) is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertip. The system is intended to be used by a single person and should not be shared.

The OneTouch Verio Sync BGMS is intended for self testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitor the effectiveness of diabetes control. The OneTouch Verio Sync BGMS should not be used for the diagnosis of or screening of diabetes or for neonatal use.

The OneTouch Verio Test Strips are for use with the OneTouch Verio Sync Blood Glucose Meter to quantitatively measure glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips.

Prescription Use _____ AND/OR Over-The-Counter Use X
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER
PAGE OF NEEDED)

Concurrence of CDRH, Office of In Vitro Devices and Radiologic Health (OIR)

Katherine Serrano

Division Sign-Off
Office of In Vitro Devices and Radiologic Health

510(k) k120708