

APR 19 2012

**Rival® PTA Dilatation Catheter**

**510(k) Summary  
21 CFR 807.92**

As required by the Safe Medical Devices Act of 1990, coded under Section 513, Part (I)(3)(A) of the Food, Drug and Cosmetic Act, a 510(k) summary upon which substantial equivalence determination is based is as follows:

**Submitter Information:**

Applicant: Bard Peripheral Vascular, Inc  
1625 West 3<sup>rd</sup> Street  
Tempe, Arizona 85281

Phone: 480-303-2662

Fax: 480-449-2546

Contact: Candace Wade, Regulatory Affairs Specialist

Date September 1, 2011

**Subject Device Name:**

Device Trade Name: Rival® PTA Dilatation Catheter

Common or Usual Name: Percutaneous Catheter (21 CFR 870.1250,  
Product Code DQY and LIT)

Classification: Class II

Classification Panel: Cardiovascular

**Predicate Devices:**

- Edwards Peripheral Dilatation Catheter (K052149; cleared September 2, 2005)
- Dorado® PTA Dilatation Catheter (K072283; cleared September 19, 2007)

**Device Description:**

The Rival® PTA Dilatation Catheter is a 0.035" guidewire compatible balloon catheter consisting of an over the wire catheter with a balloon fixed at the distal tip. The non-compliant, low profile balloon is designed to provide consistent balloon diameters and lengths. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. The hub of the catheter includes inflation and guidewire ports with female luer locks. Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use. A re-wrapping tool is also provided on the catheter shaft. A stylet is placed into the tip of the catheter to aid in rewrap/refolding of the balloon.

Attribute	Rival PTA Dilatation Catheter Product Offering
Balloon Diameter (mm)	3, 4, 5, 6, 7, 8, 9, 10
Balloon Length (cm)	2, 4, 6, 8, 10, 15
Catheter Shaft Lengths (cm)	80, 120, 135
Introducer Sheath Compatibility (compatible balloon diameters, mm)	5F: (3, 4, 5, 6) 6F: (6, 7, 8, 9) 7F: (10)

**Indications for Use of Device:**

The Rival® PTA Dilatation Catheters are intended to dilate stenoses in the peripheral arteries, treat obstructive lesions of native or synthetic A-V fistulae, and/or re-expand endoluminal stent graft elements in the aorta and iliac arteries. This device is also recommended for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

**Comparison of Indications for Use to Predicate Devices:**

The indication for use statement for the Rival® PTA Dilatation Catheter does not raise any new issues of safety and effectiveness as demonstrated through the risk analysis process based on the proposed indications for use statement as compared to the predicate devices, the Edwards Peripheral Dilation Catheter and the Dorado® PTA Dilatation Catheter. Therefore, the subject device, the Rival® PTA Dilatation Catheter, is substantially equivalent to the predicate devices.

**Technological Comparison to Predicate Devices:**

The Rival® PTA Dilatation Catheter has the following similarities to the predicate devices:

- Similar intended use (combination of both predicates)
- Similar indications for use (combination of both predicates)
- Same target population (both predicates)
- Same fundamental scientific technology (both predicates)
- Same operating principle (both predicates)
- Same packaging materials and configuration (Dorado® PTA Dilatation Catheter)
- Same sterility assurance level and method of sterilization (both predicates)

**Performance Data:**

To demonstrate substantial equivalence of the subject device, the Rival® PTA Dilatation Catheter to the predicate devices, the technological characteristics and performance criterion were evaluated. Using FDA Guidance Documents on non-clinical testing of medical devices and internal Risk Assessment procedures, the following *in vitro* tests were performed:

- Puncture Resistance
- Robustness

The results from these tests demonstrate that the technological characteristics and performance criteria of the Rival® PTA Dilatation Catheter are comparable to the predicate devices and that it can perform in a manner equivalent to devices currently on the market for the similar intended use.

**Conclusions:**

The subject device, the Rival® PTA Dilatation Catheter, met all the predetermined acceptance criteria of design verification and validation as specified by applicable standards, guidance, test protocols and/or customer inputs. The Rival® PTA Dilatation Catheter is substantially equivalent to the legally marketed predicate devices, the Edwards Peripheral Dilatation Catheter and the Dorado® PTA Dilatation Catheter.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

APR 19 2012

Bard Peripheral Vascular, Inc.  
c/o Mark Job  
Responsible Third Party Official  
Regulatory Technical Services LLC  
1394 25<sup>th</sup> Street NW  
Buffalo, MN 55313

Re: K120722  
Trade/Device Name: Rival PTA Dilatation Catheter  
Regulation Number: 21 CFR 870.1250  
Regulation Name: Percutaneous catheter  
Regulatory Class: Class II (Two)  
Product Code: LIT, DQY  
Dated: March 8, 2012  
Received: March 9, 2012

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K120722

Device Name: Rival® PTA Dilatation Catheter

Indications for Use: The Rival® PTA Dilatation Catheters are intended to dilate stenoses in the peripheral arteries, treat obstructive lesions of native or synthetic A-V fistulae, and/or re-expand endoluminal stent graft elements in the aorta and iliac arteries. This device is also recommended for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

Prescription Use X  
(Part 21 CFR 801. Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number K120722