

K121667
P. 10A2

510(k) Summary
per 21 CFR §807.92

AUG 13 2012

Submitter's Name and Address	Boston Scientific Corporation Cardiovascular, Rhythm and Vascular Division One Scimed Place Maple Grove, MN 55311 Phone: 763-494-1700 Fax: 763-494-2222		
Contact Name and Information	Margo J. Anderson Principal Regulatory Affairs Specialist Phone: 763-494-1932 Fax: 763-494-2222 E-mail: Margo.Anderson@bsci.com		
Date Prepared	19 July 2012		
Proprietary Name	NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter		
Common Name	Percutaneous Transluminal Coronary Angioplasty (PTCA) Dilatation Catheter		
Product Code	LOX		
Classification	Class II, 21 CFR Part 870.5100		
Predicate Device	NC Quantum Apex™ PTCA Dilatation Catheter	P860019 /S241	16 April 2010
Reference Devices	Emerge™ PTCA Dilatation Catheter	K113220	22 March 2012
	Quantum™ Maverick™ PTCA Dilatation Catheter	P860019 /S182	01 October 2002
Device Description	<p>The Boston Scientific NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter is a sterile, single-use, intravascular medical device. The catheter consists of a shaft with a balloon near the distal tip. The balloon is designed to provide an inflatable segment of known diameter and length at recommended pressures. There are radiopaque marker bands located under the balloon to aid in positioning the system during the procedure. Coatings are applied to the balloon and catheter to enhance insertion and withdrawal performance.</p> <p>The NC Quantum Apex Monorail PTCA Dilatation Catheter is available with balloon diameters of 2.00 mm to 5.00 mm and balloon lengths of 6 mm to 30 mm.</p>		
Intended Use of Device	The NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter is intended for dilatation of stenosis in coronary arteries or bypass grafts and for post-delivery expansion of balloon expandable stents (bare metal and drug-eluting.)		

Indications for Use	<p>The NC Quantum Apex™ Monorail™ PTCA Dilatation Catheters are indicated for the balloon catheter dilatation of the stenotic portion of a native coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. NC Quantum Apex Monorail PTCA Dilatation Catheters are also indicated for the post-delivery expansion of balloon expandable stents (bare metal and drug-eluting).</p>										
Comparison of Technological Characteristics	<p>The NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter was modified to incorporate the same device design for the hypotube, midshaft component and proximal shaft marks as referenced devices Emerge™ PTCA Dilatation Catheter (K113220, cleared March 22, 2012) and Quantum™ Maverick™ PTCA Dilatation Catheter (P860019/S182, approved October 1, 2002).</p> <p>The modified NC Quantum Apex Monorail PTCA Dilatation Catheter differs from the Boston Scientific predicate device, NC Quantum Apex™ PTCA Dilatation Catheter (P860019/S241, approved April 16, 2010), by the removal of an accessory kit pouch, flushing needle and one CLIP-IT™ hypotube clip. Otherwise the NC Quantum Apex Monorail PTCA Dilatation Catheter incorporates substantially equivalent device materials, design, packaging materials, fundamental technology, manufacturing processes, sterilization process and intended use as those featured in the predicate device.</p>										
Summary of Non-Clinical Performance Data	<p>The NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter was subjected to testing according to the requirements of <i>Guidance for Industry and FDA Staff – Class II Special Controls for Certain Percutaneous Transluminal Coronary Angioplasty (PTCA) Catheters</i>, September 8, 2010. Design verification testing was performed to demonstrate that the performance of the modified NC Quantum Apex Monorail PTCA Dilatation Catheter remains substantially equivalent to the predicate device.</p> <p>Specifically, the following tests were evaluated to support changes for the modified NC Quantum Apex Monorail PTCA Dilatation Catheter:</p> <table border="0"><tr><td>Effective Length</td><td>Balloon Deflation Time</td></tr><tr><td>Proximal Shaft Profile</td><td>Corrosion Resistance</td></tr><tr><td>Shaft and Bond Burst Pressure</td><td>Torque Strength</td></tr><tr><td>Catheter Bond Strength</td><td>Particulate Evaluation</td></tr><tr><td>Flexibility and Kink</td><td>Biocompatibility</td></tr></table>	Effective Length	Balloon Deflation Time	Proximal Shaft Profile	Corrosion Resistance	Shaft and Bond Burst Pressure	Torque Strength	Catheter Bond Strength	Particulate Evaluation	Flexibility and Kink	Biocompatibility
Effective Length	Balloon Deflation Time										
Proximal Shaft Profile	Corrosion Resistance										
Shaft and Bond Burst Pressure	Torque Strength										
Catheter Bond Strength	Particulate Evaluation										
Flexibility and Kink	Biocompatibility										
Conclusion	<p>Based on the indications for use, technological characteristics, and safety and performance testing, the NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter has been shown to be appropriate for its intended use and is considered to be substantially equivalent to the Boston Scientific predicate device, NC Quantum Apex™ PTCA Dilatation Catheter.</p>										



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room --WO66-G609
Silver Spring, MD 20993-0002

AUG 13 2012

Boston Scientific Corp.
Ms. Laura Lind
Regulatory Affairs Specialist
One Scimed Place
Maple Grove, MN 55311

Re: K121667

Trade/Device Name: NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter
Regulation Number: 21 CFR 870.5100
Regulation Name: Catheters, Transluminal Coronary Angioplasty, Percutaneous
Regulatory Class: Class II
Product Code: LOX
Dated: July 19, 2012
Received: July 20, 2012

Dear Ms. Lind:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Ms. Laura Lind

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K121667

Device Name: NC Quantum Apex™ Monorail™ PTCA Dilatation Catheter

Indications for Use:

The NC Quantum Apex Monorail PTCA Dilatation Catheters are indicated for the balloon catheter dilatation of the stenotic portion of a native coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

NC Quantum Apex Monorail PTCA Dilatation Catheters are also indicated for the post-delivery expansion of balloon expandable stents (bare metal and drug-eluting).

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER
PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K121667