

**510(K) SUMMARY
OF SAFETY AND EFFECTIVENESS**

In accordance with the Food and Drug Administration Rule to implement provisions of the Safe Medical Devices Act of 1990 and in conformance with 21 CFR 807.92, this information serves as a Summary of Safety and Effectiveness for the use of the Vital Sync™ System.

- A.1. Submitted By:** Covidien
6135 Gunbarrel Avenue
Boulder, CO 80301
- Date:** February 14, 2013
- Contact Person:** Kelsey Lee
Senior Regulatory Affairs Specialist
(303) 305-2760
- A.2. Proprietary Name:** Vital Sync™ System, Model 5000
- Common Name:** Cardiac Monitor (without alarms)
- Device Classification Regulation:** 21 CFR 870.2300 – Class II
- Device Product Code & Panel:** MWI: Monitor, Physiological, Patient (without arrhythmia detection or alarms)

Cardiovascular
- A.3. Predicate Device:** Vital Sync™ System (K093244)

A.4. Device Description

The Vital Sync™ System is being extended to offer additional device compatibility and parameter display.

The design features of the subject Vital Sync™ System are substantially equivalent to the design features of the predicate Vital Sync™ system.

A.5. Intended Use

The subject Vital Sync™ System has the same indications for use as the predicate Vital Sync™ System. The only difference is the addition of displayed parameters in the subject Vital Sync™ System and the devices that are compatible.

The Indications for use are as follows:

The Vital Sync™ System is intended for the display and recording of multiple physiological parameters of adult, pediatric and neonatal patients. It is not intended for alarm notification, nor is it intended to control any of the independent bedside devices it is connected to. A listing of supported devices and displayed parameters is attached.

WARNING: The Vital™ Sync System is not an active patient monitoring system. It is intended to supplement and not replace any part of the hospital's device monitoring

A.6. Technological Characteristics Comparison

The subject Vital Sync™ System and the legally marketed predicate Vital Sync™ System have identical indications for use and similar displayed parameters and device compatibility, have the same principles of operation, and utilize the same hardware.

The subject Vital Sync™ System includes additional device compatibility and additional parameter display.

B.1. Substantial Equivalence – Non-Clinical Evidence

Substantial equivalence was shown through driver validation and regression testing. The driver validation and regression tests were run to show that the data collected from the bedside devices was not being corrupted and would be shown accurately on the subject device. The results of the tests show that the subject Vital Sync™ System is substantially equivalent to the legally marketed predicate Vital Sync™ System.

B.2. Substantial Equivalence – Clinical Evidence

N/A – Clinical evidence was not necessary to show substantial equivalence

B.3. Substantial Equivalence – Conclusions

Substantial equivalence is shown through driver validation and regression testing. The indications for use, intended population, hardware utilized and principles of operation are identical between the subject and predicate. The subject and predicate differ in that the subject is capable of additional device compatibility and displaying additional parameters. No new questions of safety and effectiveness have been raised. From the evidence presented in the Premarket Notification, the subject device can be expected to perform at least as well as the predicate.

Displayed Parameters

Amplitude Integrated EEG Left	End Diastolic Volume Index
Amplitude Integrated EEG Right	Stat End Diastolic Volume
Anesthetic Agent	Stat End Diastolic Volume Index
Air Temperature	Exhaled Minute Volume
Air Temperature Setting	End Inspiratory Pressure
Airway Temperature	End Systolic Volume
Arterial Base Excess	End Systolic Volume Index
Arterial Bicarbonate	End Tidal CO2
Arterial pH	Environment Temperature
Arterial Temperature	Esophageal Temperature
Average Heart Rate	Exhaled Tidal Volume
Axillary Temp	Exhalation Time
Bi-level Positive Airway Pressure	Expired Positive Airway Pressure
BIS	Inspired Fraction of Oxygen
Bladder Temperature	Infusate Temp
Blood Temperature/Pulm. Artery Temperature	Flow Rate 1
Body Surface Area	Flow Rate 2
Bolus Cardiac Output	Flow Rate 3
Bolus Cardiac Index	Gastric pCO2
Brain PO2	Heart Rate
Brain Temp	Heater Output Percent
Calculated SO2	Heater Output Percent Setting
Cardiac Index	Hematocrit
Cardiac Output	Hemoglobin
Stat Cardiac Index	High Inspired Pressure Setting
Stat Cardiac Output	Humidity
Cardioplegia Line Pressure	Humidity Setting
Central Venous Pressure	Variation of Contractility Index
Cerebral Perfusion Pressure	Inspired : Expired Ratio
Cerebral Blood Flow	Inspired CO2
Ch1 rSO2	Inspired O2 Setting
Ch2 rSO2	Insp Pos Air Pressure
Ch3 rSO2	Inspiratory Pressure
Contractility Index	Inspiratory Resistance
Continuous Cardiac Index	Inspiratory Tidal Volume
Continuous Cardiac Output	Inspiratory Time
Control Temp	Intracranial Pressure
Core Temperature	Diastolic Arterial lood Pressure
Coronary Sinus Pressure	Mean Arterial Blood Pressure
Continuous Positive Airway Pressure	Systolic Arterial Blood Pressure
Delta Pressure	% Leak in Tidal Volume (Insp/Exp)
Dynamic Compliance	Left Arterial Pressure

Ejection Fraction	Left Cardiac Work
Emboli 1	Left Cardiac Work Index
Emboli 2	Left Stroke Work
Emboli 3	Left Stroke Work Index
End Diastolic Volume	Left Ventricular Ejection Time
Line Pressure	Static Compliance
Mean Airway Pressure	ST Interval
Mattress Temperature	Stat Stroke Volume Index
Mattress Temperature Setting	Stat Stroke Volume Index
Respiration Rate	Venous SO2
Ventilation Mode	Oxygenator Sweep
Myocardial Temperature	Systolic Time Ratio
Nasopharyngeal Temperature	Thoracic Fluid Index
Nitric Oxide	Tidal Volume Setting
Nitric Dioxide	Transcutaneous pCO2
Diastolic Cuff Blood Pressure	Transcutaneous pO2
Mean Cuff Blood Pressure	Tympanic Temperature
Systolic Cuff Blood Pressure	Venous pH
Oxygen Consumption	Venous Temperature
Oxygen Extraction Index	Volume Flow1
Oral Temp	Volume Flow 2
Arterial pCO2	Volume Flow 3
Arterial pO2	Water Temperature
Peak Flow	Weight
Plateau Time	CoHB
Positive Pressure Duration	MetHB
Potassium	Perfusion Index
Pressure Control	O2 Delivery
Pressure Limit	FICO2
Pressure Sensitivity	Amplitude
Pressure Support	Amplitude 2
Pulmonary Capillary Wedge Pressure	Mode
Pulmonary Arterial Diastolic Pressure	Inspired O2
Pulmonary Arterial Mean Pressure	Suppression Ratio
Pulmonary Arterial Systolic Pressure	Insp Time Set
Pulse Amplitude	Low IP Set
Pulse Rate	PIP Set
Pump Flow	Vent Rate Set
Venous pCO2	CPAP Set
Venous pO2	Dose Mode
Rate Pressure Product	Dose Value
Rectal Temperature	Flow Chk
Right Ventricular Ejection Fraction	History Cleared

Stat Right Ventricular Ejection Fraction	Infusion Mode
Arterial SO2	Label
Set Point Temperature	Primary Rate
Skin Temperature	Pri Vol Infused
Skin Temperature Setting	Vol Remaining
Pulse SO2 1	Secondary Rate
Pulse SO2 2	Sec Vol Infused
Spectral Edge Frequency	Time Remaining
Spontaneous Respiration Rate	Central Ven O2 Sat
O2 Delivery	Shunt/Tot Flow
O2 Delivery Index	Slope
Stroke Volume	Specimen Type
Stroke Volume Index	Taped At
Stroke Volume Variation	Testing Lab
Systemic Vasc Resistance	Vital Capacity
Systemic Vasc Resistance Idx	
Inter beat interval	
Heart rate variability	
Baroreflex sensitivity	
Total peripheral resistance	
Total arterial compliance	
Max. steepness of current upstroke	
Ascending aortic impedance at DIA	
Pulse Press Variation	
Systolic Press Variation	
Umbilical Artery Pressure Systolic	
Umbilical Artery Pressure Diastolic	
Umbilical Artery Pressure Mean	
Abd Girth	
Alveolar-Art	
Amp Power	
Art/Ven O2 Diff	
Body Site	
CRT	
ETT Size	
Gas Temp	
Head Circ	
Hertz	
Length	
Pain Scale	
Pathologic Source	
QP:QS	
Sample Method	

Compatible Devices

Atom Medical	Infra Warmer V505
Baxter Healthcare	AS 50
Baxter Healthcare	Colleague IP
Baxter Healthcare	FloGuard 6201
Baxter Healthcare	FloGuard 6301
Bird	VIP Gold/Sterling Ventilator
Cardiotronic	Aesculon Noninvasive Cardiac Output
Cardiotronic	ICON Hemodynamic Monitor
CDI	CDI 500 Bloodgas Monitor
Cincinnati Sub-Zero	Blanketrol II
Cincinnati Sub-Zero	Blanketrol III
Covidien	PB-840 Ventilator
Covidien	N200 Pulse Oximeter
Covidien	N295 Pulse Oximeter
Covidien	N395 Pulse Oximeter
Covidien	N595 Pulse Oximeter
Covidien	InfantStar 500 Ventilator
Covidien	InfantStar 950 Ventilator
Covidien	N600 Pulse Oximeter
Covidien	N600x Pulse Oximeter
Covidien	N85 EtCO2 Monitor
Covidien	N75 EtCO2 Monitor
Covidien	BIS VISTA (Complete Monitoring System)
Datex/Ohmeda	S/5 Patient Monitor
Draeger	Babylog 8000 Ventilator
Draeger	Babylog 8000SC Ventilator
Draeger	Evita Ventilator
Draeger	Evita 2 Ventilator
Draeger	Evita 4 Ventilator
Draeger	Infinity Patient Monitor
Draeger	SC7000
Draeger	SC8000
Draeger	SC9000XL
Edwards Life Sciences	Vigilance Hemodynamic Monitor
Edwards Life Sciences	Vigilance II Hemodynamic Monitor
GE	Dash 2000 Patient Monitor
GE	Dash 3000 Patient Monitor
GE	Dash 4000 Patient Monitor
GE	Solar 8000i Patient Monitor
GE	Solar 8000M Patient Monitor
IKARIA	INOvent NO Delivery

K123002

LiDCO	LiDCORapid Hemodynamic Monitor
Maquet/Siemens	Servo 300 Ventilator
Maquet/Siemens	Servo i Ventilator
Maquet/Siemens	Maquet Perfusion Pump
Masimo	SET Radical-7 Pulse Oximeter
Masimo	SET Radical-9 Pulse Oximeter
Masimo	SET Radical-8 Pulse Oximeter
Mennen Medical	Horizon 2000 Patient Monitor
Philips	CMS 2001 Patient Monitor
Philips	V24 Patient Monitor
Philips	V26 Patient Monitor
Philips	MP5 Patient Monitor
Philips	MP20 Patient Monitor
Philips	MP40 Patient Monitor
Philips	MP50 Patient Monitor
Philips	MP60 Patient Monitor
Philips	MP70 Patient Monitor
Philips	MP80 Patient Monitor
Philips	MP90 Patient Monitor
Philips	Capnogard EtCO2 Monitor
Somanetics	INVOS 5100B Cerebral/Somatic Oximeter
Somanetics	INVOS 5100C Cerebral/Somatic Oximeter
Sorin	Sorin S3 Perfusion Pump System
Sorin	Sorin C5 Perfusion Pump System
Sorin	Sorin S5 Perfusion Pump System
Sorin	Sorin SIII Encore Perfusion Pump System
Spectrum Medical	Spectrum M2 Flow Meter
Spectrum Medical	Spectrum M3 Flow Meter
Transonic Systems	Transonic HT 109 Flow Meter
Viasys	Viasys Avea Ventilator
Welch-Allyn	Propaq Physiomonitor



February 27, 2013

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-002

Covidien
c/o Ms. Kelsey Lee
Senior Regulatory Affairs Specialist
6135 Gunbarrel Avenue
Boulder, CO 80301

Re: K123002
Trade/Device Names: Vital Sync™ System, Model 5000
Regulatory Number: 21 CFR 870.2300
Regulation Name: Cardiac Monitor (including Cardiotachometer and Rate Alarm)
Regulatory Class: Class II (Two)
Product Code: MWI
Dated: January 31, 2013
Received: February 1, 2013

Dear Ms. Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Ms. Kelsey Lee

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Owen P. Faris -S

for Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K123002

Device Name: Vital Sync™ System, Model 5000

Indications for Use:

The Vital Sync™ System is intended for the display and recording of multiple physiological parameters of adult, pediatric and neonatal patients. It is not intended for alarm notification, nor is it intended to control any of the independent bedside devices it is connected to.

The Vital Sync™ System displayed parameters are listed on the following pages.

WARNING: The Vital Sync™ System is not an active patient monitoring system. It is intended to supplement and not to replace any part of the hospital's device monitoring.”

Prescription Use AND/OR Over-The-Counter Use
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE –
CONTINUE ON ANOTHER PAGE IF NEEDED)

CONCURRENCE OF CDRH, OFFICE OF DEVICE EVALUATION (ODE)

Owen P. Faris -S

Vital Sync™ System Displayed Parameters

Amplitude Integrated EEG Left	Exhaled Tidal Volume
Amplitude Integrated EEG Right	Exhalation Time
Anesthetic Agent	Expired Positive Airway Pressure
Air Temperature	Inspired Fraction of Oxygen
Air Temperature Setting	Infusate Temp
Airway Temperature	Flow Rate 1
Arterial Base Excess	Flow Rate 2
Arterial Bicarbonate	Flow Rate 3
Arterial pH	Gastric pCO ₂
Arterial Temperature	Heart Rate
Average Heart Rate	Heater Output Percent
Axillary Temp	Heater Output Percent Setting
Bi-Level Positive Airway Pressure	Hematocrit
BIS	Hemoglobin
Bladder Temperature	High Inspired Pressure Setting
Blood Temperature/Pulm. Artery Temperature	Humidity
Body Surface Area	Humidity Setting
Bolus Cardiac Output	Variation of Contractility Index
Bolus Cardiac Index	Inspired: Expired Ratio
Brain PO ₂	Inspired CO ₂
Brain Temp	Inspired O ₂ Setting
Calculated SO ₂	Insp Pos Air Pressure
Cardiac index	Inspiratory Pressure
Cardiac Output	Inspiratory Resistance
Stat Cardic Index	Inspiratory Tidal Volume
Stat Cardic Output	Inspiratory Time
Cardioplegia Line Pressure	Intracranial Pressure
Central Venous Pressure	Diastolic Arterial Blood Pressure
Cerebral Perfusion Pressure	Mean Arterial Blood Pressure
Cerebral Blood Flow	Systolic Arterial Blood Pressure
Ch1 rSO ₂	% Lead in Tidal Volume (Insp/Exp)
Ch2 rSO ₂	Left Arterial Pressure
Ch3 rSO ₂	Left Cardiac Work
Contractility Index	Left Cardiac Work Index
Continuous Cardiac Index	Left Stroke Work
Continuous Cardiac Output	Left Stroke Work Index
Control Temp	Left Ventricular Ejection Time
Core Temperature	Line Pressure
Coronary Sinus Pressure	Mean Airway Pressure
Continuous Positive Airway Pressure	Mattress Temperature
Delta Pressure	Mattress Temperature Setting
Dynamic Compliance	Respiration Rate
Ejection Fraction	Ventilation Mode
Emboli 1	Myocardial Temperature
Emboli 2	Nasopharyngeal Temperature
Emboli 3	Nitric Oxide
End Diastolic Volume	Nitric Dioxide
End Diastolic Volume Index	Diastolic Cuff Blood Pressure
Stat End Diastolic Volume Index	Mean Cuff Blood Pressure
Exhaled Minute	Systolic Cuff Blood Pressure
End Inspiratory Pressure	Oxygen Consumption
End Systolic Volume	Oxygen Extraction Index
End Systolic Volume Index	Oral Temp
End Tidal CO ₂	Arterial pCO ₂
Environmental Temperature	Arterial pO ₂
Esophageal Temperature	Peak Flow

Plateau Time	Flow Chk
Positive Pressure Duration	History Cleared
Potassium	Infusion Mode
Pressure Control	Label
Pressure Limit	Primary Rate
Pressure Sensitivity	Pri Vol Infused
Pressure Support	Vol Remaining
Pulmonary Capillary Wedge Pressure	Secondary Rate
Pulmonary Arterial Diastolic Pressure	Sec Vol Infused
Pulmonary Arterial Mean Pressure	Time Remaining
Pulmonary Arterial Systolic Pressure	Central Ven O2 Sat
Pulse Amplitude	O2 Delivery Index
Pulse Rate	Stroke Volume
Pump Flow	Stroke Volume Index
Venous pCO2	Stroke Volume Variation
Venous pO2	Systemic Vasc Resistance
Rate Pressure Product	Systemic Vasc Resistance Idx
Rectal Temperature	Inter Beat interval
Right Ventricular Ejection Fraction	Heart Rate Variability
Stat Right Ventricular Ejection Fraction	Baroreflex Sensitivity
Arterial SO2	Total Peripheral Resistance
Set Point Temperature	Total Arterial Compliance
Skin Temperature	Max. Steepness of Current Upstroke
Pulse SO2 1	Ascending Aortic Impedence at DIA
Pulse SO2 2	Pulse Press Variation
Spectral Edge Frequency	Systolic Press Variation
Spontaneous Respiration Rate	Umbilical Artery Systolic Pressure
Static Compliance	Umbilical Artery Diastolic Pressure
St Interval	Umbilical Artery Mean Pressure
Stat Stroke Volume Index	Abd Girth
Venous SO2	Alveolar-Art
Oxygenator Sweep	Amp Power
Systolic Time Ratio	Art/Ven O2 Diff
Thoracic Fluid Index	Body Site
Tidal Volume Setting	CRT
Transcutaneous pCO2	ETT Size
Transcutaneous pO2	Gas Temp
Tympanic Temperature	Head Circ
Venous pH	Hertz
Venous Temperature	Length
Volume Flow 1	Pain Scale
Volume Flow 2	Pathologic Source
Volume Flow 3	QP:QS
Water Temperature	Sample Method
Weight	Shunt/Tot Flow
CoHB	Slope
MetHB	Specimen Type
Perfusion Index	Taped At
O2 Delivery	Testing
PICO2	Lab
Amplitude	Vital Capacity
Amplitude 2	
Mode	
Inspired O2	
Suppression Ratio	
Insp Time Set	
Low IP Set	
PIP Set	
Vent Rate Set	
CPAP Set	
Dose Mode	
Dose Value	