

5. 510(k) SUMMARY

DATE: February 17, 2012

OWNER: Northstar Healthcare Holdings
70 Sir John Rogerson's Quay
Dublin 2, Ireland

NOV 16 2012

OFFICIAL CORRESPONDENT: Michael Riordan
Operations Manager
Telephone: 00353-21-4548255
Fax: 00353-21-4548294
Email: michael.riordan@mckesson.ie

DEVICE NAME: **Trade Name:** Textured, Blue, Latex, Powder Free Examination
Gloves with Protein Labeling Claim (50µg/dm²
Or Less of Water Soluble Protein)

Common Name: Patient Examination Gloves

Classification: Patient Examination Gloves

Class: Class I

Product Code: LYY **Regulation:** 880.6250

PREDICATE DEVICE(S):

Predicate 510(k)	Device Name	Indication	Clearance Date	Company
K062917	Latex Powder Free Examination Gloves	The latex examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	02 Nov 2006	Wear Safe SDN, BHD

DEVICE DESCRIPTION: Textured, Blue, Latex, Powder Free Examination Gloves with Protein Labeling Claim (50µg/dm² Or Less of Water Soluble Protein)

STATEMENT OF INTENDED USE:

The Latex Powder Free Examination Glove is a Disposable device intended for medical and dental purposes to be worn on the examiner's hand or finger to prevent contamination between patient and examiner.

TECHNOLOGICAL CHARACTERISTICS:

The Latex Powder Free Examination Glove is substantially equivalent to the predicate device with regard to physical characteristics, design, product features, and intended use. Both gloves are made with latex using similar manufacturing processes.

Feature	Latex Powder Free Examination Gloves K062917 Predicate	Latex Powder Free Examination Gloves (Proposed)																					
Intended Use	Intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	Same																					
Indications for Use Statement	The latex examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	The latex examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.																					
Description	Powder free, examination gloves made of latex. The gloves are provided in sizes small, medium, large and extra large.	Powder free, examination gloves made of latex and colored blue. The gloves are provided in sizes small, medium, large and extra large.																					
Presentation	Non-Sterile gloves are provided in dispenser boxes.	Same																					
Material	Latex	Same																					
Single Use	Yes	Same																					
Dimensions	Meets ASTM D3578-05	<table border="0"> <tr> <td>Length</td> <td>Small, Medium, Large, XL</td> <td>285 mm, min.</td> </tr> <tr> <td>Width</td> <td>Small</td> <td>70-90mm</td> </tr> <tr> <td></td> <td>Medium</td> <td>85-105 mm</td> </tr> <tr> <td></td> <td>Large</td> <td>101-121 mm</td> </tr> <tr> <td></td> <td>XL</td> <td>105-125 mm</td> </tr> <tr> <td>Thickness</td> <td>Finger</td> <td>0.30 mm min.</td> </tr> <tr> <td></td> <td>Palm</td> <td>0.20 mm min.</td> </tr> </table>	Length	Small, Medium, Large, XL	285 mm, min.	Width	Small	70-90mm		Medium	85-105 mm		Large	101-121 mm		XL	105-125 mm	Thickness	Finger	0.30 mm min.		Palm	0.20 mm min.
Length	Small, Medium, Large, XL	285 mm, min.																					
Width	Small	70-90mm																					
	Medium	85-105 mm																					
	Large	101-121 mm																					
	XL	105-125 mm																					
Thickness	Finger	0.30 mm min.																					
	Palm	0.20 mm min.																					
Physical Properties	Meets ASTM D3578-05	<table border="0"> <tr> <td></td> <td></td> <td colspan="2" style="text-align: center;"><u>Before aging/after aging</u></td> </tr> <tr> <td>Elongation</td> <td></td> <td>650%</td> <td>500%</td> </tr> <tr> <td>Tensile Strength</td> <td></td> <td>18MPa</td> <td>14MPa</td> </tr> </table>			<u>Before aging/after aging</u>		Elongation		650%	500%	Tensile Strength		18MPa	14MPa									
		<u>Before aging/after aging</u>																					
Elongation		650%	500%																				
Tensile Strength		18MPa	14MPa																				
Freedom from Pinholes	Meets ASTM D5151-06	Same																					
Residual Powder	Meets ASTM D6124-06	Same																					
Protein Level	Meets ASTM D5712-05	Meets ASTM D5712-10																					
Biocompatibility Tests	Passes Primary Skin Irritation in Rabbits Passes Guinea Pig Maximization	Same Same																					

**ASSESSMENT OF
NONCLINICAL DATA:**

Characteristic	Standard	Device Performance
Dimension	ASTM Standard D3578-05	Meets
Physical Properties	ASTM Standard D3578-05	Meets
Freedom from Pinholes	21 CFR 800.20; ASTM D5151-06	Meets
Powder Residual	ASTM Standard D6124-06	Meets Results generated values below 2mg of residual powder
Protein Level	ASTM Standard 5712-10	Meets Results generated values below 50 mcg/g of protein
Biocompatibility	Primary Skin Irritation in rabbits (ISO 10993-10:2010)	Gloves are non-irritating
	Dermal Sensitization in the guinea pig (ISO 10993-10:2010)	Gloves do not display any potential for sensitization

CONCLUSIONS:

The Latex Powder Free Examination Gloves meet the requirements of established standards ASTM D3578-05, ASTM D5712-10, ASTM D5151-06, ASTM D6124-06 and ISO 10993-10:2010.

Based on the comparison of intended use, design, materials and performance, the Latex Powder Free Examination Gloves are substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-002

November 16, 2012

Northstar Healthcare Holdings
C/O Mr. Ned Devine
Underwriters Laboratories, Incorporated
333 Pfingsten Road
Northbrook, Illinois 60062

Re: K123026

Trade/Device Name: Textured, Blue, Latex, Powder Free Examination Gloves with Protein Labeling Claim ($50\mu\text{g}/\text{dm}^2$ Or Less of Water Soluble Protein)

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: I

Product Code: LYY

Dated: October 25, 2012

Received: October 31, 2012

Dear Mr. Devine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Devine

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



DN: c=US, o=U.S. Government, ou=HHS,
ou=FDA, ou=People, cn=Anthony D.
Watson,
0.9.2342.19200300.100.1.1=1300092402

Anthony D. Watson, B.S., M.S., M.B.A.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K123026

Device Name: **Textured, Blue, Latex, Powder Free Examination
Gloves with Protein Labeling Claim (50µg/dm²
Or Less of Water Soluble Protein)**

Indications for Use: The examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

AND/OR

Prescription Use _____
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use X
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Digitally signed by Elizabeth F. Claverie

Date: 2012.11.16 12:42:41 -05'00'

(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K123026