

PULPDENT CORPORATION

80 Oakland Street
Watertown, MA 02472
USA

Tel: 617 926 6666
Fax: 617 926 6262
Email: Pulpdent@pulpdent.com

510(k) SUMMARY

K123265

DATE OF SUBMISSION: October 5, 2012

OWNER: *Pulpdent Corporation*

CONTACT:

Kenneth J. Berk
80 Oakland Street
PO Box 780
Watertown, MA 02472 USA

Telephone: 617-926-6666
Fax: 617-926-6262
Email: ken@pulpdent.com

DEVICE:

Trade Name: *RMGI Low Viscosity*
Classification Name: Dental cement
FDA Product Code: EMA, 21 CFR Part 872.3275

PREDICATE DEVICES:

Pulpdent Tuff-Temp 2.0
3M RelyX Luting Plus
3M Vitrebond Plus LC
GC Fuji Lining LC
GC Fuji Plus Lining Capsule

INTENDED USE:

Pulpdent RMGI Low Viscosity is a resin-modified glass ionomer preparation used by dental professionals as a liner, base or luting material in dental restorations.

DESCRIPTION:

Pulpdent RMGI Low Viscosity is a resin-modified glass ionomer preparation with both a bioactive resin matrix and bioactive glass fillers. In this context, 'bioactive' refers to the release of beneficial ions from the resin and glass fillers into the oral environment. As a resin-modified glass ionomer, *RMGI Low Viscosity* has three setting mechanisms: light cure, self-cure resin chemistry, and acid-base glass ionomer reaction. It contains no Bisphenol A, no Bis-GMA and no BPA derivatives. *RMGI Low Viscosity* is a paste-paste formulation provided in an automix syringe.

SUMMARY OF PERFORMANCE TESTING – BENCH

	<u>Self-cure</u>	<u>Light cure</u>
Flexural strength	81.0 MPa	88.0 MPa
Compressive strength	208.0 MPa	239 MPa
Flexural modulus	3.7 GPa	Not tested
Diametral tensile strength	38.0 MPa	36.0 MPa
Working time	2 min 30 seconds	Not applicable
Setting time	3-3.5 min from beginning of mix	20 seconds
Radiopacity	Equivalent to 2 mm aluminum	Equivalent to 2 mm aluminum

PULPDENT CORPORATION

80 Oakland Street
Watertown, MA 02472
USA

Tel: 617 926 6666
Fax: 617 926 6262
Email: Pulpdent@pulpdent.com

COMPARISON WITH PREDICATE PRODUCTS:

Pulpdent RMGI Low Viscosity is substantially equivalent in design, composition, performance and intended use to the predicate products:

Product	Description	Intended Use	Composition
Pulpdent RMGI Low Viscosity	Resin-modified glass ionomer in two pastes	Dental liner, base or luting agent.	Part A: Diurethane dimethacrylate and other methacrylate-based monomers and oligomers, polyacrylic acid/maleic acid copolymer, water, barium borosilicate glass, silica, reducing agents, photoinitiator. Part B: Diurethane dimethacrylate and other methacrylate-based monomers and oligomers, aluminofluorosilicate ionomer glass, silica, oxidizing agents.
Pulpdent Tuff Temp 2.0 K120784	Glass-filled resin in two parts with photo-initiator and self-cure chemistry.	Temporary crown and bridge material	Diurethane dimethacrylate and other methacrylate monomers and oligomers in two pastes.
3M RelyX Luting Plus Automix K111185	Resin-modified glass ionomer in two pastes	Dental cement / luting agent	Part A: Fluoroaluminosilicate glass, HEMA, water, opacifying agent and reducing agent Part B: Zirconia silica filler, methacrylated polycarboxylic acid, HEMA, Bis-GMA, water, potassium persulfate, photoinitiator.
3M Vitrebond Plus Glass Ionomer Liner K011200	Resin-modified glass ionomer in two parts: paste and liquid	Dental cavity liner / base and cement	Paste HEMA, BIS-GMA, fluoroaluminosilicate glass, water and initiators. Liquid Methacrylate-modified polycarboxylic acid, water, HEMA, photoinitiator.
GC Fuji Lining LC K 020117, K051427	Resin-modified glass ionomer in two parts: powder and liquid	Cavity liner/base and dental cement	Powder: Alumino-silicate glass Liquid: HEMA, polyacrylic acid, water, UDMA, photo-initiator.
GC Fuji Plus Lining Capsule K 993973	Resin-modified glass ionomer in two parts: powder and liquid	Cavity liner/base and dental cement	Powder: Aluminofluorosilicate glass. Liquid: Polyacrylic acid, HEMA, tartaric acid, dimethacrylate, water.

PULPDENT CORPORATION

80 Oakland Street
Watertown, MA 02472
USA

Tel: 617 926 6666
Fax: 617 926 6262
Email: Pulpdent@pulpdent.com

	Compressive Strength MPa	Flexural Strength MPa	Flexural Modulus GPa	Diametral Tensile strength (DTS) MPa
Pulpdent RMGI Low Viscosity (auto cure)	208	81.0	3.7	38
Pulpdent RMGI Low Viscosity (light cure)	239	88.0		36
Pulpdent Tuff Temp 2.0	205.0	84.0		---
3M-RelyX Luting Plus	110	27		15
3M-Vitrebond Plus LC	114	28		18
GC Fuji Lining LC	135	26		12
GC-Fuji Plus Lining Capsule	156	31		15

CONCLUSION:

From the above comparisons and the bench testing, it can be concluded that *Pulpdent RMGI Low Viscosity* is substantially equivalent in design, composition, performance and intended use to the predicate products listed above. The predicate products have been found substantially equivalent under the 510(k) Premarket Notification process as Class II Dental Devices under CFR 872.3275 and have been used by dental professionals for more than 20 years. A search of the relevant scientific literature shows that resin-modified glass ionomers used for cavity lining and luting have been extensively studied. See References below.

REFERENCES

Specific to Pulpdent RMGI Low Viscosity

1. Pameijer CH, et al. Biocompatibility of four experimental formulations [including RMGI Low Viscosity] in subcutaneous connective tissue of rats. Pulpdent Corporation. 2011 Nov.
2. Pameijer CH. Report on the retention of Embrace WetBond cement and a RMGI cement (Pulpdent). Pulpdent Corporation 2012 Aug.
3. Pameijer CH et al. Marginal bacterial leakage with RMGI filled restorations and castings cemented with a RMGI luting cement. University of El Salvador/Argentine Dental Association, Buenos Aires, Argentina. 2012 Jun.

PULPDENT CORPORATION

80 Oakland Street
Watertown, MA 02472
USA

Tel: 617 926 6666
Fax: 617 926 6262
Email: Pulpdent@pulpdent.com

4. Pameijer CH. Microleakage of four experimental resin-modified glass ionomer restorative materials. Pulpdent Corporation. 2011 Apr.
5. Pameijer CH, Zmener O. Histopathological evaluation of an RMGI cement, auto- and light-cured, used as a luting agent. A subhuman primate study. Pulpdent Corporation. March 2011.

General

6. Booth SE. Surface properties and the bioactivity of glass-ionomer dental cements and related materials. PhD Thesis, University of Greenwich. 2010.
7. Comisi JC. Using bioactive materials to achieve proactive dental care. Oral Health. 2011 Dec. Dec;34-46.
8. Costa CA, et al. Pulp response after application of two resin modified glass ionomer cements (RMGICs) in deep cavities of prepared human teeth. Dent. Mater. 2011 Jun;27(7):e158-70.
9. Gandolfo MG, et al. Biomimetic Remineralization of human dentin using promising innovative calcium silicate hybrid "smart" materials. Dent Mater. 2011 Nov;27(11):1055-69.
10. Goldstep F. Proactive intervention dentistry: a model for oral care through life. Compend Contin Educ Dent. 2012 Jun;33(6):394-6, 398-402; quiz 404, 416.
11. Kotsanos N, Arizos S. Evaluation of a resin modified glass ionomer serving both as indirect pulp therapy and as restorative material for primary molars. Eur Arch Paediatr Dent. 2011 Jun;12(3):170-5.
12. Loof J et al. A comparative study of the bioactivity of three materials for dental applications. Dent Mater. 2008 May;24(5):653-9. Epub 2007 Aug 28.
13. Matsuya S, et al. Structure of bioactive glass and its application to glass ionomer cement. Dent Mater J. 1999 Jun;18(2):155-66.
14. Rajesh P, Kamath MP. Application of glass ionomer cements in restorative dentistry. Indian J Dent Res. 1999 Jul-Sep;10(3):88-90.
15. Mitra SB, et al. Fluoride release and recharge behavior of a nano-filled resin-modified glass ionomer compared with that of other fluoride releasing materials. Am J Dent. 2011 Dec;24(6):372-8.
16. Mousavinasab SM, et al. Flexural strength and morphological characteristics of resin-modified glass-ionomer containing bioactive glass. J Contemp Dent Pract. 2011 Jan 1;12(1):41-6.
17. Murray PE, et al. Restorative pulpal and repair responses. JADA 2001 Apr;132:482-91.
18. Nicholson JW, Czarnecka B. The biocompatibility of resin-modified glass-ionomer cements for dentistry. Dent Mater. 2008 Dec;24(12):1702-8.
19. Peumans M, et al. Two year clinical effectiveness of a resin-modified glass-ionomer adhesive. Am J Dent. 2003 Dec;16(6):363-8.

PULPDENT CORPORATION

80 Oakland Street
Watertown, MA 02472
USA

Tel: 617 926 6666
Fax: 617 926 6262
Email: Pulpdent@pulpdent.com

20. Salata LA et al. Recent advances in the use of glass ionomers: bone substitutes. Rev Odontol Univ Sao Paulo. 1999 Apr/June;13(2):203-7.
21. Sampaio PC, et al. Effect of conventional and resin-modified glass-ionomer liner on dentin adhesive interface of Class I cavity walls after thermocycling. Oper Dent. Jul-Aug;36(4):403-12.
22. Stanley HR. Local and systemic responses to dental composites and glass ionomers. Adv Dent Res. 1992 Sep;6:55-64.
23. Tarim B, et al. Pulpal response to a resin-modified glass-ionomer material on non-exposed and exposed monkey pulps. Quintessence Int. 1998 Aug;29(8):535-42.
24. Yli-Urpo H, et al. Release of silica, calcium, phosphorus, and fluoride from glass ionomer cement containing bioactive glass. J Biomater Appl. 2004 Jul;19(1):5-20.
25. Zhou SL, et al. In vitro study of the effects of fluoride-releasing dental materials on Remineralization in an enamel erosion model. J Dent. 2012 Mar;40(3):255-63.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-002

December 7, 2012

Mr. Kenneth J. Berk
Director
Pulpdent Corporation
80 Oakland Street
Watertown, Massachusetts 02472

Re: K123265

Trade/Device Name: Pulpdent RMGI Low Viscosity
Regulation Number: 21 CFR 872.3275
Regulation Name: Dental Cement, Cavity Liner
Regulatory Class: II
Product Code: EMA, EJK
Dated: October 16, 2012
Received: October 19, 2012

Dear Mr. Berk:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Berk

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Kwame O. Ulmer

Anthony D. Watson, B.S., M.S., M.B.A.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: *Pulpdent RMGI Low Viscosity*

Indications for Use:

Pulpdent RMGI Low Viscosity is a resin-modified glass ionomer preparation used by dental professionals as a liner, base or luting material in dental restorations.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Runner DDS, MA 2012.12.07
16:19:29 -05'00'

(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K123265

Page 1 of 1