



Olive Medical Corporation
2302 S Presidents Drive
Ste. D
SLC, UT 84120
+1. 866.300.1148
www.olivemedical.com

510(k) Summary

DEC 19 2012

510(k) Submitter Olive Medical Corp.
2302 South Presidents Dr. STE D
Salt lake City, UT 84120

Contact(s) Brian Dean
Telephone (801) 953-0559
Fax (801) 823-2238

Date Prepared October 26th, 2012

Trade Name OVS1 Video System

Common Name Endoscope And/Or Accessories

Classification Name Light Source, Fiberoptic, Routine; Laparoscope, General & Plastic Surgery
Regulation Number 21 CFR 876.1500
Product Code FCW; GCJ

Predicate Devices

| | | |
|--|-----------|------------------------------------|
| Tele Pack X | Not known | Karl Storz Endoscopy- America, Inc |
| Image 1 Video Imaging System With Option | K070716 | Karl Storz Endoscopy- America, Inc |
| InnerVue Diagnostic Scope System | K072879 | Biomet Sports Medicine |
| Flexible Video Endoscope with Sheath and Video Processor | K102733 | Vision-Sciences, Inc |
| LED Light Source | K093792 | Sunoptic Technologies, Inc. |

Device Description:

The OVS1 Video System is a Camera Control Unit with integrated LED Light Source and Video Display for use in a surgical environment to view endoscopic images when interfaced with the Olive Medical TCK1 HD Camera Head, an endoscope, a light cable, and other accessory devices.

Intended Use:

Indications for Use:

The OVS1 Video System is indicated for use in diagnostic and operative endoscopic procedures to provide illumination and visualization of an interior cavity of the body through either a natural or surgical opening. The OVS1 Video System is indicated for use with a compatible Olive Medical Camera Head and other accessory devices including an endoscope, optical coupler, and light cable

This Indications for Use statement is substantially equivalent to predicate devices.

General description of the diseases or conditions that the device will diagnose, treat, prevent, cure, or mitigate:

The OVS1 Video System is intended for use in endoscopic procedures without limit to the diseases or condition being treated by the treating surgeon.

Patient Population for which the device is intended:

The OVS1 Video System is intended for use with patients undergoing endoscopic procedures for which the treating surgeon desires external video display and/or illumination.

Comparison of Technological Characteristics

The OVS1 Video System is substantially equivalent to predicate devices since the device technology and design are similar. The Indications for Use statements are similar to predicate devices, and no new issues of safety or effectiveness are raised. The minor differences between the OVS1 Video System and predicate devices have no negative effect on performance, function, or intended use of the device.

Device Comparison Tables:

| | Display Technology | Display Size | Video Output Resolution (vertical lines) | Control Unit Technology |
|--|---------------------------|---------------------|---|---|
| OVS1 Video System | Integrated LCD | 12.1" | 720 | Similar |
| Tele Pack X Karl Storz Endoscopy- America, Inc 510(k) unknown | Integrated LCD | 15" | 494* | Similar |
| Image 1 Video Imaging System With Option Karl Storz Endoscopy- America, Inc K070716 | N/A | N/A | N/A | Similar – no on-board display or Light Source |
| InnerVue Diagnostic Scope System Biomet Sports Medicine K072879 | Integrated LCD | 6.4" | 480 | Similar |
| Flexible Video Endoscope with Sheath and Video Processor Vision-Sciences, Inc K102733 | Integrated LCD | 15" | 480 | Similar |
| LED Light Source Sunoptic Technologies, Inc. K093792 | N/A | N/A | N/A | N/A |
| Lightsource or Illuminator Sunoptic Technologies, Inc. K961074 | N/A | N/A | N/A | N/A |

*Output resolution unknown. Vertical lines of input listed

Section 5: 510(k) Summary

| | Light Source Technology | Light Source Safety |
|--|--------------------------------|--|
| OVSI Video System | LED | LED Module and output identical to LED Light Source, K093792 |
| Tele Pack X Karl Storz Endoscopy- America, Inc 510(k) unknown | Metal Halide | Light Source Safety demonstrated in 510(k) submission |
| Image 1 Video Imaging System With Option Karl Storz Endoscopy- America, Inc K070716 | N/A | N/A |
| InnerVue Diagnostic Scope System Biomet Sports Medicine K072879 | Xenon Arc | Light Source Safety demonstrated in 510(k) submission |
| Flexible Video Endoscope with Sheath and Video Processor Vision-Sciences, Inc K102733 | LED | Light Source Safety demonstrated in 510(k) submission |
| LED Light Source Sunoptic Technologies, Inc. K093792 | LED | Light Source Safety demonstrated in 510(k) submission |
| Lightsource or Illuminator Sunoptic Technologies, Inc. K961074 | Xenon | Light Source Safety demonstrated in 510(k) submission |

| | Performance Comparison | Target Population and Anatomical Site | Reuse durability | Skill Level Required |
|--|--|--|-------------------------|-----------------------------|
| OVSI Video System | Similar | Similar | Similar | Similar |
| Tele Pack X Karl Storz Endoscopy- America, Inc 510(k) unknown | Equivalent | Similar | Similar | Similar |
| Image 1 Video Imaging System With Option Karl Storz Endoscopy- America, Inc K070716 | Similar, no on-board display or light source | Similar | Similar | Similar |
| InnerVue Diagnostic Scope System Biomet Sports Medicine K072879 | Similar | Similar | Similar | Similar |
| Flexible Video Endoscope with Sheath and Video Processor Vision-Sciences, Inc K102733 | Similar | Similar | Similar | Similar |
| LED Light Source Sunoptic Technologies, Inc. K093792 | Identical light output | Similar | Similar | Similar |
| Lightsource or Illuminator Sunoptic Technologies, Inc. K961074 | Similar | Similar | Similar | Similar |

| | Energy Source | Materials | Biocompatibility | Sterility |
|--|------------------------|---|-------------------------|------------------|
| OVS1 Video System | External – wall outlet | Similar | N/A | N/A |
| Tele Pack X Karl Storz Endoscopy- America, Inc 510(k) unknown | External – wall outlet | Similar | N/A | N/A |
| Image 1 Video Imaging System With Option Karl Storz Endoscopy- America, Inc K070716 | External – wall outlet | Similar - no on-board display or light source | N/A | N/A |
| InnerVue Diagnostic Scope System Biomet Sports Medicine K072879 | External – wall outlet | Similar | N/A | N/A |
| Flexible Video Endoscope with Sheath and Video Processor Vision-Sciences, Inc K102733 | External – wall outlet | Similar | N/A | N/A |
| LED Light Source Sunoptic Technologies, Inc. K093792 | External – wall outlet | Similar – no on-board display | N/A | N/A |
| Lightsource or Illuminator Sunoptic Technologies, Inc. K961074 | External – wall outlet | Similar – no on-board display | N/A | N/A |

Non-Clinical Testing:

The OVS1 Video System demonstrates substantial equivalence in safety by tested compliance with ISO 60601-1: Medical electrical equipment – Part 1: General requirements for basic safety and essential performance; and ISO 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance — Collateral standard: Electromagnetic compatibility — Requirements and tests.

Clinical Testing:

No comparison of clinical performance data was used for demonstration of substantial equivalence

Substantial Equivalence Rationale:

The OVS1 Video System fulfills criteria regarding substantial equivalence with the Predicate Devices listed for the following reasons:

1. The intended use of the OVS1 Video System is equivalent to the Tele Pack X, 510(k) unknown; Image 1 Video Imaging System, K070716; and InnerVue Diagnostic Scope System, K072879
2. The Indications for Use statements are equivalent to the indications for use of the InnerVue Diagnostic Scope System, K072879, and the Flexible Endoscope Sheath and Video Processor, K102733
3. The Technologic characteristics of the OVS1 Video System are equivalent to Predicate Devices as follows:

Section 5: 510(k) Summary

- a. Video Processing System Equivalence: Tele Pack X, 510(k) unknown; Image 1 Video Imaging System, K070716; InnerVue Diagnostic Scope System, K072879; and Flexible Endoscope Sheath and Video Processor, K102733
 - b. Light Source Similarity: Tele Pack X, 510(k) unknown; Flexible Video Endoscope with Sheath and Video Processor, K102733, Lightsource or Illuminator, K961074
 - c. Light Source Identical: LED Light Source, K093792
 - d. Display Equivalence: Tele Pack X, 510(k) unknown; InnerVue Diagnostic Scope System, K072879; and Flexible Endoscope Sheath and Video Processor, K102733
4. Performance data show no new safety concerns for the OVS1 Video System
 5. Performance of the OVS1 Video System demonstrates equivalent ability to Predicate Devices to illuminate the surgical site, process video, and display video

510(k) Summary Final Criteria Checklist

- The summary includes only information that is also covered in the body of the 510(k).
- The summary does not contain any puffery or unsubstantiated labeling claims.
- The summary does not contain any raw data, i.e., contains only summary data.
- The summary does not contain any trade secret or confidential commercial information.
- The summary does not contain any patient identification information.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

December 19, 2012

Olive Medical Corp.
% Mr. Brian Dean
VP, RA/QA
2302 South Presidents Dr. STE D
SALT LAKE CITY UT 84120

Re: K123359
Trade/Device Name: OVS1 Video System
Regulation Number: 21 CFR§ 876.1500
Regulation Name: Endoscope and accessories
Regulatory Class: II
Product Code: FCW, GCJ
Dated: October 29, 2012
Received: October 31, 2012

Dear Mr. Dean:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Herbert R. Lerner

for

Benjamin R. Fisher, Ph.D.
Director
Division of Reproductive, Gastro-Renal,
and Urological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K123359

Device Name:

OVS1 Video System

Indications For Use:

The OVS1 Video System is indicated for use in diagnostic and operative endoscopic procedures to provide illumination and visualization of an interior cavity of the body through either a natural or surgical opening. The OVS1 Video System is indicated for use with a compatible Olive Medical Camera Head and other accessory devices including an endoscope, optical coupler, and light cable

Prescription Use XX

And/Or

Over-The Counter Use _____

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Herbert P. Lerner

(Division Sign-Off)

Division of Reproductive, Gastro-Renal, and
Urological Devices

510(k) Number K123359