

K123605

510(k) (Traditional) Submission
Section 5, 510(k) Summary

510(k) Summary

This 510(k) summary is being submitted in accordance with the requirements of 21 CFR Part 807.92 for the MR II MR Compatible Ventricular Catheter.

1. Company Making the Submission:

AUG 16 2013

Name of Owner:	MR II, Inc.
Address:	5 Musick Irvine, CA 92618
Telephone:	949-900-6833
Fax:	949-900-6834
Contact:	Edward Waddell
E-mail:	ewaddell@mriinterventions.com

2. Device Name:

Common Name:	Ventricular Catheter
Proprietary Name:	MR Compatible Ventricular Catheter
Classification:	Class II
Regulation Number:	882.4100
Product Code:	HCA

3. Predicate Device:

MR II Ventricular Cannula, K102101

4. Intended Use Statement:

The SmartFlow™ Flex Ventricular Catheter is intended for injection of Cytarabine or removal of CSF from the ventricles during intracranial procedures. The device is not intended for implant. This device is intended for "single patient use only."

5. Description of Device:

The Catheter has a stepped distal tip with a 30 cm removable rigid ceramic stylet protecting the fluid lumen while providing rigidity to the distal portion of the device. The stylet is removed after insertion to the desired point. Soft tubing protects the lumen in the center portion and at the distal end where it terminates. The fluid containing central lumen is manufactured from PEEK tubing.

The Catheter will be marketed in the following sizes:

Ventricular Catheter .008 x 7.5mm tip
 Ventricular Catheter .008 x 15mm tip
 Ventricular Catheter .021 x 7.5mm tip
 Ventricular Catheter .021 x 15mm tip

Each unit will provide for an approximate 23 inch tubing extension.

6. Summary of the Technological Characteristics of the Device Compared to the Predicate Device

	MRII Ventricular Catheter	Predicate Device: MRII Ventricular Cannula	Discussion
Classification	21 CFR 882.4100	21 CFR 882.4060	Equivalent
Product Code	HCA	HCD	Equivalent
Premarket Notification	TBD	K102101	Equivalent
Intended Use	The SmartFlow™ Flex Ventricular Catheter is intended for injection of Cytarabine or removal of CSF from the ventricles during intracranial procedures. The device is not intended for implant. This device is intended for "single patient use only."	The MR Compatible Ventricular Cannula is intended for injection of Cytarabine or removal of CSF from the ventricles during intracranial procedures. The device is not intended for implant. The device is intended for "single patient use only"	Equivalent
Target Population	Any Pt's needing aspiration or injection of fluids from the brain ventricles	Any Pt's needing aspiration or injection of fluids from the brain ventricles	Equivalent
Anatomical Sites	Brain ventricle	Brain ventricle	Equivalent
Where Used	OR or MRI suite	OR or MRI Suite	Equivalent
Energy used	N/A	N/A	N/A
Human Factors	Labeling indicates size and length	Labeling indicates size and length	Equivalent
	Can be manipulated with gloved hand	Can be manipulated with gloved hand	Equivalent
Design	Designed to be placed through a prepared opening through the skull and dura into the brain ventricle	Designed to be placed through a prepared opening through the skull and dura into the brain ventricle	Equivalent
	Rigid and straight section to enter the brain, rigid section can be removed after insertion	Rigid and straight section to enter the brain	Equivalent, except the rigid section can be removed after insertion of the catheter

	MRI Ventricular Catheter	Predicate Device: MRI Ventricular Cannula	Discussion
	Hole at distal end for fluid movement	Hole at distal end for fluid movement	Equivalent
	Length of rigid section: 10.5" (30 cm)	Length of rigid section: 10.5" (30 cm)	Equivalent
	Tip design identical to that of the Ventricular Cannula	Tip design as described in K102101	Equivalent
	Sufficiently rigid to pass through brain tissue without additional support	Sufficiently rigid to pass through brain tissue without additional support	Equivalent
	Contains a channel through which fluids can be removed (aspiration) or placed (injection) into the ventricle	Contains a channel through which fluids can be removed (aspiration) or placed (injection) into the ventricle	Equivalent. Flow rates are higher for the catheter due to the shorter tubing length.
	Flow rate of: 4.0 ml/hr (.008" I.D.) to 56 ml/hr (.021" I.D.) at 0.7 PSI	Flow rate of: 0.6 ml/hr (.008" I.D.) to 34 ml/hr (.021" I.D.) at 0.7 PSI	
	Rigid body: Ceramic	Rigid body: Ceramic	Equivalent
	Through lumen PEEK	Through lumen: Polymer covered silica	Equivalent, both meet biocompatibility requirements
	MRI Safe	MRI Safe	Equivalent
Materials	Proximal Connector: Female Luer adapter connector, external to kit as described in IFU	Proximal Connector: Female Luer connector	Equivalent, both can be connected to a syringe. Catheter requires a user supplied common component to make the connection.
Biocompatibility	Tissue contact tested per ISO 10993: Biological Evaluation of Medical Devices	Tissue contact tested per ISO 10993: Biological Evaluation of Medical Devices	Equivalent
Cytotoxicity	Acceptable	Acceptable	Equivalent
Systemic Toxicity, SC & SO extracts	Acceptable	Acceptable	Equivalent
Pyrogen Material Mediated	Acceptable	Acceptable	Equivalent
Intracutaneous Study, SC & SO extracts	Acceptable	Acceptable	Equivalent

ISO Maximization SC & SO extracts	Acceptable	Acceptable	Equivalent
Sterility	Yes per ANSI/AAMI/ISO 11137-2: Sterilization of health care products -- Radiation	Yes per ANSI/AAMI/ISO 11137-2: Sterilization of health care products -- Radiation	Equivalent
Electrical Safety	N/A	N/A	Equivalent

7. Summary of difference from the predicate

The MR11 ventricular catheter differs from that of the cannula in four principal ways.

- a) The fluid delivery tubing is composed of PEEK rather than the silica of the predicate
- b) The fluid delivery tubing is shorter (23 inch nominal) versus 4 and 10 feet of the predicate.
- c) The MR11 ventricular cannula provides for a bone anchor, for subsequent infusions over a period of 24 hours.
- d) The rigid component of the catheter can be removed after insertion so that the catheter is flexible along its whole length. The cannula's rigid section is integrated into the assembly.

8. Testing:

Testing to applicable standards has been completed with acceptable outcomes. The following testing has been performed:

- Sterilization and Shelf Life
- Biocompatibility: Both the ventricular catheter and the bone anchor were tested under conditions of Good Laboratory Practices in the following tests with acceptable results.

Study	Result	Conclusion
Cytotoxicity Study, MEM Elution	"The test article extract showed no evidence of causing cell lysis or toxicity. The test article extract met the requirements of the test since the grade was less than a grade of 2 (mild reactivity)."	Non-cytotoxic
ISO Systemic Toxicity, SC and SO Extracts	"There was no mortality or evidence of systemic toxicity from the extracts. The test article extracts met the requirements of the study."	Non-toxic

Pyrogen Study – Material Mediated	“The total rise of rabbit temperatures during the 3 hour observation period was within acceptable USP limits. The test article was judged as nonpyrogenic.”	Non-pyrogenic
ISO Intracutaneous Study, SC and SO Extracts	“The test article met the requirements of the test since the difference between each test extract overall mean score and corresponding control overall mean score was 0.0 and 0.2 (0.4 bone anchor) for the SC and SO test extracts, respectively.”	No significant erythema or edema
ISO Maximization Sensitization Study, SC and SO Extracts	“The test article extracts showed no evidence of causing delayed dermal contact sensitization in the guinea pig. The test article was not considered a sensitizer in the guinea pig maximization test.”	Non-sensitizer
Drug/Device Compatibility Testing	HPLC using USP Monograph assay for cytarabine concentration	No leaching of materials or uptake of cytarabine

- Performance Testing – Bench

Test	Methodology	Conclusions
Flow Rate Testing	Set flow rate on infusion pump to specification rate. Measure amount of fluid collected over a set period of time. Measure pressure for reference.	Specified flow rates of 0.5 mL/hr and 0.3 mL/hr were achieved with the infusion pump.
Flow Rate Testing, constant pressure	Set flow rate on infusion pump so that measured pressure is 0.7 psi. Measure amount of fluid collected over a set period of time.	Flow rates at 0.7 psi for the device were greater than the predicate as expected. .008" ID Predicate 0.6 mL/hr .008" ID Device 4.0 mL/hr .021" ID Predicate 34 mL/hr .021" ID Device 56 mL/hr

510(k) (Traditional) Submission
Section 5, 510(k) Summary

High Pressure Flow Testing	Use hand held syringe to apply very high pressure via hand force. Measure amount of fluid collected over a set period of time.	High pressure flow is equivalent to or better than the predicate. The main difference comes from the shorter tubing length of the device and different applied pressure since hand force was used.
Aspiration	Use hand held syringe and apply vacuum to fixed amount of fluid. Record time to aspirate fluid.	Results were better than (higher aspiration rates achievable) the predicate. The main difference comes from the shorter tubing length of the device.
Leak/Burst Testing	Block the cannula tip. Use a hand held syringe or manual pump to apply 70 psi minimum pressure to the Catheter. Watch for leaks anywhere along the length.	Results were equivalent to the predicate device. The 70 psi specification is the same. All samples withstood 70 psi internal pressure without any leaks.
Distal Tip Compressive Strength	Apply axial compressive load to distal tip on 50A durometer material. Examine tip under magnification to check for damage.	Results were equivalent to the predicate, which has the same specification. All samples met the requirement without damage to the tip.
Distal Tip Lateral Load	Apply side load to distal tip section to a minimum force limit.	Results were substantially equivalent to the predicate, which had the same specification. All samples met the specification.
Ceramic Stylet Lateral Load	Hold Catheter with Stylet approximately 6 cm from the distal end (5 cm from the applied load). Apply lateral load and measure force to breakage.	Device met the specification, which is the same specification as the predicate. Predicate device has a fixed integrated rigid tube, not a removable internal stylet.
Catheter Axial Tensile Load	Apply a tensile load to the catheter. Pull to failure. Catheter is held at the outer capillary and at the PEEK tubing. Catheter was then held at the outer capillary and the Pebax tubing.	Results were equivalent to the predicate. The specification is the same as the predicate. All devices met the specification.
Catheter Insertion into Bone Anchor	Insert the device into the bone anchor. Observe for any interference or difficulty.	The device was inserted into the Bone Anchor without difficulty.
Catheter Retention in Bone Anchor	Apply pull force to the specified limit to determine if	All samples were retained to the specification.

510(k) (Traditional) Submission
Section 5, 510(k) Summary

	the Catheter remains in place in the Bone Anchor	
Stylet Removal from Bone Anchor	Once the catheter is tightened in the Bone Anchor, measure the pull force required to remove the Stylet from the Catheter.	The device met the specification of maximum pull force for Stylet removal. The predicate device has an integrated ceramic tube and does not have removable Stylet.
Catheter Tissue Insertion	Pass the device through bovine brain tissue to a 4 inch depth. Extract the device and check for plugging or damage.	Equivalent to the predicate device. There was no tissue in or on the device after insertion and removal.
Backflow	Insert the device into brain-simulating gel. Use an infusion pump at a set flow rate to drive dyed fluid through the device. Examine the gel for any dye that runs back up the sides of the device.	Equivalent to predicate device. There was no observed backflow.
Tip Deflection at Max Pressure	Plug a catheter tip. Mark the position of the tip. Apply 70 psi internal pressure to the catheter. Mark the position of the tip with pressure applied. Determine the difference (if any) between the marks.	Equivalent to the predicate device. There was no movement of the tip at applied internal pressure of 70 psi.
Tip Deflection at Max Aspiration	Plug a catheter tip. Mark the position of the tip. Apply 2 psi vacuum to the catheter. Mark the position of the tip with pressure applied. Determine the difference (if any) between the marks.	Equivalent to the predicate device. There was no movement of the tip at applied vacuum of 2.0 psi.
Bend Radius Withstand	Wrap the device tubing for ½ turn around a specified diameter rod. Examine the tubing for kinks or other damage.	Equivalent to the predicate device. There was no kinking or damage to the tubing at the specified bend radius.
Syringe Dropping Test	Connect the device to a syringe. Drop the syringe with the device anchored in place. After the syringe drop, check the device for damage and operation.	Equivalent to the predicate device. There was no damage sustained by the device after the syringe was dropped with the device connected.

510(k) (Traditional) Submission
Section 5, 510(k) Summary

Tip Deflection after Stylet Removal	Insert Catheter with Stylet into brain simulating gel in a clear beaker. Mark the spot on the beaker where the tip touches. Lock the catheter in place. Remove the Stylet and observe for tip deflection away from the mark.	There was no tip movement upon Stylet removal. Device met specification.
Bone Anchor Insertion Force	Insert Bone Anchor into drilled hole and measure the amount of downward force required from beginning to end of insertion.	All Bone Anchors were all inserted fully with a downward force below the requirement.
Bone Anchor Insertion Torque	Insert Bone Anchor into drilled hole and measure the amount of torque required from beginning to end of insertion.	All Bone Anchors were all inserted fully with a torque below the requirement.
Bone Anchor Side Load Force	Insert a Bone Anchor into a drilled hole. Apply a side load to the minimum specified level to four points around the bone anchor diameter. Inspect for damage or yielding.	All Bone Anchors all withstood the minimum load with no damage or yielding.
Bone Anchor Side Impact Force	Insert a Bone Anchor into a drilled hole. Apply a lateral impact load on four different points around the bone anchor diameter by dropping a known weight from a specified height to obtain an applied energy. Inspect for damage or yielding.	All Bone Anchors withstood the minimum energy without damage or yielding.
Bone Anchor Retention Force	Insert a Bone Anchor into a drilled hole. Tighten a catheter into the Bone Anchor and apply a tensile load to the catheter to the minimum specified limit.	All Bone Anchors remained in place under the applied tensile load.
Bone Anchor-Driver Detachment Force	Insert a Bone Anchor into a drilled hole using the Driver. Measure the tensile force required to remove the Driver from the Bone Anchor.	All pairs of Bone Anchors and Drivers were below the maximum limit for detachment force.

Bench testing included design verification testing and comparison testing with the predicate MRII Cannula, demonstrating the MR Compatible Ventricular Catheter functions as intended and is substantially equivalent to the legally marketed device.

9. Rx or OTC:

The MR Compatible Ventricular Catheter is an Rx prescription device per 21 CFR Part 801, Subpart D.

10. Substantial Equivalence:

The MR Compatible Ventricular Catheter is as safe and effective as the predicate MRII Cannula. The MR Compatible Ventricular Catheter has the same intended uses and identical indications, technological characteristics, and principles of operation as its predicate device. The minor technological differences between the MR Compatible Ventricular Catheter and its predicate devices raise no new issues of safety or effectiveness. Performance data demonstrate that the MR Compatible Ventricular Catheter is as safe and effective as the MRII Cannula. Thus the MR Compatible Ventricular Catheter is substantially equivalent.

MRII, Inc.



E. F. Waddell
Director of Regulatory Affairs

Date: 8/16/13



August 16, 2013

MRI Intervention Inc.
Edward Waddell
5 Musick
Irvine, CA 92618

Re: K123605
Trade/Device Name: SmartFlow Flex Ventricular Catheter
Regulation Number: 21 CFR 882.4100
Regulation Name: Catheter, Ventricular
Regulatory Class: Class II
Product Code: HCA
Dated: July 16, 2013
Received: July 17, 2013

Dear Mr. Waddell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Joyce M. Whang -S

for Victor Krauthamer, Ph.D.
Acting Director
Division of Neurological and Physical
Medicine Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number: **K123605**

Device Name: MRI Intervention, Inc. SmartFlow™ Flex Ventricular Catheter

Indications For Use:

The SmartFlow™ Flex Ventricular Catheter is intended for injection of Cytarabine or removal of CSF from the ventricles during intracranial procedures. The device is not intended for implant. This device is intended for "single patient use only."

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

<p>Joyce M. Whang -S</p> <hr/> <p>(Division Sign Off) Division of Neurological and Physical Medicine Devices (DNPMD)</p> <p>510(k) Number <u>K123605</u></p>
