

K123843

510(k) Summary  
Page 1 of 7  
15-Oct-2012

JUN 5 2013



**510(k) Summary**

**Prepared By:** Strohl Medical Technologies, Inc.  
126 Prospect Street  
Norwell, MA 02061

**Telephone:** (888) 374-8877

**FAX:** (617) 249-2120

**Contact Person:** Heather Strohl

**Date Summary prepared:** October 15, 2012

**Name of the Device:** NeuroEPG System

**Common Name:** Evoked Response System, Nerve Stimulator/Monitor

**Classification Name:** Electrical Evoked Response Stimulator (per CFR882.1870)  
Auditory Evoked Response Stimulator (per CFR 882.1900)

**Patient Population:** 18 years of age and older

**Environment of Use:** Hospital, clinic, EEG/EP technologist's, surgeon's, or  
physician's office, operating room

**Contraindications:** None

**Predicate Device(s):** Intelligent Hearing System *SmartEP* (K070608)

**Device Description:**

The NeuroEPG System is an evoked response testing and diagnostic device that is capable of eliciting, acquiring, and measuring auditory and somatosensory evoked potentials.

The NeuroEPG System is intended to objectively record evoked responses from patients 18 years of age and older, upon the presentation of sensory stimuli. The product is indicated for use as a

**510(k) Summary**

Page 2 of 7

15-Oct-2012

diagnostic aid and adjunctive tool in sensory related disorders (i.e., auditory, somatosensory) and in surgical procedures for intraoperative monitoring.

**Indications for Use:**

The NeuroEPG System is an evoked response testing and diagnostic device that is capable of eliciting, acquiring and measuring auditory and somatosensory evoked potentials.

The intended use of the NeuroEPG System is to objectively record evoked responses from patients 18 years of age and older, upon the presentation of sensory stimuli. The product is indicated for use as a diagnostic aid and adjunctive tool in sensory related disorders (i.e., auditory, somatosensory) and in surgical procedures for intraoperative monitoring.

The NeuroEPG System is intended to be used by trained personnel in a hospital, clinic, EEG/EP technologist's, surgeon's, or physician's office, operating room, or other appropriate setting.

The anatomical sites of contact for auditory evoked potential (AEP) testing are the patient's ear canal (with the contact object being a sound delivery ear bud), and the patient's head and possibly other body sites (with the contact object being electrodes that are capable of measuring bio-potentials).

The anatomical sites of contact for somatosensory evoked potential (SEP) testing are the patient's upper/lower limbs (with the contact object being skin-surface electrodes) and the patient's head and possibly other body sites (with the contact object being electrodes that are capable of measuring bio-potentials).

**Table of Comparison and Differences vs. Predicates**

The modifications described in this table are only concerned with the Auditory Evoked Potentials (AEP) and the Somatosensory Evoked Potential (SEP) testing aspects of the predicate device.

510(k) Summary  
Page 3 of 7  
15-Oct-2012

Parameter	Predicate Device <i>SmartEP (K070608)</i>	<i>The NeuroEPG System</i>
Intended Use	<p><u>SEP:</u> Stimulate, record, and process somatosensory evoked potentials</p> <p><u>AEP:</u> Stimulate, record, and process auditory evoked potentials</p>	Same
Indications for Use	The recording and analysis of physiological data necessary for the diagnosis of somatosensory and auditory related disorders.	Same
Target Population	All Ages	18 years and older
Design	External box housing circuitry connected to CPU via a USB connection	Same
Sterility	Non-sterile	Same
Anatomical Sites	<p><u>SEP:</u> Upper/lower limbs and Head</p> <p><u>AEP:</u> Head</p>	Same
Energy Delivery	<p><u>SEP:</u> Stimulation of upper or lower limbs with surface electrical signals</p> <p><u>AEP:</u> Stimulation of ears with auditory stimulus</p>	Same
Where Used	Clinical Setting	Same
Safety	Conforms to IEC 60601-1	Same

510(k) Summary  
Page 4 of 7  
15-Oct-2012

Parameter	Predicate Device <i>SmartEP (K070608)</i>	<i>The NeuroEPG System</i>
Patient Isolation	Type BF (IEC 60601-1) Fiber Optic Signal Link	Same
<b>Somatosensory Stimuli</b>		
Types	Constant Current or Voltage	Constant Current
Mode	<u>User selectable:</u> Single, Dual, or Stimulus Train	Same
Shape	Mono or Biphasic pulses	Same
Repetition Rate	0.1 - 100 Hz	Same
Phase/Polarity	Positive or Negative	Same
Duration /Pulse Width	10 - 1000 $\mu$ s	Same
Stimulus Intensity Levels	<u>Current:</u> 0 - 100 mA <u>Voltage:</u> 0 - 400V (Continuous adjustable level with user selectable maximum range into a 4000 Ohms load)	<u>Current:</u> 0 – 25 mA <u>Voltage:</u> 0 – 50V (Continuous adjustable level with user selectable maximum range into a 2000 Ohms load)
<b>Auditory Stimuli</b>		
Types	Clicks, Pure Tones, Multifrequency Stimuli	Same
Duration	25-5000 $\mu$ s	Same
Envelopes	Linear, Blackman, Gaussian, Hanning, Rectangular, Triangular, Trapezoidal, Exact Blackman, Cosine, Cosine Squared, Cosine Cubed	Same
Intensity	0-125 dB SPL	Same

510(k) Summary  
Page 5 of 7  
15-Oct-2012

Parameter	Predicate Device <i>SmartEP (K070608)</i>	<i>The NeuroEPG System</i>
Repetition Rate	1-100 Hz	Same
Test Frequencies	500-16,000 Hz	Same
Presentation	Monaural or Binaural	Same
Masking	White Noise Programmable	Same
Transducers	Insert Earphones, Bone Vibrator, Headphones, Sound Field, Ear Probe	Ear buds only
<b>Measurement Parameters</b>		
Analysis Window	Variable (up to 2.5 msec)	Same
Artifact Rejection Threshold	User Selectable	25 $\mu$ V
Measured Values	Response Level ( $\mu$ V) Noise Level ( $\mu$ V) Signal to Noise Ratio (dB) Response Latency (msec) Frequency (Hz)	Same plus left versus right comparison
<b>Computer Requirements</b>		
Computer Type	Personal Computer	Same
Operating System	Microsoft Windows 2000 or XP	Same
Interface Connection	USB (Universal Serial Bus)	Same

**Technological Characteristics:**

The NeuroEPG System hardware is based on modifications to the predicate device. Modifications were made by Intelligent Hearing Systems, the OEM of the predicate.

**510(k) Summary**

Page 6 of 7

15-Oct-2012

The NeuroEPG hardware is very similar in electronics design to the predicate device hardware, except that the electronics hardware has been repackaged into a single, smaller, stand-alone unit. Unlike the predicate devices, there is no separate pre-amplifier unit connected to the patient in the NeuroEPG System. Instead, all of the pre-amplifier electronics have been embedded inside of the NeuroEPG Data Acquisition Unit. The same patient isolation methods are used.

The software changes within the data acquisition unit are minor and do not alter the safety and effectiveness of the device in any way. The software incorporates a reduced range of testing and analysis parameters so a broader set of personnel can perform an evoked response test. Strohl Medical has additional software that operates on the data from the data acquisition unit.

**Safety and Effectiveness:**

The NeuroEPG System utilizes the same design principles, circuit designs, and operating principles as are used in the predicate device.

The NeuroEPG System meets the requirements of IEC 60601-1 and IEC 60601-1-2 for electrical safety and electromagnetic compatibility respectively. It also meets the requirement of IEC 60601-2-40 for electromyographs and evoked response equipment. The disposable (PIK) meets the requirements of ANSI/AAMI EC12 for Disposable ECG Electrodes and ANSI/AAMI EC53 for ECG Cables and Leadwires as well as clause 56.3c of IEC 60601-1 (1998) [21CFR898].

**Substantial Equivalence**

The stimulation and response mechanism are identical to the predicate. The fundamental scientific technology of the NeuroEPG System is identical to the predicate.

The differences that exist between the devices are insignificant in the terms of safety or effectiveness.

**Indications** – Equivalent to the predicate

**Technology** – The technology is identical

**Operating specifications** – Equivalent

**Materials** – The patient contact materials are identical to those in other legally marketed devices.

**Environment of Use** – Identical

**Patient Population** – Restricted to adults

**Differences:**

There are no significant differences between the proposed device and the predicate device.

**510(k) Summary**

Page 7 of 7

15-Oct-2012

**Comparative Performance and Specifications**

The NeuroEPG System is substantially equivalent to the SmartEP device marketed by Intelligent Hearing systems with FDA 510(k) clearance number K070608.



June 5, 2013

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-002

Strohl Medical Technologies  
c/o Mr. Mark Job  
Regulatory Technology Services LLC  
1394 25<sup>th</sup> Street, NW  
Buffalo, MN 55313

Re: K123843  
Trade/Device Name: NeuroEPG System  
Regulation Number: 21 CFR 882.1870  
Regulation Name: Evoked Response Electrical Stimulator  
Regulatory Class: Class II  
Product Code: GWF, GWJ, GXY  
Dated: May 14, 2013  
Received: May 15, 2013

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Joyce M. Whang -S**

for Victor Krauthamer, Ph.D.

---

Acting Director

Division of Neurological

and Physical Medicine Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K123843

Device Name: NeuroEPG System

### Indications For Use:

The NeuroEPG System is an evoked response testing and diagnostic device that is capable of eliciting, acquiring and measuring auditory and somatosensory evoked potentials.

The intended use of the NeuroEPG System is to objectively record evoked responses from patients 18 years of age and older, upon the presentation of sensory stimuli. The product is indicated for use as a diagnostic aid and adjunctive tool in sensory related disorders (i.e., auditory, somatosensory) and in surgical procedures for intraoperative monitoring.

The NeuroEPG System is intended to be used by trained personnel in a hospital, clinic, EEG/EP technologist's, surgeon's, or physician's office, operating room, or other appropriate setting.

The anatomical sites of contact for auditory evoked potential (AEP) testing are the patient's ear canal (with the contact object being a sound delivery ear bud), and the patient's head and possibly other body sites (with the contact object being electrodes that are capable of measuring bio-potentials).

The anatomical sites of contact for somatosensory evoked potential (SEP) testing are the patient's upper/lower limbs (with the contact object being skin-surface electrodes) and the patient's head and possibly other body sites (with the contact object being electrodes that are capable of measuring bio-potentials).

Prescription Use X AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21-CFR-801-Subpart-D) (21 CFR 801-Subpart C)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Joyce M. Whang -S**

(Division Sign Off)  
Division of Neurological and Physical Medicine  
Devices (DNPMD)

510(k) Number K123843