

K123912
Page 1 of 5

Kobold Applicator Kits

510(k) Summary

This summary of 510(k) safety and effectiveness information is furnished in accordance with requirements detailed in 21 CFR 807.92.

MAY 14 2013

1.
The assigned 510(k) number is K-TBD

Submitter's Identification:
Spencer Fillmore
2410 S Sumner Lane
Greenacres, WA, 99016

Correspondence:

Christina Bernstein
BB Medical Surgical, Inc.
2670 Leavenworth Street
San Francisco, CA 94133
Tel: 415-450-0515
Email tina@bbmedicalsurgical.com

Date of submission 07 December 2012

2.
Device name:
1. Fletcher type Applicator Set
Proprietary name: Kobold Fletcher Type GYN T&O Applicator Set™

2. Henschke type Applicator Set
Proprietary name: Henschke Type GYN T&O Applicator Set™

- A. Regulation Section 892.5700
- B. Classification: Class II
- C. Product Code: JAQ
- D. Panel: Radiology

3.
Intended Use:

1. The CT/MRI Compatible HDR Fletcher Type Applicator is indicated for use in any case where high dose rate (HDR) radiation treatment of the cervix and uterus is accepted by up-to-date clinical guidelines.

2. The CT/MRI Compatible HDR Henschke Type Applicator is indicated for use in any case where high dose rate (HDR) radiation treatment of the cervix and uterus is accepted by up-to-date clinical guidelines.

4.

Device Description:

1. The CT/MRI Compatible HDR Fletcher Type Applicator Set™ design is based on the conventional modified Fletcher GYN applicator with a click fit connector interface. This applicator is particularly useful in providing an HDR radiation boost to the cervix and parametrical tissues in conjunction with external beam radiation therapy. The set consists of precision manufactured intrauterine tandems and a pair of interlocking colpostats with Fletcher type unshielded ovoids ranging in size from 1cm to 3cm. All components are designed to allow for customization of the treated geometry. The distance between the ovoids is adjustable with a lock knob. Four different tandems are provided: one straight, one with a slight curvature of 15 degrees, one with medium curvature of 30 degrees, and one with a maximal curvature of 45 degrees.

2. The CT/MRI Compatible HDR Henschke Type Applicator Set™ design is based on the conventional modified Henschke GYN applicator with a click fit connector interface. This applicator is particularly useful in providing an HDR radiation boost to the cervix and parametrical tissues in conjunction with external beam radiation therapy. The set consists of precision manufactured intrauterine tandems and a pair of interlocking colpostats with Henschke type unshielded ovoids ranging in size from 2cm to 3cm. All components are designed to allow for customization of the treated geometry. The distance between the ovoids is adjustable with a lock knob. Four different tandems are provided: one straight, one with a slight curvature of 15 degrees, one with medium curvature of 30 degrees, and one with a maximal curvature of 45 degrees.

5.

Substantial Equivalence Information:

Information presented supports substantial equivalence of the Kobold Fletcher and Henschke Type GYN T&O Applicator Set™ to the predicate device. Each of the two proposed devices has the same indications for use, has similar nature of body contact, is similar in shape and design, has the same fundamental technology and uses the same sterilization procedures. Please see predicate device brochure in Appendix C.

- A. Predicate device names: Mick® Ring Type GYN Applicator Set
- B. Predicate device K numbers: 011657
- C. Comparison with predicate:

Kobold Applicator Kits

SUBSTANTIAL EQUIVALENCE TABLE	KOBOLD LLC	KOBOLD LLC	MICK RADIO-NUCLEAR INSTRUMENTS, INC.®
K-Number	TBD	TBD	K-011657
Device Description	Applicator Set, Fletcher Type	Applicator Set, Henschke Type	Applicator Set, Ring Type
Indications for Use	The CT/MRI Compatible HDR Fletcher Type Applicator is indicated for use in any case where high dose rate (HDR) radiation treatment of the cervix and uterus is accepted by up-to-date clinical guidelines.	The CT/MRI Compatible HDR Henschke Type Applicator is indicated for use in any case where high dose rate (HDR) radiation treatment of the cervix and uterus is accepted by up-to-date clinical guidelines.	The Mick Radio-Nuclear Instruments, Inc. HDR Tandem/Ring Applicator with Rectal Retractor is intended for use in Brachytherapy. The delivery of intra-cavitary radiation therapy requires not only proper visualization and localization of the applicator within the treatment volume, but precise dosimetry and a stable delivery system from which treatment can be administered. The Mick Radio-Nuclear HDR Tandem/Ring Applicator with Rectal Retractor meets these requirements.
Afterloader Compatibility	GammaMed®. VariSource® with click-fit connector.	GammaMed®. VariSource® with click-fit connector.	GammaMed®, VariSource® with click-fit connector.
Design Geometry	Flexible	Flexible	Flexible
Materials	Titanium, Ultem®, Stainless Steel, Silicone, High Temperature Vinyl	Titanium, Ultem®, Stainless Steel, Silicone, High Temperature Vinyl	Titanium, Stainless Steel, Polysulphone, Silicone
Packing	Tray	Tray	Cassette
Sterility	Provided Non-Sterile	Provided Non-Sterile	Provided Non-Sterile
Sterilization Methods, Sterilization Conditions.	Steam 18 min at 273°F (134°C), 44psi (3 bar) Dry 15 min; no Sterrad® compatibility documented	Steam 18 min at 273°F (134°C), 44psi (3 bar) Dry 15 min; no Sterrad® compatibility documented	Steam 20 min at 132°C Dry 15 min; no Sterrad® compatibility documented

Kobold Applicator Kits

Biocompatibility	Documented	Documented	Not Documented
Anatomical Sites	Uterus, Cervix	Uterus, Cervix	Cervix
Environmental Compatibility	CT/MRI Compatible	CT/MRI Compatible	CT/MRI Compatible

6.

Test Principle, Performance Characteristics:

FDA has not established special controls or performance standards for this device.

7.

Bench Top Testing

Extensive testing in accordance with known standards is documented by the manufacturer. The standard is a specification for the materials used in surgical implant applications, predicate brachytherapy devices and manufacturer's acceptance procedures are detailed in Appendix A. Mechanical testing of the finished devices is additionally described in Appendix A.

8.

Conclusions

Both versions of Kobold's Applicator Set are similar in intended use and technological characteristics to predicate devices reviewed. Each of the two versions of the device is similar with respect to indications for use and physical characteristics to predicate devices in terms of section 510(k) substantial equivalency.

Contraindications: As per clinical guidelines and standard clinical practice.

Warnings and Precautions: The precautions are provided in the device labeling for the Kobold Fletcher Style GYN T&O Applicator Set™ and Kobold Henschke Style GYN T&O Applicator Set™. There is no warning associated with this type of device.

9.

Summary

Description	Comparison with Predicate Device
Biocompatibility	Safe as Predicate Device

Kobold Applicator Kits

Performance Characteristics	Substantially equivalent
Intended Use	Substantially equivalent
Performance Tests	Not Required

Each version of the device 1. Kobold Fletcher Type GYN T&O Applicator Set™ and 2. Kobold Henschke Style GYN T&O Applicator Set™, based on the information submitted in this 510(k) application has been demonstrated to be substantially equivalent to the predicate devices.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

Kobold, LLC
% Ms. Christina Bernstein
Regulatory Director
2670 Leavenworth Street
SAN FRANCISCO CA 94133

May 14, 2013

Re: K123912

Trade/Device Name: Fletcher Type GYN T&O Applicator Set™
Regulation Number: 21 CFR 892.5700
Regulation Name: Remote Controlled Radionuclide Applicator System
Regulatory Class: Class II
Product Code: JAQ
Dated: April 19, 2013
Received: May 2, 2013

Dear Ms. Bernstein:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for

Janine M. Morris
Director, Division of Radiological Health
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K123912

Device Name: Fletcher Type GYN T&O Applicator Set™

Indications for Use:

The CT/MRI Compatible HDR Fletcher Type Applicator is indicated for use in any case where high dose rate (HDR) radiation treatment of the cervix and uterus is accepted by up-to-date clinical guidelines.

Prescription Use X

AND/OR


Over-The-Counter Use

(Part 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of *In Vitro* Diagnostics and Radiological Health (OIR)



(Division Sign Off)

Division of Radiological Health

Office of *In Vitro* Diagnostic and Radiological Health

510(k) **K123912**

Page 1 of 1