
Nit-Occlud® PDA

Instructions for Use
For USA only!

www.pfmmedical.com

Explanation of the symbols on label and packaging.

	Reference number
	Lot number
	Serial number
	Read Instructions for Use carefully
	Protect from direct sunlight
	Store in a dry place
	Expiry Date
	For single use only
	Do not resterilise
	Sterilized by Ethylene Oxide
	Do not use if packaging is damaged
	"MR conditional" = safe use of MR diagnostics under certain conditions
	Manufacturer
	Does not contain rubber latex components
	Diethylhexylphthalate (DEHP) free
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Instructions for Use

Nit-Occlud® PDA

Product Description

Nit-Occlud® PDA is a system for transcatheter occlusion of Patent Ductus Arteriosus (PDA) with spiral coils. The system consists of the following parts:

- **Nit-Occlud® PDA**




Figure 1 Nit-Occlud PDA

The spiral coil is mounted in a straightened fashion on a flexible delivery system including a disposable handle to which it is connected by means of a patented detachment mechanism.

Nit-Occlud® PDA coils are available as Flexible and Medium type.

The flexible and medium types are pre-loaded into the transportation sheath. For insertion the transportation sheath must be connected to the implantation catheter.

The Nit-Occlud® PDA coil has a cone in cone configuration which results from the fact that the proximal windings of the coil are wound in the reverse direction (see Figure 2.).




Figure 2 Nit-Occlud PDA spiral coils
(D=Distal diameter, P=Proximal diameter, Lc=Length configured).

- **Implantation catheter**

The implantation catheter is equipped with a marker ring at its distal tip for better orientation during fluoroscopy.

Indication for Use

The Nit-Occlud® PDA coil is a permanently implanted prosthesis indicated for percutaneous, transcatheter closure of small to moderate size patent ductus arteriosus with a minimum angiographic diameter less than 4 mm.

Contraindications

Medical conditions that exclude implantation of a Nit-Occlud® PDA coil include:

- Endocarditis, endarteritis or active infection at the time of the implantation
- Patients with a body weight < 5 kg
- Pulmonary hypertension (calculated PVR greater than 5 Wood Units)
- Thrombus in a blood vessel through which access to the PDA must be obtained
- Thrombus in the vicinity of the implantation site.

General Information and Warnings

These instructions for use and the information on the packaging should be read carefully before each use. The PDA coil system should be used only by physicians trained in interventional occlusion techniques.

WARNING:

- ⚠ **Do not use the product if the packaging has been opened, or is damaged, if you are not sure that it is sterile, or if the expiry date has passed.**
- ⚠ **Each product is packed separately, and is delivered in an EO-sterilised and non-pyrogenic condition. It is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization of single-use devices may result in degraded performance or a loss of functionality. Reuse of single-use devices may result in exposure to pathogens such as viruses, bacteria, fungi, or prions.**
- ⚠ **The product must be stored in dry conditions. Do not expose the packaged products to direct sunlight.**
- ⚠ **Retrieval devices should be available during implant procedures for interventional retrieval of the coil if required.**
- ⚠ **Care must be taken not to damage the coil or to dislodge it from the delivery system while unpacking or inserting it into the implantation catheter.**
- ⚠ **Since the delivery system has ferromagnetic properties, implantation must not be carried out in an MR environment.**
- ⚠ **The coil should not be removed from the delivery system. It should not be used with another delivery system since this may alter characteristics of configuration and detachability.**
- ⚠ **A detached coil should not be remounted on the core wire of the delivery system.**
- ⚠ **The configured coil should not be pulled through heart valves or ventricular chambers.**
- ⚠ **The implantation catheter is not suitable for application of contrast medium. It must not be connected to high pressure injectors.**
- ⚠ **The Nit-Occlud® PDA coil consists of a nickel-titanium alloy, which is generally considered safe. In non-clinical testing, nickel has been**

shown to be released from the device in very small amounts. Patients who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to **seek medical assistance** immediately if they suspect they are experiencing an allergic reaction. **Symptoms may include** difficulty in breathing or **swelling** of the face or throat. **While data are currently limited, it is possible that** some patients may develop an allergy to nickel if this device is implanted.

- ⌚ This product contains chemicals known to the State of California to cause cancer, birth defects, or reproductive harm.
- ⌚ **NOTE: Federal Law (USA) restricts this device to use by a physician.**

Product Identification

Each product label has peel-off labels, to allow the product to be identified precisely. These can be used for the patient file and the patient ID card.

MRI Compatibility:

The Nit-Occlud® PDA coil was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Non-clinical testing demonstrated that the Nit-Occlud® PDA coil is MR conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less
- The maximum whole-body averaged specific absorption rate (SAR) shall be limited to 2.0 W/kg (normal operating mode only) for 15 minutes of scanning.

MRI-Related Heating

- In non-clinical testing, the Nit-Occlud® PDA coil produced the following temperature rise during MRI performed for 15-min in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: **Highest temperature change +1.6°C.**

Therefore, the MRI-related heating experiments for the Nit-Occlud® PDA coil at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

Artifact Information

- MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Nit-Occlud® PDA coil. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	369 mm ²	118 mm ²	647 mm ²	739 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Potential Adverse Events

- Air embolism
- Allergic reaction to drug/contrast
- Apnea
- Arrhythmia requiring medical treatment or pacing
- Arteriovenous Fistula
- Bacterial Endocarditis
- Blood loss requiring transfusion
- Chest Pain
- Damage to the tricuspid or pulmonary valves
- Death
- Embolization of the occluder, requiring percutaneous or surgical intervention
- Endarteritis
- False aneurysm of the femoral artery
- Fever
- Headache/migraine
- Heart failure
- Hemolysis after implantation of the occluder
- Hypertension
- Hypotension or shock
- Infection
- Myocardial infarction
- Occluder fracture or damage
- Perforation of the heart or blood vessels
- Stenosis of the left pulmonary artery or descending thoracic aorta
- Stroke/TIA
- Thromboembolism (cerebral or pulmonary)
- Valvular Regurgitation
- Vessel damage at the site of groin puncture (loss of pulse, hematoma etc.).

Precautionary Measures

- An angiogram must be performed prior to implantation for measuring length and diameter of the PDA.
- The implantation catheter must be flushed with heparinized saline solution prior to introduction and during the procedure, especially after angiography.
- The pfm medical implantation catheter is specifically designed for the delivery system. Other catheters should not be used to implant the device.

- Contrast media should not be injected through the implantation catheter.
- The coil should not be pulled back into the implantation catheter using strong force.
- Administration of 50 units of heparin per kg body weight is recommended after femoral sheaths are placed.
- Antibiotic coverage before (1 dose) and after implantation (2 doses) is recommended in order to prevent infection during the implant procedure. Antibiotic prophylaxis should be performed to prevent infective endocarditis during first 6 months after coil implantation.
- A suitable lateral aortogram should be performed for measurement of PDA dimensions (see Fig. PDA Measurements):

Directions for Use




Figure 3 PDA Measurements

D1 = Narrowest diameter

D2 = Aortic Ampulla diameter

L3 = PDA length

D = Distal diameter

P = Proximal diameter

Lc = Length configured

When defining the parameters of your measurements, please consider that the anatomy of the PDA may differ among patients. According to Krichenko et al. there are five different types of PDAs (see device performance by PDA anatomy in clinical study section below).




Figure 4 Types of PDAs According to: Krichenko A. Et al. 1989, AmJ Cardiol; 63:877-880

Coil Selection

According to the measurements, the ductus type and the following recommendations, an appropriate coil should be selected:

- The distal coil diameter D should be no more than 2 mm larger than D2.
- The distal coil diameter D should be at least 3 to 4 mm larger than D1.
- Length of the configured coil Lc (see product label) should be not longer than L3.

D1	D2	Device	D1	D2	Device
----	----	--------	----	----	--------

1mm	≤ 3mm	4x4	3mm	≤ 7mm	7x6
1mm	4mm	5x4	3mm	8-9mm	9x6
1mm	≥ 5mm	6x5	3mm	≥ 9mm	9x6 or 11x6
2mm	≤ 5mm	6x5	< 4mm	9mm	11x6
2mm	6-7mm	7x6	< 4mm	10-11mm	11x6
2mm	≥ 8mm	9x6	< 4 mm	≥ 12mm	11x6

Table 1: Selection of Nit-Occlud PDA coil (according to angiographic PDA dimensions)

Coil Implantation Sequence


Step: 01

- Unpack the Nit-Occlud® PDA consisting of coil with disposable handle and implantation catheter under sterile conditions.

WARNING:

- ⚠ **Do not pull on the delivery system. If the coil is withdrawn into the Y connector, there is a danger that the system can no longer be loaded. Check all screw connections. Some screw joints may have been loosened by the sterilization process.**
- Flush the system carefully through the side access of the Y connector with heparinized saline solution, and ensure that there is no air remaining anywhere in the system.
- Check the coil position inside the transparent transportation sheath. The coil should be inside the sheath. When it is in this position, it is essential not to pull on the delivery system. If the coil is not positioned inside the transportation sheath, or shows visible signs of damage, it must be replaced with a new coil.

Step: 02 Coil Introduction



MPA: Main Pulmonary artery

Figure 5 Coil Introduction

- Using a soft guide wire, advance the implantation catheter from the right femoral vein through the right heart, across the PDA into the descending thoracic aorta.

- Remove the guide wire from the implantation catheter and flush the catheter with a heparinized saline solution.
- Attach the luer lock connector of the transportation sheath to the implantation catheter.
- Open the hemostatic valve of the Y connector. The coil is now free for advancement.

WARNING:

- ☞ **Do not pull on the delivery system in this position!**
- Advance the coil into the implantation catheter.

Step: 03 Coil Configuration




Figure 6 Coil Configuration

- Under fluoroscopic control advance the coil carefully through the implantation catheter into the aorta. This is done by moving the delivery system forward while arresting the implantation catheter.

WARNING:

- ☞ **Ensure that the pigtail aortography catheter does not become entangled with the loops of the coil.**
- Advance the coil until the first marker M1 is positioned close to the Y connector. At this position, all but one loop is configured outside the implantation catheter.

Step: 04 Coil Positioning




Figure 7 Coil Positioning

- Retract the entire system (implantation catheter, delivery system) under fluoroscopic control until the configured coil is positioned in the ampulla of the ductus (close hemostatic valve of Y connector or fix implantation catheter against delivery system).

NOTE: For longer ductus types, coil configuration inside the ductus ampulla is recommended. Here, 2–3 windings of the coil must first be configured in the aorta. Then the entire system is pulled into the ductus ampulla for further configuration of the coil.

Step: 05 Final Coil Adjustment




Figure 8 Final Coil Adjustment

- Open the hemostatic valve of the Y connector.
- Configure the last 1 or 2 loops on the pulmonary side of the ductus by simultaneously pulling back the implantation catheter (with your left hand) and pushing the delivery system (with your right hand). Advance the delivery system until the second marker M2 is close to the Y connector. At this position the coil is outside the catheter.
- Perform an aortogram to confirm that the coil is in the correct position.

NOTE: If the position or size of the coil is not satisfactory, it should be repositioned or exchanged at this point.

Repositioning

- To reposition the coil, pull it back into the implantation catheter by pulling the delivery system.

WARNING:

- ☞ **Close the gap between delivery system and coil before you retrieve the coil into the implantation catheter.**

To do so, hold the handle with one hand and move the delivery system gently forward while holding it between 2 fingers of the same hand. This movement closes the gap between coil and delivery system and should be done under fluoroscopic control. Once the gap is closed the coil can be pulled smoothly into the implantation catheter.

WARNING:

- ☞ **If a strong resistance is encountered while pulling the delivery system into the catheter, do not pull the system very hard because you risk a premature coil release of the coil.**

- To reposition the implantation catheter, the coil should be pulled back into the transportation sheath, carefully and under visual control, until the tip of the coil is in line with the marker at the distal end of transportation sheath. Fix the coil position by closing the Y connector.

WARNING:

- ☞ **If the coil is pulled back too far, there is the risk that it may not be possible to reload it into the delivery system.**

- Then flush the implantation catheter with heparinized saline solution and repeat the procedure from Step 02.

Step: 06 Coil Release




Figure 9 Coil Release

- When the coil is properly positioned, it should be released. The rotation screw should lie directly against the pusher ball. If there is any gap between the two, it must be closed.

WARNING:

⚠ **Final release should only be performed if the coil is properly positioned in the PDA. Otherwise, the coil must either be retrieved and repositioned or replaced by an appropriate substitute. Before the coil is finally released, proper position of the coil should be confirmed by angiography.**

- Remove the safety clip from the handle.
- Turn the rotation screw under fluoroscopy clockwise until the coil is released.
Note that depending on the coil type between 8-15 rotations are needed to release the coil. You will feel an increase in resistance immediately prior to release.
- Remove the delivery system and implantation catheter.
- Perform a final aortogram about 10 minutes later to document position of the coil and PDA occlusion.
- Remove the aortography catheter.

WARNING:

⚠ **Ensure that the catheter does not touch the coil**

Technical Complications and how to avoid them

Complications may be avoided or ameliorated by the following:

- Use of venous access for implants.
- Flushing all components with heparinized saline.
- Keeping delivery system and catheter straight, avoiding loops and curves on the catheterization table.
- Coil recapture in the pulmonary artery or aorta, avoiding pulling the exposed coil across heart valves or through the right ventricle.

Failure of detachment:

Complications may arise if the coil is not released successfully. “Sticking” may occur if positioning of the coil is very time-consuming and/or if the delivery system does not protrude far enough from the end of the implantation catheter. The distal part of the delivery system must be placed outside of the implantation catheter immediately before release. If, however, the coil “sticks”, the device must be retrieved and exchanged. As done with all interventional instruments, prior to implantation of the exchanged device the catheter should be flushed thoroughly to prevent coagulation, thus avoiding elevated friction or “sticking” of the system.

Device Embolization/ Premature release

The coil may embolize into the pulmonary artery if the aortic cone is too small or the coil fit is too loose. This can be prevented by accurate measurement of the PDA dimensions and choice of an appropriate coil for the PDA. Correct calibration of the angiographic measurement is a very important factor.

In case of coil embolization, interventional retrieval should be performed using a snare or a biptome.

In the event that the coil embolizes and interventional retrieval is unsuccessful, surgical retrieval should be considered.

Protrusion/ Obstruction

Protrusion of the windings into the aorta and/ or into the pulmonary artery may cause blood flow disturbances or vessel stenosis. This is avoidable by choosing an appropriate coil with a configured length (Lc) equal to or less than the PDA length (L3). The recoil force of the coil tends to return the device to its original configuration whenever possible and will retract the windings into the ampulla and/ or against the vessel wall. Pulmonary artery protrusion may be avoided by correct coil positioning during implant.

Late complications

Delayed complications such as migration or protrusion with a significant blood flow disturbance may require surgical removal of the coil.

Clinical Studies

Study description

A prospective, non-randomized, multi-center, single-arm Study and a continuing access study were performed using the same protocols at 15 centers in the United States of America to assess the safety and effectiveness of the Flex and Medium Nit-Occlud® PDA coil for occlusion of Patent Ductus Arteriosus (PDA) with minimum angiographic diameter of less than 4 mm. The primary effectiveness endpoints were echocardiographic and clinical closure rates at 12 months. The primary safety endpoint was the serious adverse event rate at 12 months. The endpoint rates were compared to an Objective Performance Criteria as follows:

- Echocardiographic closure (absence of detectable residual PDA flow on echocardiogram) greater than 85% at 12 months
- Clinical closure (absence of heart murmur) greater than 95% at 12 months
- Serious adverse event rate of less than 1% at 12 months

The following criteria were considered for their inclusion:

Inclusion criteria	Exclusion criteria
<ul style="list-style-type: none"> • PDA with 4 mm or smaller minimum diameter by color Doppler 	<ul style="list-style-type: none"> • Associated cardiac anomalies requiring surgery

<ul style="list-style-type: none"> • Patent weight \geq 5 Kg, Age 6 months to 21 years (Patients older than 21 years may have device implanted and be included in a study registry.) • Previous treatment by surgery or Nit-Occlud device with residual PDA noted at least 6 months after the procedure 	<ul style="list-style-type: none"> • Known bleeding or blood clotting disorders • Ongoing febrile illness • Pregnancy • Pulmonary hypertension/increased pulmonary vascular resistance (>5 Wood Units) • Known hypersensitivity to contrast medium
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Table 2 Inclusion/Exclusion Criteria

Study results

A total of 378 patients were enrolled and 357 patients were evaluated for safety and effectiveness. The patient's mean age was 4.26 years (range 0.5 to 21.9 years); the mean weight was 18.1 kg (range 4.7 to 109.0 kg), a total of 68.1% of the enrolled patients were female. Of the 357 evaluable patients, 347 had successful implantation of the device (technical success).

Principal safety and effectiveness results are presented in Table 3 below:

	OPC Rates	Nit-Occlud Patients	Percent	95% Lower Bound	95% Upper Bound
Technical Success at Implantation	95% ²	347/357	97.2%	95.6%	
Clinical Closure at 12 Month Follow-up	95% ¹	308/314	98.1%	96.7%	
Echocardiographic Closure at 12 Month Follow-Up	85% ¹	299/309	96.8%	95.0%	
Mortality at 12 Months	0% ¹	0	0.0%		0.95%
Serious Adverse Events at 12 Months	1% ¹	0	0%		0.95%
Total Device and Procedure Related Adverse Events at 12 Months	6%	15/316*	4.7%		7.21%
		14/316**	4.4%		6.84%
Composite Success at 12 Months	80% ³	294/309	95.1%	93.0%	

¹ Objective Performance Criteria (OPC) specified by the Multiorganization Advisory Panel to (FDA) Appendix (XII)

² Inferred from technical success rate of Gianturco coil technical success cited in Multiorganization Advisory Panel to FDA report (Appendix XII)

³ Defined in IDE protocol but not defined by the Multiorganization Advisory Panel report

* Numerator is number of events; denominator is number with 12 mos fu + 2 with AE before 12 months

** Numerator is number of person; denominator is number with 12 mos fu + 2 with AE before 12 months

Table 3 Principal Safety and Effectiveness Results

Refer to Table 4 below for procedural and fluoroscopy times by device size and type.

Catalog #	Device Size Distal x Proximal Diameter	Device Type	Number of Implants	Mean Procedure Duration [min.]	Median Procedure Duration [min.]	Mean Fluoroscopy Time [min.]	Median Fluoroscopy Time [min.]
145044	4 x 4 mm	Flex	38	68.6	66.0	17.2	14.0
145054	5 x 4 mm	Flex	27	77.8	72.0	19.6	17.0
145065	6 x 5 mm	Flex	57	91.5	82.0	19.8	18.5
145076	7 x 6 mm	Medium	110	83.3	73.5	17.0	15.0
145096	9 x 6 mm	Medium	97	92.0	79.0	18.8	16.0
145116	11 x 6 mm	Medium	26	93.0	85.0	25.5	23.5

Table 4 Procedure and Fluoroscopy Times by Nit-Occlud Device

Differing Technical Failure Rates were observed based on Angiographic Classification of the PDA on the lateral aortogram and are summarized in the Table 5 below.

Classification	N(% of Total)	Technical Failure Rate	n/N (%)
Conical (A)	267 (74.8%)	4/267	(1.5%)
Short (B)	17 (4.8%)	3/17	(17.6%)
Tubular (C)	5 (1.4%)	1/5	(20%)
Complex (D)	18 (5.0%)	1/18	(11.1%)
Elongated (E)	50 (14.0%)	1/50	(2%)
TOTAL	357 (100%)	10/357	(2.8%)

Table 5 Technical Failure rate by Angiographic Classification (See Figure above)

Study Adverse events were defined as follows:

Serious Adverse Events:

- Procedural or device related events which were life-threatening, required surgery to correct, resulted in hospitalization or prolonged hospital stay, caused long-term disability, or resulted in genetic damage or birth defect.

Major Adverse Events:

- Procedural or device related events which were not life-threatening, required interventional (catheter based) and /or medical treatment to correct up to one year follow-up evaluation but were resolved without surgical intervention.

Minor Adverse Events:

- Procedural or device related events which were not life-threatening, and were resolved without intervention or with a brief specific non-surgical intervention up to one year follow-up evaluation.

The combined studies safety results were the following:

- Mortality at 12 months: 0.0% (0/314)
- Serious Adverse Events at 12 months (device related): 0.0% (0/314)
- Serious Adverse Events at 12 months (procedure related): 0.0% (0/314)
- Total AEs (Serious, Major, and Minor) at 12 months or last follow up (related to the procedure or the device): 4.7% (15/316^{*})

The 15 Adverse Events are further described in Table 6 below:

DSMB Adjudication	Category	No of Events
Major Device Related	Device embolization	2
	Device Retrieval/Removal	2
	Obstruction of descending aorta	1
Minor Device Related	Possible Thrombus	1
Major Procedure Related	Decreased Pulse in Right Foot	1
	Reaction to anesthesia	2
Minor Procedure Related	Reaction to anesthesia	1
	Vascular access site complication	1
	Other Adverse Event	2
	Nausea	1
	Fever	1

Table 6 Adverse Events

Procedural success, effectiveness and safety results were comparable to or better than predefined objective performance criteria.*

Disposal after Use

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance

* Patients with 12 month follow up and those with an adverse event at any time

* Multiorganization Advisory Panel to FDA for Pediatric Cardiovascular Devices. Proposed Standards for Clinical Evaluation of Patent Ductus Arteriosus Occlusion Devices. Catheter Cardiovasc Interv 2000; 51:293-296.t.

with recognised medical procedure, and in compliance with the relevant legal regulations and local ordinances.

Warranty

pfm medical warrants that this medical device is free from defects in both materials and workmanship. The above warranties are in lieu of all other warranties, either expressed or implied, including any warranty of merchantability or fitness for a particular purpose. Suitability for use of the medical device for any surgical procedure shall be determined by the user. pfm medical shall not be liable for incidental or consequential damages of any kind.



pfmmedical

pfm medical ag
Wankelstr. 60, 50996 Köln, Germany
T +49 (0)2236 9641-0, F -20
www.pfmmedical.com
service@pfmmedical.com

GA034/Rev00-130801

Nit-Occlud[®] PDA

A Patient's Guide to Non-Surgical Closure of
the Patent Ductus Arteriosus (PDA)

Using the Nit-Occlud[®] PDA

Germany:

pfm medical ag
Wankelstrasse 60
50996 Cologne GERMANY
Tel: +49 (0)2236 9641-10
www.pfmmedical.com

pfmmedical
Quality and Experience

USA/Canada:

pfm medical, Inc.
1815 Aston Avenue Suite 106
Carlsbad, California 92008
Tel: + 001 760 758 8749
Fax: +001 760 758 1167

English Version

This brochure is intended to provide you with general information to discuss with your doctor. It is not intended to provide medical care or treatment. You should consult with your doctor regarding the diagnosis or treatment of your medical condition.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Warranty

pfm medical warrants that this medical device is free from defects in both materials and workmanship. The above warranties are in lieu of all other warranties, either expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Suitability for use of the medical device for any surgical procedure shall be determined by the user. pfm medical shall not be liable for incidental or consequential damages of any kind.

Inside This Book	Page
About "PDA" ("Patent Ductus Arteriosus")	5
All About Your Heart	5
Why We Have a PDA	6
Why You Need to Visit the Hospital: Closing the PDA	7
Purpose of the Nit-Occlud® PDA Device	7
Description of the Nit-Occlud® PDA Device	9
When the Device Should Not Be Used (Contraindications).....	10
Risks and Benefits of the Procedure	11
What Are the Risks?	11
Other Possible Risks	12
What Are the Benefits	13
Other Ways to Handle a PDA (Alternatives to the Procedure).....	13
What to Expect Before, During, and After the Procedure	13
Before the Procedure	13
During the Procedure	14
After the Procedure: REST!	17
Patient ID ("Identification") Card	18
Follow-up Visits With Your Doctor	19
When to Call the Doctor.....	19
Additional Information and Help	20
Questions or Concerns?	21

List of Pictures (Figures)

Page

Picture 1: A healthy heart.....	5
Picture 2: Heart with a "PDA" ("Patent Ductus Arteriosus").....	6
Picture 3: Nit-Occlud® PDA catheter tip and coil system for closing ("occluding") the PDA.....	9
Picture 4: Nit-Occlud® PDA system: coil, catheter, delivery system, and handle	9
Picture 5: How the catheter moves through the body	15
Picture 6: The catheter feeds through the heart and PDA into the aorta	15
Picture 7: Bigger picture of Picture 6. The catheter feeds through the heart and PDA into the aorta.....	16
Picture 8: The coil feeds through the catheter	16
Picture 9: The coil is pulled back and fitted to the PDA.....	16
Picture 10: The coil leaves the catheter toward the PDA	17
Picture 11: The PDA is closed	17


About "PDA" ("Patent Ductus Arteriosus")

Your heart has a small opening that allows blood to flow in a different way than a normal heart. It is called a "patent ductus arteriosus," or "PDA." This book will help you understand PDA, your procedure, and how to stay healthy afterward. **Always feel free to ask your doctor any questions!**

All About Your Heart

To better understand PDA, it is helpful to understand the way blood moves—or circulates—through and around your body. Your heart is a pump that pumps blood throughout your body. This allows your body—including all the organs, muscles, and tissues inside—to get the oxygen it needs to function.

There are two chambers on the right side and two on the left. The two chambers in the top of your heart are called the left atrium and the right atrium (**Picture 1**). The atriums receive and collect blood. The two chambers in the bottom are called "ventricles" (**Picture 1**).



Picture 1: A healthy heart.

The right ventricle pumps blood to your lungs to get fresh oxygen that your body needs to live. Inside your body, blood is carried by vessels. Vessels include veins and arteries. Some of the larger vessels that enter and leave your heart are the aorta, pulmonary


artery, and pulmonary vein. To get oxygen from the lungs, the right side of the heart pumps blood to the lungs by the pulmonary artery. It is the only artery in the body that carries blood without oxygen.

First blood is in the right side of the heart. Then it is pumped through the lungs, where it receives oxygen. After this, blood moves to the left side of your heart. There, it is pumped to your body by the left ventricle. The artery carrying blood to your body is called the aorta.

Questions? Ask your doctor!

Why We Have a PDA

During development in the womb, blood flows from your pulmonary artery to your aorta through the ductus arteriosus. (Picture 2)



Picture 2: Heart with a PDA ("Patent Ductus Arteriosus").

At birth the ductus arteriosus closes, allowing the pulmonary arteries to get oxygen from the lungs and the aorta to carry it to the rest of the body. If the ductus arteriosus remains open or "patent" after birth, it's called a "Patent Ductus Arteriosus," or "P-D-A."

No one knows why the ductus arteriosus doesn't close in some people. Being born before complete development (premature) may be one reason. PDA seems to be more

common in girls than boys. It also may be more common in kids with other heart problems called congenital heart defects.

The PDA can sometimes be detected (with a doctor's stethoscope) by the presence of a heart murmur. Other times doctors find out about a PDA after doing other tests on the heart.

Questions? Ask your doctor!

Why You Need to Visit the Hospital: Closing the PDA

Some people with PDAs that are not closed feel very tired. Others have trouble breathing or are at increased risk for infection. Still others don't grow and mature naturally. Some people have no signs of PDA at all.

If the PDA is not closed, blood can move in abnormal ways called "shunts." A shunt can carry blood from the left side of your heart to the right side through the PDA. This can overload the arteries supplying blood to your lungs, leading to "high blood pressure" in and around your lungs. In the body, high blood pressure in the main artery to the lungs—the pulmonary artery—can be very harmful. It can also lead to other problems such as inflammation or infection of the pulmonary artery or movement ("shunting") of blood from the right side of your heart to the left. Blood with too little oxygen from the right side of your heart can reach your aorta through the PDA and can cause your body to receive blood without enough oxygen.

Questions? Ask your doctor!

Purpose of the Nit-Occlud® PDA Device

To close (or occlude) your PDA, the doctor will place (or implant) a small coil in the duct connecting your pulmonary artery and aorta called the "Nit-Occlud® PDA." During the procedure, the Nit-Occlud® PDA will be passed into the PDA near your heart by a plastic tube


called a "catheter." To implant the device and close your PDA, the doctor will make a small opening in an artery and vein in your groin. Then the catheter can be passed toward your heart, and the coil moves within the catheter. Once the coil is released, it closes or occludes the PDA. The Nit-Occlud[®] PDA will stay there after your procedure and your body will grow around it. The coil will keep the PDA closed.

This kind of procedure is called "percutaneous," which means "beneath the skin" (through an artery and vein).


Questions? Ask your doctor!

Description of the Nit-Occlud® PDA Device

The Nit-Occlud PDA device is a small spring-like "coil" (**Picture 3**) shaped like a cone. The coil is made of nickel and titanium ("Nitinol").



Picture 3: Nit-Occlud® PDA catheter tip and coil system for closing ("occluding") the PDA.



Picture 4: Nit-Occlud® PDA system: coil, catheter, delivery system, and handle.

The coil can be pulled out of shape to pass through the catheter into your PDA. In this way, the metal has "memory." When it is pushed out of the catheter, it remembers its original coil shape and returns to that shape. The coil is designed to fit snugly inside the PDA. The Nit-Occlud® PDA coil is used to close small or medium PDAs—openings of

about 4 millimeters across ("in diameter") or less.

The Nit-Occlud[®] PDA includes the catheter and coil. A gold ring is fitted to the tip of the catheter so that the doctor can see it on X-rays. The coil is attached to a delivery system. The coil is pushed through the catheter using the delivery system. Once the coil has been placed in the PDA, it can be released from the delivery system using the handle (**Picture 4**).

Questions? Ask your doctor!



When the Device Should Not Be Used (Contraindications)

The Nit-Occlud[®] PDA should **NOT** be implanted if you have:

- An infection.

- A blood clot near where the Nit-Occlud[®] PDA coil will be placed. A clot is a clumping of blood cells called a "thrombus." A blood clot can block or "occlude" one of your veins or arteries and make the procedure difficult to perform (among other problems).

- A body weight of less than 11 pounds.

- High blood pressure in the pulmonary artery and smaller arteries in the lung.

- A very large PDA. Your doctor should not use the device to close PDAs that are more than 4 millimeters across (in "diameter").

- A blood clot in an artery that the Nit-Occlud[®] PDA coil needs to pass through.

- **Questions? Ask your doctor!**

Risks and Benefits of the Procedure

What Are the Risks?

All medical procedures involve risks. In the closing of PDAs using the Nit-Occlud[®] PDA, some risks are due to the procedure, and others are due to the device.

In a clinical study, a total of 357 patients at 15 hospitals in the United States received the Nit-Occlud[®] PDA to close their PDAs. These patients were ages 6 months to 21 years. The patients had PDAs with diameters of 4 millimeters or less, and the patients also weighed at least 11 pounds.

In the study, the following side effects ("adverse events") were seen. Most of these had no lasting effects.

In 3 patients out of 357 receiving the Nit-Occlud[®] PDA:

- The coil moved away from the PDA and entered a vein or artery. Patients experienced nausea and vomiting.

In 2 patients out of 357 receiving the Nit-Occlud[®] PDA:

- A blood clot formed in an artery.

In 1 patient each out of 357 receiving the Nit-Occlud[®] PDA:

- The device got trapped in a blood vessel.
- The aorta became more narrow than usual.
- A blood clot formed on the coil.
- A rash formed after tape was placed on a patient.
- The patient had trouble breathing.
- The patient got a fever (high body temperature).
- The patient had a seizure.

Other Possible Risks

Some harmful effects of the procedure not seen in the study mentioned above that can happen with this type of procedure include the following. (They may occur in less than 1 to 9 patients out of every 1,000):

- Abnormal connection between blood vessels
- Air bubbles in the blood, leading to blood vessels closing.
- Allergic drug reaction
- Allergic dye reaction
- Bruises in the region where the catheter was inserted.
- Changes in how the heart beats— arrhythmia.
- Chest pain
- Collection of blood outside the vessel (hematoma)
- Death
- Harming one of the valves in the heart.
- Harming "red blood cells," which carry oxygen. This is called "hemolysis."
- Harming the heart or vessels where the catheter was placed
- Headache/migraine
- Heart attack
- Heavy bleeding
- High or low blood pressure.
- Infection or inflammation of blood vessels or lining of the heart and valves.
- Piercing of a vessel or heart (perforation)
- Temporary absence of breathing (apnea)
- Temporary lack of oxygen to the brain (Stroke/ Transient Ischemic Attack)
- Weakening heart ("cardiac insufficiency").

What Are the Benefits

The main benefit of having the Nit-Occlud[®] PDA is that you most likely won't need open-heart surgery. Your "percutaneous" procedure is less "invasive" than surgery. This usually means:

- A shorter procedure: usually 1 to 2 hours instead of 2 or more hours for surgery.
- A shorter stay in the hospital: usually 1 to 3 days instead of 4 to 10 days for surgery.
- No scar on your chest. Very little scarring or no scarring at all (in your groin).

Other Ways to Handle a PDA (Alternatives to the Procedure)

- Surgical Closure of PDA. The chest is opened surgically. This allows the surgeon to see the whole heart and large vessels. The surgeon sews up the PDA using stitches to keep the PDA closed permanently. Sometimes this kind of surgery is needed when the PDA is large.
- Medication to help close the PDA or treat symptoms of the PDA
- No treatment: a small PDA may close on its own.

What to Expect Before, During, and After the Procedure

Before the Procedure

Before the procedure, the doctor will look at your heart with a test called an "echo" (echocardiogram), where pictures of your heart and vessels will be taken. The doctor will also perform an "ECG" ("electrocardiogram"). For this test, small pads (circles) made of plastic with metal endings ("leads") are placed on or around your chest and other parts of your body. The ECG helps your doctor make sure that your heart is beating in the right way.

The doctor will also do a "physical examination" and blood test to make sure you don't have an infection.

Just before your procedure, you will be given anesthesia through an "IV" ("intravenous") line. You may also have a "local" anesthetic rubbed on your thigh or groin area where the doctor will insert the catheters for your procedure. This is to make the area numb so you don't feel discomfort during the procedure. Some patients feel temporary stinging or burning from the anesthetic.


The procedure will take place in the "cath lab" ("catheterization laboratory") where hundreds or even thousands of procedures like yours are performed every year. The doctor will perform an "angiogram", which is a special X-ray that allows doctors to see inside your blood vessels. Doing this special test allows your doctor to get a clear picture of your PDA.

During the Procedure

The procedure itself takes 1 or 2 hours. During the procedure, cameras will move above your chest area to help guide the doctors. They will also use an ECG to monitor your heart. You may feel some discomfort in your groin area when the catheters go in. Some people feel extra heart beats from the catheter when it is inserted.

Step 1.

The doctor puts the two small plastic tubes (catheters) into a vein and artery in your groin. Then the doctor moves one of the catheters up your body toward your heart and vessels (pulmonary artery and aorta; **Picture 5**). In the other catheter, the doctor gives contrast medium. This liquid helps the doctor to see your PDA on an angiogram.




Picture 5: How the catheter moves through the body.


Step 2:

The Nit-Occlud[®] PDA device stays inside the catheter. The doctor passes the catheter (with the device inside) through an artery and vein. As he or she does, the coil gets closer and closer to your PDA.

First the coil pushes into the chambers on the right side of the heart until it reaches the pulmonary artery. The doctor uses the special X-ray (angiogram) to make sure that everything is in the correct place. When sure, the doctor feeds the catheter through the PDA and into the aorta (**Pictures 6 and 7**).




Picture 6: The catheter feeds through the heart and PDA into the aorta.



Picture 7: Bigger picture of Picture 6. The catheter feeds through the heart and PDA into the aorta.

Step 3:


Next the doctor pushes the Nit-Occlud® PDA coil out through the catheter. Once out of the catheter, the coil returns to its spiral shape (**Picture 8**).



Picture 8: Coil fed through catheter.


Step 4:

The coil is still attached to the catheter. Next, the doctor places the coil into the PDA (**Picture 9**). The coil is designed to fit snugly inside the PDA.




Picture 9: The coil is pulled back and fitted to the PDA.

The doctor releases the coil (**Picture 10**). X-rays (angiograms) guide the doctor. The coil is now inside the PDA.



Picture 10: The coil leaves the catheter toward the PDA.

Now, or soon after, the doctor takes out all catheters from your heart and body. The Nit-Occlud® PDA coil stays near your heart. The PDA is now closed (**Picture 11**). Over time, your body will grow around and over it.



Picture 11: The PDA is closed.



After the Procedure: REST!!

The doctor will take out the catheters. There may be two small marks in your groin where the catheters were inserted, and you may have some soreness or tenderness in your groin area. These will go away and you should have no scars or very small ones.

A special large adhesive bandage will be placed on your groin called a "pressure bandage." This helps to prevent bleeding at the place where the catheters were inserted. You should be able to go home (or be discharged) after 1 to 3 days in the hospital. Before discharge, a chest X-ray and an echocardiogram may be done to confirm that your device is in the correct place and your heart is OK. If the device is working, it will stay in your body for the rest of your life.

The doctor will go over your medicines and activities with you and your family or caregivers. You may need to take medication so you do not get an infection that can be serious. Your doctor will let you know if you need to take them and, if so, for how long.

During the first 4 to 6 weeks after returning home, strenuous activity such as heavy lifting should be avoided, Activities such as jumping, running, lifting, biking, or gymnastics should also be avoided. After 4 to 6 weeks, most patients can return to normal activity.

Questions? Ask your doctor!

Patient ID ("Identification") Card

As a patient with a Nit-Occlud[®] PDA device, it is important to carry a Patient ID Card to identify you as having an implanted device. Show it whenever you see a doctor or dentist for treatment. The Patient ID Card has information about safety including safe ways to do MRIs in patients with a Nit-Occlud[®] PDA. You or your mother, father, guardian, or caregiver should register the MRI information in your Patient ID Card with MedicAlert. MedicAlert provides important information to healthcare workers. This helps emergency responders, hospital staff, and others make informed decisions about your care.

Contact:

Mail: MedicAlert Foundation International

2323 Colorado Avenue

Turlock, CA 95382

Phone: (800) 432 5378 (toll-free)

Emergency: (800) 625 3780 (toll-free)

Fax: (800) 863 3429 (toll-free) or (209) 669 2495

Web: www.medicalert.org

Follow-up Visits With Your Doctor

You will probably go back and see your doctor a few times in the first year after getting your Nit-Occlud[®] PDA. It is important for you and your family or caregivers to make and keep appointments with your doctor. The doctor may wish to do an echocardiogram to make sure everything is OK.

When to Call the Doctor

If you feel any of the feelings below, you, your family, or caregiver should contact your family doctor right away:

- Pain, numbness, or a feeling of coldness or weakness in your legs or feet.
- Any pain in your back, chest (or rib cage), stomach, or groin region.
- A stronger and faster heartbeat than usual ("palpitations").
- Trouble breathing ("breathlessness").
- Fainting ("passing out"/"syncope") or dizziness.

Additional Information and Help

You, your family, or caregiver can find useful information on the Internet. Feel free to visit the following websites.

- The Adult Congenital Heart Association provides information and support for adults and older kids ("adolescents"). They provide information about heart problems that people are born with ("congenital heart disease") to patients, their families, and healthcare professionals. Click on: <http://www.achaheart.org/>
- The Congenital Heart Information Network is a non-profit group that provides support and information on congenital heart disease. They work to make people more aware of these problems. Click on: <http://tchin.org/>
- Little Hearts Inc. is a non-profit group that provides support services for families affected by heart disease. Click on: <http://www.littlehearts.org/>
- The Children's Heart Society is a non-profit group including parents and extended families of children with heart disease. Click on: <http://www.childrensheart.org/>

(The respective operators are solely responsible for the contents of the websites listed above.)

For More Information, Contact pfm medical:



Germany:

pfm medical ag
Wankelstrasse 60
50996 Cologne GERMANY
Tel: +49 (0)2236 9641-10
www.pfmmedical.com

USA/Canada:

pfm medical, Inc.
1815 Aston Avenue Suite 106
Carlsbad, California 92008
Tel: 760 758 8749
Fax: 760 758 1167