

K130266

*Intuitive Surgical, Inc.*

**Traditional 510(k)**  
***EndoWrist® One Vessel Sealer***

**510(k) Summary**  
**[As Required by 21 CFR 807.92(c)]**

**510(k) Owner:** Intuitive Surgical, Inc.  
1266 Kifer Road  
Sunnyvale, CA 94086

**Official Contact:** Crystal Ong  
Sr. Regulatory Affairs Specialist  
408-523-8636 (phone)  
408-523-8907 (fax)  
[crystal.ong@intusurg.com](mailto:crystal.ong@intusurg.com)

AUG 29 2013

**Date Summary Prepared:** February 1, 2013

**Trade Name:** *EndoWrist® One Vessel Sealer*

**Common Name:** System, Surgical, Computer Controlled Instrument

**Product Code:** NAY

**Classification:** Endoscope and Accessories, 21 CFR 876.1500

**Predicate Devices:**

- Intuitive Surgical's *EndoWrist® One Vessel Sealer*(K110639)
- Ethicon Endo-Surgery's EnSeal PTC (K071728)

**Device Description:**

The *EndoWrist One Vessel Sealer* is an electro-surgical sealing and cutting instrument to be used in conjunction with the *da Vinci Si Surgical System* and the ERBE VIO 300 D electro-surgical generator. The outside diameter (O.D.) of the instrument shaft is 8.5 mm and the working length is 38 cm. The instrument is provided sterile and is a single-use, disposable device. The distal end has jaws with bipolar electrodes for sealing vessels and contains a cutting blade that extends through the jaws to transect sealed vessels and other tissues. The sealing and cutting functions are actuated by the *da Vinci Si Surgical System* foot pedals.

The ERBE VIO 300 D generator provides the high frequency (radio frequency) electrical current for tissue sealing.

**Intended Use:**

The *EndoWrist® One Vessel Sealer* is a bipolar electro-surgical instrument for use with the *da Vinci Si* Surgical System and the ERBE VIO 300 D electro-surgical generator. It is intended for **grasping and blunt dissection of tissue and for bipolar coagulation and mechanical transection of vessels up to 7 mm in diameter and tissue bundles that fit in the jaws of the instrument.** The *EndoWrist One Vessel Sealer* has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures, and should not be used for these procedures.

**Technological Characteristics:**

The subject *EndoWrist One Vessel Sealer* is identical to the predicate *EndoWrist One Vessel Sealer* device in terms of technological characteristics. The proposed indications for use are a subset of the indications cleared for the Ethicon-EndoSurgery predicate device.

**Performance Data:**

Design validation and animal testing demonstrates that the subject device is substantially equivalent to the predicate devices and that the design output meets the design input requirements. The differences do not raise different questions of safety or effectiveness as compared to the predicate devices.

**Summary:**

Based on the intended use, indications for use, technological characteristics, and performance data, the Intuitive Surgical *EndoWrist One Vessel Sealer* is substantially equivalent to the predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-002

August 29, 2013

Intuitive Surgical, Inc.  
Crystal Ong  
Sr. Regulatory Affairs Specialist  
1266 Kifer Road  
Sunnyvale, California 94086

Re: K130266  
Trade/Device Name: EndoWrist One Vessel Sealer  
Regulation Number: 21 CFR 876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: Class II  
Product Code: NAY  
Dated: July 16, 2013  
Received: July 17, 2013

Dear Ms. Ong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

FOR **Peter D. Rumm -S**

Mark N. Melkerson  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number if known: K130266

Device Name: *EndoWrist® One Vessel Sealer*

**INDICATIONS FOR USE:**

The *EndoWrist® One Vessel Sealer* is a bipolar electro-surgical instrument for use with the *da Vinci Si* Surgical System and the ERBE VIO 300 D electro-surgical generator. It is intended for grasping and blunt dissection of tissue and for bipolar coagulation and mechanical transection of vessels up to 7 mm in diameter and tissue bundles that fit in the jaws of the instrument. The *EndoWrist One Vessel Sealer* has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures, and should not be used for these procedures.

Prescription Use  X  AND/OR Over-the-Counter Use \_\_\_\_\_

(Per 21 CFR 801 Subpart D) (Per 21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Joshua C. Nipper -S**

(Division Sign-off)

Division of Surgical Devices

510(k) Number  K130266