



JAN 13 2014

## 510(k) Summary

As required by 21 CFR 807.92, this "510(k) Summary" provides a basis for the substantial equivalence determination of the device listed below.

### General Information

**Date Prepared:** February 26, 2013  
**Submitter Name:** Xeridium Medical Devices  
[formerly Manufacturing and Research, Inc.]  
**Submitter Address:** 4700 S. Overland Dr.  
Tucson, Az 85714  
**Contact Person:** Michelle Lott, RAC  
**Phone Number:** (520) 882-7794  
**Fax Number:** (520) 882-6894

### Device Information

**Trade Name:** Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube  
**Common Name:** Gastrostomy Tube Kit (low profile)  
**Classification Name:** Tubes, gastrointestinal (and accessories)  
[21 CFR 876.5980, Product Code KNT]  
**Classification Panel:** Gastroenterology/Urology

### Predicate Device Information

**Trade Name:** EndoVive™ Low Profile Replacement Gastrostomy Tube (K083685)  
**Common Name:** Gastrostomy Tube Kit (low profile)  
**Classification Name:** Tubes, gastrointestinal (and accessories)  
[21 CFR 876.5980, Product Code KNT]  
**Classification Panel:** Gastroenterology/Urology


### Device Description

The Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube is constructed primarily of silicone. The low profile device consists of a hub which is overmolded onto a previously extruded bi-lumen shaft. At the proximal end, a silicone balloon is bonded to the shaft using a silicone RTV. The tip of the shaft is filled with radiopaque barium RTV. In use, the balloon end of the catheter is inserted through a gastrointestinal stoma tract and into the stomach lumen. The balloon is then inflated with 5mL sterile or distilled water to secure the device in place. Accessories for the Low Profile device include two extension set configurations, an inflation syringe, a catheter tip syringe, gauzes and a stoma measuring device to facilitate device placement and use.

---

Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube

Xeridium  
Special 510(k) Application  
Appendix B  
Page 1 of 3



### Intended Use/Indications for Use

Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube is indicated in adult and pediatric populations for use in percutaneous placement of an enteral feeding tube for feeding and/or administration of medication in conjunction with an established gastrointestinal stoma tract. The replacement tube may also be used for gastric decompression.

### Technological Characteristics

The table below provides a general summary of the technological characteristics of the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube compared to the predicate device.

Technological Characteristics		
	Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube (Submitted device)	EndoVive™ Low Profile Replacement Gastrostomy Tube (Predicate Device) K083685
Product Range (FR size)	12 - 24 FR	12 - 28 FR
Product Range (Length)	0.8 – 5.0 cm	1.0 – 5.5 cm
Product Range (Balloon Size)	5 cc for all FR sizes	4 cc for 12 FR device 6 cc for 14 – 28 FR devices
Hub	Same as K083685	Silicone with 2 ports: one for feeding and the other for balloon inflation
Inflation valve	Same as K083685	Thermoplastic housing with silicone valve
Balloon	Same materials as K083685, different balloon profile	Silicone balloon
Reflux (feeding) valve	Thermoplastic housing with silicone valve – different profile as K083685	Thermoplastic housing with silicone valve
Shaft	Silicone 65 duro	Silicone 60 duro
Barium stripe	Barium sulfate	Barium sulfate with blue colorant
Kit Contents	Inflation syringe, feeding syringe, gauze, Bolus Extension Set, Continuous Extension Set	Inflation syringe, feeding syringe, gauze, lubricating jelly, Bolus Feeding Set, Continuous Feeding Set
Packaging System/Materials	Same as K083685	Tyvek Pouch, Polybag, IFU, Corrugate Case & Shipper

Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube

Xeridien  
Special 510(k) Application  
Appendix B  
Page 2 of 3

### Performance Characteristics

The table below provides a general summary of the performance testing that was conducted for the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube.

### Test Method Conclusion

Test	Method	Conclusion
-Flow rate -Balloon Burst Volume -Balloon Volume Maintenance -Balloon Concentricity -Balloon Over Shaft Size -Balloon Integrity -Balloon Integrity in simulated gastric fluid	All in accordance to ASTM F2528-06 (Procedures A through G)	Candidate device met all performance requirements. Where acceptance criteria were not defined in the standard, the candidate device performed equivalent to or better than the predicate device.

In all instances the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube performed as intended.

### Substantial Equivalence – Comparison to Predicate Device

Among the information and summary tables presented in the 510(k) submission to support substantial equivalence of the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube to the legally marketed predicate device are: 1.) device description, 2.) indications for use, 3.) device comparison tables, 4.) material information, 5.) bench test results, and 6.) labeling. In particular, bench testing demonstrates that there was no difference in the performance, safety, or effectiveness between the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube and the EndoVive™ Low Profile Replacement Gastrostomy Tube.

The Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube is substantially equivalent to the EndoVive™ Low Profile Replacement Gastrostomy Tube (K083685).

The subject device and predicate device have the same intended use and similar indications, technological characteristics and principles of operation. The minor differences in technological characteristics have been identified and do not present any new issues of safety or effectiveness. Thus, the Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube is substantially equivalent to the EndoVive™ Low Profile Replacement Gastrostomy Tube (K083685).

---

**Entuit™ Thrive LP Balloon Retention Gastrostomy Feeding Tube**

Xeridiem  
Special 510(k) Application  
Appendix B  
Page 3 of 3



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

January 13, 2014

Xeridiam Medical Devices  
Michelle Lott  
RA/QA Director, RAC  
4700 South Overland Drive  
Tucson, AZ 85714

Re: K130674  
Trade/Device Name: Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube  
Regulation Number: 21 CFR§ 876.5980  
Regulation Name: Gastrointestinal tube and accessories  
Regulatory Class: II  
Product Code: KNT  
Dated: December 3, 2013  
Received: December 6, 2013

Dear Michelle Lott,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Benjamin  Fisher -S

Benjamin R. Fisher, Ph.D.  
Director  
Division of Reproductive, Gastro-Renal,  
and Urological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



4700 South Overland Drive • Tucson, Arizona 85714 • 520-882-7794 • www.xeridien.com

## Indications for Use

510(k) Number (if known): K130674

Device Name: Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube

### Indications For Use:

The Entuit™ Thrive Balloon Retention Gastrostomy Feeding Tube is indicated in adult and pediatric populations for use in percutaneous placement of an Enteral feeding tube for feeding and/or administration of medication in conjunction with an established GI stoma tract. The replacement tube may also be used for gastric decompression.

Prescription Use  AND/OR Over-The-Counter Use   
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Benjamin R. Fisher -S  
2014.01.13 15:39:18 -05'00'

Page 1 of 1