

## 510(k) Summary

JUL 12 2013

### 510(k) Owner

Yamahachi Dental Manufacturing Co.  
54-1 Ochigara, Nishiura-Cho  
Gamagori City, Aichi Pref., Japan 443-0105

Phone: 081-533-577121  
Fax: 081-533-5771764  
Contact: Lawrence Greenwald (yamahachi1@yahoo.com)

### Date Prepared

April 1, 2013

### Consultant

Denterprise International, Inc.  
110 E. Granada Blvd., Suite 207  
Ormond Beach, FL 32176

Phone: 386-672-0450  
Fax: 386-672-4402  
Primary Contact: Claude Berthoin, President (claude@510kfda.com)  
Secondary Contact: John Gillespy, VP Regulatory Affairs (john@510kfda)  
Holly Layne, Assistant (holly@510kfda)

### Device Classification

510(k) Name: Yamahachi Denture Base Resins  
Trade/Model Names: Basis, Basis HI, Basis Twin Cure, Basis Flow, Re-Fine Bright  
Common Name: Denture Base Polymer  
Classification Name: Denture Relining, Repairing, or Rebasing Resin  
Regulation Number: 872.3760  
Product Code: EBI  
Regulatory Class: II

### Predicates

K102874 Veracril, Veracril HI, Ez-Cryl (New Stetic, Columbia)  
K102640 Vertex Self-Curing, Vertex Castavaria (Vertex Dental, The Netherlands)

**Intended Use**

Yamahachi Denture Base Resins is a system of heat- and self-cure acrylic polymers intended for fabrication or repair of the denture base.

**Device Description**

Yamahachi Denture Base Resins is supplied in powder and liquid form. The powder is primarily a polymer of polymethyl methacrylate (PMMA) beads with small quantities of initiator and color pigments. The liquid is primarily the monomer methyl methacrylate (MMA) with small quantities of a cross-linking agent and activator.

To fabricate the denture base, the powder and liquid materials are mixed together and stirred to create a dough state that is packed or poured into a mold, saddle, or core. The resin then cures either through heat application (water bath, rapid immersion, or microwave oven) or a self-curing process (pressure vessel or via quick-setting).

While heat-cure resins are generally used to fabricate denture bases, self-cure acrylics are most often indicated for repair and relining of the dentures. The subject device includes three heat-cure resins (Basis, Basis HI, and Basis Twin Cure) and two self-cure resins (Basis Flow and Re-Fine Bright).

**Device Comparison To Predicates**

A comparison of technological characteristics between the subject and predicate devices follows in Tables 5A and 5B...

**Table 5A -- Comparison Of Heat-Cure Acrylics**

	<b>Subject Devices</b>	<b>Predicate Devices</b>	<b>Differences</b>
<b>Device Names</b>	<b>Basis, Basis HI, Basis Twin Cure</b>	<b>Veracril, Veracril HI, EZ-cryl</b>	NA
<b>Manufacturer</b>	Yamahachi Dental (Japan)	New Stetic (Columbia)	NA
<b>510(k)</b>	Not assigned yet	K102874 (Bundle of 5 Heat-Curing & 3 Self-Curing)	NA
<b>Classification &amp; Product Code</b>	872.3760; EBI	872.3760; EBI	No difference in Device Description or Intended Use
<b>Device Description</b>	Traditional heat-cure & microwave-cure acrylic resin for total or partial denture base and for removable prosthesis	Denture system consisting of monomer and polymer powder and liquid components (traditional heat-cured & microwave-cured acrylic resin for total or partial denture base and for removable prosthesis)	

<b>Intended Use</b>	Repair or fabrication of denture base	Repair or fabrication of denture base	
<b>Product State</b>	Polymer powder & monomer liquid	Polymer powder & monomer liquid	
<b>Mixing Ratio (Powder:Liquid)</b>	100g:43ml (Basis & Basis HI); 100g:40ml (Basis Twin Cure)	2:1 by weight or 3:1 by volume	The only difference in Method Of Use is determination of mixing ratio. This difference is related to characteristics of polymer particle size, but does not affect performance of device nor raise issues of safety or effectiveness, a conclusion supported by meeting requirements of standard ISO 20795-1.
<b>Mixing Time</b>	30 sec	30 sec	
<b>Pressing Technique</b>	Pack dough in plaster model, apply pressure, trim excess, apply final pressure (10 - 15 minute packing time)	Pack dough in flask, apply pressure, trim excess, apply final pressure	
<b>Application Time</b>	30 minutes	Approx 10 minutes	
<b>Polymerization (Curing) Method</b>	<u>Water Bath</u> : Immerse flask in water & slowly raise to boiling over 30 min, boil 30 - 40 min, air cool 30 min (Basis & Basis HI); <u>Rapid Curing</u> : Immerse flask in boiling water 30 min (Basis TC); <u>Microwave</u> : Put flask in MWO 500W for 3 min, air cool 30 min (Basis TC)	<u>Water Bath</u> : 90 min at 73° C, 30 min boiling, 30 min air cool, 15 min water cool (Veracril devices); <u>Microwave</u> : 10 min in 3 stages in MWO, air cool 30 min, 15 min water cool (EZ-cryl only)	
<b>Components (See legend after table)</b>	<u>Powder</u> : PMMA, MMA (Basis HI), initiator, pigments; <u>Liquid</u> : MMA, crosslinker, activator (Basis TC)	<u>Powder</u> : PMMA, pigments; <u>Liquid</u> : MMA, crosslinker	The only difference in Product Components are minor variances in percentages of ingredients, which do not affect product performance, safety, or effectiveness, as supported by meeting the standards of ISO 20795-1.
<b>Standards of Conformity</b>	ISP 9001:2008; ISO 13485:2003; MDD (93/42/EEC); ISO 14971	ISP 9001:2008; ISO 13485:2003	All Physical Properties within specification
<b>Biocompatibility</b>	ISO 10993-1, ISO 7405	ISO 10993-1, ISO 7405	
<b>Physical Properties</b>	ISO 20795-1	ANSI/ADA 12:2002; ISO 1567:1999	
<b>Flexural Strength (65 MPa Minimum)</b>	Basis: 84.3 MPa; Basis HI: 87 MPa; Basis Twin Cure: 81.0 MPa	Veracril: 70.8 MPa; Veracril HI: 88.1 MPa; Ez-cryl: 70.5 MPa	
<b>Flexural Modulus (2,000 MPa Min)</b>	Basis: 2,067 MPa; Basis HI: 2,178 MPa; Basis Twin Cure: 2,299 MPa	Veracril: 5,300 MPa; Veracril HI: 5,804 MPa; Ez-cryl: 5,700 MPa	
<b>Impact-Resistance (Min 1.9 MPa m<sup>1/2</sup>)</b>	Basis/Basis TC: Not Applicable; Basis HI: > 2.26 MPa m <sup>1/2</sup>	Veracril/Ez-cryl: Not Applicable; Veracril High Impact: 3.1 MPa m <sup>1/2</sup>	

<b>Residual Monomer (Maximum 2.2%)</b>	Basis: 0.4%; Basis HI: 0.7%; Basis Twin Cure: 0.2%	Veracril: 0.98%; Veracril HI: 1.88%; Ez-cryl: 0.80%	
<b>Sorption (Max 32 ug/mm3)</b>	Basis: 10.7 ug/mm <sup>3</sup> ; Basis HI: 22.8 ug/mm <sup>3</sup> ; Basis Twin Cure: 22.4 ug/mm <sup>3</sup>	Veracril: 18.1 ug/mm <sup>3</sup> ; Veracril HI: 14.5 ug/mm <sup>3</sup> ; Ez-cryl: 19.1 ug/mm <sup>3</sup>	
<b>Solubility (Max 1.6 ug/mm3)</b>	Basis: 0.6 ug/mm <sup>3</sup> ; Basis HI: 0.4 ug/mm <sup>3</sup> ; Basis Twin Cure: 0.2 ug/mm <sup>3</sup>	Veracril: 0.8 ug/mm <sup>3</sup> ; Veracril HI: 0.9 ug/mm <sup>3</sup> ; Ez-cryl: 0.72 ug/mm <sup>3</sup>	
<b>Classification (ISO 20795-1:2008)</b>	Basis & Basis HI: Type 1 Class 1; <u>Basis Twin Cure</u> : Type 5	<u>Veracril &amp; Veracril HI</u> : Type 1 Class 1; <u>EZ-cryl</u> : Type 5	No Classification difference

MMA - Methyl-methacrylate

PMMA - Poly-methyl-methacrylate

PEMA - 2-(N-pyrolyl) ethyl methacrylate

**Table 5B -- Comparison Of Self-Cure Acrylics**

	<b>Subject Devices</b>	<b>Predicate Devices</b>	<b>Differences</b>
<b>Device Names</b>	<b>Basis Flow, Re-Fine Bright</b>	<b>Vertex Self Curing, Vertex Castavaria</b>	NA
<b>Manufacturer</b>	Yamahachi Dental (Japan)	Vertex-Dental (Netherlands)	NA
<b>510(k)</b>	Not assigned yet	K102640 (Bundle of 3 Self-Curing Resins)	NA
<b>Classification &amp; Product Code</b>	872.3760; EBI	872.3760; EBI	No difference in Device Description or Intended Use
<b>Device Description</b>	<u>Basis Flow</u> : Self-curing denture base material intended as a pouring and repair acrylic. <u>Re-fine Bright</u> : Multifunctional self-polymerizing denture base material for repair and relining of full & partial dentures	<u>Vertex Self Curing</u> : Self-curing denture base material intended as a pouring and repair acrylic. <u>Vertex Castavaria</u> : Multifunctional self-polymerizing denture base material for repair and relining of full & partial dentures	
<b>Intended Use</b>	Repair and relining of full & partial dentures	Fabrication, repair and relining of full & partial dentures	

<b>Product State</b>	Polymer powder & monomer liquid	Polymer powder & monomer liquid	The only difference in Method Of Use is curing method for Re-Fine Bright, which is not pressure cured. This difference does not affect performance of device, safety, or effectiveness, as supported by meeting requirements of standard ISO 20795-1.	
<b>Mixing Ratio (Powder:Liquid)</b>	100 g : 50 - 60 ml	100 g : 60 ml		
<b>Mixing Time</b>	5 - 30 sec	20 sec		
<b>Dough Time</b>	Max 3 min	Max 4.5 - 8 min		
<b>Working Time</b>	Max 4 min	Max 5 - 8 min		
<b>Polymerization (Curing) Method</b>	<u>Basis Flow</u> : 30 min at 55° C in 0.2 atm pressure vessel, air cool 30 min; <u>Re-Fine Bright</u> : Let hard polymerize	<u>Vertex Self Curing</u> : 10 min at 55° C in 2.5 bar pressure vessel; <u>Vertex Castavaria</u> : 30 min at 55° C in 2.5 bar pressure vessel		
<b>Components (See legend before table)</b>	<u>Powder</u> : PMMA, PEMA, activators, pigments; <u>Liquid (Basis Flow)</u> : MMA, crosslinker, activator; <u>Liquid (Re-fine Bright)</u> : MMA, crosslinker, activators	<u>Powder</u> : PMMA 99.1%, Inhibitor (dibenzoyl peroxide) <1%; <u>Liquid</u> : MMA >95%, Crosslinker <5%, Accelerators <1%, UV Absorber <<1%	The only difference in Product Components are minor variances in percentages of ingredients, which do not affect product performance, safety, or effectiveness, as supported by meeting the standards of ISO 20795-1.	
<b>Standards of Conformity</b>	ISP 9001:2008; ISO 13485:2003; MDD (93/42/EEC); ISO 14971	ISO 1567, ISO 179-1, ASTM F 895-84	All Physical Properties within specification	
<b>Biocompatibility</b>	ISO 10993-1, ISO 7405	ISO 7405		
<b>Physical Properties</b>	ISO 20795-1	ISO 20795-1		
<b>Flexural Strength (60 MPa Minimum)</b>	Basis Flow: 80.1 MPa; Re-Fine Bright: 73.8 MPa	Vertex Self Curing: 68 MPa; Vertex Castavaria: 79MPa		
<b>Flexural Modulus (1,500 MPa Min)</b>	Basis Flow: 1,657 MPa; Re-Fine Bright: 1,529 MPa	Vertex Self Curing: 2,028 MPa; Vertex Castavaria: 2,316 MPa		
<b>Residual Monomer (Maximum 4.5%)</b>	Basis Flow: 4.2%; Re-Fine Bright: 3.3%	Vertex Self Curing: 3.76%; Vertex Castavaria: 3.91%		
<b>Sorption (Max 32 ug/mm<sup>3</sup>)</b>	Basis Flow: 18.8 ug/mm <sup>3</sup> ; Re-Fine Bright: 15.8 ug/mm <sup>3</sup>	Vertex Self Cure: 20.3 ug/mm <sup>3</sup> ; Vertex Castavaria: 23.2 ug/mm <sup>3</sup>		
<b>Solubility (Max 8.0 ug/mm<sup>3</sup>)</b>	Basis Flow: 1.8 ug/mm <sup>3</sup> ; Re-Fine Bright: 2.3 ug/mm <sup>3</sup>	Vertex Self Cure: 1.8 ug/mm <sup>3</sup> ; Vertex Castavaria: 1.8 ug/mm <sup>3</sup>		
<b>Classification (ISO 20795-1:2008)</b>	Basis Flow: Type 2 Class 2; Re-Fine Bright: Type 2 Class 1	Vertex SC: Type 2 Class 2; Vertex Castavaria: Type 2 Class 1		No Classification difference

### **Non-Clinical Tests**

Bench tests were performed on the subject device in conformity with ISO 20795-1. These tests confirmed the different resins (Basis, Basis HI, Basis Twin Cure, Basis Flow, and Re-Fine Bright) all met the performance criteria established by that standard.

### **Substantial Equivalence Discussion**

As noted above...

- Product Description and Intended Use are the same for the subject and predicate devices.
- The only difference in Method Of Use is determination of mixing ratio. This difference is related to characteristics of polymer particle size, but does not affect performance of device nor raise issues of safety or effectiveness, a conclusion supported by meeting requirements of standard ISO 20795-1.
- The only difference in Product Components are minor variances in percentages of ingredients, which do not affect product performance, safety, or effectiveness, as supported by meeting the standards of ISO 20795-1.
- The Physical Properties of both the subject and predicate devices conform to the standards of ISO 20795-1.
- Classification of the dental acrylic resins is the same between subject and predicate devices.

Yamahachi Denture Base Resins has the same intended use as the predicates. The subject device also has minor differences in technological characteristics that could not affect safety or effectiveness. Bench tests show conformance with performance standards in ISO 20795-1:2008 for all five models included in this submission (Basis, Basis HI, Basis Twin Cure, Basis Flow, and Re-fine Bright).

In conclusion, Yamahachi Denture Base Resins warrants a finding of substantial equivalence to legally marketed devices from New Stetic (Columbia) and Vertex Dental (The Netherlands) and of proper clearance for premarketing activities in the United States.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 12, 2013

Yamahachi Dental Manufacturing Co.  
C/O Mr. Claude Berthoin, President  
Denterprise International, Inc.  
110 E. Granada Boulevard, Suite 207  
Ormond Beach, FL 32176

Re: K131036  
Trade/Device Name: Yamahachi Denture Base Resins  
Regulation Number: 21 CFR 872.3760  
Regulation Name: Denture Relining, Repairing, or Repairing Resin  
Regulatory Class: II  
Product Code: EBI  
Dated: April 26, 2013  
Received: April 29, 2013

Dear Mr. Berthoin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary S. Runner -S

Kwame Ulmer M.S.  
Acting Division Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications For Use Statement

**Applicant:** Yamahachi Dental Manufacturing Co.

**510(k) Number (if known):** K131036

**Device Name:** Yamahachi Denture Base Resins

### Indications For Use:

Yamahachi Denture Base Resins is intended for fabrication or repair of the denture base.

*All devices are sold by or on the order of a physician. They are not for use by the general public or over-the-counter.*

Prescription Use X  
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use \_\_\_\_  
(21 CFR Part 801 Subpart C)

---

(Please Do Not Write Below This Line – Continue On Another Page If Needed)

Concurrence of CDRH, Office of Device Evaluation (ODE).

---

Division Sign-Off  
Office of Device Evaluation

510(k) \_\_\_\_\_

---

(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K131036