

K131366

## 510(k) Summary

---

### Introduction

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

### 1. Submitter Name, Address, Contact

Roche Diagnostic Corporation  
9115 Hague Rd.  
Indianapolis, IN 46250  
(317) 521-4793  
Contact Person: Nate Carrington  
Date Prepared: February 12, 2013

OCT 11 2013

### 2. Device Name

Proprietary names: ACCU-CHEK Aviva Expert System  
ACCU-CHEK Aviva Plus Test Strip  
ACCU-CHEK Botus Advisor  
ACCU-CHEK Aviva Expert Meter

Classification name: Glucose dehydrogenase, glucose test system (21 C.F.R. § 862.1345); Class II

Classification name: Drug dosing calculation (21 C.F.R. § 868.1890); Class II

LFR, Glucose Dehydrogenase  
NDC, Drug Dosing Calculator

### 3. Predicate Device

ACCU-CHEK Aviva Combo meter, cleared as a component of the ACCU-CHEK Combo System in #k111353

The ACCU-CHEK Aviva Combo blood glucose system is intended for self-testing outside the body (in vitro diagnostic use) by people with diabetes at home as an aid to monitoring the effectiveness of diabetes control. The ACCU-CHEK Aviva Combo meter is also indicated for the management of diabetes by calculating an insulin dose or carbohydrate intake based on user-entered data.

#### **4. Device Description**

The ACCU-CHEK Aviva Expert System consists of the following:

- ACCU-CHEK Aviva Expert meter
- ACCU-CHEK Bolus Advisor (a component of the Aviva Expert meter)
- ACCU-CHEK Aviva Plus test strips (k101299)
- ACCU-CHEK Aviva control solutions (k043474)

The ACCU-CHEK Aviva Expert system is a blood glucose monitoring system that makes use of the ACCU-CHEK Aviva Plus test strips (cleared under k101299) and the ACCU-CHEK Aviva control solutions (cleared under k043474). The industrial design of the ACCU-CHEK Aviva Expert meter is nearly identical to that of the ACCU-CHEK Aviva Combo meter (cleared under k111353); the only difference between the two devices is that the ACCU-CHEK Aviva Expert meter does not have Bluetooth wireless capability and therefore cannot be connected to an insulin pump, whereas the ACCU-CHEK Aviva Combo meter does have Bluetooth capability and can be connected to an insulin pump.

The ACCU-CHEK Aviva Expert system provides the user with the ability to measure capillary blood glucose levels when a sample of capillary blood is applied to the test strip. The meter also provides an optional insulin bolus calculator (the ACCU-CHEK Bolus Advisor) designed for use by individuals with diabetes who require insulin. This feature is optional in that a user can simply obtain a blood glucose value through capillary blood testing and does not need to use the insulin bolus calculator portion of the system if it is not desired. The insulin bolus calculator algorithm is identical to the bolus calculator algorithm that was cleared during the Aviva Combo 510(k) submission (k111353). For the ACCU-CHEK Aviva Expert system, this bolus calculator is meant to be used by patients with diabetes on multiple daily insulin injection (MDI) therapy. In order to calculate the appropriate bolus of insulin, the ACCU-CHEK Bolus Advisor takes the measured bG, the target bG, the carbohydrate intake, the insulin-to-carbohydrate ratio, the insulin sensitivity, health events (such as exercise), the time of day, and the active insulin into account. Before using the ACCU-CHEK Aviva Expert system, a physician or healthcare professional must provide the patient-specific target blood glucose, insulin-to-carbohydrate ration, and insulin sensitivity parameters.

#### **5. Intended Use**

The ACCU-CHEK Aviva Expert System is indicated as an aid in the treatment of insulin-requiring diabetes. The ACCU-CHEK Aviva Expert System consists of the ACCU-CHEK Aviva Expert Meter, ACCU-CHEK Aviva Plus test strips, ACCU-CHEK Aviva control solutions, and ACCU-CHEK Bolus Advisor. The ACCU-CHEK Aviva Expert System is intended to facilitate the optimization of glycemic control in patients who are trained in multiple daily insulin injection therapy and are under the supervision of healthcare professionals experienced in managing insulin treated patients.

The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended to be used for the quantitative measurement of glucose in fresh capillary whole blood samples drawn from the fingertips. The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended for self-testing outside the body (in vitro diagnostic use) by people with diabetes. The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended to be used by a single person and should not be shared. The ACCU-CHEK Aviva Expert blood glucose monitoring system should not be used for the diagnosis or screening of diabetes or for neonatal use. Alternative site testing should NOT be used with the ACCU-CHEK Aviva Expert blood glucose monitoring system. The ACCU-CHEK Aviva Expert System is intended for prescription home use only.

The ACCU-CHEK Aviva Expert meter is also indicated for the calculation of an insulin dose or carbohydrate intake based on user-entered data. The ACCU-CHEK Bolus Advisor, as a component of the Accu-Chek Aviva Expert meter, is intended for use in providing insulin dose recommendations in response to blood glucose, health events, and carbohydrate input. The ACCU-CHEK Bolus Advisor is intended to provide direction for insulin adjustment within the scope of a pre-planned treatment program from a healthcare professional. Before its use, a physician or healthcare professional must prescribe the ACCU-CHEK Aviva Expert System and provide the patient-specific target blood glucose, insulin-to-carbohydrate ratio, and insulin sensitivity parameters to be programmed into the ACCU-CHEK Bolus Advisor. Once programmed, a patient must consult with his/her physician or healthcare professional before making any changes to these ACCU-CHEK Bolus Advisor settings.

## **6. Substantial Equivalence**

The ACCU-CHEK Aviva Expert System is substantially equivalent to the ACCU-CHEK Aviva Combo System. Below is a table that provides a comparison between the ACCU-CHEK Aviva Expert System and its predicate device.

**Similarities Table**

<b>System Feature/Claim</b>	<b>Detail</b>
Test Strip	<b>Identical:</b> Both systems make use of the Aviva Plus test strip.
Test Strip Production Processes and Lot-Release Criteria	<b>Identical:</b> Both systems make use of the Aviva Plus test strip.
bG Measurement Algorithm	<b>Identical:</b> The fundamental scientific technology for the measurement of blood glucose has not changed from the predicate.
Meter Main Integrated Circuit Board	<b>Identical:</b> The firmware components and layout are identical between the two meters.
Meter Display Module, Display Frame, and Button Module	<b>Identical:</b> The same display unit, display frame, and button unit are used by the two meters.
Meter Housing and Battery Contacts	<b>Identical:</b> For both meters, the plastic parts are produced by the same molding tool using the same material, and the battery contacts are the same.
Meter Production Environment	<b>Identical:</b> Both meters are produced on one production line/process.
Underdose Detection and Meter Failsafes	<b>Identical:</b> The fundamental scientific technology for the measurement of blood glucose has not changed from the predicate.
Integrity Check for Strip	<b>Identical:</b> Early in the measurement sequence, the meter measures the resistance of the gold on the un-dosed strip to assure that it has been properly inserted and that the quality is not compromised. The meter measures the background conductivity and electrical current prior to dosing to assure that the reagent quality is not compromised or that the strip was not prematurely dosed.
Insulin Bolus Calculator Algorithm	<b>Identical:</b> Both systems use the exact same bolus calculator, the ACCU-CHEK Bolus Advisor, for insulin dosing calculations.

**Similarities Table (continued)**

<b>System Feature/Claim</b>	<b>ACCU-CHEK Aviva Expert Meter with Aviva Plus Test Strip</b>	<b>ACCU-CHEK Combo Meter with Aviva Plus Test Strip Predicate (k111353)</b>
Indications for Use	Quantitative measurement of glucose (sugar) in capillary blood from the finger tip, indicated for diabetes management by calculating an insulin dose or carbohydrate intake	Quantitative measurement of glucose (sugar) in capillary blood from the finger tip, indicated for diabetes management by calculating an insulin dose or carbohydrate intake
Test Principle	Amperometric Detection	Amperometric Detection
Enzyme	Mut. Q-GDH	Mut. Q-GDH
Sample Hematocrit	10 to 65%	10 to 65%
Maximum Altitude	10,000 feet	10,000 feet
Measuring Range	20 – 600 mg/dL	20 – 600 mg/dL
Sample Volume	0.6 $\mu$ L	0.6 $\mu$ L
Test Time	5 seconds	5 seconds
Operating Temperature and Relative Humidity	14 to 38°C (57 to 100°F) 10 to 80% r.h.	14 to 38°C (57 to 100°F) 10-80% r.h.
Coding	Code key insertion; lot-specific code key provided with each box of Aviva Plus test strips	Code key insertion; lot-specific code key provided with each box of Aviva Plus test strips
Precision	For response targets below 75 mg/dL, the SD is $\leq$ 5.0 mg/dL, and for response targets $\geq$ 75 mg/dL, the CV is $\leq$ 5.0%.	For response targets below 75 mg/dL, the SD is $\leq$ 5.0 mg/dL, and for response targets $\geq$ 75 mg/dL, the CV is $\leq$ 5.0%.
Double Dosing	No	No
Alternate Site Testing	No	No

**Similarities Table (continued)**

<b>System Feature/Claim</b>	<b>ACCU-CHEK Aviva Expert Meter with Aviva Plus Test Strip</b>	<b>ACCU-CHEK Combo Meter with Aviva Plus Test Strip Predicate (k111353)</b>
Closed and Open Vial Shelf Life Stability	18 months	18 months
Control Solutions	Aqueous, 2 levels, uses ACCU-CHEK Aviva Control Solutions	Aqueous, 2 levels, uses ACCU-CHEK Aviva Control Solutions
Primary Packaging	Standard flip top vial	Standard flip top vial
Handling	Automatic on/off with strip insertion or by pressing button	Automatic on/off with strip insertion or by pressing button
Dimensions	3.7 x 2.1 x 1 in LWH; approximately 3.6 oz. with batteries inserted	3.7 x 2.1 x 1 in LWH; approximately 3.6 oz. with batteries inserted
Data Transfer	Infrared	Infrared
Date Reminders	Yes	Yes
bG Test Reminders	Yes	Yes
Alarm Clock Reminders	Yes	Yes
Target bG Levels	Yes	Yes
Health Events	Yes	Yes
Electronic Diary	Yes	Yes
Limitations of Procedure	Galactose >15 mg/dL will cause overestimation of blood glucose results.	Galactose >15 mg/dL will cause overestimation of blood glucose results.
	Lipemic Samples >1800 mg/dL	Lipemic Samples >1800 mg/dL
	Intravenous administration of ascorbic acid which results in blood concentrations of ascorbic acid >3 mg/dL will cause overestimation of blood glucose results	Intravenous administration of ascorbic acid which results in blood concentrations of ascorbic acid >3 mg/dL will cause overestimation of blood glucose results

**Differences Table**

<b>System Feature/Claim</b>	<b>ACCU-CHEK Aviva Expert Meter with Aviva Plus Test Strip</b>	<b>ACCU-CHEK Combo Meter with Aviva Plus Test Strip Predicate (k111353)</b>
RF Wireless Capability	No	Yes; additional components to satellite board for Bluetooth communication with ACCU-CHEK Spirit Combo insulin infusion pumps
Measurement Units of Insulin Bolus Result Calculations	0.5 units (appropriate for syringe/pen administration of insulin)	0.1 units (appropriate for pump administration of insulin)
Basal Insulin*	Diary feature for tracking basal insulin	Basal insulin is controlled via pump
User Group	Diabetes patients treated with multiple daily insulin injection (MDI) therapy	Diabetes patients treated with insulin pump therapy or multiple daily insulin injection (MDI) therapy

\* The basal insulin values that are recorded do not influence the bolus advice.

**7. Data demonstrating substantial equivalence**

Performance testing on the ACCU-CHEK Aviva Expert System demonstrated that the device meets the performance requirements for its intended use. The data demonstrate that the system is substantially equivalent to the predicate device.

Below is the method comparison data for the system:

Results for glucose concentrations less than 75 mg/dL

Within $\pm 5$ mg/dL	Within $\pm 10$ mg/dL	Within $\pm 15$ mg/dL
41/48 (85.4%)	48/48 (100%)	48/48 (100%)

Results for glucose concentrations greater than or equal to 75 mg/dL

Within $\pm 5\%$	Within $\pm 10\%$	Within $\pm 15\%$	Within $\pm 20\%$
147/252 (58.3%)	222/252 (88.1%)	246/252 (97.6%)	250/252 (99.2%)

Below is the repeatability (within lot) precision for the system:

<b>Blood</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
N	100	100	100	100	100
Mean [mg/dL]	42.1	84.5	137.8	208.2	345.0
SD [mg/dL]	1.2	2.2	3.3	5.6	7.9
CV [%]	2.9	2.6	2.4	2.7	2.3

Below is the reproducibility (intermediate or day-to-day) precision for the system:

<b>Control solutions</b>	<b>Low</b>	<b>Mid</b>	<b>High</b>
N	100	100	100
Mean [mg/dL]	45.1	117.6	303.0
SD [mg/dL]	1.1	2.4	5.1
CV [%]	2.4	2.0	1.7



October 11, 2013

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Roche Diagnostics Corporation  
C/O Nathan A. Carrington, Ph.D.  
9115 Hague Rd  
INDIANAPOLIS IN 46250-0457

Re: K131366

Trade/Device Name: ACCU-CHEK Aviva Expert Blood Glucose Monitoring System  
Regulation Number: 21 CFR 862.1345  
Regulation Name: Glucose test system  
Regulatory Class: II  
Product Code: LFR, NDC  
Dated: August 26, 2013  
Received: August 29, 2013

Dear Dr. Carrington:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2—Mr. Carrington

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Carol C. Benson -S for

Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use Form

510(k) Number (if known): k131366

Device Name: ACCU-CHEK Aviva Expert Blood Glucose Monitoring System

### Indications for Use:

The ACCU-CHEK Aviva Expert System is indicated as an aid in the treatment of insulin-requiring diabetes. The ACCU-CHEK Aviva Expert System consists of the ACCU-CHEK Aviva Expert Meter, ACCU-CHEK Aviva Plus test strips, ACCU-CHEK Aviva control solutions, and ACCU-CHEK Bolus Advisor. The ACCU-CHEK Aviva Expert System is intended to facilitate the optimization of glycemic control in patients who are trained in multiple daily insulin injection therapy and are under the supervision of healthcare professionals experienced in managing insulin treated patients.

The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended to be used for the quantitative measurement of glucose in fresh capillary whole blood samples drawn from the fingertips. The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended for self-testing outside the body (in vitro diagnostic use) by people with diabetes. The ACCU-CHEK Aviva Expert blood glucose monitoring system is intended to be used by a single person and should not be shared. The ACCU-CHEK Aviva Expert blood glucose monitoring system should not be used for the diagnosis or screening of diabetes or for neonatal use. Alternative site testing should NOT be used with the ACCU-CHEK Aviva Expert blood glucose monitoring system. The ACCU-CHEK Aviva Expert System is intended for prescription home use only.

The ACCU-CHEK Aviva Expert meter is also indicated for the calculation of an insulin dose or carbohydrate intake based on user-entered data. The ACCU-CHEK Bolus Advisor, as a component of the Accu-Chek Aviva Expert meter, is intended for use in providing insulin dose recommendations in response to blood glucose, health events, and carbohydrate input. The ACCU-CHEK Bolus Advisor is intended to provide direction for insulin adjustment within the scope of a pre-planned treatment program from a healthcare professional. Before its use, a physician or healthcare professional must prescribe the ACCU-CHEK Aviva Expert System and provide the patient-specific target blood glucose, insulin-to-carbohydrate ratio, and insulin sensitivity parameters to be programmed into the ACCU-CHEK Bolus Advisor. Once programmed, a patient must consult with his/her physician or healthcare professional before making any changes to these ACCU-CHEK Bolus Advisor settings.

Prescription Use XX AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

---

Concurrence of CDRH, Office of In Vitro Devices and Radiologic Health (OIR)

Stayce Beck

Division Sign-Off  
Office of In Vitro Devices and Radiologic Health

510(k) \_\_\_\_\_

Page 1 of \_\_\_\_\_