

JAN 13 2014

510(k) Summary

Date: January 8, 2014

Trade Name: CO₂/O₂ Nasal Cannula

Common Name: Oxygen Delivery, Carbon Dioxide Sampling Nasal Cannula

Classification Name: Carbon Dioxide Gas Analyzer

Regulation: Class II per 21CFR 868.1400

Product Code: CCK

Sponsor: Southmedic Inc.
50 Alliance Blvd.
Barrie, ON
L4M 5K3
Canada
(705) 720-1902

Contact: Tish Anger

Phone: (705)720-1902

Fax: (705)728-9537

LEAGALLY MARKETED PREDICATE DEVICES

This premarket notification will demonstrate that Southmedic's CO₂/O₂ Sampling Nasal Cannula is substantially equivalent to the devices listed below:

- Salter Labs #4706 CO₂/O₂ Nasal Cannula (K892406)
- Hudson RCI BiFlow 1850 (K961150)
- UnoMedical/Hospitak 368E (K915228)

DEVICE DESCRIPTION

Southmedic's CO₂/O₂ Nasal Cannula incorporates sampling nares that channel expired carbon dioxide to a capnograph. A port located between the two nares is used to deliver oxygen to the patient as needed from an oxygen source. Its unique design delivers O₂ to both the nose and mouth area. This device when in use involves surface contact to the patient's intact skin. There

is also gas pathway contact between patient and tubing/nasal cannula. Contact is considered limited exposure (up to 24 hours).

INTENDED USE

This device is intended to provide a means for sampling end tidal carbon dioxide with the option to deliver supplemental O₂ therapy to patients for up to 24 hours.

TECHNOLOGICAL CHARACTERISTICS

Southmedic's CO₂/O₂ Nasal Cannula and the predicate devices have similar technological characteristics. Specifically Southmedic's CO₂/O₂ Nasal Cannula and the predicate devices are all designed to draw expired end-tidal CO₂ from the patient through the nasal cannula and tubing to a capnograph, and administer medical grade USP O₂ to patients from an oxygen supply.

Southmedic's CO₂/O₂ Nasal Cannula's cannula has been designed with slight differences from the predicate devices. All predicate cannulae have two separate systems built into the part to allow for the delivery of oxygen and sampling of end-tidal carbon dioxide. Each predicate takes a different approach. Exhaled carbon dioxide is collected through the nasal prongs and a single port was created central to the nasal prongs for O₂ delivery in Southmedic's model. Modification to the outer portion helps direct O₂ down towards the mouth as well as the nose.

SUBSTANTIAL EQUIVALENCE SUMMARY

Characteristic	Southmedic's CO₂ Sampling Nasal Cannula	Salter Labs #4706 CO₂/O₂ Nasal Cannula	Hudson RCI BiFlow 1850	UnoMedical/Hospitak 368E
510(k)	K131410	K892406	K961150	K915228
Intended Use	This device is intended to provide a means for sampling ET _{CO₂} with the option to deliver supplemental O ₂ therapy to patients.	Same	Same	Same
Prescription	Yes	Yes	Yes	Yes
Technological Characteristics	Provides a means to deliver exhaled CO ₂ to a capnograph via	Same	Same	Same

Characteristic	Southmedic's CO₂ Sampling Nasal Cannula	Salter Labs #4706 CO₂/O₂ Nasal Cannula	Hudson RCI BiFlow 1850	UnoMedical/Hospitak 368E
	tubing, provides a means to deliver continuous USP grade medicinal O ₂ as necessary.			
Design	Collects CO ₂ from nasal prongs and delivers O ₂ through a single port between nasal prongs	Delivers oxygen from one nasal prong and samples CO ₂ from the other.	Delivers oxygen and samples CO ₂ from both nasal prongs.	Same as Southmedic's, however, O ₂ delivery through two ports rather than one.
Material	Phthalate free PVC Not made with natural rubber latex.	PVC containing Phthalates Not made with natural rubber latex.	PVC containing Phthalates Not made with natural rubber latex.	PVC containing Phthalates Not made with natural rubber latex.
Energy Used or Delivered	Not Applicable	Not Applicable	Not Applicable	Not Applicable
Manufacturing Process	Injection molding, extrusion, assembly	Injection molding, extrusion, assembly	Injection molding, extrusion, assembly	Injection molding, extrusion, assembly
Performance	Same	Same	Same	Same
Labelling	Same	Same	Same	Same

TESTING

In order to support the specific safety of Southmedic's CO₂/O₂ Nasal Cannula, as well as confirm the suitability of the materials used in the manufacture of this product, extensive biocompatibility studies were contracted to an independent laboratory. The following table summarizes the biocompatibility testing conducted by Nelson Labs:

Biocompatibility Testing	
Test/Method:	Agar Overlay – Test for Cytotoxicity
Objective:	To determine the cytotoxicity of diffusible components from materials.
Acceptance Criteria	The United States Pharmacopeia & National Formulary (USP 87) states that the test article meets criteria if the reactivity grade is not greater than 2 or a mild reactivity. The ANSI/AAMI/ISO 10993-5 standard states the achievement of a numerical grade greater than 2 is considered a cytotoxic effect.
Results:	Three samples were investigated obtaining an average score of 1. Results were compared to both positive and negative controls (average scores of 4 and 0 respectively).
Discussion:	Southmedic's CO ₂ /O ₂ Nasal Cannula is composed of materials that do not have a cytotoxic effect. This supports the safety of the device as well as the suitability of the materials used in its production.
Test/Method:	ISO 10993 Part 10 Guinea Pig Buehler Sensitization Test
Acceptance/Evaluation Criteria	Test results were based on incidence and severity of the sensitization reaction. Incidence was defined as the percentage of animals exhibiting a sensitization reaction at each observation time point (24 and 48 hours) post challenge. Severity was calculated by dividing the sum of the scores of one (1) or greater in the test group generally indicated sensitization, provided grades of less than one (1) were observed in the control animals. If severity of one (1) or greater was noted on the controls, then the reaction of the test animals exceeded the most severe control reaction to be considered due to sensitization. In the final analysis of data consideration was given to overall patterns, intensity, duration and character of reactions of the test as compared to the control animal.
Results:	No sensitization reaction was observed in either the control or test group.
Discussion	Southmedic's CO ₂ /O ₂ Nasal Cannula is composed of materials that do not have a sensitization effect. This supports the safety of the device as well as the suitability of the materials used in its production.

Test/Method:	ISO 10993 Part 10 – Primary Skin Irritation Test in Rabbits										
Acceptance/Evaluation Criteria:	<p>After the 72 hour scoring, all erythema plus edema scores generated during the 24±2, 48±2, and 72±2 hour observations were totaled separately for each test sample and control for each animal. The resulting totals for each animal were divided by 6 (two tests/observation site, three time points) to determine separate test article and control primary irritation scores for each animal. The primary irritation index is characterized by the number (score) and description (response category) given in below:</p> <table border="1"> <thead> <tr> <th>Response Category</th> <th>Primary Irritation Index (PII)</th> </tr> </thead> <tbody> <tr> <td>Negligible</td> <td>>0 to 0.4</td> </tr> <tr> <td>Slight</td> <td>0.5 to 1.9</td> </tr> <tr> <td>Moderate</td> <td>2 to 4.9</td> </tr> <tr> <td>Severe</td> <td>5 to 8</td> </tr> </tbody> </table>	Response Category	Primary Irritation Index (PII)	Negligible	>0 to 0.4	Slight	0.5 to 1.9	Moderate	2 to 4.9	Severe	5 to 8
Response Category	Primary Irritation Index (PII)										
Negligible	>0 to 0.4										
Slight	0.5 to 1.9										
Moderate	2 to 4.9										
Severe	5 to 8										
Results:	The test subject received a Primary Irritation Index Score of 0 in all areas.										
Discussion:	Since no skin irritation was observed in the test subjects, this test demonstrates that Southmedic's CO ₂ /O ₂ Nasal Cannula is composed of materials that are a non-irritant to skin. This supports the safety of the device as well as the suitability of the materials used in its production.										

Risk evaluation performed by Southmedic determined that because of the well characterized safe and effective use of nasal cannulas for oxygen delivery and end tidal CO₂ sampling, animal or clinical testing were not necessary to support this application. In order to confirm effectiveness of Southmedic's CO₂/O₂ Nasal Cannula for its intended use, a comparative performance bench test against predicate devices was contracted from an independent laboratory. The following table summarizes the performance testing conducted by Piper Laboratory:

Study Attribute	Description
Test subjects	Southmedic's CO ₂ /O ₂ Nasal Cannula was tested and results compared against predicate devices (see above)
Objective	To measure the end tidal CO ₂ values of four oxygen delivering nasal cannulae under simulated patient conditions
Acceptance Criteria	End tidal CO ₂ values for Southmedic Nasal Cannula shall not be statistically different, or shall vary less from the true end tidal CO ₂ value than the predicate products
Apparatus	Harvard Respiratory pump and mannequin head to simulate patient head. Nasal cannula connected to oxygen supply line, CO ₂ sensing line and CO ₂ detector (capnograph)

Simulated Respiratory Settings	<p>1. Adult, Respiratory Rate 8, Tidal Volume 800</p> <p>2. Adult, Respiratory Rate 16, Tidal Volume 600</p> <p>3. Pediatric, Respiratory Rate 30, Tidal Volume 300</p> <p>All conditions had I:E ratio of 1:1 (inspiration:expiration)</p>
Deviations	No deviations were noted during the study. All equipment met predetermined operation and calibration testing before and after testing
Oxygen delivery for study	Source flows of 0, 1, 3, and 5 liters per minute for each setting was used. The system was allowed at least 3 minutes with oxygen deliver prior to sampling CO ₂ to allow the system to equilibrate.
Measurements	Each test was sampled three times at each combination of settings for a total of 144 tests (4 samples x 4 tests per sample x 3 respiratory settings x 4 oxygen flow rate settings = 144 tests total)
Results	<p>The end-tidal CO₂ values for the Southmedic CO₂/O₂ Nasal Cannula was not statistically different, or varied less from the true end tidal CO₂ value than predicate products.</p> <p>A comparison of variances of the end tidal CO₂ measurements from the actual value for all four devices indicates that the Southmedic nasal cannula has the least variance from true value when compared to the three predicate devices within 95% confidence intervals.</p>
Discussion	<p>The Southmedic CO₂/O₂ Nasal Cannula met the predetermined acceptance criteria of the test.</p> <p>Both Southmedic's CO₂/O₂ Nasal Cannula and predicate devices are intended to provide a means for sampling ETCO₂ with the option to deliver supplemental O₂ therapy to patients. This performance test demonstrated that Southmedic's CO₂/O₂ Nasal Cannula performs equally to, or better than, the predicate devices in measuring end tidal CO₂. This supports the substantial equivalence of Southmedic's CO₂/O₂ Nasal Cannula to the predicate products.</p>

CONCLUSIONS

In summary, Southmedic's CO₂/O₂ Nasal Cannula path through the 510(k) regulatory analysis is as follows:

- The Southmedic CO₂/O₂ Nasal Cannula and the predicate devices are all intended for use as a means for sampling ETCO₂ with the option to deliver supplemental O₂ therapy to patients.
- The Southmedic CO₂/O₂ Nasal Cannula has similar key technological characteristics as the predicate devices. All three are manufactured by extrusion, injection molding, and assembly. Designs are based on the same principle of one tube providing oxygen from a source and the other delivering expired CO₂ to a capnograph.
- In order to verify the safety of the materials used in the manufacture of Southmedic's CO₂/O₂ Nasal Cannula, extensive biocompatibility testing was conducted. The testing demonstrated that materials used have no cytotoxic effect. The testing also confirmed

that the materials used are a non-irritant and do not have a sensitization effect. These results help support the safety of Southmedic's CO₂/O₂ Nasal Cannula

- The use of nasal cannulas for sampling CO₂ and delivering oxygen is well established in the medical community. In order to address any difference between Southmedic's CO₂/O₂ Nasal Cannula and the predicate devices in their design, additional studies were performed. In order to demonstrate the effectiveness of the new design, a comparative bench test was conducted on Southmedic's behalf. The results demonstrated that Southmedic's CO₂/O₂ Nasal Cannula varied less or equal to the true values when compared to the predicate device. This supports the effectiveness of the device.

Therefore, Southmedic's CO₂/O₂ Nasal Cannula meets the criteria for substantial equivalence to Salter Lab's CO₂/O₂ Nasal Cannula, Hudson RCI Biflow Cannula, and the UnoMedical/Hospitak Nasal Cannula.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

January 13, 2014

Southmedic Incorporated
Ms. Tish Anger
Vice President of Quality and Regulatory Assurance
50 Alliance Blvd.
Barrie, Ontario
Canada L4M 5K3

Re: K131410
Trade/Device Name: CO2/O2 Nasal Cannula
Regulation Number: 21 CFR 868.1400
Regulation Name: Carbon Dioxide Gas Analyzer
Regulatory Class: II
Product Code: CCK
Dated: December 12, 2013
Received: December 13, 2013

Dear Ms. Anger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Tejashri Purohit-Sheth, M.D.
Clinical Deputy Director

DAGRID

FOR

Erin I. Keith, M.S.
Acting Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number: K131410
Device Name: CO₂/O₂ Nasal Cannula
Intended Use: This device is intended to provide a means for sampling end tidal carbon dioxide with the option to deliver supplemental O₂ therapy to patients for up to 24 hours.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Nayan J. Patel -S
2014.01.13 09:56:03 -05'00'