

K131517

Bayer HealthCare



## 510(k) Summary

**Submitter:**

MEDRAD, Inc.  
Bayer Medical Care, Inc.  
One Medrad Drive  
Indianola, PA 15051

AUG 02 2013

**Contact Person:**

Lisa A. Ewing  
Principal Regulatory Affairs Specialist  
Phone: (412) 406-3780  
Fax: (412) 406-4052  
Email: [lisa.ewing@bayer.com](mailto:lisa.ewing@bayer.com)

MEDRAD, INC.

One Medrad Drive  
Indianola, PA 15051  
U.S.A.

(412) 767-2400

[www.ri.bayer.com](http://www.ri.bayer.com)

**Date Prepared:**

May 24, 2013

**Device Trade Name:**

MEDRAD Twist & Go High Pressure  
Connector Tubing

**Common Name:**

Connector Tubing / Tubing, Fluid Delivery

**Classification Name:**

Tubing, Fluid Delivery [21 CFR 880.5440]

**Product Code:**

FPK

**Classification:**

Class II

**Predicate Device(s):**

The subject device is equivalent to the  
following device:  
MEDRAD High Pressure Connecting Tube  
K810924, April 17, 1981

**Device Description:**

The MEDRAD Twist & Go High Pressure  
Connector Tubing is a sterile, single-use  
only, disposable connector tube that is used  
to transfer contrast media and common  
flushing solutions from an automated  
injector syringe to an administration set  
during angiographic procedures.



Automated injection systems with which the MEDRAD Twist & Go High Pressure Connector Tubing is intended to be used include the MEDRAD Mark 7 Arterion Injection System and equivalent injectors.

**Intended Use:**

The Twist & Go High Pressure Connector Tubing is intended to be used in the delivery of contrast media and common flushing solutions. The device is indicated for single use only with MEDRAD injectors.

**Performance Data:**

Bench and laboratory testing were performed to support a determination of substantial equivalence to the predicate devices. Results from the testing provide assurance that the proposed device conforms to the requirements for its intended use. This included the following testing:

- Performance
- Package Integrity
- System level testing with a MEDRAD automated injector



Comparison of Features in MEDRAD High Pressure Connecting Tube (Predicate) and MEDRAD Twist & Go High Pressure Connector Tubing (Proposed)

Item		Predicate Device: MEDRAD High Pressure Connecting Tube (K810924)	Proposed Device: MEDRAD Twist & Go High Pressure Connector Tube
Labeling	Intended Use	The device, intended to be used by individuals with adequate training and experience in diagnostic imaging studies, are to be used to deliver contrast media during angiographic procedures.	The Twist & Go High Pressure Connector Tubing is intended to be used in the delivery of contrast media and common flushing solutions with Twist & Go Syringes. The device is indicated for single use only with MEDRAD injectors.
	Description	Polycarbonate hub insert molded to a PVC tube	Polycarbonate hub UV bonded to a flexible polyurethane nylon tube
Construction	Syringe connection	Clear Polycarbonate Female Luer	Clear Polycarbonate Twist & Go Fasturn Nut
	Administration Set Connection	Clear Polycarbonate standard male luer	Clear Polycarbonate rotating male luer
	Tubing Material	Clear Polyvinyl Chloride	Clear Polyurethane and nylon
	Tubing Volumes	1.00 – 4.78 ml	1.22 – 4.33 ml
	Tubing Lengths	25-122 cm	50-150 cm
	Outer Diameter	.478 cm (.188")	.368 cm (.145")
	Inner Diameter	.224 cm (.088")	.183 cm (.072")
	Adhesive	None	UV adhesive
Packaging	Type	Pouch	Pouch
	Material	Tyvek and clear polymer film	Tyvek and clear polymer film
	Shelf Life	3 years	3 years
Biological	Biocompatibility	Compliant to applicable sections of ISO/AAMI 10993-1:2009	Compliant to applicable sections of ISO/AAMI 10993-1:2009
	Pyrogenicity	Non-pyrogenic	Non-pyrogenic fluid path
	Latex content	Not made with natural rubber latex	Not made with natural rubber latex
	Sterilization Type	Ethylene Oxide	Ethylene Oxide
	Sterilization Assurance Level (SAL)	10 <sup>-6</sup>	10 <sup>-6</sup>
Performance	Pressure Rating	1000 psi	1200 psi



**Conclusion:**

MEDRAD considers the MEDRAD Twist & Go High Pressure Connector Tubing to be substantially equivalent to the predicate device listed above. This conclusion is based upon the devices' similarities in functional design, materials, indications for use, fundamental scientific technology, and principles of operation.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

August 2, 2013

MEDRAD, Incorporated  
Bayer Medical Care, Incorporated  
C/O Ms. Lisa A. Ewing  
Principal Regulatory Affairs Specialist  
One Medrad Drive  
INDIANOLA PA 15051

Re: K131517  
Trade/Device Name: MEDRAD Twist & Go High Pressure Connector Tubing  
Regulation Number: 21 CFR 880.5440  
Regulation Name: Tubing, Fluid Delivery  
Regulatory Class: II  
Product Code: FPK  
Dated: July 8, 2013  
Received: July 9, 2013

Dear Ms. Ewing:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Mary S. Runner -S**

Kwame Ulmer, M.S.  
Acting Division Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number: K131517  
Device Name: MEDRAD Twist & Go High Pressure Connector Tubing

### Indications for Use:

The Twist & Go High Pressure Connector Tubing is intended to be used in the delivery of contrast media and common flushing solutions. The device is indicated for single use only with MEDRAD injectors.

Prescription Use X  
(21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



Richard C. Chapman

2013.07.31 10:31:46

-04'00'

(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

Page 1 of 1

510(k) Number: K131517