

K131770

FEB 26 2014

## 510(k) Summary

**Date of Submission:** Feb 25, 2014

**Owner:**

**Name:** Shenzhen ANT Hi-Tech Industrial Co., Ltd  
**Address:** Jinhui Ave.18, Pingshan New District, Shenzhen, Guangdong, 518122, China.  
No.46, Science and Technology Road, Yuquan Industrial Park,  
Fenggang Town, Dongguan City, 523690 P.R. China

**Submission Contact Person:**

**Name:** Ms. Lynn. Fu  
**eMail:** lynn.fu@antmed.com      lynn.fu@rendermed.com  
**Tel:** 86-755-26903311      **Fax:** 86-755-26903355  
**Company:** Shenzhen ANT Hi-Tech Industrial Co., Ltd  
**Address:** Jinhui Ave.18, Pingshan New District, Shenzhen, Guangdong, 518122, China.

**Proposed Device:**

**Trade name:** Pressure Connecting Tube  
**Common name:** High pressure tube, High Pressure Line, Pressure Connecting Line  
**Classification name:** Diagnostic intravascular catheter (21CFR, 870.1200, Product Code DQO)

**Classification:**

**Class:** II  
**Product Code:** DQO  
**Regulation Number:** 870.1200

**Predicate Device:**

K 071196  
Disposable High Pressure Injection Lines with and with out rotating adapters  
Coeur, Inc

**Intended Use:**

The Pressure Connecting Tube is indicated for use during angiography procedures as a connecting line for the injection of liquid, such as radiopaque dye, saline, or other diagnostic fluids.

**Device Description:**

The Pressure Connecting Tube is a sterile, single use device which can withstand injection

pressures to 350~1200psi for use with power injectors to inject contrast media into implanted ports designed for use with power injectors.

This device is designed, like other legally marked devices, for one end to connect to the fluid source (such as an angiographic syringe) and the other end to connect to the catheter. The fluid such as radiopaque dye, saline, or other diagnostic fluid is then injected from the syringe, through the connecting tube, into the catheter. The contrast, saline, or other diagnostic fluid is then injected from the syringe through the High Pressure Tube, into the catheter.

According to the different pressure maintenance, the material of tube is PVC or PU, and the appearance is clear or braided.

**Technological Characteristics**

The design and technological characteristics of the Pressure Connecting Tube is substantially equivalent to the predicate devices except the following:

The predicate device pressure rating is 500 to 1200psi while the proposed device can withstand injection pressure from 350 to 1200psi.

The difference does not have a significant impact on the safety or effectiveness of the device.

**Summary of Technological Characteristic compared to legally market device:**

<b>Characteristic</b>	<b>Proposed device</b>	<b>Legally market device (K071196)</b>
<b>Intended Use</b>	This Pressure Connecting Tube is for use during angiography procedures as a connecting line for the injection of liquid, such as radiopaque dye, saline, or other diagnostic fluids.	The legally market device is for use during coronary angiography procedure as a connecting line for the injection of radiopaque dye, saline, or other diagnostic fluids.
<b>Material</b>	Plastic polymers	same
<b>Flexible</b>	YES	same
<b>Connector Tube</b>	Multiple lengths of braided or unbraided Clear or tinted	same
<b>Components</b>	With or without luer connector, luer lock connector. Protection cap.	Combinations of male and female luers (including rotating male luers), with and without dust caps(protection cap.
<b>Packaging</b>	Sealed Tyvek-lidded package.	same
<b>Sterility</b>	EtO	same
<b>Disposable</b>	YES	same
<b>Pressure</b>	350 to 1200psi	500 to 1200psi

**Non-Clinical Testing Data:**

The proposed device including three pressure sustaining types: 350psi, 600psi and 1200psi. The most complicated structure of each pressure type is selected to conduct the performance test. The results of the testing verified that the proposed device is capable of the intended use. The specific test items are listed below:

**Summary of the non-Clinical Testing**

<b>Items</b>	<b>Reference Standards</b>
<b>Packaging</b>	ISO 11607-1 Packaging for terminally sterilized medical devices —Part 1: Requirements for materials, sterile barrier systems and packaging systems
	<ul style="list-style-type: none"> <li>• <b>Seal strength</b> ASTM F88 Standard Test Method for Seal Strength of Flexible Barrier Materials</li> <li>• <b>Integrity</b> ASTM F1929 Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration</li> <li>• <b>Vacuum leak</b> ASTM D3078 Standard Test Method for Determination of Leaks in Flexible Packaging by Bubble Emission</li> <li>• <b>Accelerated aging</b> ASTM F1980-07 Standard Guide For Accelerated Aging Of Sterile Barrier Systems For Medical Devices. (Sterility)</li> </ul>
<b>Performance</b>	ISO 594-1 Conical Fittings With A 6% (Luer) Taper For Syringes, Needles And Certain Other Medical Equipment - Part 1: General Requirements.
	ISO 594-2 Conical Fittings With A 6% (Luer) Taper For Syringes, Needles And Certain Other Medical Equipment - Part 2: Lock Fittings.
	Pressure Maintenance: Manufacturer's standard. The product should fulfill the requirements of pressure maintenance as claimed.
<b>Sterilization</b>	<p>ISO 11135-1 Sterilization Of Health Care Products - Ethylene Oxide - Part 1: Requirements For Development, Validation And Routine Control Of A Sterilization Process For Medical Devices.</p> <ul style="list-style-type: none"> <li>• <b>EO/ECH Residuals</b> ISO 10993-7 Biological Evaluation Of Medical Devices - Part 7: Ethylene Oxide Sterilization Residuals.</li> <li>• <b>Bacterial Endotoxin</b> USP Chapter&lt;85&gt; Bacterial Endotoxins Test.</li> <li>• <b>Sterility</b> USP Chapter&lt;71&gt; Sterility Tests.</li> </ul>

**Biocompatibility of materials:**

All the biocompatibilities of materials meet the requirements of ISO 10993 as listed below:

**Summary of the Biocompatibility Testing**

<b>Items</b>	<b>Reference Standards</b>
<b>Cytotoxicity</b>	ISO 10993-5 Biological Evaluation Of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity
<b>Sensitization</b>	ISO 10993-10 Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.
<b>Irritation</b>	
<b>System toxicity (acute)</b>	ISO 10993-11 Biological Evaluation Of Medical Devices Part 11: Tests For Systemic Toxicity

**Clinical Tests: NA.**

**Conclusion:**

The predicate device has the same classification, intended use, technology designing and sterile specification. The differences between the proposed device and the predicated device will not affect the safety and efficiency.



February 26, 2014

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Shenzhen ANT Hi-Tech Industrial Co., Ltd  
c/o Ms. Lynn Fu  
Jinhui Ave. 18  
Pingshan New District  
Shenzhen, Guangdong, 518122  
China

Re: K131770

Trade/Device Name: Pressure Connecting Tube  
Regulation Number: 21 CFR 870.1200  
Regulation Name: Diagnostic Intravascular Catheter  
Regulatory Class: Class II  
Product Code: DQO  
Dated: January 5, 2014  
Received: January 17, 2014

Dear Ms. Fu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for  
Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indication for Use

510(k) Number: K131770

Device Name: Pressure Connecting Tube

### Indications for Use:

The Pressure Connecting Tube is for use during angiography procedures as a connecting line for the injection of liquid, such as radiopaque dye, saline, or other diagnostic fluids.

Prescription Use   
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



Page 1 of 1