



510(k) Premarket Notification Submission

JAN 22 2014

**510(k) Summary**

In accordance with 21 CFR 807.92 the following summary of information is provided:

**Date:** 18 June, 2013  
**Submitter:** Beijing Honkon Technologies Co., Ltd.  
 Address: No.3 Building, No.11 Yard, Kangding Street, BDA, 100176, Beijing, P.R.China

**Primary Contact Person:** Mike Gu  
 Address: 7<sup>th</sup> floor, Jingui Business Building, 982 Congyun Road, Baiyun District, 510420, Guangzhou, China  
 OSMUNDA Medical Device Consulting Co., Ltd.  
 Tel: (+86) 20 6232 1333  
 Fax: (+86) 208633 0253

**Secondary Contact Person:** Zhao Li  
 Management representative  
 Beijing Honkon Technologies Co., Ltd.  
 Tel: (+86) 10-56370050  
 Fax: (+86) 10-56370076-610

**Device:** Trade Name: CO<sub>2</sub> Laser (10600nm)  
Common/Usual Name: CO<sub>2</sub> Laser  
Regulatory Number: 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology

**Product Code:** GEX

**Predicate Device(s):** K080496

**Device Description:** The CO<sub>2</sub> Laser consists of semiconductor solid-state lasers, switching power supply, laser power, treatment handle, key switch, cooling system and accessories. The CO<sub>2</sub> Laser produces a beam of coherent infrared light with 10600nm wavelength.

The CO<sub>2</sub> laser comprises five models: YILIYA-10600AH, YILIYA-10600AL, YILIYA-10600CH, Aeslight-10600EH, Aeslight-Ultra peel; their technical parameters and functions are the same, the only difference is their appearances and sizes.

**Intended Use:** The CO<sub>2</sub> Laser (10600nm) is indicated for use in dermatological procedures requiring ablation, resurfacing and coagulation of soft tissue.

Confidential and Privileged. This document contains confidential and privileged trade secrets and other information of Beijing Honkon Technologies Co., Ltd. and as such may not be disclosed to others not employed by Beijing Honkon Technologies Co., Ltd. All rights reserved

**HONKON**  
北京宏强富瑞技术有限公司

510(k) Premarket Notification Submission

**Technology:** The CO<sub>2</sub> laser therapy device emits a highly focused laser beam. Absorbed by water of the soft tissue, the CO<sub>2</sub> laser incorporating the advanced technologies can immediately ablate discrete columns of tissue without charring. This device can be used in multiple specialties, such as Oral Surgery, Plastic Surgery, Dermatology.

**Determination  
of Substantial  
Equivalence:**

<b>Specification</b>	<b>Predicate Device K080496</b>	<b>Proposed CO<sub>2</sub> Laser</b>
<i>Manufacturer</i>	Lutronic Corporation	Honkon
<i>Intended Use</i>	It is indicated for use in dermatological procedures requiring ablation (removal), resurfacing and coagulation of soft tissue.	The CO <sub>2</sub> Laser is indicated for use in dermatological procedures requiring ablation (removal), resurfacing and coagulation of soft tissue.
<i>Energy output</i>	2-240mj	2-200mj
<i>Laser transfer method</i>	Articulated Arm with Handpiece	Articulated Arm with Handpiece
<i>CO<sub>2</sub> RF Module Maximum Power</i>	Maximum 30w at continuous wave	Maximum 30w at continuous wave
<i>Wavelength (nm)</i>	10600	10600
<i>Pulse duration (ms)</i>	5-100	5-100
<i>Frequency (Hz)</i>	10-200Hz	10-200Hz
<i>Patient contacting materials</i>	Aluminum	Aluminum
<i>Compatibility with environment and other devices</i>	Comply with the IEC 60601-1-2	Comply with the IEC 60601-1-2
<i>Electrical Safety</i>	Comply with the IEC 60601-1	Comply with the IEC 60601-1
<i>Radiation Safety</i>	IEC 60825-1; IEC 60601-2-22	IEC 60825-1; IEC 60601-2-22
<i>Biocompatibility</i>	ISO 10993-5; ISO 10993-10	ISO 10993-5; ISO 10993-10

Confidential and Privileged. This document contains confidential and privileged trade secrets and other information of Beijing Honkon Technologies Co., Ltd. and as such may not be disclosed to others not employed by Beijing Honkon Technologies Co., Ltd. All rights reserved

**HONKON**  
北京宏强富瑞技术有限公司

510(k) Premarket Notification Submission

<i>Other: Device Specific Guidance Requirements for Comparison</i>	Guidance on the Content and Organization of a Premarket Notification for a Medical Laser; 21 CFR 1040.10	Guidance on the Content and Organization of a Premarket Notification for a Medical Laser; 21 CFR 1040.10
--	---	---

Conclusion: The CO<sub>2</sub> Laser and its application comply with standards as detailed in section 9, 11 and 17 of this premarket notification. Therefore, Beijing Honkon states that the non-clinical tests determined that the CO<sub>2</sub> Laser to be as safe, as effective and performance is substantially equivalent to the predicate device(s).



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Osmunda Medical Device Consulting Co., Ltd.  
Mike Gu  
Regulatory Affairs Manager  
7th floor, Jingui Business Building, 982 Congyun Road,  
Baiyun District, 510420, Guangzhou, China

January 22, 2014

Re: K131837

Trade/Device Name: CO2 Laser  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Regulatory Class: II  
Product Code: GEX  
Dated: January 15, 2014  
Received: January 17, 2014

Dear Mr. Gu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Joshua  Xipper -S

Dr. Binita Ashar  
Acting Director  
FOR Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Premarket Notification Submission

510(k) Number: K131837

Device Name: CO2 Laser (10600nm)

Indications for Use:

The CO2 Laser (10600nm) is indicated for use in dermatology procedures requiring ablation, resurfacing and coagulation of soft tissue.

Prescription Use  AND/OR  
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use   
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

---

Neil R. Ogden -S  
2014.01.21 13:10:40  
-05'00"

(Division Sign-Off) for BSA  
Division of Surgical Devices  
510(k) Number     K131837