

**510(k) SUMMARY****Name, Address, Phone and Fax Number of Applicant**

Teleflex Medical, Incorporated  
2917 Weck Drive  
Research Triangle Park, NC 27709 USA  
Phone: 919-433-8050  
Fax: 919-433-4996

**Contact Person**

Angela Bouse  
Senior Regulatory Affairs Specialist

**Date Prepared**

December 11, 2013

**Device Name**

Trade Name: ISO-Gard® Mask  
Classification Name: Apparatus, gas-scavenging  
Product Code: CBN  
Regulation Number: 868.5430  
Classification: II  
Classification Panel: Anesthesiology

**Predicate Device**

This submission demonstrates substantial equivalence to the predicate device  
ISO-GARD ClearAir Mask - K123176

**Device Description**

The ISO-Gard Mask system is an oxygen delivery mask that actively scavenges waste anesthetic gases (WAGS) exhaled by patients recovering from surgery in the Post-Anesthetic Care Unit (PACU). Vacuum/suction for scavenging of WAGS is provided by the institution's regulated vacuum source. The proposed device allows for the delivery of supplemental / therapeutic oxygen to patients to aid in their recovery while reducing the amount of patient expelled waste anesthetic agents released to the work environment of the healthcare workers. The mask can be used with or without suction / vacuum to function as a standard oxygen mask with an ET/CO<sub>2</sub> monitoring port.

**Indications for Use**

The ISO-Gard® Mask is intended to be used to scavenge waste anesthetic gases from patients during recovery from general anesthesia and to provide supplemental oxygen.

The ISO-Gard® Mask helps to reduce the amount of anesthetic agents released to the work environment of the healthcare worker.

**Patient Population**

Patients recovering from general anesthesia in the PACU.

**Environments of use**

The environment of use is – Post-operative Care Units (PACU) in hospital, sub-acute facilities.

**Contraindications**

None

**Substantial Equivalence**

The proposed device is substantially equivalent to the predicate devices:

<b>Comparative Characteristics</b>	<b>Predicate K123176 ISO-GARD® ClearAir™ Mask</b>	<b>Proposed ISO-Gard® Mask</b>
<b>Classification Name</b>	Apparatus, gas scavenging	Same
<b>Product Code / CFR</b>	CBN 868.5430  Secondary CCK – Gaseous-Phase Carbon Dioxide Gas Analyzer 868.1400	Same
<b>Indications for Use</b>	The ISO-GARD® ClearAir™ Mask is intended to be used to scavenge waste anesthetic gases from patients during recovery from general anesthesia and to provide supplemental oxygen.  The ISO-GARD® ClearAir™ Mask helps to reduce the amount of anesthetic agents released to the work environment of the healthcare worker.	Same
<b>Trade Name</b>	ISO-GARD® ClearAir™ Mask	ISO-Gard® Mask
<b>Environment of Use</b>	Hospital, sub-acute facilities PACU	Same
<b>Patient Population</b>	Patients recovering from general anesthesia	Same

	and may need supplemental oxygen Adults	
<b>Contraindications</b>	None	Same
<b>Basic Components</b>	Mask Oxygen delivery tubing Vacuum (scavenging) tubing Mask Manifold controlling oxygen delivery and scavenging	Same
<b>Design, Features, and Specifications</b>		
<b>Mask</b>	Flexible oxygen mask with sealing foam	Same
<b>Method to hold mask on patient for seal</b>	Elastic band / strap	Same
<b>Tubing to deliver oxygen</b>	Standard oxygen tubing	Same
<b>Connects to ETCO2 monitor</b>	Yes	Same
<b>Connector to sampling line</b>	Standard female luer lock	Same
<b>Method of separating gas flows</b>	Divided manifold for separating vacuum and oxygen delivery and then a separate sampling port	Same
<b>Safety features</b>		
<b>Excess negative pressure</b>	Contains entrainment valves if the negative pressure from vacuum is too great Valves are one-way flapper/diaphragm valves that open with minimal negative pressure or flow	Same
<b>Excess positive pressure</b>	Contains entrainment valves if patient's inhalation is greater than the supply of the oxygen	Same
<b>Method to assist in sealing</b>	Foam pad around bridge of nose to assist in sealing of the mask	Same
<b>Method to separate oxygen delivery from scavenging</b>	Mask manifold body is a divided adapter which has an oxygen inlet and a scavenging outlet	Same
<b>Oxygen source</b>	Wall oxygen	Same
<b>Vacuum source</b>	Wall vacuum	Same
<b>Port for sampling end tidal CO2</b>	Port connector on exhalation side of Mask Manifold adapter	Same
<b>Typical oxygen delivered flow rates</b>	Up to 10 lpm	Same
<b>Oxygen at various Oxygen flow rates and Vacuum setting</b>	Delivered oxygen equal or greater than oxygen concentration mask	Same

<b>Mask sizes</b>	Adult	Same
<b>Performance Standards</b>	None	Same
<b>Shelf Life</b>	No shelf life	1 year shelf life
<b>Patient Contacting Materials</b>		
<b>Mask</b>	PVC	Same
<b>Star-Lumen Oxygen Tube</b>	PVC	Same
<b>Connector, Star-Lumen Oxygen Tube</b>	PVC	Same
<b>Gasketing Foam w/Adhesive</b>	Natural ester foam with acrylic pressure sensitive adhesive	Same
<b>Tethered Cap</b>	Thermoplastic Elastomer	Same
<b>One-way inhalation valves</b>	Polyisoprene	Same
<b>One-way Valve Body</b>	Polystyrene, Trans Blue	Same
<b>Oxygen Delivery Port Adaptor</b>	Polypropylene	Same
<b>White Elastic Strap</b>	White Polyester/ Polyisoprene	Same
<b>Mask Manifold</b>	Polystyrene	Same
<b>Suction/Exhalation Port</b>	Polystyrene	Same
<b>Oxygen Port Concentrator</b>	PVC	Same

**Comparison to Predicate Device**

The proposed ISO-Gard Mask is substantially equivalent to the predicate device with respect to indications for use, technology, materials, and construction. The proposed change is to add a one year shelf life to the ISO-Gard Mask labeling.

- **Indications for Use –**  
The indications for use are identical to the predicate.
- **Technology and construction -**  
The proposed device design, drawings, components, accessories, materials, packaging and product configurations remain unchanged. The proposed change is to add a one year shelf life to the ISO-Gard Mask labeling.
- **Environment of use –**  
The environment of use is identical to the predicate.

- **Patient Population -**  
The patient population is identical to the predicate.
- **Materials -**  
All patient contacting materials are identical to the predicate.
- **Performance Testing**

**Nonclinical Performance Testing Summary**

<b>Test</b>	<b>Test Objective</b>	<b>Acceptance Criteria</b>
Oxygen Delivery	To evaluate the oxygen delivery performance at variable oxygen flow rates and vacuum levels at standard Tidal Volumes of 500 ml without the use of N <sub>2</sub> O	The delivered oxygen percentage using the ISO-Gard Mask must be equal to or greater than a standard medium concentration oxygen mask for all vacuum settings
Scavenging	To evaluate the scavenging performance at variable oxygen flow rates and vacuum levels at standard Tidal Volumes of 500 ml with N <sub>2</sub> O	N <sub>2</sub> O levels must be lower than with a standard medium concentration oxygen mask
ETCO <sub>2</sub>	To evaluate the ETCO <sub>2</sub> performance in simulated conditions	The traces/waveforms during testing must be distinct and generated consistently

**Conclusion**

The ISO-Gard Mask has the same indications for use, technological characteristics, and constructions as the predicate. Performance test results demonstrate that the proposed device does not raise new questions of safety and effectiveness and because an acceptance criteria has been met, the device can be found substantially equivalent.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

January 9, 2014

Teleflex Medical, Inc.  
Angela Bouse  
Senior Regulatory Affairs Specialist  
2917 Weck Drive  
Research Triangle Park, NC 27709

Re: K132729

Trade Name: ISO-Gard® Mask  
Regulation Number: 21 CFR 868.5430  
Regulation Name: Gas-scavenging apparatus  
Regulatory Class: Class II  
Product Code: CBN  
Dated: December 11, 2013  
Received: December 12, 2013

Dear Ms. Bouse:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
Tejasri Purohit-Sheth, M.D.  
Clinical Deputy Director  
DAGRID

FOR

Erin Keith, M.S.  
Acting Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Indications for Use Statement**

Page 1 of 1

**510(k) Number:** \_\_\_\_\_ (To be assigned)

**Device Name:** ISO-Gard® Mask

**Indications for Use:**

The ISO-Gard® Mask is intended to be used to scavenge waste anesthetic gases from patients during recovery from general anesthesia and to provide supplemental oxygen.

The ISO-Gard® Mask helps to reduce the amount of anesthetic agents released to the work environment of the healthcare worker.

**Prescription Use XX**  
(Part 21 CFR 801 Subpart D)

or

**Over-the-counter use** \_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Nayan J. Patel, S  
2014.01.07 15:25:50 -05'00'

K132729