

**510(k) SUMMARY**

JUN 04 2014

**A. Name, Address, Phone and Fax Number of Applicant**

Teleflex Medical, Incorporated  
2917 Weck Drive  
Research Triangle Park, NC 27709 USA  
Phone: 919-433-8050  
Fax: 919-433-4996

**B. Contact Person**

Angela Bouse  
Senior Regulatory Affairs Specialist

**C. Date Prepared**

June 3, 2014

**D. Device Name**

Trade Name:	SOFTECH® Plus ETCO <sub>2</sub> Cannula
Classification Name:	Analyzer, Gas, Carbon-Dioxide, Gaseous-Phase
Product Code:	CCK
Additional Product Code:	CAT
Regulation Number:	868.1400
Classification:	II
Classification Panel:	Anesthesiology

**E. Predicate Device**

This submission demonstrates substantial equivalence to the predicate device Double Lumen Oxygen Delivery Gas Sampling Nasal Cannula, Hudson RCI Bi-Flo – K961150.

**F. Device Description**

The SOFTECH Plus ETCO<sub>2</sub> Cannula is a non-sterile disposable, single patient use device that acts as an adjunct to oxygen therapy with its primary function providing a means to deliver low flow oxygen, while sampling part of the patients exhaled gas. The SOFTECH Plus ETCO<sub>2</sub> Cannula has a split nare blank with oxygen delivery through one nasal prongs while allowing sampling of the patient's exhaled gas from the corresponding nasal prongs.

**G. Indications for Use**

The Hudson RCI SOFTECH® Plus ETCO<sub>2</sub> Nasal Cannula is an adjunct to oxygen therapy with its primary function being that of delivering low flow oxygen to a patient while providing a means to sample expired gas. It is intended for use in patients requiring oxygen therapy to improve blood oxygen levels while monitoring expired gas to determine ventilatory rate.

Patient Population: Infant, Pediatric, Adult

**H. Technological Characteristics and Material Comparison to the predicate**

The proposed SOFTECH Plus ETCO<sub>2</sub> Cannula is substantially equivalent to the predicate device listed above in that the indications for use, the intended use, and fundamental scientific technology remain unchanged. **Table 1** summarizes the differences between the proposed and predicate devices.

**Table 1 - Differences Between the Proposed and Predicate Devices**

Features	Predicate K961150 Double Lumen Oxygen Delivery Gas Sampling Nasal Cannula, Hudson RCI Bi-Flo	Proposed SOFTECH® Plus ETCO <sub>2</sub> Cannula	Performance Testing
Nasal Prongs (Nares)	Dual-channel nasal prong that allows oxygen delivery and gas sampling from both nares	Split nasal prong that directs oxygen flow into one nare and samples expired gas from the other nare	ETCO <sub>2</sub> Performance Testing with Simultaneous Oxygen Delivery
Oxygen supply tubing and CO <sub>2</sub> sampling line length	Oxygen supply tubing: 7 ft CO <sub>2</sub> sampling line: 7 in	Oxygen supply tubing: 7 ft and 14 ft CO <sub>2</sub> sampling line: 2 in, 7 ft and 14 ft	ETCO <sub>2</sub> Performance Testing with Simultaneous Oxygen Delivery
<b>Materials</b>			
Nares	Polyvinylchloride	Polyvinylchloride	Biocompatibility Testing: Cytotoxicity, Sensitization, Irritation, Genotoxicity, Implantation
Bolo Slide	Polycarbonate or Polyvinylchloride	Polyvinylchloride	
Tubing Clip	None	Polypropylene	
Lariat Tubing	Polyvinylchloride	Polyvinylchloride	
Oxygen Supply Tubing	Polyvinylchloride	Polyvinylchloride	
CO <sub>2</sub> Sampling Line	Polyvinylchloride	Polyvinylchloride	
Oxygen Connector	Polyvinylchloride	Polyvinylchloride	
Tubing Connector	Polyvinylchloride	Polyvinylchloride	
Male and Female Luer Lock Connectors	Acrylic	Polycarbonate	
22mm Oxygen	None	Polypropylene	

<b>Tubing Adaptor</b>			
<b>Adhesive</b>	Cyclohexanone	Cyclohexanone	
<b>PVC</b>	PVC with DEHP	Non-DEHP PVC	DEHP Testing

### I. Performance Data

The proposed device was tested for ETCO<sub>2</sub> sampling at various oxygen flow rates, biocompatibility, and DEHP testing. The test results demonstrate that the device is substantially equivalent to the predicate device.

The tests performed are summarized in **Table 2** below.

**Table 2 – Performance Testing Summary**

<b>General Description</b>
ETCO <sub>2</sub> performance in simulated conditions at various oxygen flow rates
Biocompatibility testing as per ISO 10993-1: cytotoxicity, sensitization, irritation, genotoxicity and implantation testing
DEHP Testing

### J. Conclusion

The device data and test results demonstrate that the device is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

June 4, 2014

Teleflex Medical, Inc.  
Angela Bouse  
Senior Regulatory Affairs Specialist  
2917 Weck Drive  
Research Triangle Park, NC 27709

Re: K132946  
Trade/Device Name: SOFTECH® Plus ETCO2 Cannula  
Regulation Number: 21 CFR 868.1400  
Regulation Name: Carbon Dioxide Gas Analyzer  
Regulatory Class: II  
Product Code: CCK, CAT  
Dated: May 02, 2014  
Received: May 05, 2014

Dear Ms. Bouse:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary  Bunner -S

Erin I. Keith, M.S.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K132946

Device Name  
Hudson RCI SOFTECH® Plus ETCO2 Nasal Cannula

Indications for Use (Describe)

The Hudson RCI SOFTECH® Plus ETCO2 Nasal Cannula is an adjunct to oxygen therapy with its primary function being that of delivering low flow oxygen to a patient while providing a means to sample expired gas. It is intended for use in patients requiring oxygen therapy to improve blood oxygen levels while monitoring expired gas to determine ventilatory rate.

Patient Population: Infant, Pediatric, Adult

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Todd D. Courtney -S

2014.06.03 16:33:20 -04:00



This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*