



## 510(k) SUMMARY

NOV 19 2013

## V.A.C.® Therapy Wound Dressings

<b>Submitter Information [21 CFR 807.929(a)(1)]</b>	
<b>Name</b>	KCI USA, Inc. (Kinetic Concepts, Inc.)
<b>Address</b>	6203 Farinon Drive San Antonio, TX 78249
<b>Phone number</b>	210-515-4059
<b>Fax number</b>	210-255-6727
<b>Establishment Registration Number</b>	1625774
<b>Name of contact person</b>	Melanie Avila
<b>Date prepared</b>	10/23/2013
<b>Name of the device [21 CFR 807.92(a)(2)]</b>	
<b>Trade or proprietary name</b>	V.A.C.® Negative Pressure Wound Therapy System
<b>Common or usual name</b>	Negative Pressure Wound Therapy System
<b>Classification name</b>	Negative Pressure Wound Therapy Powered Suction Pump (and components)
<b>Regulation</b>	878.4780
<b>Product Code(s)</b>	OMP
<b>Legally marketed device(s) to which equivalence is claimed [21 CFR 807.92(a)(3)]</b>	V.A.C. Therapy Wound Dressings cleared for use under multiple 510(k)s for the KCI V.A.C. Therapy Negative Pressure Wound Therapy Systems. The most recent 510(k) was K120033.
<b>Device description [21 CFR 807.92(a)(4)]</b>	V.A.C. Therapy Wound Dressings for use with KCI's V.A.C. NPWT System.
<b>Indications for use [21 CFR 807.92(a)(5)]</b>	<p>The ActiV.A.C., InfoV.A.C., V.A.C. ATS, V.A.C. Freedom, V.A.C. Via, and V.A.C. Simplicity Negative Pressure Wound Therapy Systems are integrated wound management systems for use in acute, extended and home care settings.</p> <p>When used on open wounds, they are intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. Open wound types include: chronic, acute, traumatic, subacute and dehisced wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.</p> <p>When used on closed surgical incisions, they are also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudates via the application of negative pressure wound therapy.</p>



## 510(k) SUMMARY

### V.A.C.® Therapy Wound Dressings

Summary of the technological characteristics of the device compared to the predicate device [21 CFR 807.92(a)(6)]		
Characteristic	V.A.C.® Therapy Wound Dressing System (Modified Device)	V.A.C.® Therapy Wound Dressing System (Predicate)
Indicated Wound Type	Same as predicate	Chronic, acute, traumatic, subacute and dehisced wounds; partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts
Dressing	Same as predicate	Multiple dressing components
Performance Data [21 CFR 807.92(b)]		
Summary of non-clinical tests conducted for determination of substantial equivalence [21 CFR 807.92(b)(1)]		
<p>The V.A.C.® Dressing System was evaluated to ensure conformance to the design specifications. The following tests were conducted:</p> <ul style="list-style-type: none"> <li>• Peel adhesion testing</li> <li>• Moisture vapor transmission rate (MVTR) testing</li> <li>• Biocompatibility testing according to ISO10993-1</li> <li>• Study conducted on healthy human volunteers, under design validation (210 CFR 820.30(g)), to ensure that the proposed modified device meets user requirements</li> </ul>		
Summary of clinical tests conducted for determination of substantial equivalence or of clinical information. [21 CFR 807.92(b)(2)]		
<p>There was no clinical study conducted, however, a study was conducted on healthy human volunteers, under design validation (210 CFR 820.30(g)), to ensure that the proposed modified device meets user requirements as there is no reliable method for measuring these requirements on the bench.</p>		
Conclusions drawn [21 CFR 807.92(b)(3)]		
<p>Testing indicates that the modified V.A.C.® Therapy Wound Dressing System is substantially equivalent in terms of both indications for use and fundamental scientific technology to the predicate product.</p>		



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

KCI USA, Incorporated  
Ms. Melanie Avila  
Regulatory Affairs Project Manager  
6203 Fairmon Drive  
San Antonio, Texas 78249

November 19, 2013

Re: K133276

Trade/Device Name: V.A.C.<sup>®</sup> Negative Pressure Wound Therapy System  
Regulation Number: 21 CFR 878.4780  
Regulation Name: Powered suction pump  
Regulatory Class: Class II  
Product Code: OMP  
Dated: October 23, 2013  
Received: October 24, 2013

Dear Ms. Avila:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 – Ms. Melanie Avila

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Joshua C. Nipper -S**

FOR

Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

### INDICATIONS FOR USE

510(k) Number (if known): K133276

Device Name: V.A.C.® Drape

Indications for Use:

The ActV.A.C., InfoV.A.C., V.A.C. ATS, V.A.C. Freedom, V.A.C. Via, and V.A.C. Simplicity Negative Pressure Wound Therapy Systems are integrated wound management systems for use in acute, extended and home care settings.

When used on open wounds, they are intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. Open wound types include: chronic, acute, traumatic, subacute and dehisced wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.

When used on closed surgical incisions, they are also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudates via the application of negative pressure wound therapy.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page   1   of   1  

(Posted November 13, 2003)

# Jiyoung Dang -S