

**510(k) Summary**

per 21 CFR §807.92

**DEC 20 2013**

<b>Submitter's Name and Address</b>	Boston Scientific Corporation One Scimed Place Maple Grove, MN 55311 Phone: 763-255-0877 Fax: 763-494-2222
<b>Contact Name and Information</b>	Maylin Truesdell Senior Regulatory Affairs Specialist Phone: 763-255-0877 Fax: 763-494-2222 e-mail: maylin.truesdell@bsci.com
<b>Date Prepared</b>	November 19, 2013
<b>Proprietary Name</b>	Rotablator™ Rotational Atherectomy System with the Peripheral RotaLink™ Plus
<b>Common Name</b>	Rotational Angioplasty System or Rotational Atherectomy System
<b>Product Code</b>	MCW – Catheter, Peripheral, Atherectomy
<b>Classification</b>	Class II, 21 CFR Part 870.4875 – Intraluminal Artery Stripper
<b>Predicate Devices</b>	Rotablator® Rotational Angioplasty System K901206, September 14, 1990  Rotablator™ Rotational Atherectomy System with the Peripheral RotaLink™ Plus K121774, September 13, 2012
<b>Device Description</b>	The Rotablator Rotational Atherectomy System with the Peripheral RotaLink Plus consists of an Advancer pre-connected to a Catheter. The advancer functions as a guide for the sliding elements that control burr advancement and the catheter portion of the device guides the burr through the vasculature to the treatment site. The Peripheral RotaLink Plus devices are provided sterile and non-pyrogenic. It is intended for one procedure use only. The Rotablator Rotational Atherectomy System with the Peripheral RotaLink Plus currently offers burr sizes 1.5, 1.75, 2.00, 2.15, 2.25, 2.38 and 2.5mm, with 1.25mm being introduced in this submission.
<b>Intended Use</b>	The Rotablator Rotational Atherectomy System with the Peripheral RotaLink Plus is intended to ablate occlusive material and restore luminal patency in the peripheral vasculature.
<b>Indications for Use</b>	The Rotablator Rotational Atherectomy System is intended for percutaneous use in peripheral vessels in patients with occlusive atherosclerotic disease who are acceptable candidates for endovascular procedures.

**Comparison of Technological Characteristics**

The Rotablator Rotational Atherectomy System with the Peripheral RotaLink Plus incorporates substantially equivalent device materials and design, packaging materials and design, fundamental technology, manufacturing processes, sterilization process and intended use as those featured in the Boston Scientific predicate devices, Peripheral Rotablator Rotational Angioplasty System with the RotaLink Exchangeable Catheter.

**Performance Data**

The following in-vitro performance tests were completed on the Peripheral RotaLink Plus:

Strain Relief	Burr Cutting Ability
Operational Speeds	Tensile Strength
Stall Torque	Brake Engagement
Infusate Temperature Generation	Infusate Flow Rate
Catheter Advancement	Lumen Patency
Component Interface Compatibility	Functional Life

The following biocompatibility and chemical characterization tests were completed on the Peripheral RotaLink Plus:

Natural Rubber Latex	Intracutaneous Reactivity Test (Irritation)
Hemolysis Assay: Extract Method	Acute Systemic Injection Test
Hemolysis Assay: Direct Contact Method	Materials Mediated Rabbit Pyrogen Test
Complement Activation C3a and SC5b-9 Assay	USP Physicochemical Test for Plastics
Partial Thromboplastin Time (PTT)	In vitro Cytotoxicity Test: MEM Elution
In vitro Hemocompatibility Assay	FTIR Analysis
Guinea Pig Maximization Sensitization Test: Method for Biomaterial Extracts	

**Conclusion**

Based on the indications for use, technological characteristics, and safety and performance testing, the Peripheral Rotablator Atherectomy System with the Peripheral RotaLink Plus has been shown to be appropriate for its intended use and is considered to be substantially equivalent to the Peripheral Rotablator Rotational Angioplasty System as submitted in K901206 and K121774.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

December 20, 2013

Boston Scientific Corporation  
Ms. Maylin Truesdell  
Senior Regulatory Affairs Specialist  
One Scimed Place  
Maple Grove, MN 55311-1566

Re: K133566

Trade/Device Name: Rotablator Rotational Atherectomy System with the Peripheral RotaLink Plus  
Regulation Number: 21 CFR 870.4875  
Regulation Name: Intraluminal Artery Stripper  
Regulatory Class: Class II  
Product Code: MCW  
Dated: November 19, 2013  
Received: November 20, 2013

Dear Ms. Truesdell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. However, we remind you that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

The seal of the U.S. Food and Drug Administration (FDA) is positioned behind the signature. It features a central shield with a caduceus, surrounded by the words "U.S. FOOD AND DRUG ADMINISTRATION" and "DEPARTMENT OF HEALTH AND HUMAN SERVICES".  
**Kenneth J. Cavanaugh -S**  
for  
Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known):           K133566          

Device Name: Rotablator™ Rotational Atherectomy System with the Peripheral RotaLink™ Plus

### Indications for Use:

The Rotablator Rotational Atherectomy System is intended for percutaneous use in peripheral vessels in patients with occlusive atherosclerotic disease who are acceptable candidates for endovascular procedures.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

 Kenneth J. Cavanaugh -S