

Section 8 – 510(k) Summary or 510(k) Statement

I. General Information

Submitter: Alma Lasers, Ltd,
Halamish St. POB 3021
Caesarea Industrial Park, Caesarea, Israel 38900

Contact Person: Kathy Maynor
Consultant
352-586-3113 (cell)

Summary Preparation Date: Dec 30, 2013

II. Names

Device Names: The Alma Diode Tabletop Laser

Primary Classification Names: Surgical Powered Light Instrument,

III. Predicate Devices

K #	Predicate Device
K100058	Quanta Diode Laser Family

IV. Product Description

The Alma Diode Tabletop Laser is comprised of the following major components:

1. The main console unit
2. Footswitch.
3. Accessories

V. Indications for Use

The Alma diode tabletop laser is intended for use in dermatologic and general surgical procedures.

VI. Summary of Technical Characteristics

Table 1: Salient Characteristics of the 1470nm module and the Predicate Devices

	K13 Alma diode tabletop laser	K100558 Quanta Diode Laser Family
Parameter		
Product Code & Regulation No.	GEX 21 CFR 878.4810	GEX 21 CFR 878.4810
Laser Wavelength [nm]	1470	1470
Max power	15W	15W
Light/Laser Source	Diode	Diode
Laser Delivery	Optical Fiber	Optical Fiber
Operation Mode	Continuous wave, single pulse, pulsed	Continuous wave, single pulse, pulsed
Pulse Duration	10-990ms	3ms – 2.5s
Bare fiber size	200, 300, 320, 400, 600, 800, 1000	200, 300, 320, 400, 600, 800, 1000
User Interface	LCD touch screen	LCD touch screen
Aiming beam	635nm	650nm
Electrical Requirements	100-240, V AC 50-60 Hz, 6.3 A,	100-240, V AC 50-60 Hz, 6.3 A, single phase,
Indications for Use	The Alma 1470nm diode tabletop laser is indicated for use in endovenous occlusion of the greater saphenous vein in patients with superficial vein reflux. The Alma 1470nm diode tabletop laser is further indicated for laser assisted lipolysis.	The QUANTA Diode Laser System is intended for delivery of laser light to soft tissue in the contact and non contact mode during surgical procedures including via endoscopes. The QUANTA Diode Laser System is generally indicated for use in endovenous occlusion of the greater saphenous vein in Patients with Superficial Vein Reflux. The QUANTA1470 Diode Laser is further indicated for laser assisted lipolysis.

Table 2: Salient Characteristics of the 980nm module and the Predicate Devices

	K13 Alma diode tabletop laser	K100558 Quanta Diode Laser Family
Parameter		
Product Code & Regulation No.	GEX 21 CFR 878.4810	GEX 21 CFR 878.4810
Laser Wavelength [nm]	980nm	980nm
Max power	30W	30W
Light/Laser Source	Diode	Diode
Laser Delivery	Optical Fiber	Optical Fiber
Operation Mode	Continuous wave, single pulse, pulsed	Continuous wave, single pulse, pulsed
Pulse Duration	10-990ms	3ms – 2.5s
Bare fiber size	200, 300, 320, 400, 600, 800, 1000	200, 300, 320, 400, 600, 800, 1000
User Interface	LCD touch screen	LCD touch screen
Aiming beam	635nm	650nm
Electrical Requirements	100-240, V AC 50-60 Hz, 6.3 A.	100-240, V AC 50-60 Hz, 6.3 A, single phase.
Indications for Use	The Alma 980nm diode tabletop laser is indicated for use in endovenous occlusion of the greater saphenous vein in patients with superficial vein reflux. The Alma 980nm diode tabletop laser is further indicated for laser assisted lipolysis,	The QUANTA Diode Laser System is intended for delivery of laser light to soft tissue in the contact and non contact mode during surgical procedures including via endoscopes. The QUANTA Diode Laser System is generally indicated for use in endovenous occlusion of the greater saphenous vein in Patients with Superficial Vein Reflux. The QUANTA980 Diode Laser is further indicated for laser assisted lipolysis. (This is a subset of the cleared indications for this product).

Table 3: Salient Characteristics of 810nm diode module and the predicate devices

	K13 Alma diode tabletop laser	K100558 Quanta Diode Laser Family
Parameter		
Product Code & Regulation No.	GEX 21 CFR 878.4810	GEX 21 CFR 878.4810
Laser Wavelength [nm]	810nm	808nm
Max power	30W	30W
Light/Laser Source	Diode	Diode
Laser Delivery	Optical Fiber	Optical Fiber
Operation Mode	Continuous wave, single pulse, pulsed	Continuous wave, single pulse, pulsed
Pulse Duration	10-990ms	3ms – 2.5s
Bare fiber size	200, 300, 320, 400, 600, 800, 1000	200, 300, 320, 400, 600, 800, 1000
User Interface	LCD touch screen	LCD touch screen
Aiming beam	635nm	650nm
Electrical Requirements	100-240, V AC 50-60 Hz, 6.3 A.	100-240, V AC 50-60 Hz, 6.3 A, single phase,
Indications for Use	The Alma 810nm diode tabletop laser is indicated for endoluminal or endovenous laser surgery for saphenous incompetent veins.	The QUANTA Diode Laser System is intended for delivery of laser light to soft tissue in the contact and non contact mode during surgical procedures including via endoscopes. The QUANTA 808 Diode Laser System is generally indicated for endoluminal or endovenous laser surgery for saphenous incompetent veins. (This is a subset of the cleared indications for this product).

VII. Safety and Effectiveness Information

The review of the indications for use and technical characteristics provided demonstrates that the Alma diode tabletop laser is substantially equivalent to the predicate devices.

The Diode Tabletop Laser was tested by a certified laboratory according to:

IEC 60601-1: 1988+A1:1991+A2:1995: Medical Electrical Equipment Part 1: General Requirements for basic safety and essential performance

IEC 60601-1-2: Medical electrical equipment: Part 1-2: General requirements for safety- Collateral Standard: Electromagnetic compatibility - Requirements and tests (2001 + A1(4))

IEC 60601-2-22:1995 Medical electrical equipment- Part 2-22: Particular requirements for the safety of diagnostic and therapeutic laser equipment 1995.

IEC 60825-1:2007 (2nd edition) Safety of Laser Products-Part 1: Equipment Classification and Requirements

VIII. Conclusion

The Alma diode tabletop laser was found to be substantially equivalent to the predicate devices.

The Alma diode tabletop laser shares the same or similar indications for use, similar design features, and functional features with, and thus is substantially equivalent to the predicate devices.



July 18, 2014

Alma Lasers Ltd.
% Ms. Kathy Maynor
Regulatory Consultant for Alma LTD
26 Rebecca Court
Homosassa, Florida 34446

Re: K140005
Trade/Device Name: The Alma Diode Tabletop Laser
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general
and plastic surgery and in dermatology
Regulatory Class: Class II
Product Code: GEX
Dated: April 29, 2014
Received: June 20, 2014

Dear Ms. Maynor:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

David Krause -S

for **Binita S. Ashar, M.D., M.B.A., F.A.C.S.**
Director
Division of Surgical Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K140005

Device Name
The Alma Diode Tabletop Laser

Indications for Use (Describe)

Intended Use

The Alma Diode Tabletop Laser is intended for use in dermatologic and general surgical procedures.

Indications for Use

The Alma Diode Tabletop Laser includes three possible diode laser modules depending on the customer order.

Diode Laser Modules:

The indications for use for the 810nm Alma Diode Tabletop Laser include:

- The Alma 810nm diode tabletop laser is indicated for endoluminal or endovenous laser surgery for saphenous incompetent veins.

The indications for use for the 980 nm Alma Diode Tabletop Laser include:

-The Alma 980nm diode tabletop laser is indicated for use in endovenous occlusion of the greater saphenous vein in patients with superficial vein reflux. The Alma 980nm diode tabletop laser is further indicated for laser assisted lipolysis.

The indications for use for the 1470nm Alma Diode Tabletop Laser include:

-The Alma 1470nm diode tabletop laser is indicated for use in endovenous occlusion of the greater saphenous vein in patients with superficial vein reflux. The Alma 1470nm diode tabletop laser is further indicated for laser assisted lipolysis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Neil R Ogden -S

2014.07.17 16:13:51 -04'00'