



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

November 12, 2014

Mr. Patrick Lim  
Manager  
HIOSSSEN, Incorporated  
85 Ben Fairless Drive  
Fairless Hills, PA 19030

Re: K140934  
Trade/Device Name: HIOSSSEN Implant System  
Regulation Number: 21 CFR 872.3640  
Regulation Name: Endosseous Dental Implant  
Regulatory Class: II  
Product Code: DZE, NHA  
Dated: October 6, 2014  
Received: October 10, 2014

Dear Mr. Lim:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Susan Runna DDS, MA". The signature is written in a cursive style. Behind the signature, there is a faint, semi-transparent watermark of the FDA logo.

Erin I. Keith, M.S.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure





**Hiossen Inc.**

85 Ben Fairless Dr. Fairless Hills, PA 19030  
Tel : 1-888-678-0001 / Fax : 1-267-759-7004  
www.hiossen.com

**K140934**

**510(k) Summary**

This summary of 510(k) substantial equivalence information is being submitted in accordance with requirements of 21 CFR Part 807.92.

Date : November 15, 2013

1. Company and Correspondent making the submission:

- Submitter's Name : HiOSSEN Inc.
- Address : 85 Ben Fairless Dr.  
Fairless Hills PA 19030
- Telephone No. 888 678 0001
- Contact : Mr. Patrick Lim

2. Device :

- Trade or (Proprietary) Name : HIOSSEN Implant System
- Common or usual name : Dental Implant
- Classification Name : Endosseous Dental Implant  
21CFR872.3640  
Class II  
DZE, NHA

3. Predicate Device :

- HT III SA FIXTURE SYSTEM, HIOSSEN Inc, K101096
- ETII SA Fixture System, HiOSSEN Inc., K123471
- ETIII SA Ultra Wide System, HiOSSEN Inc., K103537
- TS Implant System, OSSTEM CO., LTD., K121585
- HSII Short Fixture System, OSSTEM CO., LTD., K083633
- ET SS Implant System, OSSTEM CO., LTD., K120847
- 3i OSSEOTITE® Certain® Dental Implants, Implant Innovations, Inc., K063341

4. Description :

- 1) The HIOSSEN Implant System is similar to other commercially available products based on the intended use, the technology used, the claims, the material composition employed and performance characteristics.
- 2) The HIOSSEN Implant System include length 18mm implant  
Length 18mm implant has already predicated in K063341, 3i OSSEOTITE® Certain® Dental Implants( Length 7mm~20mm)

**Hiossen Inc.**

85 Ben Fairless Dr. Fairless Hills, PA 19030  
 Tel : 1-888-678-0001 / Fax : 1-267-759-7004  
 www.hiossen.com

3) ETIII SA Fixture

<b>Device Description</b>	Intended to be surgically placed in the bone of the upper or lower jaw arches. Fixture is supplied sterile.		
<b>Material</b>	Pure Titanium (ASTM F 67)		
<b>Surface</b>	SA surface treatment		
	<b>Diameter (mm)</b>	<b>Length (mm)</b>	<b>added/modified</b>
<b>Dimension</b>	3.77	8.7	Diameter is modified
	3.75	10.2, 11.7, 13.2, 15.2, 18.2	Length 18.2mm is added
	4.25	7.2, 8.7, 10.2, 11.7, 13.2, 15.2, 18.2	Length 18.2mm is added
	4.65	7.2	None
	4.63	8.7	Diameter is modified
	4.6	10.2, 11.7, 13.2, 15.2, 18.2	Length 18.2mm is added
	5.1	6.2, 7.2	Length 6.2mm is added
	5.08	8.7	None
	5.05	10.2, 11.7, 13.2, 15.2	None

4) ETIII SA Ultra-Wide Fixture

<b>Device Description</b>	Intended to be surgically placed in the bone of the upper or lower jaw arches. Fixture is supplied sterile.		
<b>Material</b>	Pure Titanium (ASTM F 67)		
<b>Surface</b>	SA surface treatment		
	<b>Diameter (mm)</b>	<b>Length (mm)</b>	<b>added/modified</b>
<b>Dimension</b>	5.95	6.2, 9.7	Diameter is modified/ Length 6.2mm is added
	5.92	11.2, 12.7	Diameter is modified
	6	7.2, 8.2	None
	6.8	6.2, 7.2, 8.2, 9.7, 11.2, 12.7	Diameter is modified/ Length 6.2mm is added

5) ETII SA Fixture

<b>Device Description</b>	Intended to be surgically placed in the bone of the upper or lower jaw arches. Fixture is supplied sterile.		
<b>Material</b>	Pure Titanium (ASTM F 67)		
<b>Surface</b>	SA surface treatment		
	<b>Diameter (mm)</b>	<b>Length (mm)</b>	<b>added/modified</b>
<b>Dimension</b>	3.5	8.7, 10.2, 11.7, 13.2, 15.2, 18.2	Length 18.2mm is added
	4.2	7.2, 8.7, 10.2, 11.7, 13.2, 15.2, 18.2	Length 18.2mm is added
	4.45	7.2, 8.7, 10.2, 11.7, 13.2, 15.2	None
	5.0	6.2	None
		4.9	7.2, 8.7, 10.2, 11.7, 13.2, 15.2

**Hiossen Inc.**

85 Ben Fairless Dr. Fairless Hills, PA 19030  
 Tel : 1-888-678-0001 / Fax : 1-267-759-7004  
 www.hiossen.com

6) Cover Screw

<b>Device Description</b>	Used to protect the exposed platform of the implant during healing period.
<b>Material</b>	Pure Titanium (ASTM F 67)
<b>Surface</b>	Anodizing
<b>diameters</b>	3.03, 3.58, 3.25, 3.4, 3.75, 3.9
<b>lengths</b>	5.25, 5.9, 6.25, 6.85, 6.9, 7.5
<b>Change content as compared to the predicate device</b>	Products code, Art(#)'s are changed & Design Change of screw bottom

7) Healing Abutment

<b>Device Description</b>	Used to make a soft tissue in a funnel shape before setting up prosthetics and removing cover screw after osseointegration.
<b>Material</b>	Pure Titanium (ASTM F 67)
<b>Surface</b>	None
<b>diameters</b>	4.3, 4.8, 5.3, 6.3, 7.3
<b>lengths</b>	7.5, 8.5, 9.5, 11.5, 12.5
<b>Change content as compared to the predicate device</b>	Addition of ø4.5mm abutment. Shape of middle part is changed

**- Substantial Equivalence Matrix**

	The HIOSSSEN Implant System	Predicate devices	
		HTIII SA Fixture / ETIII SA Ultra Wide System	ETII SA Fixture System
<b>510(K) No.</b>	-	K101096 / K103537	K123471
<b>Manufacturer</b>	HIOSSSEN Inc.	HIOSSSEN Inc.	HIOSSSEN Inc.

**Hiossen Inc.**  
 85 Ben Fairless Dr. Fairless Hills, PA 19030  
 Tel : 1-888-678-0001 / Fax : 1-267-759-7004  
 www.hiossen.com

<p><b>Design</b></p>			
<p><b>Intended Use</b></p>	<p>The HIOSSSEN Implant System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading. ETIII SA Ultra-Wide Fixture is intended to be used in the molar region.</p>	<p>The HTIII Fixture System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading. Ultra wide Fixture System is intended to be used in the molar region.</p>	<p>The ETII Fixture System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading.</p>
<p><b>Structure</b></p>	<ul style="list-style-type: none"> <li>- Internal Hex-connected</li> <li>- Submerged Fixture</li> <li>- Tapered body shape &amp; Straight body shape</li> </ul>	<ul style="list-style-type: none"> <li>- Internal Hex-connected</li> <li>- Submerged Fixture</li> <li>- Tapered body shape &amp; Straight body shape</li> </ul>	<ul style="list-style-type: none"> <li>- Internal Hex-connected</li> <li>- Submerged Fixture</li> <li>- Straight body shape</li> </ul>
<p><b>Diameter (D)</b></p>	<p>3.5, 3.75, 3.77, 4.2, 4.25, 4.45, 4.6, 4.25, 4.63, 4.65, 4.9, 5.0, 5.05, 5.08, 5.1, 5.92, 5.95, 6, 6.8</p>	<p>3.75, 4.25, 4.6, 5.05, 5.08, 5.1, 6, 5.92, 5.95, 6.8, 6.82</p>	<p>3.5, 4.2, 4.45, 4.9, 5.0</p>
<p><b>Length (mm)</b></p>	<p>6.2, 7.2, 8.7, 10.2, 11.7, 13.2, 15.2, 18.2</p>	<p>7.2, 8.7, 10.2, 11.7, 13.2, 15.2</p>	<p>6.2, 7.2, 8.7, 10.2, 11.7, 13.2, 15.2</p>
<p><b>Material of Fixture</b></p>	<p>Pure Titanium Grade 4 (ASTM F67)</p>	<p>Pure Titanium Grade 4 (ASTM F67)</p>	<p>Pure Titanium Grade 4 (ASTM F67)</p>
<p><b>Surface</b></p>	<p>SA</p>	<p>SA</p>	<p>SA</p>
<p><b>Sterilization</b></p>	<p>Radiation Sterile</p>	<p>Radiation Sterile</p>	<p>Radiation Sterile</p>
<p><b>Shelf life</b></p>	<p>8 Years</p>	<p>5 Years</p>	<p>5 Years</p>

**Hiossen Inc.**

85 Ben Fairless Dr. Fairless Hills, PA 19030  
Tel : 1-888-678-0001 / Fax : 1-267-759-7004  
www.hiossen.com

<b>S E</b>	<p>The HIOSSSEN Implant System is revised product that Dimension and code(art#) are changed and added from predicate product, ETII SA Fixture, HTIII SA Fixture and ETIII SA Ultra Wide Fixture therefore there is no difference about material, indication for use and design from predicate as above</p> <p>And the subject devices and the predicate devices encompass the same range of physical dimensions except length 18mm and characteristics, including implant diameter and surface treatment therefore The HIOSSSEN Implant System is substantially equivalent to the predicate devices</p>
------------	---

5. Indication for use :

The HIOSSSEN Implant System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading. ETIII SA Ultra-Wide Fixture is intended to be used in the molar region.

6. Review :

The HIOSSSEN Implant System has same material, indication for use, similar design and technological characteristics as the predicate device.

7. Summary of nonclinical testing

Non-clinical testing data submitted, referenced, or relied upon to demonstrate substantial equivalence includes engineering analysis, dimensional analysis, static and dynamic compression-bending testing according to “Guidance for industry and FDA staff Class II Special Controls Guidance Document Root-form Endosseous Dental Implants and Endosseous Dental Abutment” with the worst case scenario. Fatigue testing in air demonstrated the subject device to be equivalent to the tested predicate because shape and minimum diameter of subject devices are the same with predicate.

8. Summary of clinical testing

No clinical studies are submitted

9. Conclusions

In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807, and based on the information provided in this premarket notification HiOSSEN Inc. concludes that the The HIOSSSEN Implant System is substantially equivalent to the predicate devices as described herein