



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

August 20, 2014

Prismatik Dentalcraft, Incorporated  
Mr. Armin Zehtabchi  
Senior RA Specialist  
2212 Dupont Drive, Suite IJK  
Irvine California 92612

Re: K141788

Trade/Device Name: Obsidian™ Milling Blocks  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain powder for clinical use  
Regulatory Class: II  
Product Code: EIH  
Dated: July 22, 2014  
Received: July 23, 2014

Dear Mr. Zehtabchi:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary S. Runner -  
S 

Erin I. Keith, M.S.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)

TBD

K 141788

Device Name

Obsidian™ Milling Blocks

Indications for Use (Describe)

The Obsidian™ Milling Blocks is used to fabricate ceramic dental prostheses in the nature of crowns and bridges for posterior and anterior applications using CAD/CAM methods.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)


This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASstaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



**007\_510 (K) Summary-807.92(c)**

This 510 (k) summary is being submitted in accordance with the requirements of SMDA of 1990 and 21 CFR 807.92.

**A. SUBMITTER INFORMATION**

Company Name: Prismatik Dentalcraft, Inc.

Company Address: 2212 Dupont Dr., Suite IJK,  
Irvine, CA 92612

Company Phone: 949-225-1269

Company FAX: 949-553-0924

Facility Registration Number: 3005477956

Primary Contact Person: Armin Zehtabchi, (949) 225-1234  
Senior RA Specialist

Secondary Contact Person: Marilyn Pourazar, (949) 225-1269  
Senior Director, RA/QA

Date Summary Prepared: July 1, 2014

**B. DEVICE IDENTIFICATION**

Trade/Proprietary Name: Obsidian™ Milling Blocks

21 CFR Reference: 21 CFR 872.6660


21 CFR Common Name: Porcelain powder for clinical use

Classification: Class II, EIH

Panel: Dental

**C. IDENTIFICATION OF PREDICATE DEVICE**

Trade/Proprietary Name: Obsidian™ Ceramic Blocks-K100781



**D. PROPOSED DEVICE DESCRIPTION**


The Obsidian™ Milling Blocks are a lithium silicate ceramic to be supplied with or without an attached mandrel to be milled using CAD/CAM methods. The product can produce a variety of monolithic restorations with great esthetics, life-like translucency, and high strength including full-contour crowns, inlays, onlays, veneers partial crowns and substructures due to its excellent machining properties. The milling blocks will be available in the commonly used VITA Classical and Chromascop Bleach shades.

**E. INDICATIONS FOR USE**

The Obsidian™ Milling Blocks is used to fabricate ceramic dental prostheses in the nature of crowns and bridges for posterior and anterior applications using CAD/CAM methods.

**F. DETERMINATION OF SUBSTANTIAL EQUIVALENCE**

The comparison table below outlines and provides the similarities and the substantial equivalency of the predicate Obsidian™ Ceramic Blocks-K100781 (cleared by Prismatik Dentalcraft, Inc. on 6/21/2010) and the proposed Obsidian™ Milling Blocks, and Prismatik believes that the comparative data presented in the preceding paragraphs, demonstrate that proposed Obsidian™ Milling Blocks is essentially the same as currently marketed devices for the same indications for use, and supports our claim of substantial equivalence to predicate Class II devices under the classification of porcelain powder for clinical use (21 CFR 872.6660) that have previously been found to be substantially equivalent, and that any differences between the proposed device Obsidian™ Milling Blocks and the predicate device do not introduce any new issues of safety or effectiveness. Both the proposed device and the predicate device consist of general glass ceramic material and have the same intended use.




**Table 1 – Comparison between Predicate and Proposed Device**

Attributes	Predicate Device Obsidian™ Ceramic Blocks 510(k)-K100781	Proposed Device Obsidian™ Milling Blocks	Similarities and Difference Between the Predicate and the Proposed Device
Indications for Use	This device is used to fabricate ceramic dental prostheses in the nature of crowns and bridges for posterior and anterior applications using CAD/CAM or hot-press methods.	This device is used to fabricate ceramic dental prostheses in the nature of crowns and bridges for posterior and anterior applications using CAD/CAM methods.	Same
Composition	The average composition of the proposed device is provided. Refer to section 011.	The average composition of the proposed device is provided. Refer to section 011.	Additional materials were added. Refer to section 011.
Shades	A1, A2, A3, A35, B1, B2, B3, C1, C2, C3, D2, D3, BL1 and BL4	A1, A2, A3, A35, B1, B2, B3, C1, C2, D2, D3, BL1 and BL4	Same
Flexural Strength	>300 MPa (meeting ISO 6872 requirements)	>300 MPa (meeting ISO 6872 requirements)	Same
Chemical Solubility	< 100µg/cm <sup>2</sup> (meeting ISO 6872 requirements)	< 100µg/cm <sup>2</sup> (meeting ISO 6872 requirements)	Same
Freedom from Extraneous Material	Shall be free from extraneous materials when assessed by visual inspection (meeting ISO 6872 requirements)	Shall be free from extraneous materials when assessed by visual inspection (meeting ISO 6872 requirements)	Same
Radioactivity	Activity concentration of uranium <sup>238</sup> less than 1.0Bq g <sup>-1</sup> (meeting ISO 6872 requirements)	Activity concentration of uranium <sup>238</sup> less than 1.0Bq g <sup>-1</sup> (meeting ISO 6872 requirements)	Same
Coefficient of Thermal Expansion (25-500°C)	12.2+/-0.5 x 10 <sup>-6</sup> °C (meeting ISO 6872 requirements)	12.2+/-0.5 x 10 <sup>-6</sup> °C (meeting ISO 6872 requirements)	Same
Packaging	Boxes of 5 blocks for Milling Blocks	Boxes of 5 blocks for Milling Blocks	Same
Biocompatibility	Non-toxic and biocompatible (Meeting the ISO 10993-5 and 10993-10 Requirements)	Non-toxic and biocompatible (Meeting the ISO 10993-5 and 10993-10 Requirements)	Same
Accessories	Veneering/Add-On Powders, Paste Stains and Glaze, Low Fusing Ceramic Spray Glaze	Veneering/Add-On Powders, Paste Stains and Glaze, Low Fusing Ceramic Spray Glaze	Same

**G. SUMMARY OF NON-CLINICAL TESTING**

To meet the ISO 6872 requirements, various non-clinical and applicable tests were performed in accordance with the Design Verification Plan including a Risk Analysis addressing the impact of changes to the device. The test results for the Flexural Strength, Chemical Solubility, Freedom from Extraneous Material, Radioactivity, and Coefficient Thermal Expansion indicate the Obsidian™ Milling Blocks is comparable to the predicate device.



In addition, the proposed device, Obsidian™ Milling Blocks, have been tested for Cytotoxicity (ISO 10993-5), Sensitization (ISO 10993-10), and Irritation (ISO 10993-10) to meet the biocompatibility requirement and the reports are as follow:

- The Cytotoxicity Report shows that there was no reaction on any of the cells.
- The Sensitization Report shows that there was no reaction on the tested subject.
- The Irritation Report shows that there was no erythema or edema on the test subject.

Test Description	Results
Cytotoxicity Study using the IX MEM extraction method at 37°C	Pass
ISO Intracutaneous Study, Extract 0.9% sodium chloride USP solution (SC) and sesame oil, NF (SO)	Pass
Irritation and Skin Sensitization Study, Extract 0.9% sodium chloride USP and sesame oil, NF (SO)	Pass

Copies of the tests results are attached. Refer to **Attachment A**.

**H. CONCLUSION**

The proposed device, Obsidian™ Milling Blocks has the same performance specifications, fundamental scientific technology and intended use as that of the predicate device, Obsidian™ Ceramic Blocks. The changes to the device do not raise any new questions regarding safety or efficacy. The performance data and a declaration of conformity with design controls support a determination of continuing substantial equivalence of the proposed device to the predicate device.