



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

SIEMENS HEALTHCARE DIAGNOSTICS INC.  
ALEX HSU  
REGULATORY AND CLINICAL AFFAIRS SPECIALIST  
511 BENEDICT AVE.  
TARRYTOWN NY 10591

April 22, 2015

Re: K141999  
Trade/Device Name: ADVIA Centaur TSH  
Regulation Number: 21 CFR 862.1690  
Regulation Name: Thyroid stimulating hormone test system  
Regulatory Class: II  
Product Code: JLW  
Dated: March 25, 2015  
Received: March 26, 2015

Dear Alex Hsu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
**Katherine Serrano -S**

For : Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
k141999

Device Name  
ADVIA Centaur TSH

### Indications for Use (Describe)

For in vitro diagnostic use in the quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) in serum using the ADVIA Centaur, ADVIA Centaur XP and ADVIA Centaur XPT systems. Measurements of thyroid stimulating hormone produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# 510(k) Summary of Safety and Effectiveness

---

This 510(k) Summary of Safety and Effectiveness is being submitted in accordance with the requirements of 21 CFR 807.92 and the Safe Medical Device Act of 1990.

The assigned 510(k) Number is: k141999

## 1. Date Prepared

April 21, 2015

## 2. Applicant Information

*Contact:* Alex Hsu  
Regulatory and Clinical Affairs Specialist

*Address:* Siemens Healthcare Diagnostics Inc.  
511 Benedict Avenue  
Tarrytown, NY 10591-5097

*Phone:* 914-524-3097

*Fax:* 914-524-3579

*Email:* alex.hsu@siemens.com

## 3. Regulatory Information

**Table 1. Regulatory Information for ADVIA Centaur TSH Assay**

<b>Trade Name</b>	ADVIA Centaur <sup>®</sup> TSH
<b>Model Numbers</b>	08700387 (5-pack); 04911359 (1-pack)
<b>Common Name</b>	Radioimmunoassay, thyroid-stimulating hormone
<b>Classification Name</b>	Thyroid stimulating hormone test system
<b>FDA Classification</b>	Class II
<b>Review Panel</b>	Clinical Chemistry (75)
<b>Product Code</b>	JLW
<b>Regulation Number</b>	862.1690

## 4. Predicate Device Information

The ADVIA Centaur TSH assay was originally cleared by the FDA on 04/30/1991 (k910981) as the ACS TSH Immunoassay.

## 5. Description of Device Modifications

No changes were made to the ADVIA Centaur TSH assay reagents, calibrators or master curve value assignment in order to run on the ADVIA Centaur XPT.

Previously, the lower limit of detection was based on analytical sensitivity, whereas it is now based on Limit of Quantitation (LoQ). As result, the lower limit of detection has been revised from 0.01  $\mu$ IU/mL to 0.05  $\mu$ IU/mL. Accordingly, the analytical measuring range was also modified from 0.01–150  $\mu$ IU/mL to 0.05–150  $\mu$ IU/mL.

## 510(k) Summary of Safety and Effectiveness

Table 2. List of Assay Modifications

Item	ADVIA Centaur TSH (Unmodified Predicate Device)	ADVIA Centaur TSH (Modified Candidate Device)
<b>Platforms</b>	ADVIA Centaur ADVIA Centaur XP	ADVIA Centaur ADVIA Centaur XP ADVIA Centaur XPT
<b>Intended Use</b>	For <i>in vitro</i> diagnostic use in the quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) in serum using the ADVIA Centaur and ADVIA Centaur XP systems.	For <i>in vitro</i> diagnostic use in the quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) in serum using the ADVIA Centaur, ADVIA Centaur XP and ADVIA Centaur XPT systems. Measurements of thyroid stimulating hormone produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders.
<b>Lower Limit of Detection</b>	0.01 $\mu$ IU/mL Based on Analytical Sensitivity	0.05 $\mu$ IU/mL Based on Limit of Quantitation
<b>Expected Values</b>	Includes reference ranges for: Euthyroid Hyperthyroid Hypothyroid	Includes reference ranges for: Euthyroid

Table 3. List of Instrument Modifications

Item	ADVIA Centaur XP	ADVIA Centaur XPT
<b>User Interface CPU</b>	Sun Sparc based CPU running Solaris (UNIX based ) OS, with additional Intel based Application PC (APC) for QC, online documentation running on Windows XP	Single Intel Quad processor based PC, with a new User Interface application with integrated APC applications running on Windows 7
<b>Real Time Control CPU</b>	Sun Sparc based CPU running Solaris (UNIX based ) OS for instrument control and data collection and analysis;	Real Time application rewritten to run on a RoHS compliant ARM 9 based CPU running Nucleus OS;
<b>Microcontrollers</b>	Multiple distributed real-time Microcontrollers	Same (now RoHS compliant)
<b>QC Software</b>	ADVIA QC application providing Stored control results, Levy-Jennings plotting, and statistics, integrated on the Application PC (APC) within the product	ADVIA QC application now integrated into the UI application
<b>Display Monitor</b>	19" LCD Touch Screen color monitor with Graphical User Interface;	22" LCD Touch Screen Color monitor supporting a resolution of 1680 x 1050 with Graphical User Interface

## 510(k) Summary of Safety and Effectiveness

---

**Table 3. List of Instrument Modifications**

<b>Item</b>	<b>ADVIA Centaur XP</b>	<b>ADVIA Centaur XPT</b>
<b>External Printer</b>	No network capable printer	Support for new high speed and networked printers
<b>Remote Diagnostics</b>	External Modem for Remote Diagnostics Interface and application Software for remote diagnostics over the Internet via Server hosted on a separate Application PC (APC) within the product;	Same functionality integrated into the UI application
<b>Barcode Reader</b>	Stationary and handheld barcode scanners for identification of patient samples  Multiple barcode formats supported including 128, 2 of 5, Code39, Codabar;	Same plus added support for 2D barcodes
<b>Data Archival</b>	Data Management, instrument data can be archived to floppy disks or CD	System supports DVDs & memory sticks
<b>Cleaning Procedures</b>	Monthly cleaning procedures	Monthly cleaning has been eliminated
<b>Mounting of Reagent Compartment Refrigeration Hardware</b>	Reagent compartment mounted at the left side of the instrument.  Thermo-electric devices for refrigeration oriented in various directions.	Same physical location of reagent compartment.  Thermo-electric devices for refrigeration oriented in uniform direction.

## 510(k) Summary of Safety and Effectiveness

### 6. Comparison of Similarities and Differences between the Predicate Device and the Candidate Device

Table 4. Similarities/Differences: Unmodified and Modified ADVIA Centaur TSH Assays

Item	ADVIA Centaur TSH (Unmodified Predicate Device)	ADVIA Centaur TSH (Modified Candidate Device)
<b>Instrument Platforms</b>	ADVIA Centaur ADVIA Centaur XP	ADVIA Centaur ADVIA Centaur XP ADVIA Centaur XPT
<b>Intended Use</b>	For <i>in vitro</i> diagnostic use in the quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) in serum using the ADVIA Centaur and ADVIA Centaur XP systems.	For <i>in vitro</i> diagnostic use in the quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) in serum using the ADVIA Centaur, ADVIA Centaur XP and ADVIA Centaur XPT systems. Measurements of thyroid stimulating hormone produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders.
<b>Methodology</b>	Two-site sandwich immunoassay using direct chemiluminometric technology	Same
<b>Reagents</b>	ReadyPack Primary Reagent Pack contains both Solid Phase and Lite Reagent in separate wells	Same
<b>Lite Reagent</b>	Monoclonal mouse anti-TSH antibody labeled with acridinium ester	Same
<b>Solid Phase</b>	Polyclonal sheep anti-TSH antibody covalently coupled to paramagnetic particles	Same
<b>Specimen Type</b>	Serum	Same
<b>Sample Volume</b>	200 µL	Same
<b>Calibration</b>	2-point calibration using Calibrator B	Same
<b>Lower Limit of Detection</b>	0.01 µIU/mL Based on Analytical Sensitivity	0.05 µIU/mL Based on Limit of Quantitation
<b>Expected Values</b>	Includes reference ranges for: Euthyroid Hyperthyroid Hypothyroid	Includes reference ranges for: Euthyroid

## 510(k) Summary of Safety and Effectiveness

Table 5. Similarities/Differences: ADVIA Centaur XP and ADVIA Centaur XPT Instruments

Item	ADVIA Centaur XP	ADVIA Centaur XPT
<b>Principles of Operation</b>	Chemiluminescence using magnetic-particle Solid Phase and chemiluminescent label (acridinium ester) Lite Reagent	Same
<b>Optical System</b>	Photo Multiplier Tube (PMT) used in photon counting mode	Same
<b>Temperature Control</b>	Reactions are controlled at 37°C Reagents stored at 4°C to 8°C	Same
<b>Cleaning Procedures</b>	Monthly cleaning procedures	Monthly cleaning has been eliminated
<b>Mounting of Reagent Compartment Refrigeration Hardware</b>	Reagent compartment mounted at the left side of the instrument. Thermo-electric devices for refrigeration oriented in various directions.	Same physical location of reagent compartment. Thermo-electric devices for refrigeration oriented in uniform direction.
<b>Test Processing</b>	Random Access and Batch; Cuvettes are incubated in a circular, insulated track (Incubation Ring) that advances the cuvette at 15 second intervals and incubates the cuvette at 37°C. The incubation ring moves the cuvettes from the sample probe to the ancillary and reagent probes.	Same
<b>Assay Protocols</b>	Assay specific parameters contained in Test Definitions (TDefs) for each assay.  7.5 min incubation, single step; or 20 min incubation, single step; or 7.5 min / 20 min incubation, 2-step; or 20 min / 20 min incubation, 2-step	Same
<b>Specimens</b>	Serum or plasma; Sample cups or primary tubes; Dilutions allowed on a per-assay basis; Capability of diluting samples requiring pretreatment	Same
<b>Disposables</b>	Reaction cuvettes; Sample Pipette Tips	Same
<b>Calibration</b>	2 point user run calibration; 6 to 10 point stored calibration for each reagent; Calibrators checked with barcode; Calibrator lot numbers stored and displayed	Same

## 510(k) Summary of Safety and Effectiveness

Table 5. Similarities/Differences: ADVIA Centaur XP and ADVIA Centaur XPT Instruments

Item	ADVIA Centaur XP	ADVIA Centaur XPT
<b>Throughput</b>	120 to 240 tests/hr	Same
<b>Time to First Result</b>	15 min, 30 min, 60 min depending upon assay protocol	Same
<b>Dimensions</b>	Floor Model, 60H x 42D x 58L 1200 lbs	Same
<b>User Interface CPU</b>	Sun Sparc based CPU running Solaris (UNIX based ) OS, with additional Intel based Application PC (APC) for QC, online documentation running on Windows XP	Single Intel Quad processor based PC, with a new User Interface application with integrated APC applications running on Windows 7
<b>Real Time Control CPU</b>	Sun Sparc based CPU running Solaris (UNIX based ) OS for instrument control and data collection and analysis;	Real Time application rewritten to run on a RoHS compliant ARM 9 based CPU running Nucleus OS;
<b>Microcontrollers</b>	Multiple distributed real-time Microcontrollers	Same (now RoHS compliant)
<b>QC Software</b>	ADVIA QC application providing Stored control results, Levy-Jennings plotting, and statistics, integrated on the Application PC (APC) within the product	ADVIA QC application now integrated into the UI application
<b>Display Monitor</b>	19" LCD Touch Screen color monitor with Graphical User Interface;	22" LCD Touch Screen Color monitor supporting a resolution of 1680 x 1050 with Graphical User Interface
<b>External Printer</b>	No network capable printer	Support for new high speed and networked printers
<b>Remote Diagnostics</b>	External Modem for Remote Diagnostics Interface and application Software for remote diagnostics over the Internet via Server hosted on a separate Application PC (APC) within the product;	Same functionality integrated into the UI application
<b>Barcode Reader</b>	Stationary and handheld barcode scanners for identification of patient samples  Multiple barcode formats supported including 128, 2 of 5, Code39, Codabar;	Same plus added support for 2D barcodes
<b>Data Archival</b>	Data Management, instrument data can be archived to floppy disks or CD	System supports DVDs & memory sticks

### 7. Summary of Design Control Activities

Design control activities, as outlined in 21CFR 820.30, were completed for the ADVIA Centaur XPT system.

A risk analysis (Failure Modes and Event Analysis) was undertaken to assess risks of using the device. This evaluation followed the Siemens Healthcare Diagnostics procedure for risk analysis, which is based on ISO 14971:2007, Medical devices – Application of risk management to medical devices. The risk analysis concluded that all identified risks were properly mitigated and no unacceptable risks are present.

The performance of the ADVIA Centaur TSH assay using the ADVIA Centaur XPT system was verified to ensure equivalent performance when used on the predicate ADVIA Centaur XP system. All verification testing met pre-determined acceptance criteria. Therefore, the introduction of the ADVIA Centaur XPT system does not negatively impact the performance, safety or effectiveness of the ADVIA Centaur TSH assay.

### 8. Conclusions

The performance of the ADVIA Centaur TSH assay on the ADVIA Centaur XPT system is substantially equivalent to the ADVIA Centaur TSH assay running on the currently-marketed predicate ADVIA Centaur XP system.

The ADVIA Centaur XPT system has the same operating principles, assay performance characteristics and intended use as the predicate device, the ADVIA Centaur XP system. The results of performance testing and verification activities demonstrate that the design modifications to the ADVIA Centaur XPT do not impact its safety or effectiveness and do not alter its performance claims.

Furthermore, there have been no changes to the intended use of the ADVIA Centaur TSH assay, other than to include the ADVIA Centaur XPT, as described in the labeling, or the fundamental scientific technology of the device.