



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

November 13, 2014

Prismatik Dentalcraft, INC
Brandon Shepard
Regulatory Affairs & Quality Assurance Specialist
2212 Dupont Dr. Suite P
Irvine, California 92612

Re: K142115

Trade/Device Name: Inclusive® Titanium Abutments, compatible with Zimmer Screw-Vent, Biomet 3i Certain, and Nobel Biocare NobelReplace Implants
Regulation Number: 21 CFR 872.3630
Regulation Name: Endosseous dental implant abutment
Regulatory Class: II
Product Code: NHA
Dated: October 10, 2014
Received: October 14, 2014

Dear Mr. Shepard:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Susan Runna DDS, MA". The signature is written in a cursive style. Behind the signature, there is a faint, semi-transparent watermark of the FDA logo.

Erin I. Keith, M.S.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number: K142115

Device Name: Inclusive® Titanium Abutments, compatible with Zimmer Screw-Vent, Biomet 3i Certain, and Nobel Biocare NobelReplace Implants

Indications for Use:

Inclusive® Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

Inclusive® Titanium Abutments are compatible with:

- Biomet: 3i Certain internal hex implants in 3.4, 4.1, 5.0, 6.0 mm sizes
- Nobel Biocare: NobelReplace straight and tapered internal connection implants in NP, RP, WP, 6.0 mm sizes
- Zimmer: Screw-Vent and Tapered Screw-Vent internal hex implants in 3.5, 4.5, 5.7 mm sizes

Prescription Use 1
(Part 21 CFR 801 Subpart D)

Concurrence of CDRH, Office of Device Evaluation (ODE)



005

510(k) Summary

[As Required by 21 CFR 807.92]

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements 21 CFR 807.92.

A. SUBMITTER INFORMATION

Company Name: PRISMATIK DENTALCRAFT, INC.

Company Address: 2212 Dupont Dr., Suite P
Irvine, CA 92612

Company Phone / Fax: (949) 225-1269 / (978) 313-0850

Contact Person: Primary Contact:

- Brandon Shepard, (949) 225-1243

 Secondary Contact:

- Marilyn Pourazar, (949) 225-1269

Date Summary Prepared: August 1, 2014

B. DEVICE IDENTIFICATION

Trade/Proprietary Name: Inclusive[®] Titanium Abutments
compatible with: Zimmer Screw-Vent, Biomet 3i
Certain, and Nobel Biocare NobelReplace Implants

Common Name: Endosseous Dental Implant Abutment

Regulation Number: 872.3630

Product Code: NHA

Device Class: 2

Review Panel: Dental

C. IDENTIFICATION OF PREDICATE DEVICE

Trade/Proprietary Name:

- Inclusive™ Abutment for Zimmer, 3i and Nobel Biocare Implants (K073217)
- NobelProcera Ti Abutment (K091756)
- BIOMET 3i Dental Abutments and Restorative Components (K072642)
- Zimmer Ti Prepable Abutment (K092403).

D. DEVICE DESCRIPTION

Inclusive Titanium Abutments are endosseous implant abutments which are placed into the dental implant to provide support for a prosthetic restoration. The abutment is placed over the implant shoulder and is fastened into the implant with a screw. Abutments are made of titanium grade 23 (Ti-6AL-4V ELI) and meet ASTM F-136 Standard.

E. INDICATIONS FOR USE

Inclusive® Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

Inclusive® Titanium Abutments are compatible with:

- Biomet: 3i Certain internal hex implants in 3.4, 4.1, 5.0, 6.0 mm sizes
- Nobel Biocare: NobelReplace straight and tapered internal connection implants in NP, RP, WP, 6.0 mm sizes
- Zimmer: Screw-Vent and Tapered Screw-Vent internal hex implants in 3.5, 4.5, 5.7 mm sizes

F. NON-CLINICAL TESTING

Non-clinical test data was used to evaluate the device's safety and effectiveness, and determine substantial equivalence with predicate devices.

Non-clinical testing was performed in accordance with FDA Guidance "*Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments*" and it consisted of Reliability Calculation, Fatigue Testing and Static Load Failure Testing of finished assembled implant/abutment systems.

In addition, sterilization validation information and a recommended sterilization method based on ANSI/AAMI ST79 and ISO 17665-1 is provided in the Information for Use.

The testing performed demonstrated implant to abutment compatibility and established that the predicate device is as safe, as effective, and performs as well as the predicate device(s).

G. SUBSTANTIAL EQUIVALENCE

Inclusive[®] Titanium Abutments, compatible with Zimmer Screw-Vent, Biomet 3i Certain and Nobel Biocare NobelReplace Implants are substantially equivalent to the Inclusive[™] Abutment for Zimmer, 3i and Nobel Biocare Implants (K073217) and the OEM predicates, (K091756, K072642, and K092403), identified in Section C above. They are substantially equivalent in intended use, materials, design and performance.

(See Comparison Tables below).

Comparison of Predicate and Proposed Devices: Zimmer Screw-Vent Platform

	Predicate Device (1)	Predicate Device (2)	Proposed Device	
	Ti Prepable Abutment (K092403)	Inclusive Abutment for Zimmer, 3i, and Nobel Biocare Implants (K073217)	Inclusive Titanium Abutments compatible with: NobelReplace Implants	Similarities and Differences
Manufacturer	Zimmer Dental	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Hex Dimensions: 2.4mm Across Flats 3.0mm Across Flats	Cylindrical Diameter: 9.4mm Hex Dimensions: 2.4mm Across Flats 3.0mm Across Flats	Cylindrical Diameter: 9.4mm Hex Dimensions: 2.4mm Across Flats 3.0mm Across Flats	Same
Abutment Screw Manufacturer	OEM (<i>Zimmer</i>)	OEM (<i>Zimmer</i>)	Inclusive (<i>Prismatik</i>)	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length 8.2mm; 1-72 UNF-2A Thread	Length 8.2mm; 1-72 UNF-2A Thread	Length 8.2mm; 1-72 UNF-2A Thread	Same
Indications for Use	The Ti Prepable Abutments are designed for use as a terminal or intermediate abutment for cement retained prosthesis. It allows for preparation for the crown to be attached. It can be used for a single or multiple-unit restoration. The abutment is intended to be prepared and placed with patient specific margins.	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different wording in the Indications for Use
Platform Compatibility	Zimmer Dental Screw-Vent 3.5, 4.5, 5.7mm	Zimmer Dental Screw-Vent 3.5, 4.5, 5.7mm	Zimmer Dental Screw-Vent 3.5, 4.5, 5.7mm	Same
Connection	Internal Hex	Internal Hex	Internal Hex	Same
Design/Construction	Machined	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Flat	Flat	Flat	Same
Screw Seat	Flat	Flat	Flat	Same
Material	Titanium Alloy	Titanium Alloy	Titanium Alloy	Same

Comparison of Predicate and Proposed Devices: Biomet 3i Certain Platform

	Predicate Device (1)	Predicate Device (2)	Proposed Device	Similarities and Differences
	Biomet 3i Dental Abutments & Restorative Components (K072642)	Inclusive Abutment for Zimmer, 3i, and Nobel Biocare Implants (K073217)	Inclusive Titanium Abutments compatible with: Biomet 3i Certain Implants	
Manufacturer	Biomet	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Hex Dimensions: 2.2mm Across Flats 2.7mm Across Flats	Cylindrical Diameter: 9.4mm Hex Dimensions: 2.2mm Across Flats 2.7mm Across Flats	Cylindrical Diameter: 9.4mm Hex Dimensions: 2.2mm Across Flats 2.7mm Across Flats	Same
Abutment Screw Manufacturer	OEM (<i>Biomet 3i</i>)	OEM (<i>Biomet 3i</i>)	Inclusive (<i>Prismatik</i>)	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length: 8.4mm; M1.6 Thread	Length: 8.4mm; M1.6 Thread	Length: 8.4mm; M1.6 Thread	Same
Indications for Use	Biomet 3i Dental Abutments are intended for use as accessories to endosseous dental implant to support a prosthetic device in a partially or completely edentulous patient. A dental abutment is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be screw retained or cement retained.	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different Indications for Use statement.
Platform Compatibility	Biomet 3i Certain 3.4, 4.1, 5.0, 6.0mm	Biomet 3i Certain 3.4, 4.1, 5.0, 6.0mm	Biomet 3i Certain 3.4, 4.1, 5.0, 6.0mm	Same
Connection	Internal Hex	Internal Hex	Internal Hex	Same
Design/Construction	Machined	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°-15°	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Flat	Flat	Flat	Same
Screw Seat	Flat	Flat	Flat	Same
Material	Titanium Alloy	Titanium Alloy	Titanium Alloy	Same

Comparison of Predicate and Proposed Devices: Nobel Biocare NobelReplace Platform

	Predicate Device (1)	Predicate Device (2)	Proposed Device	
	NobelProcera Ti Abutment (K091756)	Inclusive Abutment for Zimmer, 3i, and Nobel Biocare Implants (K073217)	Inclusive Titanium Abutments compatible with: Nobel Biocare NobelReplace Implants	Similarities and Differences
Manufacturer	Nobel Biocare	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Connection Length: 2.9mm	Cylindrical Diameter: 9.4mm Connection Length: 2.9mm	Cylindrical Diameter: 9.4mm Connection Length: 2.9mm	Same
Abutment Screw Manufacturer	OEM (Nobel Biocare)	OEM (Nobel Biocare)	Inclusive (Prismatik)	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length 8.3mm; M1.8 Thread	Length 8.3mm; M1.8 Thread	Length 8.3mm; M1.8 Thread	Same
Indications for Use	The NobelProcera Ti Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different Indications for Use statement.
Platform Compatibility	Nobel Biocare NobelReplace NP, RP, WP, 6.0	Nobel Biocare NobelReplace NP, RP, WP, 6.0	Nobel Biocare NobelReplace NP, RP, WP, 6.0	Same
Connection	Internal Trilobe	Internal Trilobe	Internal Trilobe	Same
Design/Construction	Machined	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°-30°	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Flat	Flat	Flat	Same
Screw Seat	Flat	Flat	Flat	Same
Material	Titanium Alloy	Titanium Alloy	Titanium Alloy	Same