



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

Reliance Medical Systems, LLC
% Meredith L. May, MS, RAC
Empirical Consulting, LLC
4628 Northpark Drive
Colorado Springs, Colorado 80918

April 2, 2015

Re: K142217

Trade/Device Name: Reliance Interspinous Plate System
Regulation Number: 21 CFR 888.3050
Regulation Name: Spinal interlaminar fixation orthosis
Regulatory Class: Class II
Product Code: PEK
Dated: February 19, 2015
Received: February 20, 2015

Dear Ms. May:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mark N. Melkerson -S

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K142217

Device Name

Reliance Interspinous Plate System

Indications for Use (Describe)

The Reliance Interspinous Plate System is a supplemental fixation device intended as an adjunct to fusion of the thoracolumbar (T1- S1) spine by fixation to the spinous processes in skeletally mature patients, and is indicated for the following:

- degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies)
- spondylolisthesis
- trauma (i.e., fracture or dislocation)
- tumors

This device is not intended for stand-alone use.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

5. 510(K) SUMMARY

Submitter's Name:	Reliance Medical Systems, LLC
Submitter's Address:	545 West 500 South Suite 100 Bountiful, UT 84010
Submitter's Telephone:	(801) 295-3280
Contact Person:	Meredith L. May, MS, RAC Empirical Consulting, LLC 719.337.7579
Date Summary was Prepared:	31 March 2015
Trade or Proprietary Name:	Reliance Interspinous Plate System
Common or Usual Name:	Spinous Process Plate
Classification:	Class II per 21 CFR §888.3050
Product Code:	PEK
Classification Panel:	Orthopedic (87)

DESCRIPTION OF THE DEVICE SUBJECT TO PREMARKET NOTIFICATION:

The Reliance Interspinous Plate System consists of an interspinous process spacer-plate and process lock plate assembly for posterior fixation of the spine in order to achieve fusion. The interspinous process spacer-plate and process lock plate are available in multiple sizes to accommodate various patient anatomies. The insert features a lateral plate with teeth to interface with the bone of the spinous processes. The process lock plate also features teeth which contact the bone of the spinous processes. The process lock plate is passed through the insert such that, in their final position, interspinous process spacer-plate and process lock plate surround the spinous processes on both sides, and fixation is achieved via compression of the two components onto the spinous processes.

INDICATIONS FOR USE

The Reliance Interspinous Plate System is a supplemental fixation device intended as an adjunct to fusion of the thoracolumbar (T1- S1) spine by fixation to the spinous processes in skeletally mature patients, and is indicated for the following:

- Degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies)
- spondylolisthesis
- trauma (i.e., fracture or dislocation)
- tumors

This device is not intended for stand-alone use.

The indications for use for the Reliance Medical Interspinous Plate is similar to that of the predicates listed in Table 5-1 below.

TECHNOLOGICAL CHARACTERISTICS

The interspinous process spacer-plate and process lock plate are manufactured from medical grade Titanium (Ti 6Al-4V) per ASTM F136. The implants are provided non-sterile with instructions for sterilization. The interspinous plates are designed in interspinous heights of 18mm, total heights of 35-50mm, and widths of 25-30.

The subject and predicate devices have nearly identical technological characteristics and the minor differences do not raise any new issues of safety and effectiveness. Specifically the following characteristics are identical between the subject and predicates:

- Indications for Use
- Materials of manufacture
- Structural support mechanism
- Sterilization

Table 5-1 Predicate Devices

510k Number	Trade or Proprietary or Model Name	Manufacturer	Predicate Type
K133363	Interspinous Posterior Fixation System	InterBRIDGE	Additional
K073278, K131238	Affix Spinous Process Plate System	NuVasive®	Additional
K123093	StabiLink™ MIS Spinal Fixation System	Southern Spine	Additional
K100354	PrimaLOK™ Interspinous Fusion System	OsteoMed	Primary
K122509	Spinous Process Fixation Plate	VertiFlex®	Additional

PERFORMANCE DATA

The Reliance Interspinous Plate System has been tested in the following test modes:

- Static axial compression per ASTM F1717-13
- Static torsion per ASTM F1717-13
- Static axial pull-out per ASTM F1717-13
- Dynamic axial compression per ASTM F1717-13
- Static plate dissociation

The results of this non-clinical testing show that the strength of the Reliance Medical Interspinous Plate is sufficient for its intended use and is substantially equivalent to legally marketed predicate devices.

CONCLUSION

The overall technology characteristics and mechanical performance data lead to the conclusion that the Reliance Interspinous Plate System is substantially equivalent to the predicate device.