



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

February 20, 2015

Beijing Choice Electronic Technology Co., Ltd.
Mr. Lei Chen
Quality Director
North Building 3F, No.9 Shuangyuan Road
Badachu Hi-tech Zone, Shijingshan District
Beijing 100041, China

Re: K142888

Trade/Device Name: Fingertip Pulse Oximeter MD300C318T2
Regulation Number: 21 CFR 820.2700
Regulation Name: Oximeter
Regulatory Class: Class II
Product Code: DQA
Dated: January 20, 2015
Received: January 22, 2015

Dear Mr. Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Tejashri Purohit-Sheth, M.D.

Tejashri Purohit-Sheth, M.D.
Clinical Deputy Director
DAGRID/ODE/CDRH FOR

Erin I. Keith, M.S.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Section II Indications for Use Statement

Indications for Use

510(k) Number (if known): K142888

Device Name: Fingertip Pulse Oximeter MD300C318T2

Indications for Use:

The MD300C318T2 is intended for spot checking in measuring and displaying functional arterial hemoglobin (SpO₂) and pulse rate. It is intended for adult, adolescent, child and infant users in hospitals, hospital facilities and home healthcare environments.

Prescription Use √ AND/OR Over-The-Counter Use _____
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Section III 510(k) Summary

This summary of 510(k) is being submitted in accordance with the requirements of 21 CFR 807.92.

There is no prior submission for the device.

3.1 Submitter Information

- **Manufacturer Name:**
Establishment Registration Number: 3005569927
Beijing Choice Electronic Technology Co., Ltd.
Room 320, West Building 4, No.83 FuXing Road
Beijing 100039, P.R.China
- **Contact Person:**
Mr.Lei Chen
Beijing Choice Electronic Technology Co., Ltd.
North Building 3F, No. 9 Shuangyuan Road,
Badachu Hi-tech Zone, Shijingshan District
Beijing China 100041
Phone: +86-10-88798300 Ext 6020
Fax: 215-4052545
Email: cc@choicemmed.com
- **Date prepared:** Nov.06, 2014

3.2 Proposed Device Information

Device Common Name: Pulse Oximeter
Device Trade/Proprietary Name: Fingertip Pulse Oximeter
Model: MD300C318T2
Classification Name: Oximeter
Regulation Number: 870.2700
Product Code: DQA
Class: II
Panel: Anesthesiology

3.3 Predicate Device

510(k) Number: K092620
Common Name: Pulse Oximeter
Device Trade/Proprietary Name: Fingertip Pulse Oximeter

Premarket Notification 510(k) Submission—Section III 510(k) Summary

Model: MD300C318

Classification Name: Oximeter

Product Code: DQA

Regulation Number: 870.2700

Device Class: II

Panel: Anesthesiology

Manufacturer: Beijing Choice Electronic Technology Co., Ltd.

3.4 Device Description

The proposed device of Fingertip Pulse Oximeter MD300C318T2 is a battery powered fingertip device, with the Durable and Deluxe design. The new kind of fingertip pulse oximeter integrated with Bluetooth wireless technology is designed to monitor SpO₂ and pulse rate and connect with telemedicine.

The proposed device consists of power supply module, detector and emitter LED, signal collection and process module, display module, Bluetooth[®] module, voice indicator module, USB module, user interface.

The Pulse Oximeter works by applying a sensor to a pulsating arteriolar vascular bed. The sensor contains a dual light source and photo detector. The one wavelength of light source is 660nm, which is red light; the other is 940nm, which is ultra red light. Skin, bone, tissue, and venous vessels normally absorb a constant amount of light over time. The photodetector in finger sensor collects and converts the light into electronic signal which is proportional to the light intensity. The arteriolar bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases. The ratio of light absorbed at systole and diastole is translated into an oxygen saturation measurement. This measurement is referred to as SpO₂.

The power sources of the proposed device are 2 AAA alkaline batteries. The proposed device has low battery voltage indicator function and the proposed device will automatically power off when there is no signal for longer than 8 seconds.

The proposed device is not for life-supporting or life-sustaining, not for implant. The device or transducers are not sterile and the transducer is reusable and does not need sterilization or re-sterilization. The device is for prescription. The device does not contain drug or biological products.

The device is software-driven and the software validation is provided in Section of Software.

3.5 Comparison list of the technological characteristics

Table 3-1 Performance Specification Comparison

Comparison Elements	Proposed Device	Predicate Device
Product Name	Fingertip Pulse Oximeter	Fingertip Pulse Oximeter(K092620)
Model	MD300C318T2	MD300C318
Regulation No.	21 CFR 870.2700	21 CFR 870.2700
Classification	II	II
Classification Name	Oximeter	Oximeter
Product Code	DQA	DQA
Indications for Use	The MD300C318T2 is intended for spot checking in measuring and displaying functional arterial hemoglobin (SpO ₂) and pulse rate. It is intended for adult, adolescent, child and infant users in hospitals, hospital facilities and home healthcare environments.	The MD300C318 Fingertip Pulse Oximeter is indicated for continuous use or spot checking in measuring and displaying functional arterial oxygen saturation (SpO ₂) and Pulse Rate of patients in hospitals and home care. It is intended for adult and pediatric patients on finger between 0.3-1.0 inch (0.8 - 2.5 cm) thick.
Comparison Statement	The proposed device has the same indications for use and classification as the predicate device.	
Components	The proposed device consists of power supply module, detector and emitter LED, signal collection and process module, display module, Bluetooth [®] module, voice	The predicated device consists of detector and emitter LED, signal amplify unit, CPU, data display unit, power unit and with Bluetooth communication which is

Premarket Notification 510(k) Submission—Section III 510(k) Summary

		indicator module, USB module, user interface.	integrated into a telemedicine system or other health data collection system through the wireless collection.
Design Principle		Principle of the oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (RHb) and Oxyhemoglobin (HbO ₂) in red and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm glow and 940nm near infrared light) can be focused onto a human nail tip through a clamping finger-type sensor. A measured signal obtained by a photosensitive element, will be shown on the oximeter's display through process in electronic circuits and microprocessor.	Principle of the oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (RHb) and Oxyhemoglobin (HbO ₂) in red and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm glow and 940nm near infrared light) can be focused onto a human nail tip through a clamping finger-type sensor. A measured signal obtained by a photosensitive element, will be shown on the oximeter's display through process in electronic circuits and microprocessor.
Measurement Wavelength	Red	660nm	660nm
	Infrared	940nm	940nm
Comparison Statement		The proposed device has the same design and similar components as the predicate	

Premarket Notification 510(k) Submission—Section III 510(k) Summary

Device specification	device.	
Display Type	OLED	OLED
Power Supply	2 AAA-size alkaline batteries	Polymer Lithium-ion Rechargeable battery
Display Data	SpO2 , PR	SpO2 , PR
SpO2 display range	0~100%	0~100%
Measurement range	70-100%	70~100%
Accuracy	70%~100%: $\pm 2\%$; 70%~80%; $\pm 2\%$ 80%~90%; $\pm 2\%$ 90%~100%; $\pm 1\%$ 0%~69% no definition	70%~100%: $\pm 3\%$; 70%~80%; $\pm 3\%$ 80%~90%; $\pm 3\%$ 90%~100%; $\pm 3\%$ 0%~69% no definition
Resolution	1%	1%
PR display range	0~235BPM	0~235BPM
PR measurement range	30-235 BPM	30-235 BPM
PR Accuracy	30~99bpm, $\pm 2\text{bpm}$; 100~235bpm, $\pm 2\%$	30~99bpm, $\pm 2\text{bpm}$; 100~235bpm, $\pm 2\%$
Resolution	1bpm	1bpm
Wireless Transmission Range	0~10m	0~10m
Antenna Type	Internal	Internal
Transmitter	Bluetooth Compliance: Version 2.1	Bluetooth Compliance: Version 2.0
Operating Frequency	2402 to 2480 MHz	2402 to 2480 MHz
Operating Temperature	5~40°C	5~40°C
Storage	-25~+70°C	-20~55°C

Premarket Notification 510(k) Submission—Section III 510(k) Summary

	Temperature		
	Relative Humidity	15%~93% in operation ≤93% in storage	20%~85% in operation, none condensing <85% in storage, none condensing
	Atmosphere pressure	86kPa~106kPa	86kPa~106kPa
Comparison Statement		The applicant device has similar device specifications as the predicate device	
Construction Materials	Battery Cover	ABS	ABS
	Fingertip Cushion	Medical Silicon gel	Medical Silicon gel
	Enclose	ABS	ABS
Comparison Statement		The contacting materials of applicant device are as same as that of the predicate device.	
Performance Testing	Performance Testing	Meet the requirements of FDA Guidance	Meet the requirements of FDA Guidance
	Clinical Test	Conformed to ISO 80601-2-61 Clinical test for device accuracy is conducted by the Yue Bei People’s hospital. The clinical test report and protocol are provided in <i>Performance Testing-Clinical</i>	Conformed to ISO9919
Electromagnetic compatibility	Electrical Safety	Conformed to IEC60601-1	Conformed to IEC60601-1

Premarket Notification 510(k) Submission—Section III 510(k) Summary

	Electromagnetic Compatibility	Conformed to IEC60601-1-2		Conformed to IEC60601-1-2	
Software		Moderate Level of Concern		Moderate Level of Concern	
		Compliance with FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices.		Compliance with FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices.	
		Risk Management in Compliance with ISO14971		Risk Management in Compliance with ISO14971	
Biocompatibility	Black Medical Silicon gel	In Vitro Cytotoxicity	No cytotoxic potential.	In Vitro Cytotoxicity	No cytotoxic potential.
		Skin Irritation Test	No evidence of causing sensitization.	Skin Irritation Test	No evidence of causing sensitization.
		Animal Skin Sensitization Test	No evidence of significant sensitization from the test extract to rabbits.	Animal Skin Sensitization Test	No evidence of significant sensitization from the test extract to rabbits.
	ABS plastic Enclosure	In Vitro Cytotoxicity	No cytotoxic potential.	In Vitro Cytotoxicity	No cytotoxic potential.
		Skin Irritation Test	No evidence of causing sensitization.	Skin Irritation Test	No evidence of causing sensitization.
		Animal Skin	No evidence of	Animal Skin	No evidence of

Premarket Notification 510(k) Submission—Section III 510(k) Summary

		Sensitization Test	significant sensitization from the test extract to rabbits.	Sensitization Test	significant sensitization from the test extract to rabbits.
Label and Labeling		Compliance with FDA guidance		Compliance with FDA guidance	

3.6 Intended use

The MD300C318T2 is intended for spot checking in measuring and displaying functional arterial hemoglobin (SpO₂) and pulse rate. It is intended for adult, adolescent, child and infant users in hospitals, hospital facilities and home healthcare environments.

3.7 Test

Non-clinical Test

The Fingertip Pulse Oximeter MD300318T2 is designed and tested and will be manufactured in accordance with both mandatory and voluntary standards, including:

IEC60601-1 Medical Electrical Equipment-Part 1: General requirements for safety.

IEC60601-1-2:2007 Medical Electrical Equipment-Part1-2: General requirements for safety, Collateral standard: Electromagnetic compatibility -Requirements and tests.

IEC60601-1-11:2010 Medical Electrical Equipment-Part 1-11: General requirements for basic safety and essential performance-Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

ISO 80601-2-61:2011 Medical Electrical Equipment- Part 2-61: Particular requirements for basic safety and essential performance of pulse oximeter equipment.

The Software Validation is in compliance with FDA Guidance to Compliance with FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices.

The compatibility of the skin-contact component material in the finished product meets the requirement of Biocompatibility. The Biological Evaluation Tests are in compliance with the standards of ISO 10993, “Biological Evaluation of Medical Devices”.

The list of non-clinical test performed on the proposed device.

No.	Test Name
1	System Performance Test
2	Shelf Life Test
3	Performance Test according to ISO 80601-2-61
4	Electromagnetic Compatibility Test According to IEC 60601-1-2
5	Electrical Safety Test According to IEC 60601-1
6	Used in the home healthcare environment test according to IEC60601-1-11
7	Irritation ,Sensitization and Cytotoxicity Test according to ISO 10993

Premarket Notification 510(k) Submission—Section III 510(k) Summary

Clinical Test

The Clinical Test reports was conducted in accordance to EN ISO 14155-1:2009, EN ISO 14155-2:2009, ISO 9919:2005, EN ISO 9919:2009, BS EN ISO 80601-2-61:2011, and the FDA Guidance Document for Pulse Oximeters.

Subjects:

After Institutional Review Board (IRB) approval, 12 healthy adult volunteer subjects (ages 21-33yr, 46-75kg, 158-183cm, with light to dark pigmentation) were included in the study conducted July.5-6, 2013 to evaluate the SpO₂ accuracy performance of the MD300C318T2 Fingertip Pulse Oximeter.

Methods:

Each system was evaluated during steady state/non-motion conditions with various levels of induced hypoxia resulting in stable oxygen saturation levels between 100% and 70% SaO₂. Arterial blood samples were drawn during simultaneous data collection from the test devices. The blood was immediately analyzed on Reference CO-Oximetry providing functional SaO₂ for the basis of the SpO₂ accuracy comparison.

Adverse events and complications:

There are no adverse events during the clinical test.

Conclusions:

The results show the MD300C318T2 Fingertip Pulse Oximeter to pass a SpO₂ accuracy specification of 2 which Beijing Choice Electronic Tech Co., Ltd. claims during steady state conditions over the range of 70-100%.

3.8 Determination of substantial equivalence

The proposed device of Fingertip Pulse Oximeter MD300C318T2 has the same classification information, same indications and intended use, same design principle, similar product design and specifications, same performance effectiveness, performance safety as the predicate device. The differences exist in storage temperature, relative humidity, Bluetooth compliance version and power supply type. These differences are slight and do not influence the effectiveness and safety of the device. According to the non-clinical and clinical test results, the proposed device is as safe as effective and performance as well as the predicate device. So the proposed device is Substantially Equivalent (SE) to the predicate device which is US legally market device.