



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 21, 2015

Dentca, Inc.  
c/o Mr. Claude Berthoin  
Denterprise International, Inc.  
110 East Granada Blvd, Suite 207  
Ormond Beach, Florida 32176

Re: K143033  
Trade/Device Name: Dentca Denture Base  
Regulation Number: 21 CFR 872.3760  
Regulation Name: Denture Relining, Repairing, or Rebasing Resin  
Regulatory Class: II  
Product Code: EBI  
Dated: June 15, 2015  
Received: June 18, 2015

Dear Mr. Berthoin,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Tina  
Kiang -S

for Erin I. Keith, M.S.

Director

Division of Anesthesiology,

General Hospital, Respiratory, Infection  
Control, and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K143033

Device Name

Dentca Denture Base

Indications for Use (Describe)

Dentca Denture Base is a light-cured resin indicated for fabrication and repair of full and partial removable dentures and baseplates. The material is an alternative to traditional heat-cured and auto polymerizing resins.

Fabrication of dental prosthetics with Dentca Denture Base requires a computer-aided design and manufacturing (CAD/CAM) system that includes the following components not part of the device: oral casting impression, digital denture base file created in an optical impression system, stereolithographic additive printer, and curing light equipment.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

K143033

## 510(k) Summary

### Submitter

Dentca, Inc.  
3608 Griffith Avenue  
Los Angeles, CA 90011

Phone: 323-232-7505  
Fax: 323-232-7504

Contact Person: Jason Lee  
Date Prepared: July 20, 2015

### Consultant

Denterprise International, Inc.  
110 East Granada Blvd., Suite 207  
Ormond Beach, FL 32176

Phone: 386-672-0450  
eFax: 855-235-7902

Contact Person: Claude Berthoin, President

### Device Classification

Trade Name: Dentca Denture Base  
Common Name: Dental Acrylic Resin  
Regulation Number: 21 CFR 872.3760  
Regulation Name: Denture relining, repairing, or rebasing resin  
Product Code: EBI  
Submission Type: 510(k)  
Regulatory Class: 2  
Medical Specialty: Dental

### Primary Predicate Device

The following primary predicate is a legally marketed, post-amendment device:

510(k) Number: K032892  
Clearance Date: September 26, 2003  
510(k) Trade Name: Trubyte Denture Base Resin System  
Actual trade names: Trubyte Triad VLC (Visible Light Cure) Denture Base System;  
Trubyte Eclipse VLC Denture Base System  
Manufacturer: Dentsply International, Inc. (York, PA, USA)  
Regulation & PC: 872.3760; EBI

## Reference Device

The following reference predicate is a legally marketed, post-amendment device:

510(k) Number: K102776  
Clearance Date: February 18, 2011  
510(k) Trade Name: e-Dent Temporary Resin and Extraoral Curing System  
Actual Trade Name: EnvisionTEC E-Dent 100  
Manufacturer: DentaMed GmbH (Friedberg, Germany)  
Regulation & PC: 872.3770; EBG

## Device Description

Dentca Denture Base is a photosensitive resin intended to fabricate removable dentures in a CAD/CAM additive printing process.

The Dentca polymer is a viscous solution of the following compounds: methacrylate-based resins, a photoinitiator that activates at 405 nm visible (blue) light, an inhibitor, and pigments. It comes in two sizes, large and small bottles (see photo at right). It is a Type 4 (light-activated) acrylic resin as classified by ANSI/ADA Specification No. 12.

The denture fabrication process begins with a traditional casting impression of the oral region in the dentist office. This impression is sent to a dental lab for conversion to digital image in an optical impression system. The denture base is then made layer-by-layer in a stereolithographic laser printer. After attachment of preformed plastic teeth, the denture is cured in a light chamber, and, lastly, sent back to the dentist for try-in and final adjustment.

Key performance specifications include:

Flexural strength...	90.2 MPa
Flexural modulus...	2,290 MPa
Residual monomer...	Not detectable
Water sorption...	14 ug/mm <sup>3</sup>
Water solubility...	1.3 ug/mm <sup>3</sup>

## Indications For Use

Dentca Denture Base is a light-cured resin indicated for fabrication and repair of full and partial removable dentures and baseplates. The material is an alternative to traditional heat-cured and auto polymerizing resins.

Fabrication of dental prosthetics with Dentca Denture Base requires a computer-aided design and manufacturing (CAD/CAM) system that includes the following components not part of the device: oral casting impression, digital denture base file created in an optical impression system, stereolithographic additive printer, and curing light equipment.

Both the subject and predicate are light-cure resins indicated for fabrication and repair of removable dental prostheses.

## Comparison Of Technological Characteristics With Predicate

The following table compares technological and other characteristics of the subject and primary predicate devices.

**Table 5 -- Comparison Of Technical Features**

	<b>K143033 Subject Device</b>	<b>K032892 Predicate Device</b>	<b>Differences</b>
<b>Device Names</b>	<b>Dentca Denture Base</b>	<b>Dentsply Trubyte Denture Base Resin System</b>	<b>NA</b>
<b>Classification &amp; Product Code</b>	872.3760; EBI	872.3760; EBI	SAME CLASSIFICATION
<b>Intended Use</b>	Fabrication and repair of removable dentures, appliances and prostheses	Fabrication and repair of removable dentures, appliances and prostheses	SAME INTENDED USE
<b>Indication For Use</b>	Dentca Denture Base is a light-cure resin indicated for fabrication and repair of full and partial removable dentures and baseplates. The material is an alternative to traditional heat-cured and auto polymerizing resins... The device requires a CAD/CAM system that includes the following components not part of the device: oral casting impression, digital denture base file created in an optical impression system, stereo lithographic additive printer, and curing light equipment.	Trubyte Denture Base Resin System is [a light-cure resin] indicated for fabrication and repair of dentures, appliances and prostheses, and for relining of denture surfaces.	BOTH DEVICES ARE LIGHT-CURE RESINS INDICATED FOR FABRICATION AND REPAIR OF REMOVABLE DENTAL PROSTHESES.

### Device Description

<b>Acrylic Resin</b>	Light-cure resin	Light-cure resin	SAME TECHNOLOGY
<b>Chemical Characterization</b>	Methacrylate-based resins with photoinitiator, inhibitor and pigments	Methacrylate-based resins with photoinitiator, inhibitor and pigments	SAME TECHNOLOGY
<b>Polymerization (Curing) Method</b>	Visible light	Visible light	SAME TECHNOLOGY
<b>Product State</b>	Pre-mixed resin	Pre-mixed resin	SAME TECHNOLOGY
<b>Fabrication of Denture Base</b>	<u>Automated 3D printing of resin</u> in multiple layers, each light-cured before adding next layer, with post curing in light chamber	<u>Manual application of resin</u> in multiple layers, each light-cured before adding next layer, with post curing in light chamber	<u>Manual vs. automated application technology</u>
<b>Post Curing</b>	Visible light-curing unit	Visible light-curing unit	SAME TECHNOLOGY
<b>Teeth Assemble</b>	Bonding	Bonding	SAME TECHNOLOGY

The above comparison shows the only technological difference between the denture polymers is method of resin application to fabricate the denture base, with the subject applied in an automated optical method and the predicate in a manual mechanical process.

A comparison of the subject and reference showed that both devices fabricate the denture base in the same manner, via automated application of a light-cure resin in an additive CAD/CAM printer.

## **Performance Data**

The following performance data were provided in support of the substantial equivalence determination.

### Biocompatibility Testing

The biocompatibility evaluation for Dentca Denture Base was conducted in accordance with the FDA Blue Book Memorandum #G95-1 and International Standard ISO 10993-1, as recognized by FDA. The battery of testing included the following tests:

- Genotoxicity
- Cytotoxicity
- Sensitization
- Irritation
- Acute Toxicity
- Material Characterization

The denture base is considered tissue contacting for a period longer than 30 days (a removable prosthesis).

### Electrical Safety and EMC

The low voltage 3D printer used to fabricate the denture base from the Dentca Denture Base light-cure resin was certified by the printer manufacturer to conform to EN 60950 (Safety of Information Technology Equipment) and the following EMC standards: EN 55022, EN 55024, EN 61000-3-2, EN 61000-3-3, and EN 60825-1.

### Software Verification and Validation Testing

The optical scanning impression system used to convert cast impressions into 3D digital images was validated by the scanner manufacturer in conformance with ISO 12836.

### Bench Testing

Dentca Denture Base was tested for conformity with the industry consensus standard ISO 20795-1.

## **Conclusions**

The subject and predicate devices have the same intended use and substantially similar technological characteristics, with the exception of the method of fabricating the denture base. Because the subject introduces a new additive printing technology as part of a CAD/CAM manufacturing process, an acceptable reference device was brought in to demonstrate prior use of this new technology.

The non-clinical data support the substantial equivalence of the subject device and demonstrate that Dentca Denture Base should perform as intended in the specified use conditions.