



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

January 26, 2015

BSN medical GmbH
% Ms. Marcia Palma
NAMSA
4050 Olson Memorial Highway, Suite 450
Minneapolis, Minnesota 55422

Re: K143151

Trade/Device Name: Cutimed[®] Sorbact[®] Wound Dressing (Cutimed[®] Sorbact[®]
Hydroactive, Cutimed[®] Sorbact[®] Hydroactive B,
Cutimed[®] Sorbact[®] Siltec Sorbact)

Regulatory Class: Unclassified

Product Code: FRO

Dated: November 1, 2014

Received: November 4, 2014

Dear Ms. Palma:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

David Krause -S

for Binita S. Ashar, M.D., M.B.A., F.A.C.S.
Director
Division of Surgical Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K143151

Device Name

Cutimed® Sorbact® Wound Dressings (Cutimed® Sorbact® Hydroactive, Cutimed® Sorbact® Hydroactive , Cutimed® Siltec Sorbact®)

Indications for Use (Describe)

Cutimed® Sorbact® Hydroactive is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Cutimed® Sorbact® Hydroactive B is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Cutimed® Siltec Sorbact® is indicated for use in the management of moderate to heavily exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Submitter:	Martin Spengler, Head of Regulatory Affairs, Legal Manufacturer Hamburg BSN medical GmbH Quickbornstraße 24, Hamburg, 20253 Germany Phone: +49 (040) 4909-7233 Fax: +49 (040) 4909-5476
Contact Person:	Marcia Palma Senior Principle Regulatory Consultant, NAMSA 4050 Olson Memorial Hwy, Suite 450 Minneapolis, MN 55422 Telephone: 763.588.9847 Fax: 763-287-3836 Email: mpalma@namsa.com
Date Prepared:	January 22, 2015
Trade Name:	Cutimed® Sorbact® Wound Dressings (Cutimed® Sorbact® Hydroactive, Cutimed® Sorbact® Hydroactive B, Cutimed® Siltec Sorbact®)
Classification:	Unclassified
Product Code:	FRO
Predicate Device(s):	The subject device is equivalent to the following device: K063059, Sorbact® Wound Dressings, Abigo Medical AB
Device Description:	Cutimed® Sorbact® Wound Dressings combine Sorbact mesh (Sorbact acetate fabric coated with dialkylcarbamoylechloride - DACC) with a highly absorbent hydropolymer matrix or superabsorbent polyurethane foam and are covered by a semi-permeable polyurethane film. A fixation border is made of a latex free hydrogel or silicone adhesive. The dressings come in three models (Cutimed® Sorbact® Hydroactive, Cutimed® Sorbact® Hydroactive B, Cutimed® Siltec Sorbact®), multiple sizes in each model and are sterile and single use.
Intended Use:	Cutimed® Sorbact® Hydroactive is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehiscenced wounds and traumatic wounds. Cutimed® Sorbact® Hydroactive B is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected

wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Cutimed® Siltec Sorbact® is indicated for use in the management of moderate to heavily exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Comparison to Predicate:	Product	Cutimed® Sorbact® Wound Dressings	Predicate: Sorbact® Wound Dressings
	Indications for Use	<p>Cutimed® Sorbact® Hydroactive is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.</p> <p>Cutimed® Sorbact® Hydroactive B is indicated for use in the management of low to moderate exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.</p> <p>Cutimed® Siltec Sorbact® is indicated for use in the management of moderate to heavily exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.</p>	<p>Sorbact® Wound Dressings are intended for use in the management of moderate to heavily exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds) and to bind hydrophobic microbes. The dressings are indicated for post-operative wounds, trauma wounds, shallow cavity wounds, fistulas, pressure ulcers, diabetic ulcers, and venous ulcers.</p>

	Device Design	The Cutimed® Sorbact® Wound Dressings are designed with an acetate fabric coated with DACC (Sorbact®) in combination with an absorbent core, which absorbs and locks exudate from the wound, and a transparent semi-permeable polyurethane film to allow vapor and oxygen to penetrate the film. It is also designed with an adhesive layer for fixation (except the Hydroactive dressing).	Same
	Wound Contact Material	Sorbact mesh (Sorbact acetate fabric coated with dialkylcarbamoylechloride – DACC)	Same
	Exudate Absorption Material	Hydropolymer matrix or polyurethane foam	Non-woven Viscose
	Transparent Semi-permeable backing	Polyurethane film	Same
	Self Adhesive Border	Yes (Except Cutimed® Sorbact® Hydroactive)	Same
	Mechanism of Action	Used for the management of clean, colonized, contaminated, or infected wounds; rapidly absorbs and locks exudate within an absorbent core; allows water vapor and gases to evaporate from the skin to the outside protecting the skin from maceration	Same
	Sizes	Multiple sizes	Same
	Wear Time	Up to 4 days	Change at least twice per week
	Sterility	Gamma radiation (SAL 10 ⁻⁶) for Cutimed® Sorbact® Hydroactive and Cutimed® Sorbact® Hydroactive B; EO (SAL 10 ⁻⁶) for Cutimed® Siltec Sorbact®	Gamma radiation (SAL 10 ⁻⁶)
	Single Use	Yes	Same
	Storage Conditions	Dry, below 25°C	Same
	Shelf Life	3 years	5 years
	Latex Free	Yes	Same
	Biocompatible	Yes	Same
Functional and Safety Testing:	To verify that device design met its functional performance and safety requirements, representative samples of the device		

underwent testing including bench testing (absorption, MVTR, retention), biocompatibility testing (cytotoxicity, irritation, sensitization, systemic toxicity, implantation, genotoxicity), packaging testing (pouch seal and transportation), sterilization validation testing, and stability testing (AA and real time).

Conclusion:

BSN medical GmbH considers the Cutimed® Sorbact® Wound Dressings to be equivalent to the predicate device listed above. This conclusion is based upon the device's similarities in intended use, design, mechanisms of action, technology and materials.