



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Medical Intelligence Medizintechnik GmbH  
% Mr. Michael Wolff  
Regulatory Affairs Manager  
Robert-Bosch-Strasse 8  
Schwabmunchen, 86830  
GERMANY

March 4, 2015

Re: K143485  
Trade/Device Name: iGUIDE 2.1  
Regulation Number: 21 CFR 892.5050  
Regulation Name: Medical charged particle radiation therapy system  
Regulatory Class: II  
Product Code: IYE  
Dated: November 4, 2014  
Received: December 8, 2014

Dear Mr. Wolff:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



For

Robert Ochs, Ph.D.  
Acting Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K143485

Device Name

iGUIDE 2.1

Indications for Use (Describe)

The intended use of the device is the control of accurate patient positioning with assistance of a 3D Tracking System in a radiotherapy environment.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration Office  
of Chief Information Officer Paperwork  
Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# Traditional 510(k) 510(k) Summary

## iGUIDE 2.1

### iGUIDE 2.1- 510(k) Summary

- |                            |  |
|----------------------------|--|
| 1. Applicant               | Medical Intelligence Medizintechnik GmbH   |
| 2. Address                 | Robert-Bosch-Straße 8<br>86830 Schwabmünchen<br>Germany  |
| 3. Contact Person          | Michael Wolff, +49 (0) 8232 9692 701   |
| 4. Preparation Date        | December 2 <sup>nd</sup> , 2014  |
| 5. Device Submitted        | iGUIDE <sup>®</sup> 2.1 System   |
| 6. Proprietary Name        | iGUIDE <sup>®</sup> 2.1 System   |
| 7. Common Name             | iGUIDE <sup>®</sup> 2.1  |
| 8. Classification Name     | Accelerator, Linear, Medical (21 CFR 892.5050, Product Code IYE)   |
| 9. Substantial Equivalence | <p>The iGUIDE<sup>®</sup> 2.1 is substantially equivalent to the following currently marketed device:</p> <p>iGUIDE<sup>®</sup> System (K072079)</p> <p>iGUIDE<sup>®</sup> 2.1 is an advancement of the iGuide<sup>®</sup> System and has similar technological characteristics as the predicate device.</p>   |
| Device Description         | <p>The iGUIDE 2.1 System is a powered radiation therapy support assembly which provides patient positioning control prior to radiotherapy treatment. iGUIDE 2.1 System consists of the iGUIDE workstation with software, a tracking system, reference frame, calibration tool, isolation transformer, KVM Extender and monitors for the workstation and terminal.</p> <p>The iGUIDE software controls the combined working of the HexaPOD evo RT System components and provides status detection via the software GUI to the user. The HexaPOD evo RT System consists of iGUIDE and the HexaPOD evo RT Couchtop, a robotic patient-positioning platform with six degrees of freedom, which is cleared under K072898.</p> |
| 10. Indications for Use    | The intended use of the device is to control accurate patient  |

# Traditional 510(k) 510(k) Summary

## iGUIDE 2.1

positioning with assistance of a 3D Tracking System in a radiotherapy environment.

### 11. Summary of Product Changes

- Compatibility with the Brainlab ExacTrac® V6.1 system
- Control of the Elekta Precise Treatment Table
- Enhanced usability and elimination of use errors through:
  - Automatic receipt of patient data
  - Use of Remote Automatic Table Movement interface; manual repositioning function is no longer supported
  - Automated data transfer from XVI to iGUIDE & Sentinel to iGUIDE

### 12. Summary of Similarities

iGUIDE® 2.1 is similar to the predicate device in terms of:

- Intended use
- product components
- technological characteristics (regarding design, material, energy source)

### 13. Comparison with predicate device

<b>Manufacturer</b>	Medical Intelligence Medizintechnik GmbH	Medical Intelligence Medizintechnik GmbH
<b>Device Name</b>	iGUIDE 2.1 System	iGUIDE System
PreMarket Notification Number	Proposed	K072079
Indications for Use	The intended use of the device is to control accurate patient positioning with assistance of a 3D Tracking System in a radiotherapy environment.	The intended use of the device is to control accurate patient positioning with assistance of a 3D Tracking System in a radiotherapy environment.
Components	<ul style="list-style-type: none"> <li>• Software</li> <li>• Tracking System</li> <li>• Workstation</li> <li>• Monitor</li> <li>• KVM Extender</li> <li>• Isolation Transformer</li> <li>• Reference Frame</li> <li>• Calibration Tool</li> </ul>	<ul style="list-style-type: none"> <li>• Software</li> <li>• Tracking System</li> <li>• Workstation</li> <li>• Monitor</li> <li>• KVM Extender</li> <li>• Isolation Transformer</li> <li>• Reference Frame</li> <li>• Calibration Tool</li> </ul>
Protection Class	I B	I B

# Traditional 510(k)

## 510(k) Summary

### iGUIDE 2.1

<b>Manufacturer</b>	Medical Intelligence Medizintechnik GmbH	Medical Intelligence Medizintechnik GmbH
<b>Device Name</b>	iGUIDE 2.1 System	iGUIDE System
Software Safety Classification	C (Major Level of Concern)	C (Major Level of Concern)
Hardware Platform	Part of the HexaPOD evo RT System - consisting of HexaPOD evo Couchtop (K072898) and iGUIDE components (see below)	Part of the HexaPOD RT System - consisting of HexaPOD Couchtop (K071870) and iGUIDE components (see below)
Position control by 3D tracking System	yes	yes
HexaPOD evo RT Couch Top as treatment couch	yes	yes
Movement Control	yes	yes
Table Movement	Remote Automatic Table Movement	Manual prepositioning function
Accuracy of 3D Tracking System	0.35mm	0.35mm
Accuracy of total system (Treatment table + tracking system)	< 1mm; CI95: <0,5mm	< 0.5 mm (translation) < 0.25° (rotation)
Receiving patient data	Automatic data transfer to iGUIDE	Manual data transfer to iGUIDE
Sending patient data	Automated data transfer from iGUIDE	Manual data transfer from iGUIDE
Enable switches	<ul style="list-style-type: none"> <li>• Interrupts directly the power of the moving component</li> <li>• In control room and treatment room</li> </ul>	<ul style="list-style-type: none"> <li>• Interrupts directly the power of the moving component</li> <li>• In control room and treatment room</li> </ul>
Interlock connection to Linac	Overrides the enable switch control	Overrides the enable switch control

#### 14. Summary of Performance Testing

The Performance Testing is based on Benchmark Testing only.

Rationale: As the iGUIDE 2.1 is a positioning device which does not perform any treatment on its own but helps to position and re-position the patient, only physical parameters characterize the system (e.g. accuracy) and give information about substantial equivalence of the modified system compared to the predicate (unmodified) system which was cleared under K072079.

The results of the verification and validation and bench testing demonstrate that the iGUIDE 2.1 System with the extended compatibility and improved usability satisfies the intended use as described above and supports the claim of substantial equivalence.