



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 23, 2015

Jedmed Instrument Company  
% Mr. Craig R. Parks  
Regulatory Affairs/QA Manager  
5416 JEDMED Court  
St. Louis, Missouri 63129

Re: K150248

Trade/Device Name: ORL Video Nasopharyngo-laryngoscope  
Regulation Number: 21 CFR 874.4760  
Regulation Name: Nasopharyngoscope (Flexible Or Rigid) And Accessories  
Regulatory Class: Class II  
Product Code: EQN  
Dated: February 3, 2015  
Received: June 25, 2015

Dear Mr. Parks:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"

(21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Eric A. Mann -S**

for Malvina B. Eydelman, M.D.

Director

Division of Ophthalmic and Ear, Nose,  
and Throat Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K150248

Device Name

ORL Video Nasopharyngo-Laryngoscope

Indications for Use (Describe)

Jedmed ORL Video Nasopharyngo-Laryngoscopes are intended to examine the larynx, nasal cavity and nasal pharynx. They are used between the upper respiratory tracts of the nasal passage and vocal cords.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



*An Employee-Owned Company*

5416 JEDMED COURT • ST. LOUIS, MO 63129 PHONE: (314) 845-3770 • FAX: (314) 845-3771

E-MAIL: [INFO@JEDMED.COM](mailto:INFO@JEDMED.COM) WEB SITE: [www.jedmed.com](http://www.jedmed.com)

## 6. 510(k) - Summary

Submitter: JEDMED Instrument Company  
5416 JEDMED Court  
St. Louis, MO 63129

Contact: Craig Parks  
Regulatory Affairs /QA Manager  
Phone: 314-845-3770  
Email: [craigp@jedmed.com](mailto:craigp@jedmed.com)

Date: 07-21-2015

Name of Device: ORL Video NasoPharyngo-Laryngoscope

Common Name: Naso-pharyngo-laryngo-flexible scope

Classification Name of Device: Laryngo, Nasopharyngoscope  
a) product code: EQN  
b) regulation number: 874.4760

Legally Marketed Device Flexible naso-pharyngo-laryngoscope to which Equivalence is claimed  
JEDMED Instrument Company K132039  
ERGOFLEX NasoPharyngolaryngoscope  
Vision Science Flexible ENT 5000 K102733

Claimed:

Description: This flexible scope is designed with the insertion tube with its bendable distal tip w/camera, the handle, the umbilical unit which has the leakage tester, and connection for light source port and power/video port The handle incorporates the control lever to bend the distal tip.

Indications for Use: Jedmed ORL Video Nasopharyngo-Laryngoscopes are intended to examine the larynx, nasal cavity and nasal pharynx. They are used between the upper respiratory tracts of the nasal passage and vocal cords.



*An Employee-Owned Company*

5416 JEDMED COURT • ST. LOUIS, MO 63129 PHONE: (314) 845-3770 • FAX: (314) 845-3771

E-MAIL: [INFO@JEDMED.COM](mailto:INFO@JEDMED.COM) WEB SITE: [www.jedmed.com](http://www.jedmed.com)

## 6. 510(k) - Summary

Similarities/differences  
to predicates:

The submitted devices are equivalent to the flexible endoscopes referred to in K132039 Ergoflex NasoPharyngo-Laryngoscope .

The fiberscopes in the submission use the same basic design and device material as predicate device. This device does have an eyepiece for visualization, and can be connected to a camera system while the submitted device incorporates the camera in the tip and does not have direct visualization

The submitted devices are equivalent to the Flexible ENT 5000 scope from Vision Sciences 510(k)# 102733. Both devices use a camera system for visualization without the eyepiece. There are only differences regarding diameters and lengths of the device insertion tube and the outer appearance and shaping of the shell parts and control elements which has no bearing on the safety or effectiveness.

All the devices have the same intended use, to examine the larynx, nasal cavity and nasal pharynx. They are used between the upper respiratory tracts of the nasal passage and the vocal cords

Testing: The ORL Video Naso-Pharyngo-Laryngoscope was tested functionally and compared to the predicate device the ErgoFlex Naso-Pharyngo-Laryngoscope 510(k)# K132039 and found to be comparable. The test included Field of View, Depth of View, Angulation of the Bendable Tip, Color Output, Clarity, and Tip Diameter. Both scopes were tested with the same test fixtures and results compared by our knowledgeable field staff.

A cleaning and high-level disinfection test was also performed on the ORL Video Naso-Pharyngo-Laryngoscope for both Cidex 2.4% Activated Dialdehyde and Cidex OPA disinfectant solutions. The test dictated that a significant level of microbial reduction must result from tests for both vegetative microorganisms and mycobacteria. Based on the test criteria, the general disinfection processing employing both Cidex 2.4% Activated Dialdehyde and Cidex OPA (performed by Geneva Laboratories per the IFU) fully met the validation acceptance requirements for the Jedmed ORL Laryngo-Nasopharyngoscope.

The submitted devices pose the same type of questions about safety and effectiveness as the compared devices. The different technological characteristics have no influence on safety or effectiveness, therefore our device is substantially equivalent to the previously cleared device.

	<b>Device under Review</b>	<b>Predicate Device K132039</b>	<b>Predicate Device K102733</b>
<b>Device Name</b>	ORL Video NasoPharyngo-Laryngoscope	ERGOFLEX NasoPharyngo-Laryngoscope JEDMED Instrument	Vision Science ENT-5000
Indication for Use	Examination of the larynx, nasal cavity, nasal pharynx and is used between the upper nasal respiratory tracts and of the nasal passage and vocal chords.	Same	Same with the addition of arthroscopic and endoscopic visualization of an interior cavity of the body thru natural or surgical openings.
<b>Viewing Method</b>	Camera System	Same	Same
<b>Working Length</b>	310mm	310mm	300mm
<b>Field of View</b>	90°	90°	90°
<b>Depth of View</b>	10mm-55mm	10-55mm	3mm-50mm
<b>Tube Diameter</b>	4mm	4mm	3mm
<b>Angulation</b>	130° Up 130° Down	130° Up 130° Down	140° Up 140° Down
<b>Video Format</b>	NTSC---PAL/NTSC	NTSC---PAL	CCD
<b>Resolution</b>	320x240 500x582 510x492	320x240 & 500x582/510x492	
<b>Product Code</b>	EQN	EQN	EOB and HRX and GCJ
<b>Light Guide Connector</b>	ACMI-Wolf-Storz-Olympus	ACMI-Wolf-Storz-Olympus	ACMI-Wolf-Storz-Olympus
<b>Direct View Eyepiece</b>	None	None	None