



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

June 22, 2015

Summit Medical Inc.
Ms. Nicole Dove
Quality Assurance/Regulatory Affairs Manager
815 Northwest Pkwy, Suite 100
St. Paul, MN 55121

Re: K150540

Trade/Device Name: Instru-Safe[®] Instrument Protection System
Regulation Number: 21 CFR 880.6850
Regulation Name: Sterilization Wrap
Regulatory Class: II
Product Code: KCT
Dated: May 20, 2015
Received: May 22, 2015

Dear Ms. Dove:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Tejashri Purohit-Sheth, M.D.

Tejashri Purohit-Sheth, M.D.
Clinical Deputy Director
DAGRID/ODE/CDRH FOR

Erin I. Keith, M.S.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K150540

Device Name

Instru-Safe® Instrument Protection System

Indications for Use (Describe)

Instru-Safe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instru-Safe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during Amsco V-PRO Low Temperature Sterilization Cycles. The Instru-Safe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap or Aesculap rigid container. The Instru-Safe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in table 1.

Amsco V-PRO Low Temperature Sterilization Cycles:

Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle
V-PRO 1	X	N/A	N/A
V-PRO 1 Plus	N/A	X	X
V-PRO maX	N/A	X	X

Summit Cassette Model

IN-8823

IN-6105

Aesculap Container Model

*JM444

*JM742

*Validated by Summit Medical for use in Amsco V-PRO Low Temperature Sterilization systems ONLY. When using the Aesculap container as a sterile barrier, the load (Summit tray and contents), should not exceed the load claims for the container in weight or load type.

Lumen size of instrumentation validated includes:

Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens
IN-8823	3 mm	400 mm	2
IN-6105	3 mm	200 mm	1
IN-2681	1 mm	64 mm	1

The worst case validated load by vent-to-volume calculation is the IN-2681 tray.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”



Indications for Use Statement

Table 1 – Device Models

Part Number	Maximum # of Instruments	Estimated Weight - Tray w/ instruments (lbs)
IN-1315	30	3.5
IN-2006	8	2
IN-2681	13	0.64
IN-2682	13	0.89
IN-2683	13	1.39
IN-2840	36	8.75
IN-2842	24	6.5
IN-2843	36	8.75
IN-2880	56	12.1
IN-2900	22	4.18
IN-2950	12	4.1
IN-3030	34	9.5
IN-4000	20	3.2
IN-4003	30	3.25
IN-4010	10	1.98
IN-5009	8	5
IN-5401	6	2
IN-5401-02	2	1
IN-5401-03	2	1
IN-5401-08	8	3.2
IN-5401-12	12	3.25
IN-6103	2	2.15
IN-6105	2	2.15
IN-6110	3	2.15
IN-6203	2	2.75
IN-6205	2	2.75
IN-6210	2	2.75
IN-6240	2	2.75
IN-6303	2	3.28
IN-6305	2	3.28
IN-6310	2	3.28
IN-6403	2	3.28
IN-6405	2	3.28
IN-6410	2	3.28



IN-6500	30	12.5
IN-7010	2	2
IN-7012	1	1.07
IN-7020	30	8.5
IN-7030	30	10
IN-7032	2	1.1
IN-7073	10	5
IN-7120	45	11.25
IN-7123	45	12
IN-7130	45	13.5
IN-7140	45	14.5
IN-7150	8	1.9
IN-7153	6	1.7
IN-7220	45	14.5
IN-7223	10	9.2
IN-7230	45	14.5
IN-7234	45	14.5
IN-7240	45	14.5
IN-7250	45	14.5
IN-7252	25	8
IN-7260	45	14.5
IN-7273	10	6
IN-7274	30	8
IN-7322	45	14.5
IN-7323	45	14.5
IN-7343	45	14.5
IN-7344	1	4
IN-7360	45	14.5
IN-7423	45	14.5
IN-7452	10	8
IN-7453	10	8
IN-7540	45	14.5
IN-7560	45	14.5
IN-7644	45	14.5
IN-7723	15	7.18
IN-7724	15	7.2
IN-7725	10	9.5
IN-7781	45	14.5



IN-7823	45	14.5
IN-7830	45	14.5
IN-7840	45	13.5
IN-7940	20	13.25
IN-8240	20	13.5
IN-8610	2	6.65
IN-8612	2	6.8
IN-8613	2	6.1
IN-8615	2	5.8
IN-8616	2	5.8
IN-8620	3	7.2
IN-8621	4	7.18
IN-8622	4	7.18
IN-8630	3	6.5
IN-8632	3	6.45
IN-8633	3	6.8
IN-8640	4	5.35
IN-8642	4	5.35
IN-8643	5	5.35
IN-8645	4	5.35
IN-8650	4	5.85
IN-8660	4	5.35
IN-8662	4	5.35
IN-8663	3	5.5
IN-8700	40	14
IN-8810	20	13.5
IN-8820	15	8.75
IN-8823	45	14
IN-8830	15	8.75
IN-8833	45	14
IN-8840	20	13.75
IN-8850	15	8.75
IN-8853	45	14
IN-8860	15	8.75
IN-8862	30	10.5
IN-8863	45	14
IN-8880	2	3.28
IN-8882	16	12.1



IN-8883	2	3.28
IN-8884	4	5.35
IN-8885	1	2.25
IN-8886	6	12.1
IN-8889	6	12.1
IN-8891-S	1	2
IN-8891-SI-12-S	1	2
IN-8891-SI-85-S	1	2
IN-8892	12	12.1
IN-8893	9	7.5
IN-8894	5	16.1
IN-8897	8	6
IN-8898	10	10.25
IN-8899	7	6.5
IN-8901	1	2.25
IN-8902-G2	22	17
IN-8903	15	13.25
IN-8904	22	17
IN-8907	7	12.5
IN-8931	1	2.4
IN-8932	9	9.5
IN-8933	3	3.75
IN-8936	6	11.5
IN-8937	16	14.5
IN-8938	8	12.5
IN-8939	10	11.6
IN-8940	5	5.18
IN-8942	11	10
IN-8943	1	2.7
IN-8944	6	4.7
IN-8945	2	5.18
IN-8946	9	6.1
IN-8980	20	9.5
IN-8982-01	17	9.5
IN-8983-01	16	9.5
IN-8984-01	15	9.5
IN-8986-S	2	6.5
IN-8987-S	2	6.5



IN-8988-S	2	6
IN-8989-S	2	6
IN-9999-160	6	12.1
IN-9999-162-S	2	5.8
IN-9999-168-S	2	5.8
IN-9999-172-S	2	5.8
IN-9999-178-S	2	5.8



510(k) Summary

Following is a summary of 510(k) safety and effectiveness information in accordance with 21 CFR 807.92

Submitter:	Summit Medical Inc. 815 Northwest Parkway, Suite 100 St. Paul, MN 55121 Tel: (651) 789-3939
ER Number:	3008719017
Contact Person:	Nicole Dove QA/RA Manager Tel: (651) 789-3921 ndove@summitmedicalusa.com
Date Prepared:	June 22, 2015
Subject Device:	<u>Trade Name(s):</u> Instru-Safe [®] Instrument Protection System <u>Classification Name:</u> Sterilization wrap containers, trays, cassettes & other accessory (21 CFR 880.6850) <u>Common Name:</u> Instrument Tray, Sterilization Tray, Sterilization Cassettes, Instrument Delivery System <u>Device Class:</u> Class II <u>Device Code:</u> KCT <u>Panel:</u> General Hospital
Predicate Device:	Tradename: Instru-Safe Instrument Protection System 510(k) Holder: Summit Medical Inc. 510(k) #: K133015
Device Description:	Summit Medical Inc. Instru-Safe Instrument Protection System are cassettes / trays used to enclose and hold surgical instruments and accessories in an organized manner during the sterilization process and subsequent storage and transportation. The cassettes / trays by themselves do not maintain sterility.



	<p>The cassettes / trays are different sizes of the same basic configuration: a rectangular base with latchable cover. The cassettes / trays have perforations to allow sterilant penetration. The cassettes / trays contain silicone inserts in the base and/or cover to hold, organize and protect the surgical instruments within the cassette / tray.</p>																																							
Intended Use:	<p>Instru-Safe[®] Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instru-Safe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during Amsco V-PRO Low Temperature Sterilization Cycles. The Instru-Safe Instrument Protection System cassettes are intended to be used in conjunction with legally marketed wrap or Aesculap rigid container. The Instru-Safe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in table 1.</p> <p>Amsco V-PRO Low Temperature Sterilization Systems</p> <table border="1"> <thead> <tr> <th>Sterilizer</th> <th>Standard Cycle</th> <th>Lumen Cycle</th> <th>Non Lumen Cycle</th> </tr> </thead> <tbody> <tr> <td>V-PRO 1</td> <td>X</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>V-PRO 1 Plus</td> <td>N/A</td> <td>X</td> <td>X</td> </tr> <tr> <td>V-PRO maX</td> <td>N/A</td> <td>X</td> <td>X</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th>Summit Cassette Model</th> <th>Aesculap Container Model</th> </tr> </thead> <tbody> <tr> <td>IN-8823</td> <td>*JM444</td> </tr> <tr> <td>IN-6105</td> <td>*JM742</td> </tr> </tbody> </table> <p>*Validated by Summit Medical for use in Amsco V-PRO Low Temperature Sterilization Systems ONLY. When using the Aesculap container as a sterile barrier, the load (Summit tray and contents), should not exceed the load claims for the container in weight or load type.</p> <p>Lumen size of instrumentation validated includes:</p> <table border="1"> <thead> <tr> <th>Summit Cassette Model</th> <th>Minimum Inside Diameter</th> <th>Maximum Length</th> <th>Number of Lumens</th> </tr> </thead> <tbody> <tr> <td>IN-8823</td> <td>3mm</td> <td>400mm</td> <td>2</td> </tr> <tr> <td>IN-6105</td> <td>3mm</td> <td>200mm</td> <td>1</td> </tr> <tr> <td>IN-2681</td> <td>1mm</td> <td>64mm</td> <td>1</td> </tr> </tbody> </table> <p>The worst case validated load by vent-to-volume calculation is the IN-2681 tray.</p> <p>The intended use of the subject device includes the Amsco V-PRO Low Temperature Sterilization Systems. Performance testing has been completed for the Amsco V-PRO Low Temperature Sterilization Systems. These new sterilization cycles do not affect safety and effectiveness of the Instru-Safe Instrument Protection System.</p>		Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle	V-PRO 1	X	N/A	N/A	V-PRO 1 Plus	N/A	X	X	V-PRO maX	N/A	X	X	Summit Cassette Model	Aesculap Container Model	IN-8823	*JM444	IN-6105	*JM742	Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens	IN-8823	3mm	400mm	2	IN-6105	3mm	200mm	1	IN-2681	1mm	64mm	1
Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle																																					
V-PRO 1	X	N/A	N/A																																					
V-PRO 1 Plus	N/A	X	X																																					
V-PRO maX	N/A	X	X																																					
Summit Cassette Model	Aesculap Container Model																																							
IN-8823	*JM444																																							
IN-6105	*JM742																																							
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens																																					
IN-8823	3mm	400mm	2																																					
IN-6105	3mm	200mm	1																																					
IN-2681	1mm	64mm	1																																					
Comparison of Characteristics to Predicate Device:	<p>Based on a comparison of the design, technology, materials, manufacturing, performance, specifications and methods of use, the Instru-Safe Instrument Protection System is equivalent to the identified 510(k) cleared predicate device.</p>																																							
Element	New Device	Predicate (K133015)																																						



<p>Intended Use</p>	<p>Instru-Safe[®] Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instru-Safe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during Amsco V-PRO Low Temperature Sterilization Cycles. The Instru-Safe Instrument Protection System cassettes are intended to be used in conjunction with legally marketed wrap or Aesculap rigid container. The Instru-Safe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in table 1. Sterilization methods and configurations</p> <p>Amsco V-PRO Low Temperature Sterilization Systems</p> <table border="1" data-bbox="414 1098 755 1220"> <thead> <tr> <th>Sterilizer</th> <th>Standard Cycle</th> <th>Lumen Cycle</th> <th>Non Lumen Cycle</th> </tr> </thead> <tbody> <tr> <td>V-PRO I</td> <td>X</td> <td>NA</td> <td>NA</td> </tr> <tr> <td>V-PRO I Plus</td> <td>NA</td> <td>X</td> <td>X</td> </tr> <tr> <td>V-PRO maX</td> <td>NA</td> <td>X</td> <td>X</td> </tr> </tbody> </table> <table border="1" data-bbox="402 1255 857 1371"> <thead> <tr> <th>Summit Cassette Model</th> <th>Aesculap Container Model</th> </tr> </thead> <tbody> <tr> <td>IN-8823</td> <td>*JM444</td> </tr> <tr> <td>IN-6105</td> <td>*JM742</td> </tr> </tbody> </table> <p>*Validated by Summit Medical for use in Amsco V-PRO Low Temperature Sterilization Systems ONLY. When using the Aesculap container as a sterile barrier, the load (Summit tray and contents), should not exceed the load claims for the container in weight or load type.</p>	Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle	V-PRO I	X	NA	NA	V-PRO I Plus	NA	X	X	V-PRO maX	NA	X	X	Summit Cassette Model	Aesculap Container Model	IN-8823	*JM444	IN-6105	*JM742	<p>Instru-Safe Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instru-Safe System cassettes are intended to allow sterilization of the enclosed medical devices during a pre-vacuum steam sterilization cycle. The Instru-Safe System cassettes are intended to be used in conjunction with legally marketed wrap or Aesculap rigid container. The Instru-Safe System cassettes are not intended on their own to maintain sterility.</p> <p>Sterilization methods and configurations</p> <ul style="list-style-type: none"> Autoclave Sterilization Parameter: <ul style="list-style-type: none"> Cycle: Pre-vacuum Temperature: 270°F (132°C) Exposure Time: 4 minutes Minimum Dry Time: 30 minutes <table border="1" data-bbox="919 1064 1334 1236"> <thead> <tr> <th>Summit Cassette Model</th> <th>Aesculap Container Model</th> </tr> </thead> <tbody> <tr> <td>IN-8823-AE</td> <td>*JN444</td> </tr> <tr> <td>IN-2880</td> <td>*JK444</td> </tr> <tr> <td>IN-6105</td> <td>*JN742</td> </tr> </tbody> </table> <p>*Validated by Summit Medical for use in steam prevacuum sterilizers ONLY operating at 270°F (132°C) for 4 minutes exposure time. Consult container instructions to ensure that contents do not exceed the sterilization container intended load claims.</p>	Summit Cassette Model	Aesculap Container Model	IN-8823-AE	*JN444	IN-2880	*JK444	IN-6105	*JN742
Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle																													
V-PRO I	X	NA	NA																													
V-PRO I Plus	NA	X	X																													
V-PRO maX	NA	X	X																													
Summit Cassette Model	Aesculap Container Model																															
IN-8823	*JM444																															
IN-6105	*JM742																															
Summit Cassette Model	Aesculap Container Model																															
IN-8823-AE	*JN444																															
IN-2880	*JK444																															
IN-6105	*JN742																															

	<p>Lumen size of instrumentation validated includes:</p> <table border="1"> <thead> <tr> <th>Summit Cassette Model</th> <th>Minimum Inside Diameter</th> <th>Maximum Length</th> <th>Number of Lumens</th> </tr> </thead> <tbody> <tr> <td>IN-8823</td> <td>3mm</td> <td>400mm</td> <td>2</td> </tr> <tr> <td>IN-6105</td> <td>3mm</td> <td>200mm</td> <td>1</td> </tr> <tr> <td>IN-2681</td> <td>1mm</td> <td>64mm</td> <td>1</td> </tr> </tbody> </table> <p>The worst case validated load by vent-to- volume calculation is the IN-2681 tray.</p>	Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens	IN-8823	3mm	400mm	2	IN-6105	3mm	200mm	1	IN-2681	1mm	64mm	1	
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens															
IN-8823	3mm	400mm	2															
IN-6105	3mm	200mm	1															
IN-2681	1mm	64mm	1															
Material Composition	No changes from predicate device	The cassette contains components made of anodized aluminum, stainless steel, blue																
Physical Properties	<p>Instru-Safe Instrument Protection System cassettes include</p> <ul style="list-style-type: none"> - perforated base - perforated cover - silicone inserts (hold-it / hold down) - Handles - Latches - Feet - Posts (optional) - Divider (optional) - Shelf (optional) 	<p>Instru-Safe Instrument Protection System cassettes include</p> <ul style="list-style-type: none"> - perforated base - perforated cover - silicone inserts (hold-it / hold down) - Handles - Latches - Feet 																
Chemical Properties	Not Applicable	Not Applicable																
Configurations / Dimensions	Various configurations / dimensions	See table located in predicate device																
Air permeance	Not Applicable	Not Applicable																

Percent of surface performances	Not Applicable	Not Applicable
Performance	New Device	Predicate
Sterilant Penetration	Amsco V-PRO Low Temperature Sterilization Systems.	Pre-Vacuum Steam Cycle: Pre-vacuum Temperature: 270°F (132°C)
Microbial Barrier Properties (Packaging Integrity)	Not Applicable	Not Applicable
Material Compatibility	No changes from predicate device	Refer to predicate device K133015
Toxicological Properties (Biocompatibility, including	MEM Elution Cytotoxicity (ISO 10993-5) - The test samples meet the USP and ISO 10993-5 requirements for this	Refer to predicate device K133015



Sterilant Residue Limits)	test. All controls were acceptable and the test considered valid. The test samples PASSED and are considered NON-CYTOTOXIC under the test conditions employed.	
Shelf Life	No Change	Reusable (5 year accelerated shelf life study)
Drying Time	Not Applicable	Autoclave Sterilization Parameter: Cycle: Pre-Vacuum Temperature: 270°F (132°C) Minimum Exposure Time: 4 minutes Minimum Dry Time: 30 minutes
Aeration Time	Not Applicable	Not Applicable
Technological Characteristics:	The technological characteristics of the subject devices are equivalent to the predicate devices. The cassettes / trays are made of standard medical grade materials and do not incorporate any new technological characteristics.	
Performance Data:	Sterilization validation testing was performed to demonstrate Instru-Safe Instrument Protection System compatibility when used in an Amsco V-PRO Low Temperature Sterilization Systems with a legally marketed wrap or Aesculap rigid container.	
Conclusion:	Based upon intended use, performance data and technical information provided in this pre-market notification, the Instru-Safe Instrument Protection System described herein is substantially equivalent to the predicate device [Instru-Safe Instrument Protection System (K133015)].	



Table 1 – Device Models

Part Number	Maximum # of Instruments	Estimated Weight - Tray w/ instruments (lbs)
IN-1315	30	3.5
IN-2006	8	2
IN-2681	13	0.64
IN-2682	13	0.89
IN-2683	13	1.39
IN-2840	36	8.75
IN-2842	24	6.5
IN-2843	36	8.75
IN-2880	56	12.1
IN-2900	22	4.18
IN-2950	12	4.1
IN-3030	34	9.5
IN-4000	20	3.2
IN-4003	30	3.25
IN-4010	10	1.98
IN-5009	8	5
IN-5401	6	2
IN-5401-02	2	1
IN-5401-03	2	1
IN-5401-08	8	3.2
IN-5401-12	12	3.25
IN-6103	2	2.15
IN-6105	2	2.15
IN-6110	3	2.15
IN-6203	2	2.75
IN-6205	2	2.75
IN-6210	2	2.75
IN-6240	2	2.75
IN-6303	2	3.28
IN-6305	2	3.28
IN-6310	2	3.28
IN-6403	2	3.28
IN-6405	2	3.28
IN-6410	2	3.28



IN-6500	30	12.5
IN-7010	2	2
IN-7012	1	1.07
IN-7020	30	8.5
IN-7030	30	10
IN-7032	2	1.1
IN-7073	10	5
IN-7120	45	11.25
IN-7123	45	12
IN-7130	45	13.5
IN-7140	45	14.5
IN-7150	8	1.9
IN-7153	6	1.7
IN-7220	45	14.5
IN-7223	10	9.2
IN-7230	45	14.5
IN-7234	45	14.5
IN-7240	45	14.5
IN-7250	45	14.5
IN-7252	25	8
IN-7260	45	14.5
IN-7273	10	6
IN-7274	30	8
IN-7322	45	14.5
IN-7323	45	14.5
IN-7343	45	14.5
IN-7344	1	4
IN-7360	45	14.5
IN-7423	45	14.5
IN-7452	10	8
IN-7453	10	8
IN-7540	45	14.5
IN-7560	45	14.5
IN-7644	45	14.5
IN-7723	15	7.18
IN-7724	15	7.2
IN-7725	10	9.5



IN-7781	45	14.5
IN-7823	45	14.5
IN-7830	45	14.5
IN-7840	45	13.5
IN-7940	20	13.25
IN-8240	20	13.5
IN-8610	2	6.65
IN-8612	2	6.8
IN-8613	2	6.1
IN-8615	2	5.8
IN-8616	2	5.8
IN-8620	3	7.2
IN-8621	4	7.18
IN-8622	4	7.18
IN-8630	3	6.5
IN-8632	3	6.45
IN-8633	3	6.8
IN-8640	4	5.35
IN-8642	4	5.35
IN-8643	5	5.35
IN-8645	4	5.35
IN-8650	4	5.85
IN-8660	4	5.35
IN-8662	4	5.35
IN-8663	3	5.5
IN-8700	40	14
IN-8810	20	13.5
IN-8820	15	8.75
IN-8823	45	14
IN-8830	15	8.75
IN-8833	45	14
IN-8840	20	13.75
IN-8850	15	8.75
IN-8853	45	14
IN-8860	15	8.75
IN-8862	30	10.5
IN-8863	45	14



IN-8880	2	3.28
IN-8882	16	12.1
IN-8883	2	3.28
IN-8884	4	5.35
IN-8885	1	2.25
IN-8886	6	12.1
IN-8889	6	12.1
IN-8891-S	1	2
IN-8891-SI-12-S	1	2
IN-8891-SI-85-S	1	2
IN-8892	12	12.1
IN-8893	9	7.5
IN-8894	5	16.1
IN-8897	8	6
IN-8898	10	10.25
IN-8899	7	6.5
IN-8901	1	2.25
IN-8902-G2	22	17
IN-8903	15	13.25
IN-8904	22	17
IN-8907	7	12.5
IN-8931	1	2.4
IN-8932	9	9.5
IN-8933	3	3.75
IN-8936	6	11.5
IN-8937	16	14.5
IN-8938	8	12.5
IN-8939	10	11.6
IN-8940	5	5.18
IN-8942	11	10
IN-8943	1	2.7
IN-8944	6	4.7
IN-8945	2	5.18
IN-8946	9	6.1
IN-8980	20	9.5
IN-8982-01	17	9.5
IN-8983-01	16	9.5



IN-8984-01	15	9.5
IN-8986-S	2	6.5
IN-8987-S	2	6.5
IN-8988-S	2	6
IN-8989-S	2	6
IN-9999-160	6	12.1
IN-9999-162-S	2	5.8
IN-9999-168-S	2	5.8
IN-9999-172-S	2	5.8
IN-9999-178-S	2	5.8