



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

October 16, 2015

Everyway Medical Instruments Co., Ltd.
Robert Tu
President And Operator Owner
3 Fl., No. 5, Lane 155, Section 3, Peishen Rd.,
Shenkeng Hsiang, Taipei Hsien, 222 TW

Re: K151479

Trade/Device Name: EverywayEV-820 OTC Pain Relief TENS
Regulation Number: 21 CFR 882.5890
Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief
Regulatory Class: Class II
Product Code: NUH
Dated: September 12, 2015
Received: September 15, 2015

Dear Mr. Tu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Michael J. Hoffmann -S

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K151479

Device Name

Everyway EV-820 OTC Pain Relief TENS

Indications for Use (Describe)

The Everyway EV-820 OTC Pain Relief TENS is intended for temporary relief of pain associated with sore and aching muscles in the low back as well as upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(K) SUMMARY

This summary of 510(k) safety and effectiveness information is being prepared in accordance with the requirements of SMDA 1990 and 21 CFR 807.92 at October 01, 2015.

The assigned 510(k) number is: K151479 .

1. **Submitter's Identifications:**

Establishment: EVERYWAY MEDICAL INSTRUMENT CO., LTD.

Address: 3Fl., No. 5, Lane 155, Sec. 3, Peishen Rd., Shenkeng Hsiang, Taipei Hsien
222, Taiwan

Registration Number: 9616877

Operations: Manufacturer

Owner/Operator: EVERYWAY MEDICAL INSTRUMENT CO., LTD.

Address : 3Fl., No. 5, Lane 155, Sec. 3, Peishen Rd., Shenkeng Hsiang, Taipei
Hsien 222, Taiwan

Contact Person: Robert Tu

Phone : 886-2-2662-0038

Fax No: 886-2-2664-5566

e-mail : tu922@ms35.hinet.net

2. **Name of the Device:**

Everyway EV-820 OTC Pain Relief TENS.

3. **Information of the 510(k) Cleared Device (Predicate Device):**

K133723: OTC Patch/ Model WL-2406.

K091757: OTC TENS For Arm & Leg Pain Relief/ Model WL-2407

4. **Classification Information:**

Trade/Device Name: Everyway EV-820 OTC Pain Relief TENS.

Regulation Number: 21 CFR 882.5890

Classification Name: Stimulator, Nerve, Transcutaneous For Pain Relief.

Regulatory Class: II

Product Code: NUH

5. **Device Description:**

The Everyway EV-820 OTC Pain Relief TENS is a dual channel transcutaneous electrical nerve stimulator used for pain relief by applying an electrical current to electrodes, which are attached on the user's skin. The output and waveform characteristic is fixed for every operation mode, only the intensity is adjustable within specified limit.

The Everyway EV-820 OTC Pain Relief TENS consists mainly of two parts: the stimulus generator, electrode. The stimulus generator generates the output current specified as the input of controller. The output port transmits the output current to the electrode, which is attached to the user's skin so as to transmit this stimulus current for pain relief.

The stimulation modes for Everyway EV-820 OTC Pain Relief TENS is pre-program modes with fixed pulse width, pulse rate, frequency. The timer and amplitude for each program are adjustable. This operation way is considered the simple design change from the pre-program modes of comparison clear model, existing 510(k) cleared devices, Low Back Pain Relief System/ Model EV-820(K110716). Every operation mode of Everyway EV-820 OTC Pain Relief TENS System has its individual stimulation operation cycle.



3Fl., No. 5, Lane 155, Sec. 3, Peishen Rd, Shen Keng Hsiang, Taipei Hsien, Taiwan,

For the device included in this submission, we use the following of our 510(K) legally marketed electrodes: K062675, "Everyway Lifecare Electrode", Wire Series/ model no. TFK5050, size 50x50mm, wire type as standard accessory.

With the combination of the main device parts, the device can be placed on the treatment locations as recommended in the user manual for temporary relief of pain associated with sore and aching muscles in the low back as well as upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities.

6. Intended Use:

The Everyway EV-820 OTC Pain Relief TENS is intended for temporary relief of pain associated with sore and aching muscles in the low back and/or upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities.

7. Comparison to the 510(k) Cleared Device (Predicate Device):

The following features are completely identical among the predicate device and our devices.

Features	510(K) Cleared Models		New Model
Model	WL-2406	WL-2407	EV-820
510(K) No.	K133723	K091757	Unknown
Prescription or OTC	OTC	OTC	OTC
Indication for use	temporary relief of pain associated with sore and aching muscles in the low back as well as upper and lower extremities (arm and/or leg) due to strain from exercise or normal Household and work activities	temporary relief of pain associated with sore and aching muscles in the upper And lower extremities (arm and/or leg) due to strain from exercise or normal Household and work activities	temporary relief of pain associated with sore and aching muscles in the low back as well as upper and lower extremities (arm and/or leg) due to strain from exercise or normal Household and work activities
FDA product code	NUH	NUH	NUH
Electrode Used	Self Adhesive Electrode (13x7 cm)/K082065	Belt Electrode and/or Self Adhesive Electrode (5 X 5 cm)/K082065	Self Adhesive Electrode (5x5 cm)/K083302

8. Significant output characteristics comparison table:

Comparison feature	510(K) Cleared Model		New Model
	WL-2406(K133723)	WL-2407(K091757)	EV-820
Net charge	0	0	0
Max. phase charge	13 uc	20.8 uc	20.8 uc
Max. current Density	0.04875 mA/cm ²	0.04992 mA/cm ²	0.0998 mA/cm ²
Max. Average current (RMSA)	500 Ω	50 mA	80 mA
	2K Ω	22.5 mA	30 mA
	10K Ω	7.5 mA	10 mA
Max. Power Density	0.001219 Watts/ cm ²	0.00200 Watts/ cm ²	0.00399 Watts/ cm ²
Burst Mode	Yes	Yes	Yes



Comparison of Unit Characteristics & Output Specification

Mode or Program Name	Predicate Device		New Device
	WL-2406	WL-2407	EV-820
510(K) Number	K133723	K091757	Unknown
Waveform (e.g., pulsed monophasic, biphasic)	Biphasic	Biphasic	Biphasic asymmetric
Shape (e.g., rectangular, spike, rectified sinusoidal)	Rectangular	Rectangular	Rectangular
Maximum Output Voltage (volts) (+/- 20 %)	25V @500Ω 45V @2KΩ 75V @10KΩ	40V @500Ω 60V @2KΩ 100V @10KΩ	40V @500Ω 80V @2KΩ 100V @10KΩ
Maximum Output Current (mA) (+/- 20 %)	50mA @500Ω 22.5mA @2KΩ 7.5mA @10KΩ	80mA @500Ω 30mA @2KΩ 10mA @10KΩ	80mA @500Ω 40mA @2KΩ 10mA @10KΩ
Duration of primary phase (usec)	260 max	260 max	260 max
Pulse Duration (usec)	8333 max	16666 max	8333 max
Frequency (Hz) [or Rate (pps)]	120 max	60 max	120 max
For multiphasic waveforms only:	Yes	Yes	Yes
	Not applicable	Not applicable	Not applicable
Power Source(s)	1.5Vx2 (AAA Size)	1.5Vx3 (AAA Size)	9V x 1 (6F22 Size)
- Method of Line current Isolation	Type BF	Type BF	Type BF
- Patient Leakage Current	---	---	---
- Normal condition (uA)	Under 0.1	Under 0.1	Under 0.1
- single Fault condition (uA)	Under 0.5	Under 0.5	Under 0.5
Average DC current through electrodes when device is on but no pulses are being applied (uA)	Not applicable	Not applicable	Not applicable
Number of Output Modes	8	8	8
Number of Output Channels:	Synchronous Output Coil	Synchronous Output Coil	Synchronous Transformer
Regulated Current or Regulated Voltage?	Voltage	Voltage	Current
Software/Firmware/Microprocessor control?	Yes	Yes	Yes
Automatic Overload Trip?	No	No	No
Automatic No-Load Trip?	Yes	Yes	Yes
Automatic Shut Off?	Yes	Yes	Yes
User Override control?	No	No	No
Indicator Display:	Yes	Yes	Yes
	Yes	Yes	Yes
	Yes	Yes	Yes
Timer Range (Minutes)	10-60	5-60	15, 30,60 and Continuous
Compliance with Voluntary Standards?	IEC 60601-2-10	IEC 60601-2-10	IEC 60601-2-10
Compliance with 21 CFR 898?	Yes	Yes	Yes
Weight (g) including battery	51.6	125.5	142
Dimensions (mm.) [W x H x D]	68x60x17.5	90x50.8x12.7	110x62x27
Housing Materials and construction	ABS	ABS	ABS
Pulse per burst	9	9	6
Burst per second	2	2	2
Burst duration	260us	260us	260us
Duty Cycle	Same for each program	Same for each program	Same for each program
Method of achieving zero net charge for net charge/pulse	Biphasic symmetric wave for each pulse	Biphasic symmetric wave for each pulse	Biphasic asymmetric wave for each pulse

9. Discussion of Non-Clinical Tests Verification Activities Performed to Determine the Safety and Performance of the devices are as the followings:

Compliance to applicable voluntary standards includes IEC 60601-2-10, as well as IEC 60601-1, and IEC 60601-1-2 requirement. In addition to the compliance of voluntary standards, the software verification has been carried out according to the FDA software guidance.

Discussion: The Compliance to applicable voluntary standards as above mentioned indicates that the new device in this submission used the same standards as that of predicate device. Therefore; we consider that the compliance of standards included in our submission is adequate for the determination of substantial equivalence.



10. Discussion of Clinical Test Validation Activities Performed to Determine the Effectiveness of Device are as the followings:

No particular Clinical Test was conducted for Everyway EV-820 OTC Pain Relief TENS.

11. Summary for the technology comparison.

Basically the Everyway EV-820 OTC Pain Relief TENS has the similar technological characteristics with the predicate device in the product design, material, energy source type, main program mode and the main output waveform...etc. There exists some difference in the detailed output parameters (mainly in the output intensity and electrode sizes). Through the detailed calculation comparison of stimulation output energy for each operation mode (in particular the output current density and power density), we found the output level in each operation mode for our Everyway EV-820 OTC Pain Relief TENS and predicate device are very close and within the acceptable range. So we believe the difference in detailed output parameters does not affect the determination of substantial equivalence.

12. Conclusions

The Everyway EV-820 OTC Pain Relief TENS has the same intended use and the similar technological characteristics as the cleared devices. Moreover, verification and validation tests contained in this submission demonstrate that the difference in the submitted models could maintain the same safety and effectiveness as that of cleared device.

In the other words, those engineering difference do not: (1) affect the intended use or (2) alter the fundamental scientific technology of the device.