



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

February 12, 2016

Clinical Laserthermia Systems AB  
% Mr. David Makanani  
OMEDtech, LLC  
1725 Signal Ridge Drive, Suite 150  
Edmond, Oklahoma 73013

Re: K151569

Trade/Device Name: Tranberg CLS Laser Fiber  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic  
surgery and in dermatology  
Regulatory Class: Class II  
Product Code: GEX  
Dated: January 4, 2016  
Received: January 7, 2016

Dear Mr. Makanani:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Jennifer R. Stevenson -S

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K151569

Device Name

TRANBERG CLS Laser fiber

Indications for Use (Describe)

The TRANBERG CLS Laser fiber is indicated for use to necrotize or coagulate soft tissue through interstitial irradiation or thermal therapy in medicine and surgery in cardiovascular thoracic surgery (excluding the heart and the vessels in the pericardial sac), dermatology, ear nose-throat surgery, gastroenterology, general surgery, gynecology, head and neck surgery, neurosurgery, plastic surgery, pulmonology, radiology, and urology at a wavelength of 1064nm.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

---

## 510(K) SUMMARY

**Date** June 1, 2015

**SUBMITTER** Lars-Erik Eriksson, CEO  
Clinical Laserthermia Systems, AB  
Scheelevagen 2  
Lund, Sweden 22381

**CONTACT PERSON** Lars-Erik Eriksson, CEO  
Clinical Laserthermia Systems, AB  
Scheelevagen 2  
Lund, Sweden 22381  
Tel: +4646152100  
Email: lee@clinicallaser.se

**DEVICE NAME**

Classification	Class II
Trade Name	TRANBERG <sup>CLS</sup>   Laser fiber
Common Name	TRANBERG <sup>CLS</sup>   Laser fiber
Classification	21 CFR 878.4810
Product Code	GEX - Powered Laser Surgical Instrument
Review Panel	General and Plastic Surgery

**PREDICATE DEVICE:** K053087, Visualase Cooled Laser Application System (VCLAS)

**INTENDED USE:** The TRANBERG<sup>CLS</sup> | Laser fiber is indicated for use to necrotize or coagulate soft tissue through interstitial irradiation or thermal therapy in medicine and surgery in cardiovascular thoracic surgery (excluding the heart and the vessels in the pericardial sac), dermatology, ear-nose-throat surgery, gastroenterology, general surgery, gynecology, head and neck surgery, neurosurgery, plastic surgery, pulmonology, radiology, and urology, at a wavelength of 1064nm.

**DEVICE DESCRIPTION:**  
The TRANBERG<sup>CLS</sup> | Laser fiber is used to transfer laser energy from the laser unit to the location for the treatment.

The laser fiber is an optical fiber with a core of 550 mic and radial diffusor. The length is 3m and it has a standard connector SMA 905 to fit the laser unit. The numerical aperture is at 0.22. The material in contact with human tissue is biocompatible.

The TRANBERG<sup>CLS</sup> Laser fiber is delivered sterile and for single use only.

**TECHNOLOGICAL CHARACTERISTIC AND SUBSTANTIAL EQUIVALENCE:**

The following table provides more detailed information regarding the basis for the determination of substantial equivalence:

Table 5-1  
Equivalence Comparison

Laser fiber		
Product name	Tranberg <sup>CLS</sup> Laser fiber	Visualase Laser fiber LDF
Manufacturer	Clinical LaserThermia Systems CLS, Sweden	BioTex Inc., US
Intended use / Indications for use.	The Tranberg <sup>CLS</sup> Laser fiber is indicated for use to necrotize or coagulate soft tissue through interstitial irradiation or thermal therapy in medicine and surgery in cardiovascular thoracic surgery (excluding the heart and the vessels in the pericardial sac), dermatology, ear-nose-throat surgery, gastroenterology, general surgery, gynecology, head and neck surgery, neurosurgery, plastic surgery, pulmonology, radiology, and urology, at a wavelength of 1064nm.	The LDF is indicated for use to necrotize or coagulate soft tissue through interstitial irradiation or thermal therapy in medicine and surgery in cardiovascular thoracic surgery (excluding the heart and the vessels in the pericardial sac), dermatology, ear-nose-throat surgery, gastroenterology, general surgery, gynecology, head and neck surgery, neurosurgery, plastic surgery, pulmonology, radiology, and urology, for wavelengths 800nm through 1064nm.
Device Regulatory Classification	Accessory to powered surgical laser instrument	Accessory to powered surgical laser instrument

	FDA 878.4810	FDA 878.4810
Product code	GEX	GEX
Device class	Accessory to powered surgical laser instrument Class 2	Accessory to powered surgical laser instrument Class 2
510(k) No		K053087
Fiber core diameter:	550 $\mu\text{m}$	400 $\mu\text{m}$ (200 – 1000)
Numerical aperture:	0.22	0.37
Fiber length:	3 m	3-12m standard
Proximal connector:	SMA 905	SMA 905
Wavelength:	1064 nm	532-1064nm
Laser operation mode:	Continuous Wave	Continuous Wave
Diffusing region length:	1 mm	7,5-30mm
Diffusing tip assembly diameter:	1.55 mm	0.6-1.4mm
Lesion Shape:	Round	Ellipsoidal / Round
Max power:	8 W for 550 $\mu\text{m}$	8 W for 400 $\mu\text{m}$
Lesion volume (refer to Report DV-2015-024)	0,8 $\text{cm}^3$	0,8 $\text{cm}^3$

#### PERFORMANCE TESTING - (NON-CLINICAL) BENCH

The Tranberg <sup>CL<sup>S</sup></sup> Laser fiber has been determined through engineering bench testing to support substantial equivalence with this device and the predicates. This testing showed the Tranberg <sup>CL<sup>S</sup></sup> Laser fiber to meet applicable ISO, IEC and FDA safety and performance standards,

Non-clinical bench performance testing completed:

- Engineering comparative temperature testing
- Biocompatibility

#### PERFORMANCE TESTING – ANIMAL/CLINICAL

There are no animal or clinical data submitted with this Notification.

#### CONCLUSION:

Based on the results of non-clinical testing, the TRANBERG<sup>CL<sup>S</sup></sup> Laser Fiber performs according to specifications, and as intended, and the comparative discussion of intended use, principle of operation, and technological characteristics, has determined that the Tranberg <sup>CL<sup>S</sup></sup> Laser fiber is substantially equivalent to the predicate device.