



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

September 23, 2015

Specialty Water Technologies
Alex Clark
Compliance Manager
1020 Industrial Drive
Orlinda, TN 37141

Re: K151637
Trade/Device Name: UPT Series Medical Reverse Osmosis Water Treatment System
Regulation Number: 21 CFR 876.5665
Regulation Name: Water purification system for hemodialysis
Regulatory Class: II
Product Code: FIP
Dated: August 21, 2015
Received: August 24, 2015

Dear Alex Clark,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,


Benjamin R. Fisher -S

Benjamin R. Fisher, Ph.D.
Director
Division of Reproductive, Gastro-Renal,
and Urological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K151637

Device Name

UPT Series Medical Reverse Osmosis Water Treatment System

Indications for Use (Describe)

UPT Series Medical Reverse Osmosis Water Treatment System is intended to be used to remove organic and inorganic contaminants from a tap water supply to dilute a dialysate concentrate for hemodialysis treatments, as well as for use for dialyzer reprocessing and dialysis equipment disinfecting and rinsing.

This System is strictly for multi patient use, such as for a dialysis clinic or hospital serving multiple patients with a Central Water Purification System and is intended for use only with conventional hemodialysis treatments. It is not intended for use with high-flux, CRRT, or any other form of dialysis.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary

510(k) Number: _____

Manufacturer: Specialty Water Technologies

Address: 1020 Industrial Drive
Orlinda, TN 37141
(615) 654-4441

Official Contact: Alex Clark
Compliance Manager

Trade Name: UPT Series Medical Reverse Osmosis Water Treatment System

Common Name: Water Purification System

Classification Name: Water Purification System for Hemodialysis

Product Code: 78 FIP

Device Class: II

Classification Reg: 876.5665

Specialty Water Technologies has provided the following information to the US Food and Drug Administration to support substantial equivalence of the UPT Series Medical Reverse Osmosis (RO) Water Purification System to other RO water purification systems currently cleared for sale in the U.S.

1. Device Description

The UPT Series Medical Reverse Osmosis Water Treatment System is an accessory device that is intended for use with hemodialysis applications and is intended to remove organic and inorganic substances and microbial contaminants from tap water used to dilute dialysis concentrate to form dialysate. The UPT Series Medical Reverse Osmosis Water Treatment System will produce water as prescribed by the Association for the Advancement of Medical Instrumentation (AAMI) Standard ANSI/AAMI/ISO 13959:2014 Water for Hemodialysis and Related Therapies provided adequate flow rate of the feed (tap) water and compliance with existing drinking water standards and proper pre-treatment.

2. Intended Use

UPT Series Medical Reverse Osmosis Water Treatment System is intended to be used to remove organic and inorganic contaminants from a tap water supply to dilute a dialysate concentrate for hemodialysis treatments, as well as for use for dialyzer reprocessing and dialysis equipment disinfecting and rinsing.

This System is strictly for multi patient use, such as for a dialysis clinic or hospital serving multiple patients with a Central Water Purification System and is intended for use only with conventional hemodialysis treatments. It is not intended for use with high-flux, CRRT, or any other form of dialysis.

3. Substantial Equivalence to Predicate Device

| Specifications | UPT Series Medical Reverse Osmosis Water Treatment System | Predicate Devices | |
|-------------------------|--|--|--|
| | | Matrix Series Digital Reverse Osmosis | Complete Water System Models MD610, 620, 630 and 640 |
| Manufacturer | Specialty Water Technologies | Better Water, Inc | Isopure Corporation |
| 510(k) Number | - | K051620 | K041163 |
| Decision Date | - | 11/18/2005 | 10/7/2004 |
| Classification | 876.5665, Class II | 876.5665, Class II | 876.5665, Class II |
| Product Code | 78 FIP | 78 FIP | 78 FIP |
| Water Contact Materials | FDA-NSF Compliant | FDA-NSF Compliant | FDA-NSF Compliant |
| Pre-Treatment Options | Based upon site specific feed water chemical analysis | Based upon site specific feed water chemical analysis | Based upon site specific feed water chemical analysis |
| | Feed water Pre-Treatment may include: | Feed water Pre-Treatment may include: | Feed water Pre-Treatment may include: |
| | Carbon Filtration | Carbon Filtration | Carbon Filtration |
| | Softeners | Softeners | Softeners |
| | Temperature Control | Temperature Control | Temperature Control |
| | Pressure Compensation | Pressure Compensation | Pressure Compensation |
| | Special Filtration | Special Filtration | Special Filtration |
| | Ultra-Violet (UV) | Ultra-Violet (UV) | Ultra-Violet (UV) |
| Operational Features | Electronic Touch Screen Interface | Electronic Touch Screen Interface | Electronic Touch Screen Interface |
| | All components mounded on an open frame for accessibility | All components mounted on an open frame for accessibility | All components mounted on an open frame for accessibility |
| | <u>Electronic Monitoring of:</u> | <u>Electronic Monitoring of:</u> | <u>Electronic Monitoring of:</u> |
| | Feed Water Pressure / Flow | Feed Water Pressure / Flow | Feed Water Pressure |
| | Product Water Pressure / Flow | Product Water Pressure / Flow | Product Water Pressure |
| | Membrane Feed Pressure | Membrane Feed Pressure | Membrane Feed Pressure |
| | Pump Pressure | Pump Pressure | Pump Pressure |
| | Reject Water Pressure / Flow | Reject Water Pressure / Flow | Reject Water Pressure / Flow |
| | Product Water Conductivity | Product Water Conductivity | Product Water Conductivity |
| | Feed Water Conductivity | Feed Water Conductivity | Feed Water Conductivity |
| Safety Features | 100% Diversion to drain above product conductivity set-point | 100% Diversion to drain above product conductivity set-point | 100% Diversion to drain above product conductivity set-point |
| | Low Supply Pressure Alarm | Low Supply Pressure Alarm | Low Supply Pressure Alarm |
| | High Product Water Conductivity Alarm | High Product Water Conductivity Alarm | High Product Water Conductivity Shut Down |
| | High Supply Water Temperature Shut Down | High Supply Water Temperature Shut Down | High Supply Water Temperature Shut Down |
| | Power Disturbance Alarm | Power Disturbance Alarm | Power Disturbance Alarm |
| | Remote Alarm (Nurses Station) | Remote Alarm (Nurses Station) | Remote Alarm (Nurses Station) |
| Post-Treatment | May include 0.5 micron or better ultra-filtration | May include 0.5 micron or better ultra-filtration | May include 0.5 micron or better ultra-filtration |
| Performance | Delivery of AAMI Standard Water | Delivery of AAMI Standard Water | Delivery of AAMI Standard Water |

4. Safety and Effectiveness

The intended use and technological characteristics of the UPT Series Medical Reverse Osmosis Water Treatment System are similar or equivalent to the predicate devices indicated. Any differences between the device and the predicate devices have no significant influence on the safety or effectiveness of the product.

5. Clinical Performance Data

Not required for determination of substantial equivalence for this type and class of device.

6. Electromagnetic Compatibility & Safety

The water system does not contain any wifi, wireless or other transmitting or receiving devices. Any electrical connections between components are wired to include any remote alarms.

The electrical control components purchased and used in the construction of the components used for this water system comply with:

- EN/IEC 61000-4-2 Electrostatic Discharges
- EN/IEC 61000-4-3 Radiated Fields and Transients / Burst
- EN/IEC 61000-4-4 Electrical Fast Transients / Burst
- EN/IEC 61000-4-8 Magnetic Field
- UL Approved / Recognized for safety

All pumps used are UL recognized (UR) and CSA Approved for safety

7. Conclusion Drawn from Clinical and Nonclinical Test Data

As compared with the predicate devices, K051620 and K041663, the UPT Series Medical Reverse Osmosis Water Treatment System will produce water that will meet the minimum water quality requirements as specified by Association for the Advancement of Medical Instrumentation (AAMI) standard ANSI/AAMI/ISO 13959:2014, Water for Hemodialysis and Related Therapies; and ANSI/AAMI/ISO 26722:2014, Water Treatment Equipment for Hemodialysis Applications and Related Therapies, when used as directed.