



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Nanjing Jusha Display Technology Co., Ltd.  
% Ma Jing  
Certification Engineer  
301 Hnazhongmen Street, 8F Block A, No. 1  
Nanjing International Service Outsourcing  
Mansion, Nanjing Jiangsu 210036  
CHINA

September 10, 2015

Re: K151861  
Trade/Device Name: JUSHA-C23C LCD Monitor  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: II  
Product Code: PGY  
Dated: August 13, 2015  
Received: August 17, 2015

Dear Ma Jing:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Robert Ochs". The signature is written in a cursive style with a prominent "R" and "O".

Robert Ochs, Ph.D.  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K151861

Device Name

JUSHA-C23C LCD Monitor

Indications for Use (Describe)

JUSHA-C23C LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

In accordance with 21 CFR 807.92 the following summary of information is provided:

Date:	May 6, 2014
Submitter:	Nanjing Jusha Display Technology Co., Ltd  Add: 8A, Block 1. Nanjing International Service Outsourcing Mansion, No. 301, Hanzhongmen street, Nanjing, China.
Contact Person:	Ma Jing  Certification Engineer  Nanjing Jusha Display Technology Co., Ltd  Tel: +86-25- 83305050  Fax: +86-25- 58783271
Device Trade Name:	JUSHA-C23C LCD Monitor
Common/Usual Name:	2MP Color LCD Monitor
Classification Name: Product Code:	System, image processing 21CFR 892.2050 PGY
Predicate Device(s):	RADIFORCE RX240; K113844
Device Description:	<p>JUSHA-C23C LCD Monitor is the display system with the high resolution (1600*1200), high luminance 600cd/m<sup>2</sup>), and 1024 simultaneous shades of gray out of a palette of 4096, 8 DICOM look up table and 3 GAMMA look up table inside. JUSHA-C23C has ambient brightness adapting and presence induction system, with these this display can automatic adjustment according to different requirements in order to achieve the best results.</p> <p>The product is consisted of the following components:</p> <ul style="list-style-type: none"><li>- 21.3inches, Color-TFT LCD Panel</li><li>- JUSHA- SMS_19inch motherboard/FR-4/REV:0.1</li><li>- JUSHA-C23C LCD Monitor software</li><li>- Power Adapter</li><li>- Data Cable.</li></ul>

	<p>The LCD Monitor is designed, tested, and will be manufactured in accordance with both mandatory and voluntary standards:</p> <ol style="list-style-type: none"> <li>1. IEC 60601-1 Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance 2005+CORR.1(2006)+CORR.2(2007)</li> <li>2. IEC 60601-1-2 Edition 3:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.</li> </ol>
Intended Use:	<p>JUSHA-C23C LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.</p>
Technology:	<p>JUSHA-C23C LCD Monitor is the display system with the high resolution monitor (2 megapixels) with electronic capabilities for evaluation of high resolution medical images, high luminance (600 cd/m<sup>2</sup>) and 1024 simultaneous shades of gray out of a palette of 4096, 8 DICOM look up table and 3 GAMMA look up table inside. In particular, JUSHA-C23C LCD Monitor contains CGA function, it is specially made by JUSHA, it can automatic identify gray and color signals, then gray area calls DICOM LUT and color area calls GAMMA2.2 LUT. JUSHA-C23C has ambient brightness adapting and presence induction system, with these this display can automatic adjustment according to different requirements in order to achieve the best results.</p>
Determination of Substantial Equivalence:	<p><u>Summary of Non-Clinical Tests:</u></p> <p>The LCD Monitor complies with voluntary standards as following:</p> <ol style="list-style-type: none"> <li>1 IEC 60601-1 Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance 2005+CORR.1(2006)+CORR.2(2007)</li> <li>2 IEC 60601-1-2 Edition 3:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.</li> </ol> <p>JUSHA-C23C is substantially equivalent to RADIFORCE RX240. JUSHA-C23C employs the maximum resolution values same as that of RADIFORCE RX240. Comparison table of the principal characteristics of 2 devices is shown in the Attachment 1.</p>

	<p>Attachment 1</p> <p>The following quality assurance measures were applied to the development of the system:</p> <ul style="list-style-type: none"> <li>• Risk Analysis</li> <li>• Requirements Reviews</li> <li>• Design Reviews</li> <li>• Raw materials verification</li> <li>• Testing on unit level (Module verification)</li> <li>• Integration testing (System verification)</li> <li>• Final acceptance testing (Validation)</li> <li>• Performance testing (Verification)</li> <li>• Safety testing (Verification)</li> </ul> <p><u>Summary of Clinical Tests:</u></p> <p>The subject of this premarket submission, LCD Monitor, did not require clinical studies to support substantial equivalence.</p> <p>The proposed device is Substantially Equivalent (SE) to the predicate device which is US legally market device. Therefore, the subject device is determined as safe and effectiveness.</p>
<p>Conclusion:</p>	<p>Nanjing Jusha Display Technology Co., Ltd Considers the JUSHA-C23C LCD Monitor to be as safe, as effective, and performance is substantially equivalent to the predicate device(s).</p>

## 12.1 Product Comparison

This comparison identifies the similarities and differences of the proposed JUSHA-C23C LCD Monitor device to the legally marketed predicate EIZO RADIFORCE RX240 LCD Monitor device to which substantial equivalency is claimed.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	RADIFORCE RX240	JUSHA-C23C LCD Monitor	
510(k) Number	K113844	/	
Display Performance/Specifications			
Screen technology	21.3" Color TFT LCD Panel	21.3" Color TFT LCD Panel	Same
Viewing angle (H, V)	Horizontal 176 °;Vertical 176 °	Horizontal 176 °;Vertical 176 °	Same
Resolution	1600 x 1200/1200 x 1600	1600 x 1200/1200 x 1600	Same
Display area	432.0 (H) x 324.0(V) mm	432.0 (H) x 324.0(V) mm	Same
Contrast Ratio	1200:1	1400:1	Due to the different panel.
Scanning frequency (H; V)	31~100 kHz;59~61Hz	52~76 kHz;59~61Hz	This difference between the two device is caused by the different no display area defined by different manufacturers , nothing to do with the display function.
Recommended Luminance	400cd/m <sup>2</sup>	400cd/m <sup>2</sup>	Same
Pixel Pitch	0.27x0.27 mm	0.27x0.27 mm	Same.
Backlight	LED	LED	Same.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	RADIFORCE RX240	JUSHA-C23C LCD Monitor	
510(k) Number	K113844	/	
DICOM LUT	10-bit (Display Port) : 1024  8-bit: 256	12-bit:4096	The JUSHA-C23C LCD Monitor uses a grayscale expansion technology to improve image display quality, the image clarity is better than the image displayed on the predicate device.
Luminance calibration	Built in calibration sensor provided	Built in calibration sensor provided	Same
Video Signal Input			
Input signals	DVI standard 1.0,  DisplayPort 1.1a	DVI standard 1.0,  DisplayPort 1.1a	Same
Input terminational	DVI-D x 1,  DisplayPort x 1	DVI-D x 1,  DisplayPort x 1	Same
Display controller	Off the shelf	Off the shelf	Same
Power Related Specification			
Power Requirement	AC 100~240V 50~60Hz	AC 100~240V 50~60Hz	Same



Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	RADIFORCE RX240	JUSHA-C23C LCD Monitor	
510(k) Number	K113844	/	
Power Consumption/Save Mode	52W/less than 1.6W	65W/less than 2.5W	The differences caused by different adapter and components used in the LCD Monitor, the JUSHA-C23C has more power consumption. This only shows the power consumption is different, nothing to do with the display function
Power Management	DVI DMPM DisplayPort 1.1a	DVI DMPM DisplayPort 1.1a	Same
Miscellaneous Features/Specifications			
USB Ports/standard	1 upstream (endpoint), 2 downstream/ Rev. 2.0	1 upstream (endpoint), 2 downstream/ Rev. 2.0	Same
Dimensions w/o stand (W x H x D)	Without stand: 376mmx505mmx98mm  With stand: 376mmx599mmx245.5 mm	Without stand: 382mm x490mm x75mm  With stand: 382mm x533mm x238mm	Different housing design due to the different panel size.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	RADIFORCE RX240	JUSHA-C23C LCD Monitor	
510(k) Number	K113844	/	
Indication for use	<p>RADIFORCE RX240 is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.</p>	<p>JUSHA-C23C LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device does not support the display of mammography images for diagnosis.</p>	Same
Applicable standard	<p>1 IEC 60601-1 Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance 2005+CORR.1(2006)+CORR.2(2007)</p> <p>2 IEC 60601-1-2 Edition 3:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.</p>	<p>1 IEC 60601-1 Medical equipment medical electrical equipment - Part 1: General requirements for basic safety and essential performance 2005+CORR.1(2006)+CORR.2(2007)</p> <p>2 IEC 60601-1-2 Edition 3:2007, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests.</p>	Same

## *PERFORMANCE DATA:*

The following performance data were provided in support of the substantial equivalence determination.

### **Bench testing:**

Bench testing was conducted to demonstrate the JUSHA-C23C meets all performance standards as follows:

- Measurement of the angular dependency of luminance response in horizontal, vertical and diagonal directions
- Measurement of the luminance non-uniformity characteristics of the display screen as specified in TGI18 guideline.
- Measurement of the chromaticity non-uniformity characteristics of the display screen as specified in TG18 guideline.
- Measurement of small-spot contrast ratio.
- Measurement of temporal response
- Performance data on luminance stability

### **Electrical safety and electromagnetic compatibility (EMC)**

Electrical safety and EMC testing were conducted on the JUSHA-C23C. The device complies with the IEC 60601-1 standard for safety and the IEC 60601-1-2 standard for EMC.

### **Animal and clinical study**

The subject of this premarket submission, JUSHA-C23C, does not require animal or clinical studies to support substantial equivalence.

## **CONCLUSIONS**

JUSHA-C23C Medical Display is substantially equivalent to the predicate device with respect to technical characteristics, performance, application and intended use. The non-clinical data support the safety of the device. The device should perform as intended in the specified use conditions. Nanjing Jusha Display Technology Co., Ltd considers the JUSHA-C23C Medical Display does not raise any new issues of safety or effectiveness.