



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

January 8, 2016

TNI Manufacturing, Inc.  
Ms. Marianne Grunwaldt  
Director, Quality Assurance & Regulatory Affairs  
4635 NW 103<sup>rd</sup> Avenue  
Sunrise, Florida 33351

Re: K152876

Trade/Device Name: Long Sheath  
Regulation Number: 21 CFR 870.1250  
Regulation Name: Percutaneous Catheter  
Regulatory Class: Class II  
Product Code: DQY  
Dated: December 3, 2015  
Received: December 7, 2015

Dear Ms. Grunwaldt:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Carlos L. Peña-S

Carlos L. Peña, PhD, MS  
Director  
Division of Neurological  
and Physical Medicine Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K152876

Device Name

Long Sheath

Indications for Use (Describe)

The Long Sheath is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasulature.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## **510(k) SUMMARY**

### **Submitter's Name and Address**

TNI Manufacturing, Inc.  
4635 NW 103<sup>rd</sup> Avenue  
Sunrise, Florida 33351  
Telephone: 1-954-742-5988  
Facsimile: 1-954-742-5989

### **Contact Information**

Marianne Grunwaldt  
Director, Quality Assurance &  
Regulatory Affairs  
InNeuroCo, Inc  
4635 NW 103<sup>rd</sup> Avenue  
Sunrise, Florida 33351  
Telephone: 1-305-495-3883  
Facsimile: 1-954-742-5989  
E-Mail: Marianne@InNeuroCo.com

### **Date Prepared**

September 25, 2015

### **Device Trade or Proprietary Name**

Long Sheath or LS

### **Device Common or Classification Name:**

Catheter, Percutaneous, 21CFR870.1250, Class II

### **Product Code:**

DQY

**Identification of the Legally Marketed Devices to which Equivalence is Being Claimed**

<b>Name of Predicate Device</b>	<b>Name of Manufacturer</b>	<b>510(k) Number</b>
NEURON MAX SYSTEM	Penumbra, Inc. Alameda, California	K111380

**Device Description**

The Long Sheath or LS is a variable stiffness catheter and has a catheter shaft reinforced with a stainless steel braid and has a radiopaque Platinum/Iridium marker band on the distal end. The catheter has a nominal outer diameter of 0.109 inches. It is available with a nominal inner diameter of 0.088 inches. They are available in three working lengths: 70 cm, 80 cm, and 90 cm. The Long Sheath has a PTFE-lined lumen, which is braid-reinforced, flexible, and has a hydrophilic coating. The LS are inserted through a guide catheter or vascular sheath, provide access to the target site and once in place, provide a reinforcing conduit for other intravascular devices. Accessories included with the device are a Tuohy-Borst hemostasis valve with an extension luer and a dilator. The LS is supplied sterile, non-pyrogenic, and intended for single use only.

**Indications for Use**

The Long Sheath is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.

**Comparison to Predicate Device**

	<b>Predicate Device Penumbra NEURON MAX SYSTEM</b>	<b>TNI Manufacturing, Inc. Long Sheath</b>
<b>510(k) Number</b>	K111380	K152876
<b>Classification</b>	21CFR870.1250, Class II	21CFR870.1250, Class II
<b>Product Code</b>	DQY	DQY

<b>Review Panel</b>	Cardiovascular	Cardiovascular
<b>Indication For Use</b>	The Neuron™ MAX System is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.	The Long Sheath is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.
<b>Components Supplied</b>	Sheath, Vessel Dilator, Hemostasis Valve,	Sheath, Vessel Dilator, Hemostasis Valve
<b>Catheter Shaft Material</b>	Nylon, Polyurethane	Polyether Block Amide, (PEBAX)
<b>Inner Liner</b>	PTFE	PTFE
<b>Hub Material</b>	Polycarbonate	Polycarbonate
<b>Strain Relief</b>	Polyolefin	Polyolefin
<b>Catheter Shaft Reinforcement</b>	Stainless Steel Braid	Stainless Steel Braid
<b>Lubricious Coating</b>	Hydrophilic Coating	Hydrophilic Coating
<b>Radiopaque Marker Band</b>	Platinum/ Iridium	Platinum/ Iridium
<b>Packaging</b>	Tyvek Pouch, polyethylene support tube, packaging card, SBS carton	Tyvek/Nylon Pouch, polyethylene support tube, packaging card, SBS carton
<b>Sterilization</b>	Ethylene Oxide	Ethylene Oxide
<b>Pyrogenicity</b>	Nonpyrogenic	Nonpyrogenic
<b>Working Lengths</b>	70, 80, 90, and 100 cm	70, 80, 90 cm
<b>Internal Diameter</b>	0.088 inches	0.088 inches
<b>Outer Diameter</b>	0.108 inches	0.109 inches
<b>Shelf Life</b>	3 years	3 years

<b>Vessel Dilator</b>	Aid in sheath introduction during procedure	Aid in sheath introduction during procedure
<b>Hemostasis Valve</b>	Minimizes blood loss	Minimizes blood loss
<b>Luer Tapered Hub</b>	Yes	Yes

### Summary of Non Clinical Data

Biocompatibility tests conducted with the Long Sheath and accessories were selected in accordance with ISO 10993-1 guidelines (Biological Evaluation of Medical Devices) for limited duration (<24 hours), external communicating devices, contacting circulating blood. Studies were conducted pursuant to 21CFR58, Good Laboratory Practices. Biocompatibility testing found the Long Sheath to be biocompatible and non-pyrogenic.

The conclusions drawn from the physical, mechanical, and performance testing of the subject LS and accessories demonstrates that the product is Substantially Equivalent to the legally marketed predicate device for its labeled indications.

### Long Sheath Performance Testing

Test	Test Method Summary	Results
Biocompatibility	Testing completed per ISO 10993-1	All units tested meet the acceptance criteria. All samples met the acceptance criteria
Radiographic Detectability	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Visual Inspection	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Working Length	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Hub compatibility	Testing completed per ISO 594-1 and ISO 594-2	All units tested must meet the acceptance criteria. All samples met the acceptance criteria

Test	Test Method Summary	Results
Outside diameter	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Leak – Air	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Leak – Liquid	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Particulates	Testing completed per USP 788	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Catheter Burst	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Tensile	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Corrosion	Testing completed per ISO 10555-1	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Packaging – Dye Leak	Testing completed per ASTM F1929-12	All units tested must meet the acceptance criteria. All samples met the acceptance criteria
Packaging – Peel	Testing completed per ASTM F88-09	All units tested must meet the acceptance criteria. All samples met the acceptance criteria