

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center - WO66-G609 Silver Spring, MD 20993-0002

May 20, 2016

Implants Microdent System S.L Mr. Jordi Clapés Donadeu Quality Director C/ Carles Buigas, I-Can Magre Santa Eulalia de Roncana, Barcelona 08187 SPAIN

Re: K153650

Trade/Device Name: Microdent Ektos Implant System Regulation Number: 21 CFR 872.3640 Regulation Name: Endosseous Dental Implant Regulatory Class: II Product Code: DZE, NHA Dated: April 11, 2016 Received: April 18, 2016

Dear Mr. Donadeu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<u>http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm</u> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours, Tina Kiang

for Erin I. Keith, M.S. Director Division of Anesthesiology, General Hospital, Respiratory, Infection Control, and Dental Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

#### DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

#### Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: January 31, 2017 See PRA Statement below.

510(k) Number (if known) K153650

Device Name Microdent Ektos Implant System

#### Indications for Use (Describe)

Microdent Ektos Implant System is indicated for surgical placement in the upper or lower jaw arches, for single-stage or two-stage surgical procedures and cemented, screw retained restorations or overdentures. Microdent Ektos Implant System is intended for immediate placement and function on single tooth and/or multiple tooth applications when good primary stability is achieved, with appropriate occlusal loading, in order to restore chewing function. Small diameter implants are indicated only for replacement of central and lateral incisors in the maxilla and mandible.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPAR	ATE PAGE IF NEEDED.
This section applies only to requirements of	of the Paperwork Reduction Act of 1995.
*DO NOT SEND YOUR COMPLETED FORM TO	THE PRA STAFF EMAIL ADDRESS BELOW.*
The burden time for this collection of information is estir time to review instructions, search existing data sources and review the collection of information. Send comment of this information collection, including suggestions for r	nated to average 79 hours per response, including the s, gather and maintain the data needed and complete ts regarding this burden estimate or any other aspect reducing this burden, to:
Department of Heal Food and Drug Adm Office of Chief Infor Paperwork Reduction PRAStaff@fda.hhs.	th and Human Services ninistration mation Officer on Act (PRA) Staff gov
"An agency may not conduct or sponsor, and a pe information unless it displays a	rson is not required to respond to, a collection of currently valid OMB number."

FORM FDA 3881 (8/14)

Page 1 of 1

PSC Publishing Services (301) 443-5740 EF



# 510(k) SUMMARY (21 CFR 807.92)

This 510(k) summary is submitted in accordance with the requirement of 21 CFR 807.92

The assigned 510(k) number is: K153650

Premarket Notification [510(k)] Summary

#### A. General Information

Submitter's Name:	Implant Microdent Systems
Address:	C/ Carles Buigas, 1 – Can Magre
	Sta Eulalia de Ronçana
	E-08187 Barcelona, Catalonia
	Spain

Telephone:	+34-902-402-420
Fax Number:	+34-94-844-7893
Contact Person:	Jordi Clapes Donadeu
Date Prepared:	May 9th, 2016

B. Device

Trade Name:	Microdent Ektos Implant System
Classification Name:	Endosseous Dental Implants
Product Code:	DZE and NHA (DZE is considered the primary product code)
Class:	П
Regulation Number:	21 CFR 872.3640

# C. Identification of Legally Marketed Predicate Device

Legacy System Dental Implants – K073033 (primary predicate) Microdent Genius Implant System – K141188 (reference predicate) Legacy Abutment System – K060063 (reference predicate)

# D. Description of the Device

Microdent Ektos Implant System is comprised of dental implants and prosthetic components. Microdent Ektos dental implants are internal connection endosseous implants machined in titanium that can be used with deferred load or immediate load techniques. The implant is designed to be inserted at crestal level. Microdent Ektos Implant is recommended to be placed at crestal bone level.



The implant has an internal hexagonal connection.

Microdent Ektos Implants are provided with blasted surface.

The implants are supplied sterile and the abutments are provided non-sterile.

• Implants are provided in the following diameters and length.

Ø platform	Ø core	Length	Reference
	3.70	8	EK3708CP
		10	EK3710CP
2 50		12	EK3712CP
5.50		14	EK3714CP
		16	EK3716CP
		18	EK3718CP
		8	EK4208CP
		10	EK4210CP
2 50	4 20	12	EK4212CP
5.50	4.20	14	EK4214CP
		16	EK4216CP
		18	EK4218CP
	4.80	8	EK4808CP
		10	EK4810CP
4.50		12	EK4812CP
4.50		14	EK4814CP
		16	EK4816CP
		18	EK4818CP
5.70	5.80	8	EK5808CP
		10	EK5810CP
		12	EK5812CP
		14	EK5814CP
		16	EK5816CP



• The system also contains Microdent Ektos abutments made of Ti-6AL 4-V-ELI alloy as follows:

The abutments Conical, Multi-function (supplied with the implant) and Mini Capitel are used for cemented and screw-retained restorations.

EKPCCP3501H	Microdent Ektos hex. Conical with flap abutment	Plat. Ø 3.50 Height 1, 3 and 5 mm
EKPCCP4501H	Microdent Ektos hex. Conical with flap abutment	Plat. Ø 4.50 Height 1, 3 and 5 mm
EKPCCP5701H	Microdent Ektos hex. Conical with flap abutment	Plat. Ø 5.70 Height 1 mm
EKPCCP3501R	Microdent Ektos Circular Conical with flap abutment	Plat. Ø 3.50 Height 1, 3 and 5 mm
EKNPCCP4501R	Microdent Ektos Circular Conical with flap abutment	Plat. Ø 4.50 Height 1, 3 and 5 mm
EKPCCP5701R	Microdent Ektos Circular Conical with flap abutment	Plat. Ø 5.70 Height 1 mm

Conical abutment with flap as indicated in the following table:

Conical abutment without flap as indicated in the following table:

EKPCSP3500H	Microdent Ektos Hex. Conical without flap abutment	Plat. Ø 3.50 Height 0 to 5 mm
EKPCSP4500H	Microdent Ektos Hex. Conical without flap abutment	Plat. Ø 4.50 Height 0 to 5 mm
EKPCSP5700H	Microdent Ektos Hex. Conical without flap abutment	Plat. Ø 5.70 Height 0 and 1 mm
EKPCSP3500R	Microdent Ektos Circular Conical without flap abutment	Plat. Ø 3.50 Height 0 to 5 mm
EKPCSP4500R	Microdent Ektos Circular Conical without flap abutment	Plat. Ø 4.50 Height 0 to 5 mm
EKPCSP5700R	Microdent Ektos Circular Conical without flap abutment	Plat. Ø 5.70 Height 0 and 1 mm



EKPC13500H	Microdent Ektos hex. immediate	Plat. Ø 3.50
	loading conical abutment	Height 0 mm
	Microdent Ektos hex. immediate	Plat. Ø 4.50
	loading conical abutment	Height 0 mm
	Microdent Ektos hex. immediate	Plat. Ø 5.70
	loading conical abutment	Height 0 mm
	Microdent Ektos Circular	Plat Ø 3 50
EKPCI3500R	immediate loading conical	Height 0 mm
	abutment	rieight o min
	Microdent Ektos Circular	Plat Ø 4 50
EKPC4500R	immediate loading conical	Height 0 mm
	abutment	
	Microdent Ektos Circular	Plat Ø 5 70
EKPCI5700R	immediate loading conical	
	abutment	

Immediate loading conical abutments as indicated in the following table:

Mini Capitel abutment as indicated in the following table:

	Microdent Ektos Circular Mini	Plat. Ø 3.50
ERCAPINSSUIR	Capitel Abutment	Height 1 to 4 mm
	Microdent Ektos Circular Mini	Plat. Ø 4.50
ERCAPIN450TR	Capitel Abutment	Height 1 to 4 mm
	Microdent Ektos Circular Mini	Plat. Ø 5.70
ERCAPING/UTR	Capitel Abutment	Height 1 to 2 mm
	Microdent Ektos Mini Capitel	<i>a</i> 4 90
UTSNFC4A	cementable coping	Ø 4.80
	Microdent Ektos Mini Capitel	Ø 5 60
	cementable coping	0 5.00
	Microdent Ektos Mini Capitel	<i>(</i> <b>1</b> 4 90
UTSNCF4A	Protective Cap	Ø 4.00
	Microdent Ektos Mini Capitel	Ø E 60
	Protective Cap	0.60



• The following Cover Screws are provided to protect the inner configuration of the implant and are supplied sterile with the implant.

EKTC35	Microdent Ektos cover screw	Plat. Ø 3.50
EKTC45	Microdent Ektos cover screw	Plat. Ø 4.50
EKTC57	Microdent Ektos cover screw	Plat. Ø 5.70

• Healing abutment to shape the soft tissue during the healing phase as indicated in the following table:

	Microdent Ektos Healing	Ø 4.50
LKF K3502	abutment	Height 2 to 5 mm
	Microdent Ektos Healing	Ø 6.50
EKFK4502	abutment	Height 2 mm
	Microdent Ektos Healing	Ø 5.50
EKPR5/02	abutment	Height 1 to 6 mm
	Microdent Ektos Healing	Ø 3.50
EKPCK3502	abutment	Height 2 to 5 mm
	Microdent Ektos Healing	Ø 4.50
ERFCR4302	abutment	Height 2 to 5 mm
EKPCR5702	Microdent Ektos Healing	Ø 5.70
	abutment	Height 2 mm

- The Retention Screws are used for securing the abutments to the implant.
- Overdenture retention consists of a titanium alloy socket that attached to a threaded post for use with titanium endosseous implants having an internal threaded socket. Both devices have a plastic component that has a shape on one and that mate into the titanium socket, while the other end with metal cap is attached to the denture.

	Microdent Ektos Osscilia	Ø 3.50
	retention abutment	Height 0 to 4 mm
	Microdent Ektos Osscilia	Ø 4.50
ERE0334500	retention abutment	Height 0 to 4 mm
EKEOSS5700	Microdent Ektos Osscilia	Ø 5.70



	retention abutment	Height 0 to 1 mm
CSUTGOSS	Metal cap with soft, middle	Titanium grade 5
	and strong Osscilia retainer.	+POM

All abutments include appropriate features and dimensions to mate with Microdent Ektos implants. No compatibility with other implant systems is claimed.

# E. Intended Use

Microdent Ektos Implant System is indicated for surgical placement in the upper or lower jaw arches, for single-stage or two-stage surgical procedures and cemented, screw retained restorations or overdentures. Microdent Ektos Implant System is intended for immediate placement and function on single tooth and/or multiple tooth applications when good primary stability is achieved, with appropriate occlusal loading, in order to restore chewing function.

Small diameter implants are indicated only for replacement of central and lateral incisors in the maxilla and mandible.

# F. Summary of Testing and Comparison to the Predicate Device

The devices are designed and manufactured in accordance with the following standards:

ISO 5832-2:1999 Implants for surgery - Metallic materials - Part 2: Unalloyed titanium
ISO 5832-3:1996 Implants for surgery -- Metallic materials -- Part 3: Wrought titanium
6-aluminium 4-vanadium alloy

ISO 14971 Second edition 2007-03-01 Medical devices - Application of risk management to medical devices

ISO 14801 Second edition 2007-11-15 Dentistry-Implants-Dynamic fatigue test for endosseous dental implants

ISO 10993-1:2009 Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process

ISO 10993-5:2009 Biological evaluation of medical devices -- Part 5: Tests for In Vitro cytotoxicity

USP 37<sup>th</sup> ed. 2014<85> Bacterial Endotoxins Test

USP 37<sup>th</sup> ed. 2014<151> Pyrogen Test

ISO 14698-1:2003 Cleanrooms and Associated Controlled Environments -Biocontamination Control - Part 1: General Principles and Methods

ISO 14644-1:1999 Cleanrooms and Associated Controlled Environments - Part 1: Classification of Air Cleanliness

ISO 14644-3:2005 Cleanrooms and associated controlled environments - Part 3: Test methods



ISO 11737-1:2006 (R)2011 Sterilization of medical devices - Microbiological methods Part 1: Determination of the population of microorganisms on product, 2ed

ISO 11737-2:2009Sterilization of medical devices -- Microbiological methods --Part 2: Tests of sterility performed in the definition, validation and maintenance of asterilization process

ISO 11607-1:2006/(R)2010 Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

ISO 11607-2:2006/(R)2010 Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

ISO 11137-1:2006/(R) 2010 Sterilization of health care products - Radiation - Part 1: Requirements for development, validation, and routine control of a sterilization process for medical devices

ISO 11137-2:2012 Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose

ISO 11137-3:2006/(R)2010 Sterilization of health care products - Radiation - Part 3: Guidance on dosimetric aspects.

ASTM F1980-07 (Reapproved 2011), Standard Guide For Accelerated Aging Of Sterile Barrier Systems For Medical Devices. (Sterility).

ANSI / AAMI ST79: 2010& A1:2010 & A2:2011 & A3:2012 & A4:2013 Comprehensive guide to steam sterilization and sterility assurance in health care facilities

# Comparison of Technological Characteristics

Characteristic	Subject Device	Own predicate	Predicate Device	SE / Comments
Device	Microdent Ektos	Microdent Genius	Legacy System	Yes (Microdent
	implant	implant ( K141188)	Dental Implants	Ektos and
			(K073033)	Legacy have the
				same hex
				connection.
				Microdent
				Genius has
				connection
				tapered with
				indexation, rest
				of characteristics
				same as
				Microdent
				Ektos)
Intended	Microdent Ektos Implant	Microdent Genius Implant	The intended use of the	Yes (Microdent



use	System is indicated for	System is indicated for	Legacy System Implants with	Ektos specified
	surgical placement in the	surgical placement in the	HA Coating is identical to the	the use of small
	upper or lower jaw arches,	upper or lower jaw arches, for	intended use of the predicate	diameter
	for single-stage or two-stage	single-stage or two-stage	abutments. These implants	implants,
	surgical procedures and	surgical procedures and	are two-piece implants for	Legacy System
	cemented, screw retained	cemented, screw retained	single-stage or two-stage	are coated
	restorations or	restorations or overdentures.	surgical procedures. The	blasted or HA
	overdentures. Microdent	Microdent Genius Implant	Legacy implants are intended	plasma)
	Ektos Implant System is	System is intended for	for use in the mandible and	r
	intended for immediate	immediate placement and	maxilla, in support of single or	
	placement and function on	function on single tooth and/or	multiple-unit cement or screw	
	single tooth and/or multiple	multiple tooth applications	receiving fixed restorations	
	tooth applications when	when good primary stability is	and for retention and support	
	good primary stability is	achieved, with appropriate	of overdentures. The implants	
	achieved, with appropriate	occlusal loading, in order to	are intended for immediate	
	occlusal loading, in order to	restore chewing function.	placement and function for	
	restore chewing function.	Small diameter implants are	support of single tooth and/or	
	Small diameter implants are	indicated only for replacement	multiple tooth restorations,	
	indicated only for	of central and lateral incisors	recognizing bone stability and	
	replacement of central and	in the maxillar and mandible.	appropriate occlusal load	
	lateral incisors in the		requirements	
	maxilla and mandible.			
Indication	Delayed or	same	same	Yes
	inmediate load			
Material	Comercially pure	same	Titanium Alloy	Yes (Different
	titanium (grade		(grade 5 as ISO	degree but with
	4 as ISO 5832-2)		5832-3)	very similar
				properties)
Design;	Threaded, root	same	same	Yes (Legacy
	form implant.			implant has
	<b>1</b>			quadruple-lead
				threads near top)
		-1		
		33	13	
		38	13	
	•		•	
Length (mm)	8, 10, 12, 14, 16	8, 10, 12, 14, 16 and	8, 10, 11.5, 13 and	Yes
	and 18 mm.	18 mm.	16 mm	
Platform	Diameter ranges:	Diameter ranges:	Diameter ranges:	Yes
	3.5 mm, 4.5 mm	3.5 mm, 4.5 mm,	3.0 mm, 3.5 mm.	
	and 5.7 mm	5.0 and 5.5 mm	45  mm and $57$	



D: /	D' (	D' (	D' (	N/
Diameter	Diameter ranges:	Diameter ranges:	Diameter ranges:	Yes
(mm)	3.7 mm, 4.2 mm,	3.5 mm, 4.5 mm,	3.2 mm, 3.7 mm,	
	4.8 mm, 5.8 mm	5.0 and 5,5 mm.	4.2 mm, 4.7 mm,	
	and 7.00 mm.		5,2 mm, 5.7 mm	
			and 7.0 mm.	
Connection	Internal	Internal connection	Internal connection	Yes with Legacy
type	connection with	tapered with	with hex.	(Hex
	hex.	indexation.		connection)
				Different
				connection from
				Genius implant
Surface	Blasting	same	SBM (Soluble Blast	Yes (Blasted
treatment	(roughness 0,82		Media) or	treatment with
	μm		HA-coated. (SBM	similar result)
	peak-to-valley).		roughness 1.5-2.3	
			μm peak-to-valley).	
Starilization	Storilo (Commo			Vac
Stermzation	irradiation)	Same	same	105
De alva ain a		Dealso and in starila	Immon alaona ta	Vag (hut Lagaar)
Packaging	miler sleeve to	vial with cover	sugmend the	includes
	implant/fixture m	viai with cover	implant/fixture me	autonder)
	implant/fixture-in	screw	unt aggembly	extender)
	ount assembly		unt assembly	
	packaged inside		packaged inside	
	all Outer Vial		with a con	
	Bealed with a cap.		Paakaging also	
	rackaging also		rackaging also	
	action of the series of the se		and a surgical	
	multi function		ovtender and	
	abutmont		tomporary coning	
Shalf Life	5 yoors	sama	sama	Vac
Mating	J years	All Miero dent		ICS Voc
Iviating	All Microdent	All Microdent	All Legacy Destorative	ICS
Components	Kestorative	Restorative	Restorative	
	Components	Components	Components	

#### Table 2: General Prosthetic Device Comparison



Characteristic	Subject Device	Own predicate	Predicate Device	SE / Comments
Device Name	Microdent Ektos	Microdent Genius	Legacy Abutment	yes
	implant abutments	implant abutments	System	
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants.).
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating diameters,	
	function on single tooth and/or	placement and function on	lead-in bevels, internal hex	
	multiple tooth applications	single tooth and/or multiple	sizes, and 1-72UNF internal	
	when good primary stability is	tooth applications when	threads, as shown in the	
	achieved, with appropriate	good primary stability is	Zimmer Dental Tapered	
	occlusal loading, in order to	achieved, with appropriate	Screw-Vent Surgical Manual.	
	restore chewing function.	occlusal loading, in order to	Implant Direct LLC will monitor	
	Small diameter implants are	restore chewing function.	the compatible implants for	
	indicated only for replacement	Small diameter implants are	modifications to ensure future	
	of central and lateral incisors	indicated only for	compatibility. In the event of	
	in the maxilla a nd mandible.	replacement of central and	any modification, Implant	
		lateral incisors in the	Direct LLC will either modify	
		maxillar and mandible.	the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Material	Titanium alloy	same	Titanium alloy.	Yes
	(Grade 5).			
Surface	Polished.	same	same	Yes
treatment				
Sterilization	Non-sterile	same	same	Yes
Packaging	Blister	same	same	Yes

Device Name	Microdent Ektos	Microdent Genius	Legacy cover screw	Yes
	cover screw	cover screw		



510K	NA	K141188	K073033	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The intended use of the	Yes (Microdent
	System is indicated for	System is indicated for	Legacy System Implants with	Ektos specified
	surgical placement in the	surgical placement in the	HA Coating is identical to the	the use of small
	upper or lower jaw arches, for	upper or lower jaw arches,	intended use of the predicate	diameter
	single-stage or two-stage	for single-stage or two-stage	abutments. These implants	implants,
	surgical procedures and	surgical procedures and	are two-piece implants for	Legacy System
	cemented, screw retained	cemented, screw retained	single-stage or two-stage	are coated
	restorations or overdentures.	restorations or	surgical procedures. The	blasted or HA
	Microdent Ektos Implant	overdentures. Microdent	Legacy implants are intended	plasma)
	System is intended for	Genius Implant System is	for use in the mandible and	1 /
	immediate placement and	intended for immediate	maxilla, in support of single or	
	function on single tooth and/or	placement and function on	multiple-unit cement or screw	
	multiple tooth applications	single tooth and/or multiple	receiving fixed restorations	
	when good primary stability is	tooth applications when	and for retention and support	
	achieved, with appropriate	good primary stability is	of overdentures. The implants	
	occlusal loading, in order to	achieved, with appropriate	are intended for immediate	
	restore chewing function.	occlusal loading, in order to	placement and function for	
	Small diameter implants are	restore chewing function.	support of single tooth and/or	
	indicated only for replacement	Small diameter implants are	multiple tooth restorations,	
	of central and lateral incisors	indicated only for	recognizing bone stability and	
	in the maxilla a nd mandible.	replacement of central and	appropriate occlusal load	
		lateral incisors in the	requirements.	
		maxillar and mandible.		
Design;	3 diameters and	One diameter	same	Yes (Microdent
	metric thread M1.8			Genius has one
				diameter and
				smaller)
Collar Height	Without heights	same	same	Yes
(mm, min	C			
-max)				
Seating	Diameter 3.5 mm,	Diameter 2.9 mm.	same	Yes (Microdent
Surface (mm)	4.5 mm and 5.7			Genius has one
	mm.			diameter and
				smaller)
Connection	Internal connection		Internal connection	Yes
type	tapered.		tapered.	
Sterilization	Sterile (Gamma	same	same	Yes.
	irradiation)			
	, , , , , , , , , , , , , , , , , , , ,			Packaged
				sterile with the
				implant



Device Name	Microdent Ektos	Microdent Genius	Legacy Healing	yes
	Healing Abutment	Healing Abutment	Collar	
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants).
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating diameters,	
	function on single tooth and/or	placement and function on	lead-in bevels, internal hex	
	multiple tooth applications	single tooth and/or multiple	sizes, and 1-72UNF internal	
	when good primary stability is	tooth applications when	threads, as shown in the	
	achieved, with appropriate	good primary stability is	Zimmer Dental Tapered	
	occlusal loading, in order to	achieved, with appropriate	Screw-Vent Surgical Manual.	
	restore chewing function.	occlusal loading, in order to	Implant Direct LLC will monitor	
	Small diameter implants are	restore chewing function.	the compatible implants for	
	indicated only for replacement	Small diameter implants are	modifications to ensure future	
	of central and lateral incisors	indicated only for	compatibility. In the event of	
	in the maxilla a nd mandible.	replacement of central and	any modification, Implant	
		lateral incisors in the	Direct LLC will either modify	
		maxillar and mandible.	the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Design;	4 diameters and	One diameter	One diameter and	Yes
	metric thread M1.8		metric thread M1.8	
Collar Height	0 mm to 4 mm	1 mm to 6 mm	3 and 5 mm	Yes
(mm, min				
-max)				
Seating	Diameter from 3.5,	Diameter from	Diameter from	Yes
Surface (mm)	4.5, 5.6 and 6.5	4.5, 5 and 5.5 mm.	3.7, 4.7, 5.7 and 6.5	
	mm.		mm.	
Connection	Internal connection	same	same	Yes





type

tapered.

Device Name	Microdent Ektos	Microdent Genius	Legacy Fixation	Yes
	Retention screw	Retention screw	screw	
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants).
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating diameters,	
	function on single tooth and/or	placement and function on	lead-in bevels, internal hex	
	multiple tooth applications	single tooth and/or multiple	sizes, and 1-72UNF internal	
	when good primary stability is	tooth applications when	threads, as shown in the	
	achieved, with appropriate	good primary stability is	Zimmer Dental Tapered	
	occlusal loading, in order to	achieved, with appropriate	Screw-Vent Surgical Manual.	
	restore chewing function.	occlusal loading, in order to	Implant Direct LLC will monitor	
	Small diameter implants are	restore chewing function.	the compatible implants for	
	indicated only for replacement	Small diameter implants are	modifications to ensure future	
	of central and lateral incisors	indicated only for	compatibility. In the event of	
	in the maxilla a nd mandible.	replacement of central and	any modification, Implant	
		lateral incisors in the	Direct LLC will either modify	
		maxillar and mandible.	the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Design;	This screw having	same	same	Yes
	a threaded fuse and			
	a head with a			
	hexagon 1.20 mm			
	flat to flat.			
Seating	This is a single	same	same	Yes
Surface (mm)	Retention screw for			
	all platforms			



Device Name	Microdent Ektos	Microdent Genius	Legacy straight	Yes
	Conical abutment	Conical abutment	abutments	
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants)
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating diameters,	
	function on single tooth and/or	placement and function on	lead-in bevels, internal hex	
	multiple tooth applications	single tooth and/or multiple	sizes, and 1-72UNF internal	
	when good primary stability is	tooth applications when	threads, as shown in the	
	achieved, with appropriate	good primary stability is	Zimmer Dental Tapered	
	occlusal loading, in order to	achieved, with appropriate	Screw-Vent Surgical Manual.	
	restore chewing function.	occlusal loading, in order to	Implant Direct LLC will monitor	
	Small diameter implants are	restore chewing function.	the compatible implants for	
	indicated only for replacement	Small diameter implants are	modifications to ensure future	
	of central and lateral incisors	indicated only for	compatibility. In the event of	
	in the maxilla a nd mandible.	replacement of central and	any modification, Implant	
		lateral incisors in the	Direct LLC will either modify	
		maxillar and mandible.	the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Design	3 diameters and	3 diameters and	4 diameters and hex	Yes, but
	hex connection	same connection	connection	Microdent
	(Microdent Ektos).	(Microdent	(Legacy).	Ektos and
		Genius).		Legacy have
				the same hex
		Abutment with		connection.
	Abutment with flap	flap and without	Abutments	Microdent
	and without flap.	flap.	Contoured and	Genius has
			Full-contoured.	connection
				tapered with



	Microdent Ektos inmediate loading abutment with slots on outer taper surface.	Microdent Genius inmediate loading abutment with slots on outer taper surface.	Design for inmediate loading is called Titanium Temporany Abutments and have not slots on taper surface.	indexation, rest of characteristics same to Microdent Ektos.
	Multifunction abutment.	Not available.	Fixture-Mount.	Yes, but Fixture-Mount has the same function than multifunction abutment using two pieces.
Collar Height (mm, min -max)	1 mm to 5 mm	1 mm to 6 mm	1, 2 and 3 mm	Yes (Legacy has not height 4 and 5 mm, but Genius yes)
Seating Surface (mm)	Diameter from 4.5, 5.6 and 6.2 mm.	Diameter from 4.5, 5 and 5.5 mm.	Diameter from 3.7, 4.7, 5.7 and 6.5 mm.	Yes
Connection type	Internal hex connection (H) or circular (R).	Internal connection tapered with indexation (H) or only tapered (R).	Internal hex connection (Engaging) or Non-Engaging.	Yes Microdent Ektos and Legacy have the same hex connection. Microdent Genius has connection tapered with indexation, rest of characteristics same to Microdent Ektos.



Device Name	Microdent Ektos	Microdent Genius	Legacy	Yes
	Mini Capitel	Mini Capitel	Screw-Receiving	
	abutment	abutment	Overdenture	
			Abutments	
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants.)
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating	
	function on single tooth and/or	placement and function on	diameters,lead-in bevels,	
	multiple tooth applications	single tooth and/or multiple	internal hex sizes, and	
	when good primary stability is	tooth applications when	1-72UNF internal threads, as	
	achieved, with appropriate	good primary stability is	shown in the Zimmer Dental	
	occlusal loading, in order to	achieved, with appropriate	Tapered Screw-Vent Surgical	
	restore chewing function.	occlusal loading, in order to	Manual.	
	Small diameter implants are	restore chewing function.	Implant Direct LLC will monitor	
	indicated only for replacement	Small diameter implants are	the compatible implants for	
	of central and lateral incisors	indicated only for	modifications to ensure future	
	in the maxilla a nd mandible.	replacement of central and	compatibility. In the event of	
		lateral incisors in the	any modification, Implant	
		maxillar and mandible.	Direct LLC will either modify	
			the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Design	One prosthetic	One prosthetic	One prosthetic	Yes
	diameter of 4,8	diameter of 4,8	diameter of 5 mm.	
	mm.	mm.		
Collar Height	1 mm to 4 mm	1 mm to 6 mm	1, 2 and 4 mm	Yes
(mm, min				



-max)				
Seating	Diameter ranges:	Diameter ranges:	Diameter ranges:	Ves
	Diameter Tanges.	Diameter Tanges.	Diameter ranges.	105
Surface (mm)	3.5 mm, 4.5 mm	3.5 mm, 4.5 mm,	3.0, 3.5 mm, 4.5	
	and 5.7 mm.	5.0 and 5,5 mm.	mm and 5.7 mm.	
Connection	Internal connection	Internal	Internal connection	Yes (These
type	tapered (R).	connection	tapered	Ektos
		tapered (R).	(Non-engaging)	abutments are
				used only for
				multiple screw
				prosthesis)
Conince	Commental 1. comine	Comontal 1	T	Ver (Leese
Copings	Cementable coping	Cementable	Iwo cementable	Yes (Legacy
	of Titanium (Grade	coping of	coping of Gold and	has more
	5).	Titanium (Grade	titanium.	components)
		5)		
		,	Plastic coping.	
			Gold&Plastic	
			Coping.	
			Titanium	
			Temporary Coping.	



Device Name	Microdent Ektos	Microdent Genius	Legacy Ball	Yes
	Osscilia retention	Osscilia retention	Abutments	
	abutment	abutment		
510K	NA	K141188	K060063	Yes
Intended use	Microdent Ektos Implant	Microdent Genius Implant	The Legacy Abutment System	Yes (Ektos
	System is indicated for	System is indicated for	is intended for use with dental	Microdent
	surgical placement in the	surgical placement in the	implants in the maxillary	specified the
	upper or lower jaw arches, for	upper or lower jaw arches,	and/or mandibular arches to	use of small
	single-stage or two-stage	for single-stage or two-stage	provide support for crowns,	diameter
	surgical procedures and	surgical procedures and	bridges, or overdentures for	implants.)
	cemented, screw retained	cemented, screw retained	edentulous or partially	
	restorations or overdentures.	restorations or	edentulous patients.	
	Microdent Ektos Implant	overdentures. Microdent	The Legacy Abutment System	
	System is intended for	Genius Implant System is	is compatible with implants	
	immediate placement and	intended for immediate	that have mating	
	function on single tooth and/or	placement and function on	diameters,lead-in bevels,	
	multiple tooth applications	single tooth and/or multiple	internal hex sizes, and	
	when good primary stability is	tooth applications when	1-72UNF internal threads, as	
	achieved, with appropriate	good primary stability is	shown in the Zimmer Dental	
	occlusal loading, in order to	achieved, with appropriate	Tapered Screw-Vent Surgical	
	restore chewing function.	occlusal loading, in order to	Manual.	
	Small diameter implants are	restore chewing function.	Implant Direct LLC will monitor	
	indicated only for replacement	Small diameter implants are	the compatible implants for	
	of central and lateral incisors	indicated only for	modifications to ensure future	
	in the maxilla a nd mandible.	replacement of central and	compatibility. In the event of	
		lateral incisors in the	any modification, Implant	
		maxillar and mandible.	Direct LLC will either modify	
			the Legacy abutment to	
			ensure compatibility, or cease	
			claiming compatibility to the	
			modified Zimmer Dental	
			Screw-Vent implants.	
Design;	Diameter ranges:	One prosthetic	Diameter ranges:	Yes
	3.5 mm, 4.5 mm	diameter of 3.5	3.5 mm, 4.5 mm	
	and 5.7 mm.	mm.	and 5.7 mm.	
	~	~ .	~	
	Compensation 17°.	Compensation	Compensation 28°.	
		1/°.		
Potention:	Three types of	Three times of	Housing and a time	Var
	retainers (with	retainers (with	Trousing and a type	105



	dimensional	dimensional	of retainer.	(Microdent
	variations) varies	variations) varies		Genius and
	the cap retainer.	the cap retainer.		Ektos have 3
			1000-81 1000-82	retainers with different
				friction and
			Ĺ	Legacy only
				one.)
Collar Height	0 mm to 4 mm	1 mm to 6 mm	2 and 4 mm.	Yes
(mm, min				
-max)				
Seating	Diameters 3.5 mm,	Diameter 3.5 mm.	Diameters 3.5 mm,	Yes
Surface (mm)	4.5 mm and 5.7		4.5 mm and 5.7	
	mm.		mm.	
Connection	Internal connection	Internal	Internal connection	Yes
type	tapered.	connection	tapered.	
		tapered.		

Microdent Ektos Implant System is substantially equivalent intended use as the identified predicates. Microdent Ektos Implant System is similar in fundamental scientific technology to the predicate devices in that they all have been designed, manufactured and tested in compliance with FDA'S Class II special controls guidance document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments.

Microdent Ektos Implant System is substantially equivalent in materials, indications and intended use, packaging, labeling and performance to the predicate devices currently marketed in the U.S.

The only differences the subject device and the predicate are slight differences in design and dimensions with Legacy and the connection type with Genius.

# Non-Clinical performance tests

The proposed devices have been subject to bench testing to determine fulfillment of design and performance requirements. Bench testing followed the recommendations provided in FDA Guidance Document – Class II Special Controls Guidance Document for Endosseous Dental Implants and Endosseous Dental Implant Abutments.

Specific performance testing on the EKTOS implant system included:

- Static and dynamic fatigue testing in accordance with ISO 14801;
- Sterilization validation has been performed following standards ISO 11137-1 and ISO 11137-2.



Test data from the predicate device Microdent GENIUS implant system was leveraged for the following performance tests based on the use of identical implant and packaging materials as well as identical manufacturing, packaging and cleaning processes for the Microdent EKTOS implant system:

- Biocompatibility evaluation in accordance with ISO 10993-1 and specifically cytotoxicity testing in accordance with ISO 10993-5;
- Package integrity testing in accordance with ASTM F-1980-07;
- Pyrogen testing in accordance with USP <151>
- Bacterial endotoxin LAL test USP <85>
- ESEM analyses to determine adequate surface finish and roughness
- EDS analyses to determine adequate surface cleaning
- TOC analyses to determine adequate surface cleaning

All specific non-clinical testing demonstrated that the Microdent EKTOS implant system fulfills the related requirements. In addition, all leveraged non-clinical testing demonstrated that materials and processes identical to those used for the Microdent EKTOS implant system fulfilled the related requirements.

#### G. Clinical Testing

No clinical testing was performed. Non-clinical testing was used to support the determination of substantial equivalence.

#### H. Conclusion of Substantial Equivalence

Based on the similarities observed and results of non-clinical testing performed, we conclude that the proposed devices are substantially equivalent to the predicate devices.