



May 26, 2016

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

Visco Vision Inc.
% Ms. Jennifer Ting, Manager
Jens Medical Consulting Ltd.
6F. No. 39, Ln 224, Jixian Rd.
Luzhou Dist. 247 New Taipei City
Taiwan R.O.C

Re: K160344

Trade/Device Name: Si-Hy (olifilcon B) Spherical Silicone Hydrogel Soft Contact Lens;
Si-Hy (olifilcon B) Multifocal Silicone Hydrogel Soft Contact
Lenses; Si-Hy (olifilcon B) Toric Silicone Hydrogel Soft Contact
Lens

Regulation Number: 21 CFR 886.5925

Regulation Name: Soft (Hydrophilic) Contact Lens

Regulatory Class: Class II

Product Code: LPL, MVN

Dated: April 20, 2016

Received: April 22, 2016

Dear Ms. Ting:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing

(21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"

(21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Kesia Alexander

for Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic and Ear,
Nose and Throat Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K160344

Device Name
Si-Hy (olifilcon B) Silicone Hydrogel Soft (Hydrophilic) Contact Lens

Indications for Use (Describe)

The Si-Hy Spherical Silicone Hydrogel Soft Contact Lenses are indicated as daily wear single use soft contact lens for the correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who exhibit refractive astigmatism of 1.00D or less where the astigmatism does not interfere with visual acuity.

The Si-Hy Toric Silicone Hydrogel Soft Contact Lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) with astigmatism in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters and astigmatic corrections are from -0.75 to -2.25 diopters.

The Si-Hy Multifocal Silicone Hydrogel Soft Contact lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) and emmetropia with presbyopia in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters with add powers from +0.75 to +2.75 diopters. The lenses may be worn by persons who exhibit astigmatism of 1.00 diopters or less where the astigmatism does not interfere with visual acuity.

Eye care practitioners may prescribe the lens for daily wear (disposable use) single use. The lenses are to be discarded upon removal. Therefore, no cleaning or disinfecting is required.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

005_510(k) Summary

Preparation Date: May 6, 2016

5.1 Establishment Information:

Name Visco Vision Inc.
Address No. 1, Xingye St., Guishan Dist., Taoyuan City 33341 Taiwan
Contact Ted TT Huang
Phone No. 886-3-3490202
Fax No. 886-3-3596868

5.2 Contact Person:

Company Jens Medical Consulting Ltd.
Name Jennifer TING
Phone No 886-2-82823192
Fax No 886-2-82867686
e-mail: jen.medical@msa.hinet.net

5.3 Device Identification:

Proprietary Name Si-Hy (olifilcon B) Silicone Hydrogel Soft (Hydrophilic) Contact Lens
Common Name Soft (hydrophilic) Contact Lenses
Classification Name Lenses, Soft Contact, Daily Wear,
(21 CFR 886.5925, Product Code LPL)
Lenses, Soft Contact, Daily Wear (Disposable),
(21 CFR 886.5925, Product Code MVN)
Classification II

5.4 Legally Marketed Equivalent Device:

Indication

Vistakon (narafilecon A), K073485, Product Code: MVN
Johnson & Johnson Vision Care Inc.

Material

Visco (olifilcon A) soft contact lens, K141348, Product Code: LPL/MVN
Visco Vision Inc.

5.5 Device Description

The Si-Hy (olofilcon B) Silicone Hydrogel Soft Contact Lens is made of silicone hydrogel material with UV blocker available as spherical lens, toric lens and multifocal lens. The composition of the lens is 53% olifilcon B and 47% water. A light blue color tinted with “reactive Blue19” is for handling visibility purpose. A benzotriazole UV absorbing monomer is used to block UV radiation. The transmittance characteristics are about 2.2% (<5%) in the UVB range of 280-315nm and 8.7% (<50%) in the UVA range of 316-380nm. It is supplied in a sterile state packaged in a buffered saline solution.

5.6 Indication for Use:

The Si-Hy **Spherical** Silicone Hydrogel Soft Contact Lenses are indicated as daily wear single use soft contact lens for the correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who exhibit refractive astigmatism of 1.00D or less where the astigmatism does not interfere with visual acuity.

The Si-Hy **Toric** Silicone Hydrogel Soft Contact Lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) with astigmatism in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters and astigmatic corrections are from -0.75 to -2.25 diopters.

The Si-Hy **Multifocal** Silicone Hydrogel Soft Contact lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) and emmetropia with presbyopia in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters with add powers from +0.75 to +2.75 diopters. The lenses may be worn by persons who exhibit astigmatism of 1.00 diopters or less where the astigmatism does not interfere with visual acuity.

Eye care practitioners may prescribe the lens for daily wear (disposable use) single use. The lenses are to be discarded upon removal. Therefore, no cleaning or disinfecting is required.

5.7 Technological characteristic

Si-Hy SPHERICAL Soft Contact Lenses characteristics:

- Diameter Range : 13 to 15 mm
- Base Curve : 8.0 to 9.2 mm
- Center Thickness : 0.08 mm for -3.00D (varies with power)
- Power : +20.00 to -20.00 D

Si-Hy TORIC Soft Contact Lenses characteristics:

- Diameter Range : 13 to 15 mm
- Base Curve : 8.0 to 9.2 mm
- Center Thickness : 0.08 mm for -3.00D (varies with power)
- Power : +20.00 to -20.00 D
- Cylinder: -0.75D ~ -2.25D
- Axis: 10° to 180° (in 10° increments)

Si-Hy MULTIFOCAL Soft Contact Lenses characteristics:

- Diameter Range : 13 to 15 mm
- Base Curve : 8.0 to 9.2 mm
- Center Thickness : 0.08 mm for -3.00D (varies with power)
- Power : +20.00 to -20.00 D
- Additional Powers: +0.75D ~ +1.25D (LOW)
+1.50D ~ +2.00D (MID)
+2.25D ~ +2.75D (HIGH)

5.8 Comparison table:

The characteristic comparison to predicate device is summarized in the following table.

Indication for Use			
Lens Design	Si-Hy (olifilcon B) soft contact lens	VISTAKON (Narafilcon A) Contact lens (K073485)	Visco (olifilcon A) soft contact lens (K141348)
Spherical	The Si-Hy Spherical Silicone Hydrogel Soft Contact Lenses are indicated as daily wear single use soft contact lens for the correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who exhibit refractive astigmatism of 1.00D or less where the astigmatism does not interfere with visual acuity.	The VISTAKON (narafilcon A) contact Lens is indicated for daily wear single use only for the optical correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who may have 1.00D or less of astigmatism.	The VISCO Soft (Hydrophilic) Contact Lenses is indicated for daily wear for the correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who exhibit refractive astigmatism of 2.00D or less where the astigmatism does not interfere with visual acuity.
Toric	The Si-Hy Toric Silicone Hydrogel Soft Contact Lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) with astigmatism in aphakic and	The VISTAKON (narafilcon A) contact Lens is indicated for daily wear single use only for the optical correction of visual acuity in phakic or aphakic	NA

	non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters and astigmatic corrections are from -0.75 to -2.25 diopters.	persons with non-diseased eyes that are hyperopic or myopic and may have 10.00D or less of astigmatism.	
Multifocal	The Si-Hy Multifocal Silicone Hydrogel Soft Contact lenses are indicated as daily wear for the correction of ametropia (myopia and hyperopia) and emmetropia with presbyopia in aphakic and non-aphakic persons with non-diseased eyes and whose powers are from -20.00 to +20.00 diopters with add powers from +0.75 to +2.75 diopters. The lenses may be worn by persons who exhibit astigmatism of 1.00 diopters or less where the astigmatism does not interfere with visual acuity.	The VISTAKON (narafilecon A) contact Lens is indicated for daily wear single use only for the optical correction of distance and near vision in presbyopic, phakic or aphakic persons with non-diseased eyes who may have 0.75D of astigmatism or less.	NA
Multifocal Toric	NA	The VISTAKON (narafilecon A) contact Lens is indicated for daily wear single use only for the optical correction of distance and near vision in presbyopic phakic or aphakic persons with none-diseased eyes who may have 10.00D of astigmatism or less.	NA
Replacement	Eye care practitioners may prescribe the lens for daily wear (disposable use) single use. The lenses are to be discarded upon removal. Therefore, no cleaning or disinfecting is required. The Si-Hy (olifilcon B) Silicone Hydrogel Soft (Hydrophilic) Contact Lens helps protect against transmission of harmful UV radiation to the cornea and into the eye.	Eye Care Professionals should prescribe the lenses for daily wear single use only (see "Wearing Schedule"). The lenses are to be discarded upon removal. Therefore, no cleaning or disinfection is required.	Eye care practitioners may prescribe the lens for either single-use disposable wear, or for frequent/planned replacement wear, with cleaning, disinfection, and scheduled replacement. When prescribed for frequent replacement, the lens may be disinfected using a chemical disinfection system only. The VISCO Soft (Hydrophilic) Contact Lenses help protect against transmission of harmful UV radiation to the cornea and into the eye.

Similarities and differences			
Item	Device	Predicate (K141348)	Predicate (K073485)
Product Name	Si-Hy (olifilcon B) Silicone Hydrogel soft contact lens	Visco (olifilcon A) soft contact lens	VISTAKON (Narafilcon A) Contact lens
Manufacturer	VISCO VISION Inc.	VISCO VISION Inc.	Johnson & Johnson Vision Care Inc.
Intended Use	Myopia, Hyperopia, astigmatism, Presbyopia	Myopia, Hyperopia	Myopia, Hyperopia, astigmatism, Presbyopia
Lens Design	Spherical, toric, or multifocal	Spherical	spheric, aspheric, toric or multifocal
Replacement Schedule	Daily Disposable (Single use)	Monthly	Daily Disposable (Single use)
Chemical composition	Silicone Hydrogel	The same	The same
Classification	Group 1 (Nonionic, Low water)	The same	The same
Water Content	47 % (<50%)	47 % (<50%)	46 % (<50%),
Oxygen Permeability (DK, 35°C)	120 (Fatt method)	150 (Fatt method)	100 (Fatt method)
Base Curve Range (mm)	8.0~9.2	8.0~9.2	7.80 ~ 10.00
Diameter (mm)	13.0~15.0	13.0~15.0	12.0 ~15.0
Center Thickness	Varies with design and power (0.08 mm at -3.00D)	Varies with design and power (0.08 mm at -3.00D)	Varies with design and power (0.085 mm at -3.00D)
Powers	-20.00D to +20.00D	-20.00D to +20.00D	The same
Refractive Index	1.410	1.410	1.410
Light Transmittance	94%	94%	> 85%
Blue handling tint	Reactive Blue19	Reactive Blue19	Reactive Blue Dye #4
Method of Manufacture	Molded	Molded	Molded

5.9 Nonclinical Tests Performed

5.9.1 Physiochemical studies were conducted according to ISO 18369 First edition 2006-08-15, Ophthalmic optics - Contact lenses (Ophthalmic). The physical,

optical and chemical properties of the lens are within established specifications for the lenses.

5.9.2 Toxicology studies report shows that the lenses are non-toxic and biocompatibility result is acceptable in ocular environment.

5.10 Clinical Studies

A three-month clinical study was conducted to demonstrate the safety and effectiveness of the Si-Hy (olofilcon B) silicone hydrogel soft contact lens by comparison with Vistakon (narafilecon A) Contact Lens when worn on a daily wear basis. At least 50 evaluable subjects were participated. Parameters measured include visual acuities, adverse reactions, symptom, problem and complaints, slit lamp findings, as well as lens wearing time. It was found that the mean VA for both lenses was similar at each visit. The visual acuity could be corrected to 0.1 (log MAR) or better. No significant slit lamp findings (i.e. Grade 3 or Grade 4) were noted. No significant differences in slit-lamp findings were found between groups.

The clinical study provided data to establish substantial equivalence with the predicate, control lens, Vistakon (narafilecon A) soft contact lens, *JOHNSON & JOHNSON VISION CARE INC.* (K073485), with respect to the safety and effectiveness.

5.11 Conclusion

Comparison to the predicate device for chemical composition, physical and optical properties, it shows that “Si-Hy silicone hydrogel soft contact lens” is as safe, as effective and performs as well as the predicate device.