



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

June 14, 2017

IMMUNALYSIS CORPORATION
MICHELLE BODIEN, ASSOCIATE DIRECTOR, REGULATORY AFFAIRS
829 TOWNE CENTER DRIVE
POMONA, CA 91767

Re: k161216

Trade/Device Name: Immunalysis SEFRIA Fentanyl Urine Enzyme Immunoassay and
Immunalysis Fentanyl Urine Calibrators

Regulation Number: 21 CFR 862.3650

Regulation Name: Opiate test system

Regulatory Class: II

Product Code: DJG, DLJ

Dated: May 30, 2017

Received: June 1, 2017

Dear Ms. Bodien:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,


Kellie B. Kelm -S

for Courtney H. Lias, Ph.D.
Director
Division of Chemistry and Toxicology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

k161216

Device Name

Immunalysis SEFRIA Fentanyl Urine Enzyme Immunoassay,
Immunalysis Fentanyl Urine Calibrators

Indications for Use (Describe)

The Immunalysis SEFRIA Fentanyl Urine Enzyme Immunoassay is an enzyme immunoassay with a cutoff of 1.0 ng/mL. The assay is intended for use in laboratories for the qualitative analysis of Fentanyl in human urine with automated clinical chemistry analyzers. This assay is calibrated against Fentanyl. This in vitro diagnostic device is for prescription use only.

The Immunalysis SEFRIA Fentanyl Urine Enzyme Immunoassay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography / Mass Spectrometry (GC-MS) or Liquid Chromatography / Mass Spectrometry (LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Immunalysis Fentanyl Urine Calibrators:

The Immunalysis Fentanyl Urine Calibrators are used as calibrators in the Immunalysis SEFRIA Fentanyl Urine Enzyme Immunoassay for the qualitative determination of Fentanyl in urine on automated clinical chemistry analyzers.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

A. General Information

1. Applicant Name: Immunalysis Corporation
829 Towne Center Drive
Pomona, CA 91767
2. Company Contact: Michelle Bodien
Associate Director, Regulatory Affairs
Phone: (858) 805-2283
Email: michelle.bodien@alere.com
3. Date prepared: June 13, 2017

B. Device Identification

1. Trade Name: Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay
Immunalysis Fentanyl Urine Calibrators
2. Common Name: Fentanyl Urine Enzyme Immunoassay
Fentanyl Urine Calibrators

C. Regulatory Information

1. Device Classification: II
2. Regulation Number: 21 CFR 862.3650 Enzyme Immunoassay, Opiates
21 CFR 862.3200 Calibrators, Drug Specific
3. Panel: Toxicology (91)
4. Product Code: DJG
DLJ
5. Predicate Device: Emit® II Plus Buprenorphine Assay
6. Predicate Company: Siemens Healthcare Diagnostics Inc.
7. Predicate K Number: K150606

D. Device Description

The assay consists of antibody/ substrate reagent and enzyme conjugate reagent. The antibody/ substrate reagent includes EA protein and rabbit antibodies to Fentanyl in PIPES buffer with sodium azide as a preservative. The enzyme conjugate reagent includes ED peptide labeled with Fentanyl and CPRG substrate in malic acid buffer with sodium azide as a preservative.

Calibrators and controls are included as part of the test system and provided separately. The Fentanyl calibrators consist of a Level 1 calibrator at 1 ng/mL, a Level 2 calibrator at 2 ng/mL, and a Level 3 calibrator at 4 ng/mL. The control set contains a LOW control at 0.5 ng/mL and a HIGH control at 1.5 ng/mL.

Automated clinical chemistry analyzers capable of maintaining a constant temperature, pipetting samples and reagents, mixing reagents, timing the reaction accurately and measuring enzymatic rates at 570nm can be used to perform the assay.

The SEFRIA™ technology is based on artificial fragments of the *E. coli* enzyme β -galactosidase. A mutant enzyme, termed Enzyme Acceptor (EA), is created by deletion of 28 amino acids in the amino-terminal region of the sequence. EA is inactive, but can combine with peptides, termed Enzyme Donors (ED's), containing the deleted sequence, to form active β -galactosidase. This process is termed complementation, and the active enzyme formed as a result can be measured by hydrolysis of a chromogenic substrate such as chlorophenolred β -D-galactopyranoside (CPRG). The ED peptides can be modified by attachment of a derivative of fentanyl, which does not interfere with the formation of active β -galactosidase. However antibodies to fentanyl bind to the ED-fentanyl conjugate, and block complementation. The assay is based on the competition of fentanyl in a urine sample with the ED-fentanyl conjugate for the fixed amount of antibody binding sites. In the absence of the free drug in the sample, the antibody binds the ED-fentanyl conjugate, resulting in inhibition of enzyme formation. As the fentanyl concentration in the sample increases, ED-fentanyl becomes available for complementation, creating a dose response relationship between fentanyl concentration in the urine and enzyme formation. The β -galactosidase activity is determined spectrophotometrically at 570 nm by the conversion of CPRG (orange) to chlorophenol red (red) and galactose.

E. Intended Use

1. The Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay is an *in vitro* diagnostic test for the qualitative analysis of Fentanyl in human urine with automated clinical chemistry analyzers. This test is an enzyme immunoassay with a cutoff of 1.0 ng/mL. The assay is intended for use in laboratories and for prescription use only. This assay is calibrated against Fentanyl.

The Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/ Mass Spectrometry (GC-MS) or Liquid Chromatography / Mass Spectrometry (LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

2. Immunalysis Fentanyl Urine Calibrators

The Immunalysis Fentanyl Urine Calibrators are used as calibrators in the Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay for the qualitative determination of Fentanyl in urine on automated clinical chemistry analyzers.

F. Comparison With Predicate

Attribute	Predicate Device K150606 Emit® II Plus Buprenorphine Assay	Candidate Device Immunalysis Fentanyl Urine EIA
Similarities		
Test System	Homogeneous enzyme immunoassay	Same
Sample Matrix	Urine	Same
User Environment	For use in laboratories	Same
Mass Spectrometry Confirmation	Required for preliminary positive analytical results	Same
Storage	2 – 8° C	Same
Materials	Antibody/substrate reagents and enzyme labeled conjugate	Same
Calibrator Form	Liquid	Same
Differences		
Intended Use	For the qualitative and semi-quantitative determination of the presence of Buprenorphine in human urine at a cutoff of 5 ng/mL	For the qualitative determination of the presence of Fentanyl in human urine at a cutoff of 1 ng/mL
Detection	Absorbance change measured spectrophotometrically at 340 nm	Absorbance change measured spectrophotometrically at 570 nm
Measured Analytes	Buprenorphine	Fentanyl
Cutoff Levels	5 ng/mL of Buprenorphine	1 ng/mL of Fentanyl
Antibody	Mouse monoclonal antibody to buprenorphine	Rabbit antibodies to Fentanyl
Reagents Form	R1 and R2: Liquid – Ready to use	EA and ED: Liquid – Ready to use
Calibrator Levels	One negative and four levels (0, 2.5, 5, 15 and 25 ng/mL)	One negative and three levels (0, 1, 2 and 4 ng/mL)

G. Performance Characteristics:

The following laboratory performance studies were performed to determine substantial equivalence of the Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay to the predicate. All studies utilized the Beckman Coulter AU 400e instrument.

1. Precision/ Cutoff Characterization/ Reproducibility – Study was performed for 10 days, 2 runs per day in replicates of 4 on drug free urine (N=80) spiked with fentanyl to concentrations of $\pm 25\%$, $\pm 50\%$, $\pm 75\%$, and $\pm 100\%$ of the cutoff (1.0 ng/mL). The spiked concentrations were confirmed by mass spectrometry (MS). The study verified that the cutoff serves as a boundary between a negative and positive interpretation of a qualitative result.

The following is a summary table precision results for the 1.0 ng/mL cutoff test data results.

Table 1 - Precision Results			
Concentration (ng/mL)	% of cutoff	# of determinations	Total Result
0	-100%	80	80 Negative
0.25	-75%	80	80 Negative
0.5	-50%	80	80 Negative
0.75	-25%	80	80 Negative
1.0	Cutoff	80	32 Negative/48 Positive
1.25	+25%	80	80 Positive
1.5	+50%	80	80 Positive
1.75	+75%	80	80 Positive
2.0	+100%	80	80 Positive

2. Specificity and Cross-Reactivity – Structurally similar compounds were spiked into drug free urine at levels that will yield a result that is equivalent to the cutoff. The study verified assay performance relative to the ability of the device to exclusively determine certain drugs.

The following is a summary table of results:

Table 3 - Structurally Related Compounds			
Compound	Concentration Tested (ng/mL)	Result	Cross-Reactivity (%)
Fentanyl	1	Positive	100.0000
Butyryl Fentanyl	0.8	Positive	125.0000
Acetyl Fentanyl	1	Positive	100.0000
Despropionyl Fentanyl	40	Positive	2.5000
Sufentanil	175	Negative	<0.5714
Haloperidol	1,250	Positive	0.0800
Pipamperone	1,500	Positive	0.0667
Risperidone	2,500	Positive	0.0400
Norfentanyl	20,000	Positive	0.0050
Trazodone	10,000	Positive	0.0100
Labetalol	15,000	Positive	0.0067
Clomipramine	45,000	Positive	0.0022
Benzylpiperazine	50,000	Positive	0.0020
Fluoxetine	60,000	Positive	0.0017
Fenfluramine	60,000	Positive	0.0017
Methamphetamine	70,000	Positive	0.0014
Amitriptyline	75,000	Positive	0.0013
Chlorpromazine	75,000	Positive	0.0013
PCP	100,000	Positive	0.0010
Pentazocine	75,000	Positive	0.0013
6-Acetyl Codeine	100,000	Negative	<0.0010
6-Acetyl Morphine	100,000	Negative	<0.0010
Buprenorphine	100,000	Negative	<0.0010
Bupropion	100,000	Negative	<0.0010
Codeine	100,000	Negative	<0.0010
Cyclobenzaprine	100,000	Positive	0.0010
Desipramine	100,000	Positive	0.0010
Diacetyl Morphine	100,000	Negative	<0.0010
Dihydrocodeine	100,000	Negative	<0.0010
Diphenhydramine	100,000	Positive	0.0010
Doxepin	100,000	Positive	0.0010
EDDP	100,000	Negative	<0.0010
EMDP	100,000	Negative	<0.0010
Hydrocodone	100,000	Negative	<0.0010
Hydromorphone	100,000	Negative	<0.0010
Imipramine	100,000	Positive	0.0010
Levorphanol	100,000	Negative	<0.0010
Meta-chlorphenyl piperazine	100,000	Positive	0.0010
Meperidine	100,000	Negative	<0.0010

Table 3 - Structurally Related Compounds			
Compound	Concentration Tested (ng/mL)	Result	Cross-Reactivity (%)
Methadone	100,000	Negative	<0.0010
Morphine	100,000	Negative	<0.0010
Morphine-3-glucuronide	100,000	Negative	<0.0010
Morphine-6-glucuronide	100,000	Negative	<0.0010
Nalorphine	100,000	Negative	<0.0010
Naloxone	100,000	Negative	<0.0010
Naltrexone	100,000	Negative	<0.0010
Norcodeine	100,000	Negative	<0.0010
Nordiazepam	100,000	Negative	<0.0010
Normorphine	100,000	Negative	<0.0010
Nortriptyline	100,000	Positive	0.0010
Oxycodone	100,000	Negative	<0.0010
Oxymorphone	100,000	Negative	<0.0010
Propoxyphene	100,000	Negative	<0.0010
Protryptiline	100,000	Negative	<0.0010
Tramadol	100,000	Negative	<0.0010
Trimethoprim	100,000	Negative	<0.0010
Trimipramine	100,000	Negative	<0.0010
Venlafaxine	100,000	Negative	<0.0010

3. Interference – Structurally unrelated compounds were evaluated by spiking the potential interferent into drug free urine containing fentanyl at $\pm 50\%$ of the cutoff. All potential interferents analyzed verified that assay performance is unaffected by externally ingested compounds. The results of this study are indicated in Table 5 below.

Table 5 - Structurally Unrelated Compounds			
Compound	Concentration Tested (ng/mL)	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
11-hydroxy-delta-9-THC	75,000	Negative	Positive
11-nor-9 carboxy THC	100,000	Negative	Positive
1S, 2R(+)-Ephedrine	100,000	Negative	Positive
7-Aminoclonazepam	100,000	Negative	Positive
7-Aminoflunitrazepam	100,000	Negative	Positive
7-Aminonitrazepam	100,000	Negative	Positive
Acetaminophen	500,000	Negative	Positive
Amobarbital	100,000	Negative	Positive
Barbital	100,000	Negative	Positive
Benzoylcegonine	100,000	Negative	Positive
Bromazepam	100,000	Negative	Positive
Butabarbital	100,000	Negative	Positive
Caffeine	100,000	Negative	Positive
Cannabidiol	100,000	Negative	Positive
Cannabinol	75,000	Negative	Positive
Carbamazepine	100,000	Negative	Positive

Table 5 - Structurally Unrelated Compounds			
Compound	Concentration Tested (ng/mL)	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
Carisoprodol	100,000	Negative	Positive
Chlordiazepoxide	100,000	Negative	Positive
cis-Tramadol	100,000	Negative	Positive
Clobazam	100,000	Negative	Positive
Clonazepam	100,000	Negative	Positive
Cotinine	100,000	Negative	Positive
Delta-9-THC	100,000	Negative	Positive
Demoxepam	100,000	Negative	Positive
Ecgonine	100,000	Negative	Positive
Ecgonine methyl ester	100,000	Negative	Positive
Ethyl beta-D-glucuronide	100,000	Negative	Positive
Flunitrazepam	100,000	Negative	Positive
Heroin	100,000	Negative	Positive
Hexobarbital	100,000	Negative	Positive
Ibuprofen	100,000	Negative	Positive
Ketamine	100,000	Negative	Positive
Lamotrigine	100,000	Negative	Positive
Lidocaine	100,000	Negative	Positive
Lorazepam Glucuronide	50,000	Negative	Positive
LSD	100,000	Negative	Positive
Mephobarbital	100,000	Negative	Positive
Methaqualone	100,000	Negative	Positive
Naproxen	100,000	Negative	Positive
Nitrazepam	100,000	Negative	Positive
Normorphine	100,000	Negative	Positive
Norpseudoephedrine	100,000	Negative	Positive
Oxazepam	100,000	Negative	Positive
Pentobarbital	100,000	Negative	Positive
Phenobarbital	100,000	Negative	Positive
Phenylephedrine	100,000	Negative	Positive
Salicylic Acid	100,000	Negative	Positive
Secobarbital	100,000	Negative	Positive
Temazepam	100,000	Negative	Positive
Phenytoin	100,000	Negative	Positive
PMA	100,000	Negative	Positive
Propranolol	100,000	Negative	Positive
(+)-MDA	75,000	Negative	Positive
4-Bromo-2,5-Dimethoxyphenethylamine	75,000	Negative	Positive
Desalkylflurazepam	75,000	Negative	Positive
Dextromethorphan	75,000	Negative	Positive
Diazepam	75,000	Negative	Positive
Flurazepam	75,000	Negative	Positive
Lorazepam	75,000	Negative	Positive
Lormetazepam	75,000	Negative	Positive
Maprotiline	75,000	Negative	Positive
Medazepam	75,000	Negative	Positive

Table 5 - Structurally Unrelated Compounds			
Compound	Concentration Tested (ng/mL)	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
Meprobamate	75,000	Negative	Positive
Methyphenidate	75,000	Negative	Positive
Midazolam	75,000	Negative	Positive
N-Desmethyltapentadol	75,000	Negative	Positive
Oxazepam glucuronide	75,000	Negative	Positive
Phentermine	75,000	Negative	Positive
Phenylpropanolamine	75,000	Negative	Positive
R,R(-)-Pseudoephedrine	75,000	Negative	Positive
Ranitidine	75,000	Negative	Positive
Ritalinic Acid	75,000	Negative	Positive
Sertraline	75,000	Negative	Positive
Theophylline	75,000	Negative	Positive
Thioridazine	75,000	Negative	Positive
Zolpidem Tartrate	75,000	Negative	Positive
Cocaine	40,000	Negative	Positive
MDEA	50,000	Negative	Positive
MDMA	50,000	Negative	Positive
S-(+) Amphetamine	50,000	Negative	Positive
Triazolam	50,000	Negative	Positive
Trifluoromethylphenyl-piperazine	40,000	Negative	Positive

4. Interference – Endogenous compounds were evaluated by spiking the potential interferent into drug free urine containing fentanyl at $\pm 50\%$ of the cutoff. All potential interferents analyzed verified that assay performance is unaffected by internally existing physiological conditions. The results of this study are indicated in Table 6 below:

Table 6 - Endogenous Compounds			
Compound	Concentration Tested (ng/mL)	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
Acetone	1.0 g/dL	Negative	Positive
Ascorbic Acid	0.56 g/dL	Negative	Positive
Bilirubin	2.0 mg/dL	Negative	Positive
Creatinine	0.5 g/dL	Negative	Positive
Ethanol	1.0 g/dL	Negative	Positive
Galactose	10 mg/dL	Negative	Positive
γ -Globulin	0.5 g/dL	Negative	Positive
Glucose	2.0 g/dL	Negative	Positive
Hemoglobin	0.5 g/dL	Negative	Positive
Human Serum Albumin	0.5 g/dL	Negative	Positive
Oxalic Acid	0.1 g/dL	Negative	Positive
Riboflavin	7.5 mg/dL	Negative	Positive
Sodium Azide	1% w/v	Negative	Positive
Sodium Chloride	6.0 g/dL	Negative	Positive
Sodium Fluoride	1% w/v	Negative	Positive
Urea	2.0 g/dL	Negative	Positive

5. Boric Acid Interference – Boric acid at a concentration of 1% w/v was evaluated by spiking the potential interferent into drug free urine containing fentanyl at $\pm 50\%$ of the cutoff. The results of this study are indicated in Table 7 below.

Table 7 – Boric Acid			
Compound	Concentration Tested (ng/mL)	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
Boric Acid	1% w/v	Negative	Positive

6. pH Interference – To evaluate potential interference from the effect of urine pH, device performance was tested using a range of urine pH values (3.0, 4.0, 5.0, 6.0, 7.0, 8.0, 9.0, 10.0 and 11.0). All test samples were prepared in drug free urine containing fentanyl at $\pm 50\%$ of the cutoff. No positive or negative interference was observed at urine pH values ranging from 3.0 to 11.0 for each test mode. The results of this study are indicated in Table 8 below.

Table 8 - Effect of pH		
pH Value	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
3.0	Negative	Positive
4.0	Negative	Positive
5.0	Negative	Positive
6.0	Negative	Positive
7.0	Negative	Positive
8.0	Negative	Positive
9.0	Negative	Positive
10.0	Negative	Positive
11.0	Negative	Positive

7. Specific Gravity Interference - To evaluate potential interference from the specific gravity of urine, device performance was tested using a range of physiologically relevant urine specific gravity values (1.000, 1.002, 1.005, 1.010, 1.015, 1.020, 1.025 and 1.030). All test samples were prepared in drug free urine containing fentanyl at $\pm 50\%$ of the cutoff. No positive or negative interference was observed at urine specific gravity values ranging from 1.000 to 1.030 for each test mode. The results of this study are indicated in Table 9 below.

Table 9 - Effect of Specific Gravity		
Specific Gravity Value	-50% Cutoff (0.5 ng/mL)	+50% Cutoff (1.5 ng/mL)
1.000	Negative	Positive
1.002	Negative	Positive
1.005	Negative	Positive
1.010	Negative	Positive
1.015	Negative	Positive
1.020	Negative	Positive
1.025	Negative	Positive
1.030	Negative	Positive

8. Linearity/ Recovery - Not applicable, this device is intended for qualitative use only

9. Method Comparison – Eighty deidentified, unaltered leftover clinical urine samples obtained from clinical testing laboratories were analyzed for fentanyl at an assay cutoff of 1.0 ng/mL with the Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay compared to results by mass spectrometry (LC-MS/MS). The instruments used were a Beckman Coulter AU 400e and an Agilent 6430 Liquid Chromatography Tandem Mass Spectrometer.

Table 12 –Performance Verified by LC/MS					
SEFRIA™ Fentanyl Urine EIA Result	Fentanyl Concentration by LC/MS				Agreement (%)
	< 0.5 ng/mL	0.5 ~ 0.9 ng/mL	1.0 ~ 1.5 ng/mL	> 1.5 ng/mL	
Positive	0	1	9	31	100% (40/40)
Negative	30	9	0	0	98% (39/40)

The following is a summary table of qualitative discordant results for the 1.0 ng/mL cutoff:

Table 13 - Discordant Result Summary			
Sample ID	In-House ID	1ng Cutoff	LC/MS Confirmation
		Test Device	
Negative 36	21378	Positive	0.9 ng/mL

10. Immunalysis Fentanyl Urine Calibrators

- Traceability – all components of the calibrators have been traced to a commercially available fentanyl solution.
- Value Assignment – Calibrators are manufactured and tested by mass spectrometry. The negative calibrator is a processed, drug free urine matrix. The negative calibrator is compared to a reference negative standard to ensure that it is free of analyte. The non-zero calibrators are prepared by spiking a known concentration of fentanyl in the negative calibrator matrix. If any of the analytes are not within the acceptable range, then the calibrator is adjusted and re-tested. Values are assigned to the calibrators once the mass spectrometry results are within the acceptable ranges.

H. Conclusion

The information provided in this pre-market notification demonstrates that the Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay is substantially equivalent to the legally marketed predicate device for its intended use.