



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

September 15, 2016

MedTech Products, Inc.  
% Ms. Kathryn Coressel  
Emergo Global Consulting, LLC  
816 Congress Avenue, Suite 1400.  
Austin, Texas 78701

Re: K161294

Trade/Device Name: Compound W Wart Removal System Dual Power  
Regulation Number: 21 CFR 878.4350  
Regulation Name: Cryosurgical unit and accessories  
Regulatory Class: Class II  
Product Code: GEH  
Dated: July 29, 2016  
Received: August 2, 2016

Dear Ms. Coressel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803); good manufacturing practice requirements as set forth

in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Christopher J. Ronk -S**

<sup>FOR</sup> Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K161294

Device Name

Compound W® Wart Removal System Dual Power

Indications for Use (Describe)

For the over-the-counter removal of common and plantar warts.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## Section 5 – 510(k) Summary

### Compound W<sup>®</sup> Wart Removal System Dual Power

K # 161294

**1. Submission Sponsor**

MedTech Products Inc.

660 White Plains Road

Tarrytown, NY, 10591

USA

Phone number: (914) 524-6836

Contact: Czarina Ochoa

Title: Regulatory Affairs Manager

**2. Submission Correspondent**

Emergo Global Consulting, LLC

816 Congress Avenue, Suite 1400

Austin, TX 78701

Cell Phone: PHONE

Office Phone: (512) 327-9997

Contact: Kathryn A Coressel, Senior Consultant, RA

Email: [project.management@emergogroup.com](mailto:project.management@emergogroup.com)

**3. Date Prepared**

September 14, 2016

#### 4. Device Identification

Trade/Proprietary Name: Compound W® Wart Removal System Dual Power  
Common/Usual Name: OTC Wart Removal System  
Classification Name: Cryosurgical unit and accessories  
Regulation Number: 878.4350  
Product Code: GEH  
Device Class: Class II  
Classification Panel: General & Plastic Surgery

#### 5. Legally Marketed Predicate Device(s)

Dr Scholl's® Dual Action Freeze Away™ Wart Remover (K052259)

#### 6. Device Description

The Compound W® Wart Removal System Dual Power is a cryosurgical system for the over-the-counter removal of common and plantar warts.

Compound W® Wart Removal System Dual Power consists of the following:

- Pressurized aerosol spray canister with a mixture of dimethyl ether and propane (K032271)
- Reusable actuator/cap that releases the cryogen onto the disposable applicators
- Disposable foam applicators
- Salicylic acid gel wart remover provided in a squeeze tube (in compliance with OTC monograph 21 CFR 358 subpart B)
- Comfort pads
- Instruction Leaflet

#### 7. Indication for Use Statement

Compound W® Wart Removal System Dual Power is intended for the over-the-counter removal of common and plantar warts.

#### 8. Substantial Equivalence Discussion

The following table compares the Compound W® Wart Removal System Dual Power to the predicate device with respect to indications for use, principles of operation, technological characteristics, and materials. The comparison of the devices provides more detailed information regarding the basis for the determination of substantial equivalence. The subject device does not raise any new issues of safety or effectiveness based on the similarities to the predicate device.

**Table 5A – Comparison of Characteristics**

<b>Manufacturer</b>	<b>Schering-Plough HealthCare Products, Inc.</b>	<b>MedTech Products, Inc.</b>	<b>Significant Differences</b>
<b>Trade Name</b>	Dr Scholl's Dual Action Freeze Away Wart Remover	Compound W® Wart Removal System Dual Power	
<b>510(k) Number</b>	K052259	TBD	Not Applicable
<b>Product Code</b>	GEH	SAME	NONE
<b>Regulation Number</b>	878.4350	SAME	NONE
<b>Regulation Name</b>	Cryosurgical unit and accessories	SAME	NONE
<b>Indications for Use</b>	Dr. Scholl's® Dual Action Freeze Away TM Wart Remover <u>is indicated for the over-the-counter removal of common and plantar warts.</u>	Compound W® Wart Removal System Dual Power <u>is indicated for the over-the-counter removal of common and plantar warts.</u>	SAME
<b>Presentation of kit</b>	Kit contains both the cryogen spray and salicylic acid components in the same kit.	SAME	NONE
<b>Sterile</b>	No components are provided as sterile.	SAME	NONE
<b>Single-Use</b>	Only the disposable applicators are single-use.	SAME	NONE
<b>Target Population</b>	Adults and children 4 or older	SAME	NONE
<b>Biocompatibility</b>	Complies with ISO 10993	SAME	NONE
<b>Cryogen spray</b>			
<b>Application</b>	Portable cryosurgical system comprised of an aerosol pressurized canister containing cryogen and an applicator that applies the cryogen to the wart.	SAME	NONE

<b>Manufacturer</b>	<b>Schering-Plough HealthCare Products, Inc.</b>	<b>MedTech Products, Inc.</b>	<b>Significant Differences</b>
<b>Trade Name</b>	Dr Scholl's Dual Action Freeze Away Wart Remover	Compound W® Wart Removal System Dual Power	
<b>Formulation</b>	A mixture of dimethyl ether and propane	SAME	NONE
<b>Freeze time</b>	Dependent upon type and size of wart	Identical treatment times as predicate	NONE
<b>Salicylic Acid</b>			
<b>Application</b>	Post cryo-treatment with a reusable applicator to apply the acid to the wart.	Application is via a squeezable tube	Difference does not raise new concerns for safety or effectiveness.
<b>Formulation</b>	A 17% /salicylic acid <b>solution</b>	A 17% /salicylic acid <b>gel</b>	Difference does not raise new concerns for safety or effectiveness.
<b>Frequency of Use</b>	One drop of the liquid can be applied to the wart 24 hours post cryo-treatment.	Apply 1 drop at a time to cover the wart. Repeat once or twice a day as needed.	Difference does not raise new concerns for safety or effectiveness.