



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Varian Medical Systems, X-ray Products – InfiMed  
% Ms. Catherine Mulcahy  
Regulatory Affairs Manager  
121 Metropolitan Drive  
LIVERPOOL NY 13088

September 6, 2016

Re: K161459

Trade/Device Name: Nexus DR™ Digital X-ray Imaging System (with PaxScan  
4336Wv4)

Regulation Number: 21 CFR 892.1680

Regulation Name: Stationary x-ray system

Regulatory Class: II

Product Code: MQB

Dated: August 5, 2016

Received: August 8, 2016

Dear Ms. Mulcahy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Robert Ochs". The signature is written in a cursive style with a clear, legible font.

Robert Ochs, Ph.D.  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K161459

Device Name

Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4)

Indications for Use (Describe)

The Varian Nexus DR™ Digital X-ray Imaging System is a high resolution digital imaging system intended to replace conventional film techniques, or existing digital systems, in multipurpose or dedicated applications specified below. The Nexus DR™ Digital X-ray Imaging System enables an operator to acquire, display, process, export images to portable media, send images over a network for long term storage and distribute hardcopy images with a laser printer. Image processing algorithms enable the operator to bring out diagnostic details difficult to see using conventional imaging techniques. Images can be stored locally for temporary storage. The major system components include an image receptor, computer, monitor and imaging software.

The Varian Nexus DR™ Digital X-ray Imaging System is intended for use in general radiographic examinations and applications (excluding fluoroscopy, angiography, and mammography).

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**Section 3: 510(k) Summary**

**Date Prepared:** 5/25/2016

**Contact Person:** Catherine Mulcahy  
Regulatory Affairs Manager

**Telephone:** 315-234-6853

**Fax:** 315-234-6801

**Submitter Name:** Varian Medical Systems, X-Ray Products - InfiMed  
121 Metropolitan Drive  
Liverpool, NY 13088

**Device Trade Name:** Nexus DR™ Digital X-ray Imaging System  
(with PaxScan 4336Wv4)

**Common Name:** Digital Radiographic System

**Regulation:** 21 CFR 892.1680

**Classification Name:** Stationary X-ray System

**Primary Product Code:** MQB

**Primary Predicate Device:**  
Stingray DR Digital Radiographic System

**510(k) Number:** K992794  
**Product Code:** MQB

**Reference Predicate Device:**  
Nexus DRF™ Digital X-ray Imaging System

**510(k) Number:** K130318  
**Product Code:** JAA, MQB

**Device Description:**

The Varian Nexus DR™ Digital X-ray Imaging System is a high resolution digital imaging system designed for digital X-ray imaging through the use of an X-ray detector. The Varian Nexus DR™ Digital X-ray Imaging System is designed to support general radiographic (excluding fluoroscopy, angiography, and mammography) procedures through a single common imaging platform.

The modified device consists of an X-ray imaging receptor, Varian PaxScan 4336Wv4, computer, monitor, and the digital imaging software.

**Intended Use:**

The Varian Nexus DR™ Digital X-ray Imaging System is a high resolution digital imaging system intended to replace conventional film techniques, or existing digital systems, in multipurpose or dedicated applications specified below. The Nexus DR™ Digital X-ray Imaging System enables an operator to acquire, display, process, export images to portable media, send images over a network for long term storage and

distribute hardcopy images with a laser printer. Image processing algorithms enable the operator to bring out diagnostic details difficult to see using conventional imaging techniques. Images can be stored locally for temporary storage. The major system components include an image receptor, computer, monitor and imaging software.

The Varian Nexus DR™ Digital X-ray Imaging System is intended for use in general radiographic examinations and applications (excluding fluoroscopy, angiography, and mammography).

**Technological Characteristics Comparison:**

The Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4) supports the same modality as the predicate devices with similar components or imaging concepts, and delivers equivalent or better image quality as the predicate devices. The comparison chart below reveals that functions performed by the predicate devices are performed by the modified device for the DR application. Therefore, the modified device is substantially equivalent to the predicate devices.

However, the modified device, Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4), has the ability to interface with the Varian PaxScan 4336Wv4 wireless detector with vTrigger capability. The Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4) is not integrated with the x-ray generator. With vTrigger, there is no intercepting of the Prep/Expose signals from the generator. The receptor is placed in a “ready” mode by the user and on the Nexus DR™ User Interface. The user operates the generator (technique and exposure switch) when the receptor indicates it is “ready.” The receptor captures the exposure once it arrives and acquires the data.

**Comparison Chart:**

<b>Feature/Item</b>	<b>Nexus DRF™ Digital X-ray Imaging System</b>	<b>Stingray DR Digital Radiographic System</b>	<b>Nexus DR™ Digital X-ray Imaging System</b>
<b>Device Type</b>	<b>Reference Predicate</b>	<b>Primary Predicate</b>	<b>Subject Device</b>
<b>510 (k) Number</b>	<b>K130318</b>	<b>K992794</b>	<b>K161459</b>
<b>Flat Panel Detector</b>	Varian PaxScan 4343R	Trixell Pixium 4600	Varian PaxScan 4336Wv4
<b>Detector Material</b>	a-Si sensor array with CsI or Gd <sub>2</sub> O <sub>2</sub> S:TB scintillator	a-Si sensor array with CsI scintillator	a-Si sensor array with CsI or Gd <sub>2</sub> O <sub>2</sub> S:TB scintillator
<b>Detector Dimensions</b>	17” x 17”	17” x 17”	17” x 14”
<b>Pixel Size</b>	139 x 139 microns	143 x 143 microns	139 x 139 microns
<b>Detector Element Matrix</b>	3072 x 3072	2981 x 3021	3072 x 2560
<b>Dynamic Range</b>	14 bits	14 bits	16 bits

<b>Uniform Density</b>	1.63	N/A	1.52
<b>Spatial Resolution</b>	3.2 lp/mm	3.51 lp/mm	3.2 lp/mm
<b>Sensitivity</b>	128 @ 1.1uGy/frame 395 @ 3.9uGy/frame 1060 @ 10uGy/frame 3143 @ 30uGy/frame <i>(14-bit reference panel)</i>	N/A	540 @ 1.1uGy/frame 1206 @ 3.9uGy/frame 4290 @ 10uGy/frame 12804 @ 30uGy/frame <i>(16-bit subject panel)</i>
<b>Signal to Noise Ratio</b>	67 @ 2.8uGy/frame 96 @ 7uGy/frame 147 @ 17uGy/frame 275 @ 50uGy/frame	N/A	73 @ 2.8uGy/frame 117 @ 7uGy/frame 174 @ 17uGy/frame 285 @ 50uGy/frame
<b>Modulation Transfer Function</b>	0.521 @ 1cycle/mm 0.206 @ 2cycles/mm 0.08 @ 3cycles/mm	N/A	0.551 @ 1cycle/mm 0.234 @ 2cycles/mm 0.099 @ 3cycles/mm
<b>Detective Quantum Efficiency</b>	0.242 @ 1cycle/mm 0.125 @ 2cycles/mm 0.04 @ 3cycles/mm	N/A	0.232 @ 1cycle/mm 0.15 @ 2cycles/mm 0.07 @ 3cycles/mm
<b>External Connectivity</b>	DICOM 3.0 Compatible	DICOM 3.0 Compatible	DICOM 3.0 Compatible
<b>Operator Console</b>	Graphical User Interface	Graphical User Interface	Graphical User Interface
<b>Image Processor</b>	Intel Based PC, nvidia GPU	Pentium Based PC	Intel CPU Based PC
<b>Image Storage</b>	Hard Drive	Hard Drive	Hard Drive
<b>Operating System</b>	Windows 7	Windows NT	Windows 10
<b>Total Image Processing Time</b>	10 seconds per image	30 seconds per image	10 seconds per image
<b>Power Requirements</b>	110/120V, 230/240V, 50/60 Hz	110/120V, 230/240V, 50/60 Hz	110/120V, 230/240V, 50/60 Hz

**Non-clinical Tests Discussion:**

Non-clinical Data submitted is consistent with FDA guidance document “Guidance for Industry and/or for FDA Reviewers/Staff and/or Compliance: Guidance for the Submission of 510(k)’s for Solid State X-ray Imaging Devices” available at the website <http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm073781.pdf>.

Validation was completed in accordance with the Validation Protocols included with this submission. Protocols were designed, executed and documented according to the Design Validation process with predetermined test methods and corresponding acceptance criteria. In conclusion, all release criteria have been met and the Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4) is as safe and effective as the predicate devices and does not raise different questions of safety and effectiveness.

**Clinical Tests Discussion:**

Clinical Data submitted is consistent with FDA guidance document “Guidance for Industry and/or for FDA Reviewers/Staff and/or Compliance: Guidance for the Submission of 510(k)’s for Solid State X-ray Imaging Devices” available at the website <http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm073781.pdf>.

**Conclusion:**

Based upon the results of Verification and Validation testing, the Nexus DR™ Digital X-ray Imaging System (with PaxScan 4336Wv4) has no new indications for use, has no significant technological differences, and is as safe and effective as, does not raise different questions of safety and effectiveness and is therefore substantially equivalent to the above listed current legally marketed predicate devices.